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**VISCOSYSTEM
AVS 500**

Operating Instructions

**VISCOSYSTEM
AVS 500**

Mode d'emploi

**VISCOSYSTEM
AVS 500**

SCHOTT

Gebrauchsanleitung	1	44
Operating Instructions	45	88
Mode d'emploi	89	132

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Important notes: Read these operating instructions carefully before using the unit. For safety reasons the unit may only be used for the purposes described in the operating instructions. This product is subject to continuous modification to adapt it to the state of the art. For this reason, despite great care, it may be possible that these operating instructions do not describe the properties of the unit in full. If you are in doubt please contact the Technical Application Department at our company.

Please also observe the operating instructions for the units to be connected.

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Operating Instructions

VISCOSYSTEM AVS 500

CONTENTS	PAGE
1 The ViscoSystem AVS 500	47
1.1 Functional principles of the VISCOSYSTEM AVS 500	47
1.2 Capillary viscometry	47
1.3 Measuring principles	48
Opto-electronic light barriers	48
Viscometers with thermistor sensors (TC-Viscometers)	48
1.4 Warnings and safety instructions	48
Declaration of Conformity	50
Technical specifications and data of the VISCOSYSTEM AVS 500	51
2 Setting-up, commissioning and transportation instructions	54
2.1 Unpacking	54
2.2 Mounting the sample transfer unit	54
2.3 Inserting the measuring amplifier and pneumatic modules	54
2.4 Hose and plug connections	54
2.5 Safety set	55
2.6 Used thermostats	55
3 Functional units of the VISCOSYSTEM AVS 500	56
3.1 The LED displays	56
3.2 Logic status displays	56
3.3 The keyboard	57
3.4 Sample transfer unit	60
3.5 Sample turntable	61
3.6 Measuring amplifier cards	61
3.7 Pneumatic modules	64
3.8 Printer, inserting the paper	65
3.9 The data memory	65
3.10 Data transmission	66
4 Measuring with the VISCOSYSTEM AVS 500	68
4.1 Overview, brief instructions	68
4.2 Determining the absolute kinematic viscosity	68
Changing the stored data	69
Measurements for the absolute kinematic viscosity	69
4.3 Determining the relative viscosity	71
Changing the stored data	71
Measurements for the relative viscosity	72
4.4 Determining the flow time	73
Changing the stored data	73
Measuring the flow through time	74
4.5 Determining the inherent viscosity	75
Changing the stored data	75
Measurements for the inherent viscosity	77
4.6 Brief instructions for dilution viscometry	79
4.7 Determining the limiting viscosity number (dilution viscometry)	79
Changing the stored data	80
Measurements for the dilution viscometry	81
4.8 Terminating the measurement	82

5	Routine checks and trouble shooting	82
5.1	Faults which are not reported by the ViscoSYSTEM AVS 500	82
5.2	Error messages reported by the ViscoSYSTEM AVS 500 by print-out or warning lamp	83
5.3	Replacing functional units	84
	Replacing a defective sample transfer unit	84
	Replacing the dosing unit	84
	Filter change and needle replacement	84
6	Maintenance and care of the ViscoSYSTEM AVS 500	85
7	Index	86

1 The ViscoSystem AVS 500

1.1 Functional principles of the VISCOSYSTEM AVS 500

The VISCOSYSTEM AVS 500 can carry out measurements of the absolute kinematic and relative viscosity independently ranging from calculation to documentation. 12 sample bottles with 100 ml or 50 ml capacity each, or 36 sample glasses with 20 ml capacity each, can be inserted. The position of the sample in the sample turntable is recognized automatically.

The sample liquid is drawn into a calibrated glass cylinder through a moving stainless steel or glass needle and then delivered from this cylinder into the viscometer.

The sequence of the measurements is programmed via the keyboard. Numerical and function LED indicators show the measurement results and the entered parameters. The results are printed on plain paper.

The pneumatic modules and the measuring amplifier cards are inserted into the rear panel of the VISCOSYSTEM AVS 500.

Two measuring amplifier cards are available, one for opto-electronic sensing and the other for TC-Viscometers. Thus black and absolutely opaque liquids can be measured as well as clear transparent ones.

Flow times up to 9999.99 s can be measured. The absolute, relative and inherent viscosity can be calculated from the time measurement, with optional inclusion of the Hagenbach correction.

In addition to the date and time of day, the position and sample number assignments can be entered and the thermostatic bath temperature can be read from a calibrated thermometer.

If you wish a determination of the limit viscosity number, a piston burette from SCHOTT-GERÄTE (preferably type no. T 110/50) has to be connected additionally.

The VISCOSYSTEM AVS 500 can also be operated with two viscometers simultaneously. This nearly doubles the sample throughput. An individual pneumatic module and an individual measuring amplifier card are required for each viscometer.

The VISCOSYSTEM AVS 500 is equipped with a data output conforming to RS-232-C for data transmission. Connections are provided for a Pt 1000 temperature sensor and for safety devices for dependable safeguarding of the measuring action sequence.

A thermostat must be used in order to achieve highest possible accuracy. The control deviation of the thermostat must be ± 0.02 K or better. Clear view (transparent) thermostats from SCHOTT-GERÄTE (e. g. type no. CT 1450/2) meet this requirement.

1.2 Capillary viscometry

Capillary viscometry is the most accurate of all methods for determining the viscosity of liquids with Newtonian flow behavior. The actual reading is a time measurement. The time taken by a defined volume of liquid to flow through a capillary tube with accurately known width and length, is measured. Conventionally, such measurements used to be made with the eye of the observer and a stopwatch. In the VISCOSYSTEM AVS 500 - as in all viscosity measuring units from SCHOTT-GERÄTE - the liquid meniscus is detected opto-electronically or by thermal conductivity measurement.

The VISCOSYSTEM AVS 500 permits viscosity measurements from 0.3 to approx. 1200 mm²/s at a sample temperature of 20...25 °C. The upper limit is determined by the combination of the hose equipment (inner diameter) / sample transfer logic and the switch positions of the measuring amplifier card (more detailed information, see chapter 3.4 „Sample transfer unit“. Different hose assemblies are available to you, coordinated to the range of viscosity.

1.3 Measuring principles

Opto-electronic light-barriers

A measuring stand AVS/S (anodized aluminium) or AVS/SK (PVDF) is required for opto-electronic sensing of the meniscus position. These precision devices ensure constant high precision of measurements made according to the principle of capillary viscometry, even when the measuring stand and the viscometer are interchanged.

The light produced in the upper part of the LED's measuring stand (in the close infrared range) is lead to the measuring plane with the help of a fiber-optical light guide cable. The beam of light passes through the viscometer and reaches at the other side a fiber-optical light guide cable, which leads the light to a photo-diode at the upper side. A precise measuring signal is produced.

Viscometers with thermistor sensors (TC-Viscometers)

TC-Viscometers use glass-sheathed thermistors sealed-in at the measuring levels. When the meniscus passes through the measuring level, the thermal equilibrium at the thermistor is changed due to the different thermal conductivities of air and sample liquid.

The thermistors of the TC-Viscometers are hermetically sealed in glass sheaths, so that the interior of these viscometers is chemically resistant with respect to all kinds of substances with the exception of concentrated alkalis, solutions containing fluoride or highly concentrated hot phosphate solutions.

TC-Viscometers are protected by German Industrial Design No. 85 04 764 and USA Patent No. 4 685 328.

1.4 Warnings and safety instructions

For safety and functional reasons, the VISCOSYSTEM AVS 500 must basically be opened by authorized persons only; e.g. work at the electric appliance must only be carried out by trained specialists. With unauthorized operation of the VISCOSYSTEM AVS 500 as well as careless or deliberate damaging, the warranty is no longer valid. For safety reasons the VISCOSYSTEM AVS 500 must only be used for purposes described in these operations instructions.

Please also read and observe the operating instructions for the viscometers and thermostats and for all other connected equipment units.

The VISCOSYSTEM AVS 500 conforms to Protection Class I. It has been built and tested according to DIN 57 411, Part 1/VDE 0411, Part 1, Protection precautions for electronic measuring equipment and has left the factory in a safe and properly functioning state. In order to preserve this state, the user must take note of and observe the instructions and warnings which are contained in this operating instructions manual.

Before switching-on, make sure that the mains voltage setting of the VISCOSYSTEM AVS 500 corresponds to the nominal value of the available mains voltage. The mains plug may be inserted only into a power outlet with grounding contact. The ground connection must not be interrupted by an extension cable without ground conductor. Any interruption of the grounding conductor inside or outside the VISCOSYSTEM AVS 500, or disconnection of the grounding conductor, can bring the VISCOSYSTEM AVS 500 into a dangerous state. Deliberate interruption of the grounding connection is forbidden.

Make sure that only fuses of the specified type and nominal current rating are used for replacements. The use of repaired fuses or short-circuiting of the fuseholder is forbidden.

If there is reason to believe that non-dangerous operation is not possible, the VISCOSYSTEM AVS 500 should be stopped and secured against unintentional operation.

Please switch off the VISCOSYSTEM AVS 500, disconnect the mains cable and call the SCHOTT-GERÄTE customer service. It must be assumed that safe operation is no longer possible:

- when the VISCOSYSTEM AVS 500 shows visible damage
- after severe stress during transportation.

The built-in safety devices, in particular the needle protection, should never be disabled.

In the case of grossly incorrect manual control, sample liquid may come into the pneumatic module. In order to prevent such damage, it is advisable to use the safety sets provided with the delivery. The fill level sensor for viscometers (Fig. 24 and 25, Pos. 11) may be used only when the ignition temperature of the sample liquid lies above 250 °C.

The safety devices (Fig. 24, 25 and 26) are connected to the VISCOSYSTEM AVS 500 via the distributor type no. VZ 6247 for safety sensors (Fig. 24, 25 and 26, Pos. 7) according to the connecting diagram (Fig. 2, 24, 25 and 26) and are to be considered as typical examples for all safety devices.

A safety set type no. VZ 6247 is required for each measuring amplifier. The weight parameter for the liquid level sensor type no. VZ 6244 for the waste bottle type no. VZ 5379 is set to 2.5 kg. This value is valid for liquids with a density of 0.8 kg/l and empty weight of approx. 0.9 kg for the waste bottle complete with bottle adapter NS 29, type no. VZ 6240.

New settings are required for other weights: Turning the milled screw to the left at the front on the liquid level sensor corresponds to a heavier weight, see symbol on the milled screw. The normal practical procedure is to fill 2 l of liquid into the bottle and then make a new setting of the switching point.

The capacitive safety sensor type no. VZ 6245 for the suction line (Fig. 24, 25 and 26, Pos. 10) is set for 5 ml of water dielectric constant = 80.

The statutory regulations for the prevention of accidents (e. g. safety goggles, safety clothing, respirators, etc.) must be observed as protection against dangers arising from the solvents and/or samples used.

The VISCOSYSTEM AVS 500 may not be operated or stored in damp conditions.

Declaration by the manufacturer

We certify herewith that the VISCOSYSTEM AVS 500 is equipped with radio frequency interference suppression in accordance with the stipulations of the Official Bulletin Amtsbl Vfg 1046/1984. Commercial distribution of this unit has been reported to the German Federal Post Office which has been given the right to check the serial production for compliance with the regulations.

All specifications in this instruction manual are guidance values which are valid at the time of printing. However, SCHOTT-GERÄTE may make changes at any time for technical or commercial reasons or in the necessity to comply with the statutory stipulations of various countries.

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Technical dates Viscosity Automatic Sampler VISCOSYSTEM AVS 500 (State June 28, 1995)

Sampler system:

Sample bottles:	50 ml or 100 ml screw- / standard taper bottles (12 pieces); 20 ml circular glasses (36 pieces)
Sample turntable:	type no. VZ 6204 for 12 pieces 100 ml screw- and standard taper bottles type no. VZ 6205 for 12 pieces 50 ml screw bottles type no. VZ 6206 for 12 pieces 50 ml standard taper bottles brown type no. VZ 6207 for 36 pieces 20 ml circular glasses

Measuring:

Principle:	meniscus detection by optoelectronic or thermal conductivity, see table 1: measuring amplifier cards and pneumatic modules
Measuring parameters:	flow through time in seconds [s] temperature in degree Celsius [°C]

Calculated parameters:

absolute kinematic viscosity, relative viscosity, intrinsic viscosity,
inherent viscosity, reduced viscosity (limiting viscosity number),
Hagenbach correction (HC), average values, standard deviation

Selected parameters:

Keyboard: <...>	decimal keyboard 0 ... 9; Clear, Enter, Function, constants, t_0 runtimes, temperature (nominal, tolerance, calibration), number of measurements, sample number, Hagenbach correction, turntable position, temperature equilibration, rinsing, statistical evaluation, date / time of day, <ml> input solvent quantity (dilution viscometry), paper feed, turntable rotation, Start, Stop / Reset
Number of measurements:	1 ... 99
Temperature equilibration time:	0 ... 99 min in steps of 1 min
Number of viscometer rinsing:	0 ... 9 with next sample (check minimum quantity of sample)
Sample identification number:	8 digits with turntable position
Memory:	measuring parameters are stored until modification (permanent memory)

Measuring ranges:

Viscosity:	0,3 ... 600 mm ² /s at a sample temperature of approx. 20 ... 25 °C
Time:	up to 9999,99 s, resolution 0,01 s
Suction pressure:	full automatic control
Viscometer:	Ubbelohde (DIN, ASTM, Micro, TC), Cannon-Fenske-Routine, Micro-Ostwald and Dilution Viscometer (see table 2): viscometer types, measuring and viscometer stands

Measuring accuracy:

Accuracy of the time measurement:	$\pm 0,01 \text{ s} \pm 1 \text{ digit}$, precision may not exceed 0.01 % The measuring uncertainty, while measuring the absolute kinematic viscosity, also depends on the uncertainty of the numerical value for the viscomete constant and on the measuring conditions, especially the measurement temperature.
--------------------------------------	--

Displays: <<...>>

Readout of measurement results:	
Flow through time:	6 digits LED display with floating decimal point, character height 14.2 mm
Temperature:	4 digits, resolution $\pm 1 \text{ digit}$,
Date / clock:	date; h, min, s
Dilution:	ml for the solvent quantity (dilution viscometry)

Viscometer constants or

t_0 -runtimes:	6 digits with floating decimal point, (input possible up to 8 digits)
LED functional displays:	<<s>>, <<min>>, <<d>>, <<h>>, <<ml>>, <<set>>, <<tol>>, <<°C>>, <<abs>>, <<rel>>, <<t>>, <<c → 0>>, <<RUN>>, <<meniscus passages>>, <<!>>, (error) keys: <°C>, <°C cal>, <n>, <N°>, <type HC>, <pos>, <temp>, <rinsing>, <Δ %>, <clock / date>, <ml>, <K/ t_0 >, <start>

Evaluations / results:

Correction: Hagenbach correction (HC) for Ubbelohde-, Cannon-Fenske-Routine, Micro-Ubbelohde- and Micro-Ostwald-Viscometer

Statistical evaluation: standard deviation, search for runaway values

Printer: 24 characters per line in 5 x 7 dot matrix format, character height = 2.5 mm

Paper: ordinary paper, width = 58 mm

Ambient temperature conditions:

Ambient temperature: 0 ... + 50 °C

Relative humidity of air: max. 85 % rel.

Electromagnetic emission: electrical shielding, radio frequency interference suppression in accordance with the stipulation of the Official Bulletin Amtsblattverfügung (Amtsbl Vfg) 1046/1984, according to VDE 0781 B

Device safety: conforms to DIN 57 411, Part 1/VDE 0411, Part 1, Protection Class I

Case:

Materials: fully shielded plastic/stainless steel/aluminium cabinet with chemically resistant two-component coating

Dimensions: 700 x 560 x 460 mm (W x H x D)

Weight: approx. 33 kg for the VISCOSYSTEM AVS 500 with measuring amplifier card and pneumatic module for 1 viscometer

Connections:

Pneumatic connections: screw connections for viscometer

Electrical connections: concentric plug-and-socket with bayonet catch for measuring stands and TC-Viscometer

Viscometer: 1 or 2 viscometer each connectable with 1 measuring amplifier card and 1 pneumatic module: types see table 2; for dilution viscometry: 1 Dilution Viscometer on measuring point 1

Temperature: Pt 1000 sensor: Resistance Thermometer Pt 1000; offset calibrated

Thermostat: signal 10 mV/°C; Reference: Type CT 1450/2 from SCHOTT-GERÄTE

Connection for measuring sensor Pt 1000 has priority over connection for the thermostat

Data interfaces:

- 1) serial interface RS-232-C:
plug connection: female 25 pin connector subminiature D
- 2) interface IEEE 488:
plug connection: female 25 pin connector series 57 Micro Ribbon®/CHAMP
- 3) interface for Piston Burettes VISCODOSER® AVS 20, TITRONIC® T 110, TITRONIC® T 200:
plug connection: male 37 pin connector subminiature D

Safety:

- 1) connection of the safety set (box) VZ 6247
 - a) with Level sensor VZ 6246 for Ubbelohde Viscometer permitted up to a viscosity of $\leq 200 \text{ mm}^2/\text{s}$ and an inflammation point of the measuring solution above 250 °C)
 - b) with capacitive safety sensor VZ 6245 for suction line, preset by manufacturer to 5 ml of water at a relative permittivity = 80
 - c) with overflow protection VZ 6244 preset by manufacturer to 2.5 kg for waste bottle VZ 5379
- 2) additional potential free controlling working contact for any desired signal device (max. 250 V~ / 10 A~)

Mains supply: Euro panel plug DIN 49 457 B with safety fuse, mains filter

Data transfer:

Interface: bidirectional serial interface according to EIA RS-232-C,

Data format: 7 bit ASCII code, 2 stop bits, 4800 Baud, no parity

Power supply:

Voltage supply: The VISCOSYSTEM AVS 500 conforms to Protection Class I
230 V~ or 115 V~, 50 ... 60 Hz

Mains fuse: instrument fuse (5 x 20 mm): for 220 V~, 0,8 A delayed
for 110 V~, 1,6 A delayed

Power consumption: 150 VA

Table 1: Measuring amplifier cards and pneumatic modules

- Measuring amplifier card VZ 6211 for meniscus detection by optoelectronic
- Measuring amplifier card VZ 6212 for meniscus detection by thermal conductivity, TC-Viscometer

Pneumatic unit type	VZ 6203 Membrane pump	VZ 6208 Jet pump	VZ 6213 Membrane pump	VZ 6214 Membrane pump
Material	Viton	PTFE	PTFE	EPDM
Liquids	aromatic and chlorinated	aggressive	neutral not aggressive	Ketones
Examples for liquids:				
Formic acid	—	+	○	—
Sulfuric acid	—	+	○	—
Dichloro acetic acid	—	+	—	—
Trifluoro acetic acid	—	+	—	—
Phenol/Dichloro benzene	+	+	+	○
m-Cresol	+	+	+	○
o-Chlorophenol	+	+	+	○
Hexa-Fuoro-Isopropanol	—	+	○	—
Cyclo-Hexanon	—	+	+	+
Aceton	—	+	+	+
Methyl-Ethyl-Keton	—	+	+	+
Dichloromethylene	+	+	○	○
Chloroform	+	+	○	○
Ethanol	+	+	+	+
Methanol	○	+	+	+
Paraffinic oils	+	+	+	○

+ = applicable; — = not applicable; ○ = conditionally applicable, in doubt it has to be tested.

Table 2: Viscometer types, measuring and viscometer stands

Viscometer		Frame	Measuring stand		Hose assembly
Type	Type no.	Type no.	Type no.		Type no.
Ubbelohde (DIN)	532.. 530..	053 92	AVS/S	AVS/SK	VZ 6270*
Ubbelohde (ASTM)	526..	053 92	AVS/S	AVS/SK	VZ 6270*
Micro-Ubbelohde	537.. 538..	053 92	AVS/S	AVS/SK	VZ 6272
Cannon-Fenske-Routine	520..	----	AVS/S-CF	AVS/SK-CF	VZ 6270*
Dilution Viscometer	531..	----	AVS/S-V	AVS/SK-V	VZ 6270 VZ 6284
Micro-Ostwald	517..	053 97	AVS/S	AVS/SK	VZ 6272
TC-Ubbelohde	547..	053 93	Bracket VZ 5931		VZ 6270* VZ 6226

* above 200 mm²/s to 600 mm²/s: type no. VZ 6271, but without liquid level sensor.

2 Setting-up, commissioning and transportation instructions

2.1 Unpacking

The VISCOSYSTEM AVS 500 has been carefully checked and packed at the factory. When unpacking, please take care to remove all parts out of the packing material. Please check that the nominal mains voltage at the intended operating location corresponds to the value marked on the type specification plate.

The individual parts are packed separately for safe transportation. The mounting instructions are contained in the following sections.

Please check that the mains voltage actually available agrees with the nominal value shown on the type specification plate (see rear of the VISCOSYSTEM AVS 500 and Fig. 4, Pos. 2).

2.2 Mounting the sample transfer unit

Remove the transport securing devices on the sample transfer unit and on the pneumatic module.

For mounting, first remove the front panel of the sample transfer unit. For this purpose, release the 4 screws at the front (Fig. 5, Pos. 1). Then push the assembly into the VISCOSYSTEM AVS 500 on the runners. The plug connector thereby inserts itself into the socket in the basic unit, with automatic mechanical self-adjustment. Then tighten the three screws shown in Fig. 6, Pos. 1 (see also Section 5.3.1).

The sample transfer unit must be taken out again for transportation of the VISCOSYSTEM AVS 500.

2.3 Inserting the measuring amplifier cards and pneumatic modules

As viewed from the rear, the order of the cards and modules (Fig. 4 and 5) from right to left is:

Measuring amplifier card no. 1

Measuring amplifier card no. 2

Pneumatic module no. 1

Pneumatic module no. 2

Remove the blank panels before mounting the measuring amplifier cards and pneumatic modules.

Fix the cards and modules with the associated screws. The VISCOSYSTEM AVS 500 is then ready for operation (Fig. 14 to 19).

Proper electrical contacts have been established when the rear panels of the cards and modules are flush.

2.4 Hose and plug connections

The viscometer type and hose assembly are selected according to table 2. The viscometer is mounted in the frame (Fig. 20 to 23).

The frame containing the viscometer is inserted into the measuring stand for optoelectronic sensing of the meniscus. The orientation must be such that the breakout in the bottom plate faces to the front and latches into the nose provided for it. Apply light pressure against the frame to make the viscometer lock into the holding spring on the measuring stand. Before making measurements, always check that the viscometer is correctly seated in the measuring stand (Fig. 20 and 21).

Establish the electrical and pneumatic connections between the measuring stand and the viscometer on the one hand and the VISCOSYSTEM AVS 500 on the other hand (Fig. 1, 2 and 24 to 26).

All hose connections must be checked carefully to make sure that they do not leak. A safety vessel each is interposed in the suction and discharge line (Fig. 2, 24 to 26, Pos. 10). When checking for leaks, it is advisable to first make a test run with water or some other non-aggressive solvent.

A measuring stand is not required for TC-Viscometers. A stainless steel frame VZ 5931, conspiring with a stainless steel stand 05393 which has been specially adapted for this type of viscometer provides the means for inserting the latter into the thermostat. The electrical connector is attached directly to the viscometer and is connected with the help of a connecting cable VZ 6226 for TC-Viscometer with the measuring amplifier card. A separate measuring amplifier card is required for TC-Viscometers. In parallel with the electrical connection, the TC-Viscometer is also connected pneumatically to the pneumatic module and to the sample transfer unit. The pneumatic modules are the same for viscometers with opto-electronic sensing and for TC-Viscometers (Fig. 15 and 17).

Connectors are provided on the rear of the VISCOSYSTEM AVS 500 (Fig. 13) for the data cables and measuring sensors:

- for mV output of the thermostat head (only from SCHOTT-GERÄTE, type no. CT 1450/2)
- for the Resistance Thermometer Pt 1000 for measuring the bath temperature
(The Pt 1000 connection has priority over the connection for the thermostat head)
- for the serial interface RS-232-C
- for the Piston Burettes, preferentially ViscoDOSER® AVS 20, TITRONIC® T 110, type no. T 110/50

2.5 Safety set

The safety set must be connected in order to ensure early detection of possible functional disturbance of the VISCOSYSTEM AVS 500. The safety set indicates the following events:

Full waste bottle, see Fig. 24, 25 and 26, Pos. 4.

The viscometer is too full, see Fig. 24 and 25, Pos. 1; **only Ubbelohde-Viscometers up to 200 mm²/s.**

Fill level sensor, see Fig. 24 and 25, Pos. 11.

Liquid in the safety vessel suction line, see Fig. 24, 25 and 26, Pos. 3.

Optionally: Safety vessel waste line (logically combined with safety vessel suction line), see Fig. 24, 25 and 26, Pos.5).

2.6 Used thermostats

In order to obtain the high accuracy of measurement results possible with the VISCOSYSTEM AVS 500, it is necessary to pay special attention to temperature equilibration and constancy. A thermostat which ensures a bath temperature constancy of ± 0.02 °C is required, e. g. a transparent thermostat type no. CT 1450/2 from SCHOTT-GERÄTE (Fig. 1 and 2).

As a rule of thumb, it may be assumed that the temperature deviation in °C multiplied by 10 gives the percentage deviation of the viscosity measurement results. For example, a temperature deviation of 0.05 degrees (x 10) would correspond to a viscosity measurement error of 0.5 %.

3 Functional units of the VISCOSYSTEM AVS 500

3.1 The LED displays

The displays 1 and 2 show the measured flow times (1.0...9999.99 s) after the liquid has run through. The measurement result remains in the display after print-out until the next measurement has been terminated (Fig. 1 and 3).

After depressing a function key, the last entered value or the current value (e. g. the time of day) is displayed, including the physical dimension (small LEDs).

The entered values are displayed as they are being entered via the decimal keyboard. To clear the displays, depress the key <C> (Clear).

The top right numerical LED displays shows exclusively the temperature (nominal temperature, temperature tolerance, calibration temperature). All temperature values are displayed in °C.

3.2 Logic status displays

<<abs>>: Calculation of the absolute kinematic viscosity [mm²/s]. The values which have been entered for the constants are used for the evaluation.

<<rel>> Calculation of the relative viscosity. The entered t_0 runtimes are used for the evaluation. For relative viscosity measurements it is very advisable to specify the approximate values for the constants in the mode <<abs>> because the dosing rate and the suck off time are selected automatically by the VISCOSYSTEM AVS 500 according to the values of the constants.

<<t>>: For pure time measurements. Only the average of the measured flow times is determined. It is very advisable to specify the approximate constants in the <<abs>> mode because the dosing rate and the suck off time are selected automatically by the VISCOSYSTEM AVS 500 according to the values of the constants.

<<c → 0>>: Measuring sequence "dilution viscometry".

LED indicators in the viscometer symbols: The corresponding LED symbols for both top and bottom measuring levels illuminate on and turn off as the liquid passes on its upward and downward travel. As soon as the upper or lower reference mark is passed through during the flowing-out of the liquid, the LED symbol, which belongs to it, expires. While sucking up the liquid, it applies to the opposite order.

LED indicators in the keys: These illuminate when the function or parameter is called-up by depressing the key. Each display (and the corresponding value) can be suppressed by depressing another key. But the previous value remains stored in memory, provided that it has been acknowledged with <Enter> at some earlier time.

Dimension indicators for the measured values and parameters: These are lighted according to their meanings: <<s>> for measuring the flow time, <<min>> for the temperature equilibration time, <<d>> and <<h>> for the date and time of day, <<ml>> for the solvent quantity in dilution viscometry, <<°C>> for the temperature.

<<Exclamation mark>>: This LED flashes in the case of a fault and blinks when the needle is put down in the sampling plate. An additional warning device, e. g. a warning buzzer, connected at the rear of the VISCOSYSTEM AVS 500 can be switched on by alarm contacts inside the VISCOSYSTEM AVS 500.

3.3 The keyboard

The corresponding numerical value appears in one of the displays when a key is depressed. When another key is depressed, the new value is displayed and removes (overwrites) the former value in the display. In each case the green LED in the last depressed key lights as orientation aid for the operator. All keys except <Stop/Reset> are locked during program running.

The functions of the individual keys are:

The numerical keyboard is used for entering:

- viscometer constants
- t_0 runtimes
- nominal temperature, calibration temperature and temperature tolerance
- number of measurements
- sample number
- viscometer type, Hagenbach correction and evaluation procedure
- turntable position of the first and last sample on the turntable
- temperature equilibration time
- manner and number of rinsings
- tolerance specification and statistical evaluation
- date and time of day
- solvent quantity for dilution viscometry

<Start>: Starts the measuring sequence, the green LED and the display <<RUN>> light. The program runs until the preprogrammed end or until the <Stop/Reset> key is depressed.

<Stop/Reset>: Stops a running measuring program immediately. All data so far measured for this sample are erased from memory.

When the <Stop/Reset> key is pressed again, the VISCO SYSTEM AVS 500 goes into ready status (the green LED in the <Start> key and <<RUN>> are not lit). Only print-out of the AVS memory takes place.

Symbol <arrow circle>: Depressing this key makes the sample turntable turn to the next sample. This key can be depressed up to 8 times to store and successively execute up to this number of sample position increments.

<P> (Printer): The printer paper is advanced by one line. Holding this key depressed continuously advances the printer paper.

<C> (Clear): Whenever it is desired to change the value of a measuring parameter, the value standing in the display is deleted on depressing the <C> key. The value in memory is changed only when the newly entered value is confirmed with <Enter>.

<Enter>: After entering a numerical value (now seen in the display window), this value is stored into memory after depressing the <Enter> key. Numerical values which have not been confirmed with <Enter> are ignored in the program and in the measuring sequence.

<F (Function): Selection for calculation of the viscosity: absolute viscosity <<abs>>, relative viscosity <<rel>>, measurement of the runout time <<t>> and limiting viscosity number <<c → 0>>.

<°C>: The nominal value temperature or the permitted temperature tolerance is displayed in °C (the smallest permitted tolerance is 0.1 °C). Pressing the key <°C> repeatedly toggles between display of the nominal value temperature <<Set>> and display of the temperature tolerance <<Tol>>. Examples of entries: <25.00> (°C) and <0.50> <K>.

When calibrating the temperature for the first time it is appropriate to enter a tolerance of 99.99 K. Overshoot of the temperature tolerance, which is indicated by the print-out "Temp. Tol." overshoot and flashing of the LED <<Exclamation mark>>, can be reset by new calibration followed by the keystroke sequence <Clear> and <Enter>. The temperature measurement is not made continuously but depends on the operating states of the VISCO SYSTEM AVS 500.

<°C cal>: For entering the exact temperature determined with a calibrated thermometer. This exact temperature value is taken over by the VISCO SYSTEM AVS 500 as calibration temperature. The calibration temperature should not differ by more than 3 °C from the working temperature. Example of an entry: <24.85> = measured exact temperature [°C] in the thermostat bath for a working temperature of 25.00 °C.

<n>: The number of individual measurements. This number is the same for all samples on the sample turntable. Up to 99 measurements for each sample can be specified. The default value for "0" is 1.

<No>: For entering the sample numbers. The entry can contain up to 8 digits. The display at the same time shows the position number of the sample turntable. The input can be made after the erasure of the display. Delete the display while pressing down the <C> (clear) key. Example of an entry: <5 12345678>, <Enter>.

The input mode for continuous position numbers, independent from each plate position, can be reached as follows: as described before, the position and the sample number have to be put in:

Example: : 0187654321, <Enter>, decimal point, display <02 00000000>, <C> (clear),

display: <<02.....>>. Further sample numbers according to the continuous position numbers can be put in with the input of the numeral run. <Enter>, <C> (clear).

If a decimal point is pressed down between "position" and "sample number", each position number belonging to the sample plate will be displayed together with the numeral run 00000000 after <C> (clear) and the desired sample number can be put in. **Example:** The plate is in position 5: Press down decimal point:

Display: <<05 00000000>>. Input: <12345678>. Display: <<05 12345678>>. Confirmation of the input with <Enter>. After <C> (clear) is pressed down, <<06.>> will be on display and the next sample number can be put in. On depressing <C> (Clear), the display shows <<6.>> and the next sample number then can be entered.

<Type HC>: For calculating the Hagenbach correction [$HC = E/(K \cdot t)^2$], e. g. according to DIN 51 562, Part 1, the viscometer type must be entered, e. g. Ubbelohde-Viscometer = No. 1.

If a Hagenbach correction is to be made in pure time measurement mode, the actual or nominal constant of the viscometer must be entered in <<abs>> mode.

The viscometer type and the evaluation mode are separated by a decimal point in this entry. Example of an entry: <1.1>.

Entry before the decimal point: 1 = Ubbelohde-Viscometer according to DIN (also TC-Viscometer)

2 = Micro-Ubbelohde-Viscometer

3 = Cannon-Fenske-Routine

4 = Dilution Viscometer

5 = Micro-Ostwald-Viscometer

6 = Ubbelohde-Viscometer ASTM

Entry after the decimal point: 0 = calculation without HC

1 = with HC

2 to 9 are for calculation of the inherent viscosity.

The concentration can be entered in 4 steps, whereby

2 = 0.25 g/100 ml

3 = 0.5 g/100 ml

4 = 0.75 g/100 ml

5 = 1.00 g/100 ml

The calculation is made **without** Hagenbach correction in the cases 2 to 5.

The numbers 6 to 9 have the same meanings, but **with** Hagenbach correction.

<Pos>: Position of the first sample to be measured and position of the last sample to be measured. If only the start position is entered followed by a decimal point and <00>, then the measurement sequence runs continuously (example of an entry: <01.00>). **Attention:** An entry such as <1.3> is not permitted; the correct form to be entered is <01.03>.

The position of the last sample to be measured can also be marked with the magnetic rod "last sample". The magnetic rod has priority.

Symbol <Temperature Equilibration>: For entry and display of the temperature equilibration time in minutes. This entry is the same for all samples on the turntable. Example of an entry: <05> (min). The Program accepts the entry <0>. Only integer minute specifications are allowed.

Symbol <Rinsing>: For entry and display of the number of rinsing cycles with the next sample or pre-selection of rinsing with fresh solvent.

Example of an entry: <1.3>.

Entry before the decimal point: 1 = rinse with next sample

Entry after the decimal point: Number of rinsing cycles (max. 9);
the entry <0> is possible and will be obeyed.

< Δ %>: For rejecting measured values which differ from the calculated average by more than the entered tolerance limit. This possibility applies only to the number of measurements $n = 3$ or 4 or 5. The printed-out average value is determined from the 3 or 4 or 5 values which lie within the tolerance limit. If the desired number of measurements within the tolerance limit has not been found after 10 measurements, the average of all measurements is calculated and the message "TIME TOLERANCE EXCEEDED" is printed. All values always appear in the protocol and the selected values are designated. If measured values have been rejected using this procedure, a standard deviation is not calculated, not even from all values including the runaway ones. Default when no entry is made: No search for runaway values, but calculation of standard deviation. The entry <0> means: No " Δ %" evaluation.

The general application of the method according to the elimination of the measurement values is measuring technical not correct. According to the application of quality directives, for example according to ISO 9000 ff, it is important to note the unacceptance of these screening tests.

This procedure of rejecting deviating measured values is technical incorrect but often helpful in routine work, especially if the first value is discrepant because of inadequate temperature equilibration time. The measuring system must always be checked if large deviations or scatter ranges are found. Otherwise it could happen in < Δ %> mode that measurements which agree by coincidence will be designated as correct results.

Symbol <Clock>: Entry and display of the date and the time of day. Entry: <YYMMDDHHMMSS>; example of an entry: <880512083500> for 12 May 1988, 8 h 35 m 0 s. The date and time are preserved even after switching off the equipment, but the date and time must be reset if the mains plug is disconnected.

<ml>: For entry and display of the solvent volume to be used for dilution viscometry. Each dilution step with the corresponding solvent volume is entered separately. The starting volume which the automatic sampler transfers from the sample bottle into the viscometer is 15 ml. The maximum total of added solvent volumes is 60.00 ml.

Example of an entry for 3 dilution steps:

```
<1>  <10.00>  <Enter>  <Clear>
<2>  <12.00>  <Enter>  <Clear>
<3>  <15.00>  <Enter>
```

The digits after the decimal point must be entered.

<V1>, <V2>: The viscometers 1 or 2 are preselected or deselected with the <Clear> key. The constant or the t_0 runtime can then be entered.

<K/ t_0 >: Means "Constant" when <<abs>> is lit or << t_0 runtime>> when <<rel>> is lit.

Example of an entry for constant: <0.1052>. Example of an entry for t_0 runtime: <123.85>.

Any number with not more than 8 digits can be entered as constant. The decimal point may lie in any position. If there are 6 digits after the decimal point, the zero before the decimal point is not displayed. The maximum value for the constant is $K = 100$. This is a practical limit set by the physical dimensioning of the viscometers. Higher values can be entered.

This key has no effect in "t" mode (time measurement). But if a Hagenbach correction is to be made in pure time measurements, the absolute or relative constant of the viscometer must be entered in <<abs>> mode.

3.4 Sample transfer unit

The sample transfer unit is used to draw the liquid from the sample bottle and to deliver it into the viscometer. The sample is drawn via a stainless steel or glass needle and a flexible PTFE hose into a calibrated cylinder. Dosing is made according to the piston burette principle. The piston is made of PTFE and designed for smallest possible dead volume. Switchover of the needle to the individual viscometers is done with a ceramic plate valve. A highly inert Kalrez gasket is located between the borosilicate glass cylinder and the PVDF part of the valve. The filling quantities for the individual viscometers are controlled via different piston strokes (Fig. 5 to 10). For specific usages (for example higher viscose samples) a special program is available to you that is included with the sample transfer unit. Here the suction velocity and the piston strokes are adapted to the demands. Especially the first stroke for the filling is selected very high. With the addition of all negative tolerances of the viscometers selected and hoses used, it may occur that this volume is too high. In this case, while adjusting the upper limit switches of the dosing unit, find a compensation (also see Fig. 8 and 10). The easily accessible dosing unit (piston, cylinder and valve) has three ports in the front: two openings for the two viscometers and a third opening.

The third opening is without any function and has to remain closed. Only steel screws may be used to attach the PTFE hoses on the front openings. Openings which are not used must be closed with steel plugs. If such a screw is missing, a corresponding error message is printed and the ViscoSYSTEM AVS 500 does not start its functioning operation. This safety precaution prevents the possibility of liquid leaking out and causing personal injury or damaging the ViscoSYSTEM AVS 500.

You can insert a filter between the needle and the hose. This prevents solid suspended particles (e. g. fluffs) in the solution from entering the viscometer. This protective measure is highly recommended when using viscometers with an approximate constant $K \leq 0,01 \text{ mm}^2/\text{s}^2$ because of the small capillary diameters.

The filter plate with pore size 60...100 μm (ISO 4793, P 100) or the filter plate with pore size 100...160 μm (ISO 4793, P 160) is fitted according to Fig. 10, (not recommended for Micro Viscometers).

Filter case		type no. VZ 6236
Filter plate	60...100 μm	type no. VZ 6237
Filter plate	100...160 μm	type no. VZ 6238

In cases of usage where the slightly higher dead volume of the filter equipment described above for glass filter is disturbing (for example Micro-Viscometers) the filter can be used for ECTFE fibrous filter discs.

Filter case		type no. VZ 6295
Clamping ring		type no. VZ 6296
Filter disk ECTFE	25 μm	type no. VZ 6297
Filter disk ECTFE	60 μm	type no. VZ 6298
Filter disk ECTFE	115 μm	type no. VZ 6299

The needle which moves up and down is surrounded by a needle protector which is dimensioned in the serially manufactured unit such that when using 100 ml bottles it is not possible to place a finger between the bottle and the needle protector. For sample bottles or vessels of a different size, the height of the needle protector must be adjusted accordingly (Fig. 8, Pos. 6). For this purpose, release the two screws, adjust the height of the needle protector and then tighten the screws again. After this adjustment the distance between the needle protector and the top edge of the bottle or vessel should be less than 5 mm. When the needle protector is in the raised position, the drive motor is disabled. Over the bottom and of the needle protection a splash protection is installed, that also functions as a stripping-off device when using 20 ml vessels in the 36 sample plate even if, despite the clamping effect of the vessel holding device, a sample vessel of the sample needle is transported upwards.

The materials used for the construction of the dosing system (glass, PTFE, PVDF, Kalrez) permit dosing of most liquids without damaging the system. When working with very aggressive liquids it is advisable to consult a table of chemical resistance values.

3.5 Sample turntable

The following sample turntables are available for operating with different sample vessels (Fig. 28 and 29):

12 sample turntable for	100 ml	screw-cap and ground stopper bottles,	type no. VZ 6204
12 sample turntable for	50 ml	screw-cap bottles,	type no. VZ 6205
12 sample turntable for	50 ml	ground stopper bottles, brown,	type no. VZ 6206
36 sample turntable for round bottom tubes, 20 ml,			type no. VZ 6207

When working with very volatile liquids it is advisable to seal the sample bottles with film (Fig. 27).

Adjacent to each opening for the sample bottles on the sample turntables there is a hole for inserting the rod magnet indicating the "last sample". In a sequence of measurements this rod magnet marks the last sample and thus sets the stop position.

The rod magnet may be moved to any other position while the measurements are running, provided that the last sample has not yet been reached. No measurement is possible when the rod magnet is located at the 1st sample. Print-out: **Measuring sequence finished**. Remedy: insert the rod magnet only after the needle has reached the bottom position! The rod magnet is recognized only together with a sample glass. Remedy: place the rod magnet only in positions containing a sample bottle.

The steel basin in which the sample turntable is placed incorporates an optical sample detector. This ensures that the needle will not move down in an empty sample position. When such empty positions are encountered, the turntable will continue to rotate until a position with a sample bottle or the "**last sample**" position comes. The position numbers are marked on the turntable. The correct sample is recognized electronically with the aid of a coding inside the VISCOsystem AVS 500.

When inserting the sample turntable, it must be turned by hand until it is keyed downward into position. This can be felt unambiguously. The automatic coding and proper function of the VISCOsystem AVS 500 is therewith ensured. While inserting the sample turntable, the needle protector must be raised slightly with the finger (Fig. 11).

In the latched in state the sample turntable must **not** be turned by hand; use the corresponding key on the front panel.

3.6 Measuring amplifier cards

Two different measuring amplifier cards are available:


Measuring amplifier card VZ 6211 for viscometers with optoelectronic sensing (Fig. 14)

Measuring amplifier card VZ 6212 for TC-Viscometers (Fig. 15)

Three switch positions are available on the measuring amplifier cards.

See figure next page!

Layout Meßverstärkerkarte
 Layout Measuring amplifier card
 Layout Carte d'amplification

The switch VZ 6203/VZ 6208 (SW 201) has to be set accordingly for the operation with the pneumatic module (see chapter 3.7 "Pneumatic Modules"). Different working parameters can be set with both switches in parallel SW 202 and SW 302 with the printings 0 and 1:

SW203	SW202	binär	Wartezeit nach Messen Delay time after measurement Attente après mesure Minimum	Wartezeit Teiler Divisor delay time Temps d'attente	Nachpumpen nach Spülen Discharge time after rinsing Vidange après rinçage
1	1	0 0	10 sec	07 = 14,3%	25 sec
1	0	0 1	25 sec	06 = 16,7%	50 sec
0	1	1 0	55 sec	04 = 25,0%	60 sec
0	0	1 1	55 sec	04 = 25,0%	185 sec

These switch positions make it possible to select the appropriate time combination for the respective viscosity range. The following cross-reference list only serves as an example to make a selection easier and does not replace the test for an optimum selection in connection with the respective usage case. The maximum viscosity range that can be reached lies at approx. 1200 mm²/s at a room temperatur of approx. 25 °C.

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This is valid for both VZ 6211 and VZ 6212.

Due to measuring technical reasons you only may use the same viscometers and also the same measuring amplifier cards in a parallel way.

During the measurement run the measuring amplifier cards automatically acquires the flow time and reports this time at the termination of the measurement to the central computer unit. The measuring amplifier card is connected electrically to the measuring stand or to the TC-Viscometer. Internally the measuring amplifier cards control the pneumatic modules.

3.7 Pneumatic modules

For the optimum adaption of the VISOsystem AVS 500 to the respective measuring task, four different pneumatic modules are at your dispositionn:

- Pneumatic module VZ 6203 operative range: aromatic and chlorinated solvents
- Pneumatic module VZ 6208 operative range: aggressive solvents (formic acid and others)
- Pneumatic module VZ 6213 operative range: higher viscose samples (mineral oils etc.)
- Pneumatic module VZ 6214 operative range: special ketones as solvents

The pneumatic modules VZ 6203 and 6214 differ from each other due to the materials used for the gasket and membranes. The pneumatic module VZ 6213 has a continuous current-motor for the operation of the membrane pump. The main parts of each pneumatic module are the pump and the valve blocks. Pneumatic connections from the viscometer are made with the pneumatic module (Fig. 18, 18a and 19).

The pump pressure regulation is automatic, depending on the viscosity of the measured good and the used viscometer. These parameters serve as a control basis for the sucking pressure of the pump.

The valves of the module are mounted on a valve block. This design increases the operational reliability. When setting-up the VISOsystem AVS 500, make sure that the measuring amplifiers and the pneumatic modules are inserted in pairs into the correct slots. Thus it is possible, for example, when using only one measuring position to make this the measuring position no. 2 leaving position no. 1 unoccupied. This must be taken into account also when connecting the viscometer to the sample transfer unit of the VISOsystem AVS 500: Connect the filling hose to "Visco 2" (see also Section 2.2).

The pneumatic module VZ 6208 uses air-jet pumps to produce a low pressure according to the Venturi-principle. That is why no valve block is necessary for this pneumatic module. The adaptation of the pump pressure is carried out as described above.

3.8 Printer, inserting the paper

The built-in printer operates with plain paper and ink ribbon, and prints 24 characters per line in 5 x 7 dot character matrix format. A roll of paper type no. VZ 6248 has been installed in the factory.

To insert a new roll of paper, take off the front of the printer, see Fig. 12. Take out the empty paper roll together with the spindle. For this purpose, pull the metal spindle in the direction of the operator. Place a new roll of paper onto the spindle and insert in reversed sequence such that the paper strip faces downward. Pass the paper between the ink ribbon and the ink ribbon holder, whereby it is convenient to take out the ink ribbon (this is easily possible). Reinsert the ink ribbon after passing the paper between the ink ribbon and the ink ribbon holder. Then pass the paper through the slot in the cover and remount the cover.

The paper is advanced automatically after each line print. Several blank lines of paper feed are generated after each result print-out, so that the record can be torn off. The paper is advanced electrically by one line each time the <P> key is depressed. The paper is advanced continuously when this key is depressed in and held.

3.9 The data memory

The VISCOsystem AVS 500 has a data memory for the most recently entered parameter values "number of measurements", "constants", " t_0 runtimes", "time tolerance", "temperature equilibration time", "bath temperature", "viscometer type", "calibration temperature", "rinsing", "date" and "time of day".

The data stored in memory is printed when the VISCOsystem AVS 500 is powered-up:

```
*****
*                                     *
*   AVS 500 MEMORY                   *
*                                     *
*****
Constant 1  =   0.03124
Constant 2  =   0.03095
t0 runtime 1 = 145.12   s
t0 runtime 2 = 144.85   s
No. of meas. =    3
Precond.time =    2     m
Hagenb.corr. DIN-Ubbeloh.
Delta %     =    0.5    %
Date        =   19.02.90
Time        =   08h52m26s
*****
```

This print-out gives an overview for the operator and the values stored in memory make it unnecessary to enter the values again when making series of measurements with the same viscometers.

The data memory can also be called up by pressing the <Stop/Reset> key repeatedly when the VISCOsystem AVS 500 is not operating (the LED in the start key is lit and <<RUN>> is not lit). The same print-out then appears as when the VISCOsystem AVS 500 is switched-on.

The print-out is always the current status according to the most recently entered data.

Note: However, the data are transferred permanently to memory only after at least one measurement has been started after data entry.

3.10 Data transmission

The RS-232-C data interface

This interface provides the following facilities:

- Transfer of results to the computer
- Transfer of the memory contents to the computer
- Changing the measuring parameters from the computer
- Remote control of the start command from the computer

Transferring the results to the computer

After completion of a measurement and print-out of the measurement results on the built-in printer of the VISOsystem AVS 500, the results are sent in the same form via the data interface to the computer (see Fig. 13, Pos 2).

Example: Measurement No. 1 in the viscometer 1 is transferred as "1/001 = 123.45 s".

Protocols after the end of a sequence of measurements are sent to the computer in the same form as described for the individual modes.

Read-out of the memory contents by the computer

In order to get a listing of the current memory contents of the VISOsystem AVS 500 from the computer, it suffices to send the two letter combination "MY" (MEMORY) **followed by the control characters "CR"** (Carriage Return) and "LF" (Line Feed) to the VISOsystem AVS 500. The line terminator "LF" is an ASCII character.

The VISOsystem AVS 500 replies by transmitting its memory contents in the following format:

```
*****
*                                     *
*   AVS 500 MEMORY                   *
*                                     *
*****
Constant 1   =   0.01004
t0 runtime 1 = 145.12   s
No. of meas. =    3
Precond.time =    5     m
Hagenb.corr. DIN-Ubbeloh.
Delta %      =    0.5   %
Date        =   19.02.90
Time        =   08h52m26s
=====
```

Technical specifications and plug pinout of the RS-232-C interface

This interface used the following format:

- 7 bit ASCII Code
- no parity bit
- 2 stop bits
- data transmission rate = 4800 Baud

These parameters are permanently programmed and cannot be changed.

Plug connector pinout:

- Pin 1 = Protection ground
- Pin 2 = Data output to computer (TxD)
- Pin 3 = Data input from computer (RxD)
- Pins 4 and 5 are jumpered at the instrument
- Pins 6 and 20 are jumpered at the instrument
- Pin 7 = Logic ground
- All other pins are not used and not connected.

Changing system parameters at the computer

The control codes for changing system parameters at the computer are listed in the following table:

Control code	Meaning	Format	Characters
MO	Mode	MO = 1..4	(1byte)
V1	Viscometer 1	V1 = 0..1	(1byte)
V2	Viscometer 2	V2 = 0..1	(1byte)
VH	Hagenbach correction	VH = 0..1	(1byte)
VA	Viscometer type	VA = 1..9	(1byte)
K1	Constant 1	K1 = 0.1002	(8 + 1 byte)
K2	Constant 2	K2 = 1.002	(8 + 1 byte)
T1	Blank value 1	T1 = 123.456	(8 + 1 byte)
T2	Blank value 2	T2 = 234.567	(8 + 1 byte)
NM	Number of meas.	NM = 01..99	(2 byte)
TT	Temp. equilib. time	TT = 00..99	(2 byte)
DP	Delta percent	DP = 0.5	(5 + 1 byte)
RO	Sequence	RO = 01.36	(5 + 1 byte)
CL	Rinsing	CL = 1.2	(3 + 1 byte)
VO	Volumes	VO = 01/01.00	(6 + 2 byte)
NR	Sample number NR	NR = 01/12345678	(10 + 1 byte)
MY	Memory print-out	MY	
ST	Start ST	ST	

The control character sequence **must be** terminated with "CR" and "LF".

The bytes specified in the brackets have the following meaning:

- 1 byte: The number of digits required after the letter sequence is an ASCII digit.
- 2 byte: two ASCII digits
- 3 + 1 byte: three ASCII digits and one special character, e. g. "."
- 5 + 1 byte: five ASCII digits and one special character, e. g. "."
- 6 + 2 byte: six ASCII digits and two special characters, e. g. "/" "
- 8 + 1 byte: eight ASCII digits and one special character, e. g. "."
- 10 + 1 byte: ten ASCII digits and one special character, e. g. "/"

Data entry errors

Wrong entries made at the computer for the AVS 500 are responded to with the message: "wrong entry".

A wrong entry has no effect on the previous memory contents.

4 Measuring with the VISCOSYSTEM AVS 500

4.1 Overview, brief instructions

- The VISCOSYSTEM AVS 500 is switched-on at the mains switch. The parameter values stored in memory are printed-out automatically ("AVS 500 Memory") and all displays light-up (segment test).
- Connect-up the viscometers and insert them into the thermostatic bath. Check that the hose connections do not leak.
- If necessary, calibrate the temperature display of the VISCOSYSTEM AVS 500.
- Correct the date and time of day if necessary.
- Select the desired function mode with the key <F>: <<abs>> for absolute measurements, <<rel>> for relative measurements, <<t>> for measuring the flow times, <<c → 0>> for intrinsic viscosity determinations.
- The respective values stored in memory can be recalled to the display by depressing one of the keys <K/t₀>, <n>, <No>, <Type HC>, <Precond. time>, <Δ %>, <Rinse>. A new value can be entered after clearing the old value from the display by depressing the key <C> (Clear). To store the new value into memory, it must be confirmed by depressing the <Enter> key.
- Insert the sample turntable with the samples.
- Enter the start position and optionally the end position on the sample turntable, or mark the end position by inserting the magnetic rod for "last sample" there.
- Depress the start key.
- The measuring sequence now runs automatically.
- At the end of the measuring sequence, remove the printed protocol.

4.2 Determining the absolute kinematic viscosity

After switching-on the VISCOSYSTEM AVS 500, the contents of the non-volatile memory are printed-out automatically, e. g.:

```
*****
*                               *
*   AVS 500 MEMORY             *
*                               *
*****
Constant 1   =   0.03124
Constant 2   =   0.03095
t0 runtime 1 = 145.12   s
t0 runtime 2 = 144.85   s
No. of meas. =    3
Precond.time =    2     m
Hagenb.corr. DIN-Ubbeloh.
Delta %      =    0.5    %
Date        =   19.02.90
Time        =   08h52m26s
=====
```

The function mode last used is displayed.

LED indicator: <<abs>>, <<rel>>, <<t>> or <<c → 0>>.

If necessary, press the function key <F> to set the mode to <<abs>>.

Changing the stored data

The stored value can be recalled into the display by depressing the respective parameter key (e. g. <n>). If the displayed value corresponds to the desired value, the next parameter can be recalled into the display by depressing the next key. This makes no change in the memory contents (see also Section 3.3).

If a displayed value is to be changed, depress the key <C> (Clear) to delete the old value in the display. A new value can then be entered. The new value is written into memory on depressing the <Enter> key.

- The date and time of day are called-up with the symbol <Clock>. The new entry must have the format <YYMMDDHHMMSS>.
- Calibration temperature: Called-up with the key <°C cal>. Enter the temperature reading of the calibrated thermometer as calibration temperature. The calibration temperature should not differ from the operating temperature by more than 3 K.
- Entering the nominal temperature and the temperature tolerance: Press the key <°C>, enter the nominal temperature, press <°C> again and then enter the temperature tolerance (not smaller than 0.1 °C).
- On pressing the key <K/t₀>, the stored constant is shown in the display. If necessary, enter the new constant on the keyboard.
- The number of measurements is displayed on pressing the key <n> (maximum number of measurements n = 99). If necessary, enter the new number of measurements on the keyboard.
- The sample number for the current position of the sample turntable is displayed on depressing the <No> key. A new entry consists of two digits for the position and up to eight digits for the sample number (see chapter 3.3 <No>).
- Depress the <Temperature equilibration time> key (thermometer symbol) to display the current value of the temperature equilibration time in minutes and to enter a new value if desired.
- The time tolerance for sorting-out runaway measured values can be entered after depressing the <Δ %> key.
- Depress the <Type HC> key to preselect the calculation of the Hagenbach correction, see also chapter 3.3.
- The next sample is used for rinsing the viscometer. Depress the key with the <Rinse> symbol to make entries for method of rinsing and the number of rinsing cycles. The entry <1.3> means that three rinsing cycles are to be executed with the next sample.

Measurements for the absolute kinematic viscosity

- Any position may be taken as starting position on the sample turntable, thus also the sample which is currently under the needle. If a different start position is desired, depress the <Pos> key and then enter the start position and optionally the last position. The last measurement of the sequence can also be marked with the magnetic rod "last sample".
- Depress the <Start> key to start the program run with the entered data.
- If a temperature equilibration time has been specified, when using optoelectronic tracing, the liquid is continually pumped-up and allowed to flow downward in the viscometer within for this time period, to establish better temperature equilibrium, without making any time measurements, whereas when using TC-Viscometers the sample in the viscometer rests for a temperature adjustment. The remaining equilibration time still left to elapse is continually shown in the display. When this time has expired, the actual measurement commences.

Print-out, evaluation:

```

1/001      =   180.32 s
1/002      =   180.38 s
1/003      =   180.42 s
=====

```

Viscometer 1

```

Table pos.  =         01
Sample no.  =   12345678
Date        =   19.02.90
Time        =  10h30m20s
Bath temp.  =    25.00 C
=====

```

```

1/001      =   180.32 s
1/002      =   180.38 s
1/003      =   180.42 s
=====

```

```

Average      =   180.373 s
HC.CORR      =     0.049 s
Av.Corr.     =   180.324 s
Std. deviat. =     0.050
=====

```

```

Constant 1   =    0.03124s
Visc (abs)   = 5.633 mm²/s
=====

```

- After completion of each measurement, the measured flow time is printed and stored in memory. After print-out of the last measured value, the average of the measured values is determined and the absolute viscosity is calculated and printed. The standard deviation "S" of the time measurements is calculated only when the <Δ %> key has not been programmed.
- Other evaluations are possible in function mode <<abs>>, provided that it can be expressed as a multiplication of the flow time by a constant number. The result is then the conventional operational value instead of the absolute viscosity! It is sensible to use the <<rel>> mode for such an evaluation in order to avoid any conflict with correct constant entry. For relative viscosity measurements it is very advisable to specify a nominal value for the constant because the VISCO SYSTEM AVS 500 automatically chooses the dosing rate and the suck off time according to the value of the constant.

4.3 Determining the relative viscosity

After switching the VISCO SYSTEM AVS 500 on, the contents of the non-volatile data memory is printed-out, e. g.:

```
*****
*                               *
*   AVS 500 MEMORY             *
*                               *
*****
Constant 1   =   0.03124
Constant 2   =   0.03095
t0 runtime 1 = 145.12   s
t0 runtime 2 = 144.85   s
No. of meas. =    3
Precond.time =    2     m
Hagenb.corr. DIN-Ubbeloh.
Delta %      =    0.5   %
Date        =  19.02.90
Time        =  08h52m26s
*****
```

The last active function mode before switch-off is displayed.

LED indicator: <<abs>>, <<rel>>, <<t>> or <<c → 0>>.

If necessary, set the mode to <<rel>> with the function key <F>.

Changing the stored data

The stored value can be recalled into the display by depressing the respective parameter key (e. g. <n>). If the displayed value corresponds to the desired value, call-up the next parameter with the next key. The old value remains unchanged, see also section 3.3.

If a displayed value is to be changed, depress the key <C> (Clear) to delete the old value in the display. A new value can then be entered. The new value is written into memory when the <Enter> key is depressed.

- The date and time of day are called-up by depressing the key with the <Clock> symbol. The new value must be entered with the format <YYMMDDHHMMSS>.
- The calibration temperature is called-up by depressing the key <°C cal>. Enter the temperature reading of the calibrated thermometer as calibration temperature. The calibration temperature should not differ from the operating temperature by more than 3 K.
- To enter the nominal temperature and the temperature tolerance: Depress the key <°C>, enter the nominal temperature, depress the key <°C> again, enter the temperature tolerance (not smaller than 0.1 °C).
- Depress the key <K/t₀₀ runtime to the display. If necessary, enter the new t₀ runtime with the numerical keyboard, t₀ runtimes up to 9999.99 s are possible.
- Depress the <n> key to see the number of measurements (maximum number of measurements n = 99). If necessary, enter the new number of measurements with the numerical keyboard.
- Depress the <No> key to display the sample number for the current position of the sample turntable. The new entry must consist of two digits for the position and up to eight digits for the sample number (see chapter 3.3 <No>).
- To enter the temperature equilibration time in minutes, depress the corresponding key <Temperature Equilibration> (thermometer symbol).
- The time tolerance beyond which values are to be considered as runaway values, can be entered after depressing the <Δ %> key.
- For calculating the Hagenbach correction "HC", the constant (or the approximate constant) of the viscometer must be entered in <<abs>> mode. This entry is not necessary if a Hagenbach correction is not required. When measuring the t₀ runtime, the Hagenbach correction is always necessary if it is taken into account in the determination of the flow through times for the sample solutions.
- The next sample is used to rinse the viscometer. Depress the key with the <Rinse> symbol to enter the method of rinsing and the number of rinsing cycles. The entry <1.3> means that three rinsing cycles are to be run with the next sample.

Measurements for the relative viscosity

- The start position can be any position of the sample turntable, also the sample which happens to be under the needle. If a different start position is desired, depress the <Pos> key and then enter the start position and optionally also the last position. The position for the last measurement can also be marked with the magnetic rod "last sample".
- Depress the <Start> key to start the program run with the entered data.
- If a temperature equilibration time has been specified, the liquid is continually pumped-up and allowed to flow downward in the viscometer within this time period, without making any time measurements, for better temperature equilibration. The remaining time is continually shown in the display. The actual measurement commences when this time has expired.

Print-out, evaluation:

```

1/001      =   182.39 s
1/002      =   182.41 s
1/003      =   182.43 s
=====
                Viscometer 1
Table pos.  =           01
Sample no.  =   12345678
Date        =   19.02.90
Time        =  10h30m20s
Bath temp.  =   25.00 °C
=====
1/001      =   182.39 s
1/002      =   182.41 s
1/003      =   182.43 s
=====
Average     =   182.410 s
HC.CORR     =     0.048 s
Av.Corr.    =   182.362 s
Std. deviat. =     0.020
=====
t0 runtime 1 =   145.12 s
Visc (rel)   =     1.257
=====

```

- After completion of each measurement, the determined flow time is printed and stored in memory. After the last measured value has been printed, the average of the measured values is determined and the relative viscosity is calculated and printed. The standard deviation "S" of the time measurements is calculated only when the <Δ %> key has not been programmed.
- Other evaluations are possible in function mode <<rel>>, provided that it can be expressed as division of the flow time by a constant number. The result is then the conventional operational value instead of the relative viscosity.

Because during the viscosity determination with capillary viscometers, the kinematic viscosity is determined fundamentally without external pressure, the calculation is only valid in case that the density of measurement liquid and comparisons liquid (e.g. pure solvent) are the same.

4.4 Determining the flow time

After switching-on the VISCO SYSTEM AVS 500, the contents of the non-volatile data memory is printed-out, e. g.:

```
*****
*                               *
*   AVS 500 MEMORY             *
*                               *
*****
Constant 1   =   0.03124
Constant 2   =   0.03095
t0 runtime 1 = 145.12   s
t0 runtime 2 = 144.85   s
No. of meas. =    3
Precond.time =    2     m
Hagenb.corr. DIN-Ubbeloh.
Delta %      =    0.5   %
Date         =   19.02.90
Time         =   08h52m26s
=====
```

The last active function mode before switch-off is displayed.

LED indicator: <<abs>>, <<rel>>, <<t>> or <<c → 0>>.

If necessary, press the function key <F> to set the mode to <<t>>.

Changing the stored data

The stored value can be recalled into the display by depressing the respective parameter key (e. g. <n>). If the displayed value corresponds to the desired value, call-up the next parameter with the next key. The old value remains unchanged (see Section 3.3).

If a displayed value is to be changed, depress the key <C> (Clear) to delete the old value in the display. A new value can then be entered. The new value is written into memory when the <Enter> key is depressed.

- The date and time of day are called-up by depressing the key with the <Clock> symbol. The new value must be entered with the format <YYMMDDHHMMSS>.
- The calibration temperature is called-up by depressing the key <°C cal>. Enter the temperature reading of the calibrated thermometer as calibration temperature. The calibration temperature should not differ from the operating temperature by more than 3 K.
- To enter the nominal temperature and the temperature tolerance: Depress the key <°C>, enter the nominal temperature, depress the key <°C> again, enter the temperature tolerance (not smaller than 0.1 °C).
- Depress the <n> key to see the number of measurements (maximum number of measurements n = 99). If necessary, enter the new number of measurements with the numerical keyboard.
- Depress the <No> key to display the sample number for the current position of the sample turntable. The new entry must consist of two digits for the position and up to eight digits for the sample number (see Chapter 3.3 <No>).
- To enter the temperature equilibration time in minutes, depress the corresponding key <Temperature Equilibration> (thermometer symbol).
- The time tolerance beyond which values are to be considered as runaway values, can be entered after depressing the <Δ %> key.
- If the Hagenbach correction "HC" is to be calculated, the actual or approximate constant of the viscometer must be entered in <<abs>> mode. This entry is not necessary if a Hagenbach correction is not desired.
- The next sample is usually used to rinse the viscometer. Depress the key with the <Rinse> symbol to enter the method of rinsing and the number of rinsing cycles. The entry <1.3> means that three rinsing cycles are to be run with the next sample.

Measuring the flow through time

- The start position can be any position of the sample turntable, also the sample which happens to be under the needle. If a different start position is desired, depress the <Pos> key and then enter the start position and optionally also the last position. The position for the last measurement can also be marked with the magnetic rod "last sample".
- Depress the <Start> key to start the program run with the entered data.
- If a temperature equilibration time has been specified, when using optoelectronic tracing, the liquid is continually pumped-up and allowed to flow downward in the viscometer within for this time period, to establish better temperature equilibrium, without making any time measurements, whereas when using TC-Viscometers the sample in the viscometer rests for a temperature adjustment. The remaining equilibration time still left to elapse is continually shown in the display. When this time has expired, the actual measurement begins.

Print-out, evaluation:

1/001	=	145.09 s
1/002	=	145.12 s
1/003	=	145.15 s
=====		
Viscometer 1		
Table pos.	=	01
Sample no.	=	12345678
Date	=	19.02.90
Time	=	10h30m20s
Bath temp.	=	25.00 °C

1/001	=	145.09 s
1/002	=	145.12 s
1/003	=	145.15 s

Average	=	145.120 s
HC.CORR	=	0.076 s
Av.Corr.	=	145.044 s
Std. deviat.	=	0.030

- After completion of each measurement, the determined flow time is printed and stored in memory. After the last measured value has been printed, the average of the measured values is determined and the relative viscosity is calculated and printed. The standard deviation "S" of the time measurements is calculated only when the <Δ %> key has not been programmed.

4.5 Determining the inherent viscosity

After switching-on the VISCO SYSTEM AVS 500, the contents of the non-volatile data memory is printed-out, e. g.:

```
*****
*                               *
*   AVS 500 MEMORY             *
*                               *
*****
Constant 1   =   0.03124
Constant 2   =   0.03095
t0 runtime 1 = 145.12   s
t0 runtime 2 = 144.85   s
No. of meas. =    3
Precond.time =    2     m
Hagenb.corr. DIN-Ubbeloh.
Delta %      =    0.5   %
Date        =  19.02.90
Time        =  08h52m26s
=====
```

The last active function mode before switch-off is displayed.

LED indicator: <<abs>>, <<rel>>, <<t>> or <<c → 0>>.

If necessary, press the function key <F> to set the mode to <rel>.

Changing the stored data

The stored value can be recalled into the display by depressing the respective parameter key (e. g. <n>). If the displayed value corresponds to the desired value, call-up the next parameter with the next key. The old value remains unchanged (see also Section 3.3).

If a displayed value is to be changed, depress the key <C> (Clear) to delete the old value in the display. A new value can then be entered. The new value is written into memory when the <Enter> key is depressed.

- The date and time of day are called-up by depressing the key with the <Clock> symbol. The new value must be entered with the format <YYMMDDHHMMSS>.
- The calibration temperature is called-up by depressing the key <°C cal>. Enter the temperature reading of the calibrated thermometer as calibration temperature. The calibration temperature should not differ from the operating temperature by more than 3 K.
- To enter the nominal temperature and the temperature tolerance: Depress the key <°C>, enter the nominal temperature, depress the key <°C> again, enter the temperature tolerance (not smaller than 0.1 °C).

- Depress the key <K/t₀> to recall the t₀ runtime to the display. If necessary, enter the new t₀ runtime with the numerical keyboard, t₀ runtimes up to 9999.99 s are possible.
- Depress the <n> key to see the number of measurements (maximum number of measurements n = 99). If necessary, enter the new number of measurements with the numerical keyboard.
- Depress the <No> key to display the sample number for the current position of the sample turntable. The new entry must consist of two digits for the position and up to eight digits for the sample number (see chapter 3.3 <No>).
- The temperature equilibration time can be entered after pressing the <temperature equilibration> key (thermometer symbol).
- The time tolerance beyond which values are to be considered as runaway values, can be entered after depressing the <Δ %> key.
- The constant (or approximate constant) of the viscometer must be entered in <<abs>> mode for calculation of the Hagenbach correction "HC". This entry is not necessary if the Hagenbach correction is not required. Always bear in mind that the Hagenbach correction is also required for determination of the t₀ runtime if it is taken into account in the calculation of the run through time of the sample solutions.
- To enter the concentration, depress the key <Type HC>. Enter the viscometer type before the decimal point (see also Section 3.3). The digit value after the decimal point specifies the concentration of the solution. The values 2 to 5 are for calculation without Hagenbach correction:

2 = 0.25 g/100 ml

3 = 0.50 g/100 ml

4 = 0.75 g/100 ml

5 = 1.00 g/100 ml

The values 6 to 9 have the same respective concentration meanings,
but for calculation with Hagenbach correction.

- The next sample is used to rinse the viscometer. Depress the key with the <Rinse> symbol to enter the method of rinsing and the number of rinsing cycles. The entry <1.3> means that three rinsing cycles are to be run with the next sample.

Measurements for the inherent viscosity

- The start position can be any position of the sample turntable, also the sample which happens to be under the needle. If a different start position is desired, depress the <Pos> key and then enter the start position and optionally also the last position. The position for the last measurement can also be marked with the magnetic rod "last sample".
- Depress the <Start> key to start the program run with the entered data.
- If a temperature equilibration time has been specified, the liquid is continually pumped-up and allowed to flow downward in the viscometer within this time period, without making any time measurements, for better temperature equilibration. The remaining time is continually shown in the display. The actual measurement commences when this time has expired.
- After completion of each measurement, the determined flow time is printed and stored in memory. After print-out of the last measured value, the average of the measured values is determined and the inherent viscosity is calculated and printed out. The standard deviation "S" of the time measurements is calculated too.
- Calculation of the inherent viscosity is based on the relative dynamic viscosity.

Viscosity determinations with capillary viscometers without applying any external pressure always yield the kinematic viscosity, so that the calculation is valid exactly only when the density of the sample solution is the same as the density of the reference solution (e. g. pure solvent).

Print-out, evaluation:

```
1/001      =    184.51 s
1/002      =    184.47 s
1/003      =    184.49 s
```

```
=====
```

Viscometer 1

```
Table pos. =          01
Sample no. =    12345678
Date       =    19.02.90
Time       =    10h30m20s
Bath temp. =    25.00 °C
```

```
-----
1/001      =    184.51 s
1/002      =    184.47 s
1/003      =    184.49 s
```

```
-----
Average    =    184.490 s
HC.CORR    =    2.080 s
Av.Corr.   =    182.410 s
Std. deviat. =    0.020
```

```
-----
t0 runtime 1 =    145.12 s
Concentrat. =    0.5 g/100 ml
Visc.(rel)  =    1.257
Visc.(spc)  =    0.257
Visc.(red,J) =    51.392
Visc.(inh)  =    35.73
```

```
=====
```

The calculation basis is:

Absolute kinematic viscosity:

$$\nu = K \cdot t$$

Relative viscosity:

$$\eta_{\text{rel}} = \frac{\eta}{\eta_0} = \frac{t}{t_0}$$

Inherent viscosity (logarithmic viscosity):

$$\eta_{\text{inh}} = \frac{\ln \cdot \eta_{\text{rel}}}{c}$$

Specific viscosity:

$$\eta_{\text{spec}} = \frac{\eta}{\eta_0} - 1 = \frac{t}{t_0} - 1$$

Reduced viscosity (viscosity number J):

$$\eta_{\text{red}} = \frac{\eta_{\text{spez}}}{c}$$

- ν = absolute kinematic viscosity
- K = viscometer constant
- η = dynamic viscosity of the solution in [mPa · s]
- η_0 = dynamic viscosity of the solvent in [mPa · s]
- t = flow time of the solution in [s]
- t_0 = flow time of the solvent in [s]
- c = concentration of the solution in [g/100 ml]

4.6 Brief instructions for dilution viscometry

- See also Section 4.1 "Overview, brief instructions".
- Measure the t_0 runtime of the pure solvent (in this respect, see also Section 4.4 "Determining the flow time").
- After selecting the mode <<rel>>, enter the t_0 runtime for the viscometer 1 (V1).
- With the key <F> select the function mode for limiting viscosity determination <<c → 0>>.
- Enter the individual dilution steps by pressing the key <ml>.
- A new value can be entered after pressing the key <C> (Clear). Press the <Enter> key to confirm.
- To transfer the entered dilution parameters to working memory, press the <ml> key after completing the entry.

Example of entry for 3 dilution steps:

Key	Numerical keypad	Terminator key
<ml>	<1 10.00>	<Enter>
<C>	<2 12.00>	<Enter>
<C>	<3 15.00>	<Enter>
<ml>		
<START>		

4.7 Determining the limiting viscosity number (dilution viscometry)

Before switching on the VISCOsystem AVS 500, check that a piston burette VISCODOSER® AVS 20, TITRONIC® T 110 or T 90 (in each case with 50 ml exchangeable unit) is connected and switched on (first filling). If this is not the case, connect the piston burette using the data cable type no. TZ 1952 to the output "Dos. [ml]" on the rear of the VISCOsystem AVS 500, see Fig. 13, Pos. 3).

The contents of the permanent memory are printed automatically after switching on the VISCOsystem AVS 500, e. g.:

```

*****
*                               *
*   AVS 500 MEMORY             *
*                               *
*****
Constant 1   =   0.03124
Constant 2   =   0.03095
t0 runtime 1 = 145.12   s
t0 runtime 2 = 144.85   s
No. of meas. =    3
Precond.time =    2     m
Hagenb.corr. DIN-Ubbeloh.
Delta %      =    0.5    %
Date        =   19.02.90
Time        =   08h52m26s
=====

```

The last active function mode before switch-off is displayed.

LED indicator: <<abs>>, <<rel>>, <<t>> or <<c → 0>>.

If necessary, set the mode to <<c → 0>> with the function key <F>.

Determination of the limiting viscosity is possible only with a viscometer in position 1.

Changing the data stored in memory

If only the number of dilutions is to be changed, i. e. the volumes remain the same but, for example, 2 instead of 3 dilutions are to be carried out, press the key <ml> to display the volumes 1 and 2 in cyclic sequence and confirm them by pressing <Enter> and pressing the <ml> key again (see also Section 3.3).

Important: Dilution measuring sequences are possible only with 2 or more (max. 10) steps.

The stored value can be recalled into the display by depressing the respective parameter key (e. g. <n>). If the displayed value corresponds to the desired value, call-up the next parameter with the next key. The old value remains unchanged.

If a displayed value is to be changed, depress the key <C> (Clear) to delete the old value in the display. A new value can then be entered. The new value is written into memory when the <Enter> key is depressed.

- The date and time of day are called-up by depressing the key with the <Clock> symbol. The new value must be entered with the format <YYMMDDHHMMSS>.
- The calibration temperature is called-up by depressing the key <°C cal>. Enter the temperature reading of the calibrated thermometer as calibration temperature. The calibration temperature should not differ from the operating temperature by more than 3 K.
- To enter the nominal temperature and the temperature tolerance: Depress the key <°C>, enter the nominal temperature, depress the key <°C> again, enter the temperature tolerance (not smaller than 0.1 °C).
- Depress the key <K/t₀> to recall the t₀ runtime to the display. If necessary, enter the new t₀ runtime with the numerical keyboard, t₀ runtimes up to 9999.99 s are possible.
- Depress the <n> key to see the number of measurements (maximum number of measurements n = 99). If necessary, enter the new number of measurements with the numerical keyboard. The default value of 1 is used if no entry has been made.
- Depress the <No> key to display the sample number for the current position of the sample turntable. The new entry must consist of two digits for the position and up to eight digits for the sample number (see chapter 3.3 <No>).
- The temperature equilibration time in minutes can be entered after pressing the <temperature equilibration> key (thermometer symbol).
- The time tolerance beyond which values are to be considered as runaway values, can be entered after depressing the <Δ %> key.
- If the Hagenbach correction "HC" is to be calculated, the actual or nominal constant of the viscometer must be entered in <<abs>> mode. This entry is not necessary if a Hagenbach correction is not desired. The Hagenbach correction must also be made when determining the t₀ runtime if it is to be taken into account when determining the flow times of the sample solutions.
- To enter the concentration, press the key <Type HC>. The viscometer type is specified before the decimal point. The entry <4.3> corresponds to a dilution viscometer based on DIN and concentration 0.5 g/100 ml (see also Section 3.3).
- The digit value after the decimal point specifies the concentration of the solution. The values 2 to 5 for calculation without Hagenbach correction are:

2 = 0.25 g/100 ml

3 = 0.50 g/100 ml

4 = 0.75 g/100 ml

5 = 1.00 g/100 ml

The values 6 to 9 have the same respective concentration meanings, but for calculation with Hagenbach correction.

- Only entry of a Dilution Viscometer according to DIN is permitted; all other entries lead to incorrect functioning of the VISCO SYSTEM AVS 500.
- The next sample is usually used to rinse the viscometer. Depress the key with the <Rinse> symbol to enter the method of rinsing and the number of rinsing cycles. The entry <1.3> means, for example, that three rinsing cycles are to be run with the next sample. For exact dosing **it is essential to rinse twice**.
- A measurement with only one dilution is not possible; the minimum required number of dilutions is 2! For changing the <ml> entry (see also Section 3.3, paragraph <ml>).

Measurements for dilution viscometry

- The start position can be any position of the sample turntable, also the sample which happens to be under the needle. If a different start position is desired, depress the <Pos> key and then enter the start position and optionally also the last position. The position for the last measurement can also be marked with the magnetic rod "last sample".
- Depress the <Start> key to start the program run with the entered data.
- A temperature equilibration time must always be chosen because adequate time must be allowed to elapse for mixing and temperature equilibrium establishment of the sample and the added solvent.
- If a temperature equilibration time has been specified, the liquid is continually pumped-up and allowed to flow downward in the viscometer within this time period, without making any time measurements, for better temperature equilibration. The remaining time is continually shown in the display. The actual measurement commences when this time has expired.

Print-out, evaluation:

1/001	= 184.51 s	1/001	= 142.67 s	1/001	= 125.25 s
1/002	= 184.47 s	1/002	= 142.71 s	1/002	= 125.28 s
1/003	= 184.49 s	1/003	= 142.63 s	1/003	= 125.31 s
=====					
Viscometer 1		Viscometer 1		Viscometer 1	
Sampler pos.	= 01	Sampler pos.	= 01	Sampler pos.	= 01
Sample no.	= 12345678	Sample no.	= 12345678	Sample no.	= 12345678
Date	= 19.02.90	Date	= 19.02.90	Date	= 19.02.90
Time	= 10h30m20s	Time	= 10h40m30s	Time	= 11h00m40s
Bath temp.	= 25.00 °C	Bath temp.	= 25.00 °C	Bath temp.	= 25.00 °C

1/001	= 184.51 s	1/001	= 142.67 s	1/001	= 125.25 s
1/002	= 184.47 s	1/002	= 142.71 s	1/002	= 125.28 s
1/003	= 184.49 s	1/003	= 142.63 s	1/003	= 125.31 s

Average	= 184.490 s	Average	= 142.670 s	Average	= 125.280 s
HC.CORR	= 2.080 s	HC.CORR	= 3.490 s	HC.CORR	= 4.490 s
AVG. corr.	= 182.410 s	AVG. corr.	= 139.180 s	AVG. corr.	= 120.790 s
Std. deviat.	= 0.020	Std. deviat.	= 0.040	Std. deviat.	= 0.030

t ₀ runtime 1	= 104.55 s	t ₀ runtime 1	= 104.55 s	t ₀ runtime1	= 104.55 s
Concentrat.	= 1.00g/100 ml	Concentrat.	= 0.50g/100ml	Concentrat.	= 0.25g/100ml
Visc.no. (J)	= 74.472	Visc.no.(J)	= 66.246	Visc.no.(J)	= 62.133
				Lim.visc.no.	= 58.020
=====					

A new protocol is printed after each dilution step. The intrinsic viscosity number (intrinsic viscosity) is calculated and printed after the last protocol.

4.8 Terminating the measurement

Various conditions - depending on the particular situation - lead to termination or abortion of the measuring program:

- Automatic termination when the program has executed for all samples with the specified number of measurements and evaluation. This is the normal case.
- <Stop/Reset>: The measurements are aborted immediately without evaluation. The state as after switch-on and after print-out of the memory contents is reestablished. The measured values are cleared. The data for the measuring parameters in the non-volatile memory (AVS 500 memory) are preserved. This method of termination becomes necessary in the case of danger, breakage of the viscometer or other kinds of disturbance.
- By response of the safety devices: These safety devices may respond, for example, when the waste vessel gets too full. The program is aborted immediately, as when the <Stop/Reset> key is depressed, but without memory print out. Instead, a corresponding error message is issued.
- Mains voltage failure: After return of the mains voltage, the state is the same as after first switch on. The memory contents are printed, the measuring parameters are preserved. The pure measured values were deleted at the time of mains voltage failure. However, the individual time measurements are never lost because the printer has printed the measured value after every time measurement. Thus "manual" evaluation is always still possible.

The VISCO SYSTEM AVS 500 must be emptied and cleaned after the measurements, particularly if there is a waiting time until the next measurements are due to be made. This is very important for all parts which have come into contact with the sample liquids.

5 Routine checks and trouble shooting

5.1 Faults which are not reported by the VISCO SYSTEM AVS 500

The mains switch (power switch) must be in setting <Off> when connecting the mains cable and measuring stand cables and when changing the measuring amplifier plug in modules or the pump insert unit, otherwise the VISCO SYSTEM AVS 500 could become damaged in a way which can be repaired only by service technicians.

Symptoms	Cause	Remedy
After connecting all cables and switching on the instrument, the green pilot lamp remains dark and no functions are accessible.	Mains voltage not connected.	Connect the mains cable and check the mains voltage.
	Main fuse blown.	Replace the main fuse in the Euro power connector (on the rear of the instrument). Instrument fuse (5 x 20 mm): for 220 V~, 0,8 A delayed for 110 V~, 1,6 A delayed
Loud noise in the pneumatic module.		Remove the transport securing device.
The transfer unit does not function.	Needle protection not in right position.	Push the needle protection to the bottom stop position.
No pump runs after switching on the instrument but "AVS Memory" is printed and a keyboard entry is possible when the printing has finished.	The pump insert unit is missing or not fitted correctly.	If the fault persists after correct push in the switched off state of the instrument, the pump insert unit is defective. Push a new pump insert unit.

5.2 Error messages reported by the VISCO SYSTEM AVS 500 by print-out or warning lamp

	Error (text printed)	Remedy
1.	Hose missing on the sample transfer unit, change and then press start key.	Connect the hose correctly. Screw the blind stopper into the opening which is not used.
2.	Sample bottle missing!	The table automatically continues to turn until a position with a sample bottle is reached.
3.	Bad entry!	Make a correct entry via the RS-232-C data interface.
4. 5.	Measuring amplifier 1 missing! Measuring amplifier 2 missing!	Check that the measuring amplifier plug in card is in correct position. If necessary, push it in properly.
6.	Measuring amplifiers missing!	Push in the measuring amplifier plug in card properly.
7.	No parameters entered for the dilution sequence.	Press the parameter key, clear the display with the <C> key and then enter the correct values.
8.	Constant 1 missing	same
9.	Constant 2 missing	same
10.	No runtime $t_0/1$	same
11.	No runtime $t_0/2$	same
12.	No type entry	same
13.	Viscometer does not match the turntable!	Please observe that with small sample glasses you cannot use viscometers which require more volume than is available.
14. 15.	Defective NTC Viscometer no. 1! Defective NTC Viscometer no. 2!	Check TC Viscometer no. 1 and measuring cable. Check TC Viscometer no. 2 and measuring cable.
16.	No report back from turntable! Please call Service!	Insert the turntable! Call SCHOTT-GERÄTE Customer Service.
17. 18.	Visc. 1 overfilled! Visc. 2 overfilled!	Check the mounting clearance of the fill level sensor in the viscometer (Fig. 24, Pos. 11). (The sensor must not dip into the viscometer supply vessel).
19. 20.	Waste bottle 1 too full! Waste bottle 2 too full!	Empty the waste bottle. Check that the viscometer frame is seated correctly in the stand.
21. 22.	Safety vessel 1 too full! Safety vessel 2 too full!	If it is clogged, please clean the aperture in the safety bottle.
23.	Temp. Tol.	Choose a proper nominal temperature (see also Section 3.3, paragraph: <°C>); if necessary choose a larger temperature tolerance (greater than 0.1 °C) and check the thermostat bath.

5.3 Replacing functional units

Caution: Carry out only when the VISCOSYSTEM AVS 500 is disconnected from the mains voltage!

Replacing a defective sample transfer unit

The complete sample transfer unit can be replaced after removing the front panel (Fig. 5):

- a) Release and unscrew the upper mounting bolt (Fig. 6, Pos. 1).
- b) Remove the two screws at the bottom left and right (Fig. 6, Pos. 1).
- c) Take the sample turntable out of the VISCOSYSTEM AVS 500 (Fig. 11) because otherwise the transfer unit could hit against it because the center of gravity lies far to the rear and the weight is considerable (large tilt torque).
- d) Take out the transfer unit to the front.
- e) Mount the new unit in reversed sequence.

Replacing the dosing unit

- a) Prise the 10 pole ribbon cable connector out of the plug by pressing on the locking device (Fig. 7).
- b) Unscrew the transfer hose at the needle holder (above the needle) (Fig. 6).
Caution: Dripping liquid can be dangerous for the user and for the VISCOSYSTEM AVS 500!
- c) Turn on the upper cog belt wheel of the spindle to set the piston in the top position such that the piston plug can slide out of the guide groove (Fig. 7).
Caution: Do not turn to the stop position.
- d) Press down the retaining spring and at the same time pull away the dosing unit to the front (Fig. 7).

Filter change and needle replacement

- a) The filter can be changed as follows after taking off the front sheet metal panel of the sample transfer unit (Fig. 5):
- b) The measuring VISCOSYSTEM AVS 500 must be in the switched off state. Thus the needle is in its top stop position (Fig. 10).
- c) The screw fitting of the needle can be released with a SW 13 wrench for replacing the needle (Fig. 10).
When replacing the needle, a glass sample needle can also be fitted (type no. VZ 5360).
- d) The filter sieve can be replaced after opening the milled nut of the needle adapter (Fig. 10).
- e) When mounting and assembling, make sure that the gaskets used are positioned correctly.

6 Maintenance and care of the VISOsystem AVS 500

To retain the function of the VISOsystem AVS 500 the following maintenance and inspection work is to be carried out.

Maintenance intervals

Normal operation:

- As a rule, the max. intervals for carrying out all maintenance work is 6 months.

Under particular stress:

- Once per quarter, the glass cylinder of the dosing module VZ 6239 has to be inspected visually in order to check for damages, and a test according to DIN 12 650, part 6, or ISO/DIS 8655, part 3 has to be carried out. These verifications are particularly necessary if solutions have been used which are likely to attack glass. e.g. alkaline, fluoride or phosphatic solutions.
- If there is suspicion that a solution is attacking the glass excessively, the maintenance intervals are to be reduced accordingly.
- When working with aggressive solvents such as sulfuric acid, dichloroacetic acid or formic acid, one has to ensure in particular that the permanently adhering wall film inside the dosing module does not form an excessive accumulation of liquid drops. In such cases the maintenance and cleaning intervals have to be shortened correspondingly.
- Once per quarter, the electrical contacts have to be checked, if the VISOsystem AVS 500 is used in premises with an occasional occurrence of corrosive matters in their atmosphere.

In case of disturbances:

- The work has to be carried out immediately if a disturbance, an error or another defect becomes obvious.

Maintenance work to be carried out

- Check the piston and cylinder of the dosing module for visible signs of mechanical damage, such as cracks, fractures and damage of the surface.
- Check the piston- and driving needle for corrosion.
- Check the hoses, screw unions and seals for signs of visible damage, contamination and leaks.
- Check the electrical plug contacts for corrosion and mechanical damage (on dosing valve and cables).
- Defect parts must be repaired or replaced with new ones. Defect glass parts must always be replaced.
- After each maintenance operation, a test for metrological reliability according to DIN 12 650, part 6 or ISO/DIS 9001, part 3 has to be carried out.

Periods without operation

- If the transfer unit of the VISOsystem AVS 500 is not used over a long period of time, the liquids contained in the system, in particular aggressive solutions, have to be drained. If the liquid is left in the system, one has to reckon with corrosion and an alteration of solutions used over the time. Since the state of the art does not provide plastic hoses which are absolutely free of diffusion occurrences, this precaution measures apply in particular to the hose-line section.

Cleaning

- The interior of the instrument may only be cleaned if it is **completely disconnected** from the electric supply. The mechanic parts can be cleaned with a rag or usual water or household cleaners that contain ethanol. A rinsing with water, or ethanol should basically be avoided in the transfer area, because humidity rests can lead to electric short circuits and therefore to disturbances.
- The housing of the ViscoSystem 500 may also be cleaned on the outside with a rag or household cleaner. The rear side has to be cleaned using dry rags. It has to be **avoided** that water gets inside the instrument.

7 Index

Absolute viscosity	56	Enter	57
Accuracy	47	Error messages	83
Advancing of printer paper by one line	57	Evaluation procedure	57
Alarm contacts	56	Evaluation, statistical	57
Ambient temperature	52	Exclamation mark	56
Arrow circle	57	Extension cable	48
ASCII code	66	Fill level sensor	49
Average of measured flow times	56	Filter case	60
AVS Memory	57	Filter sieve	84
Baud	66	Filter, changing	84
Blank panels	54	Flow time	56
Bottle adapter	49	Frame	54
C cal (°Celsius cal)	58	Function	57
°C Celsius	57	Functional units	56
Calibrating	57	Functional units, replacing	84
Capillary viscometry	47	Functions indicators	47
Changing system parameters	67	Fuses	48
Clear	57	Gaskets	84
Clear view transparent thermostat	47	German Federal Post Office	49
Commissioning	54	Glass cylinder	60
Constants	56	Glass needle	60
Control characters		Ground conductor	48
CR (Carriage Return)	66	Grounding contact	48
LF (Line Feed)	66	Hagenbach correction	47; 57; 58
Damage, visible	49	Holding spring	54
Damaging	48	Hose assembly	54
Data cables	55	Inherent viscosity	58
Data memory	65	Inserting the paper	65
Data output	47	K_t/to, constant	59
Data transmission	47; 66	Kalrez gasket	60
Date	57	Keyboard	47
Declaration by the manufacturer	49	Leaks	54
Delta %	59	LED displays	56
Density of liquids	49	Limiting viscosity number	57
Determining the viscosity	71	Liquid level sensor	49
Dielectric constant	49	Liquids	
Dilution viscometry	56; 57; 79	newtonian	47
Dimension indicator		volatile	61
measured value	56	Magnetic rod	58
parameters	56	Mains cable	49
DIN 51 562	58	Measuring	68
Display		Measuring amplifier cards	47; 54; 61
Solvent volume	59	Measuring sensors	55
Temperature equilibration	59	Measuring sequence	57
Disturbance	55	Meniscus position	48
Dosing rate	56	ml Display	59
Dosing unit, replacing	84		
Drive motor	60		

n , entering of sample numbers	58	Sensor	
Needle		capacitive	49
glass needle	60	Series of measurement	65
replacement	84	Setting up	54
stainless steel needle	60	Solvent volume	
Needle protection	49	maximum	59
Needle protector	60	Stainless steel frame	55
New calibration	57	Stainless steel needle	60
Nominal value temperature	57	Start of the measuring sequence	57
Number of measurements	57	Start position	58
Numerical keyboard	57	Status displays, logic	56
Solvent quantity	57	Stop	57
Operating instructions	48	Stop bits	66
Optical sample detector	61	Stop/Reset	57
Parity bit	66	Stopwatch	47
Particles in the solution	60	Suck off time	56
Piston burette	47; 55	Symbol	
Plain paper	65	Arrow circle	57
Plug connections	54	Clock	59
Pneumatic modules	47; 49; 54; 55; 64	Printer	57
Pos., position 1st sample	58	Rinsing	59
Print out	65	Temperature Equilibration	59
Printer	57; 65	Table of chemical resistant value	60
Pt 1000 Temperature Sensor	55	TC Viscometers	47; 48; 55
Read out of the memory contents	66	Temperature	
Regulations for prevention of accidents	49	calibration	57
Respirators	49	nominal	57
Rinsing	57	tolerance	57
Rod magnet	61	Temperature deviation	55
Roll of paper	65	Temperature equilibration time	56; 57
Routine checks	82	Temperature Sensor Pt 1000	47
RS-232-C interface	55; 66	Temperature tolerance	57
RUN	57	Terminating the measurement	82
Runout time	57	Thermistors	48
Runtime t₀	59	Thermostat	55
Safety clothing	49	Thermostat head	55
Safety devices	47; 49	Thermostatic bath	47
Safety goggles	49	Time measurement	47
Safety instructions	48	Time of day	57
Safety set	55	Tolerance specification	57
Safety devices	49	Transferring results	66
Sample bottles	61	Transmission of memory contents	66
Sample glasses	47; 61	Transmission rate	66
Sample number	57	Transport securing devices	54
Sample transfer unit	54; 55; 60	Transportation	49; 54
Sample transfer unit, changing	84	Transportation instructions	54
Sample turntable	61	Trouble shooting	82
Sensing		Turning of the turntable	57
conductivity	48	Turntable position	57
opto-electronic	47; 48; 54	Type HC, Hagenbach correction	58
		Ubbelohde	
		ASTM	58
		DIN	58
		Unpacking	54

V1, Viscometer 1	59	Viscometer, too full	55
V2, Viscometer 2	59	Viscosity	
Viscometer		absolute	56, 57
Cannon Fenske	58	inherent	47; 58; 75
Dilution Viscometer	58	kinematic	47; 56
Micro-Ostwald	58	relative	56, 57
Micro-Ubbelohde	58	Warning buzzer	56
Ubbelohde ASTM	58	Warnings instructions	48
Ubbelohde DIN	58	Waste bottle, full	49, 55
Viscometer constants	57	Weight parameter	49
Viscometer symbols	56	Wrong entry	67
Viscometer type	54; 57		

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SCHOTT-GERÄTE GmbH