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OPERATING INSTRUCTIONS

Digital pH Meter with mV-range, manual/automatic temperature-compensation (Pt 100), Dead-Stop-Connection (10 μ A) and Chart Recorder Connection (100 mV/1000 mV).

Commencing Operation

1. Connect electrode* and mains cable (see illustration). Mains voltage 220 V $\pm$ 15 %, 50 ... 60 cps.
2. Switch on unit (slider switch on rear of unit).

Calibration of pH Meter to the function of electrode

3. Set "FUNKTION" switch to position "pH man".
4. Select two buffer solutions: one with a pH value in the vicinity of the electrode zero point, e. g. pH = 6.87 (according to DIN 19 266 and NBS) or pH = 7.00, and the other with a pH value differing therefrom by about 3 pH units, e. g. pH = 4.01 (according to DIN 19 266 and NBS) or pH = 3.00.
5. Equilibrate temperature of buffer solutions and electrode at 25 °C for instance (with a tolerance of $\pm$ 0.2 °C if a reproducibility of approx. $\pm$ 0.01 pH units is desired). In case a reproducibility of $> \pm$ 0.1 pH units is required, a special temperature equilibration can usually be dispensed with.
6. Set "TEMPERATUR" knob to buffer solution temperature.
7. Remove cap from refilling port of electrode. Rinse electrode with distilled water. (see Operating Instructions for Electrodes).
8. Immerse electrode in buffer solution with pH value in vicinity of zero point (e. g. pH = 6.87 or pH = 7.00). Set digital display to value of buffer solution by turning "NULL-PUNKT" (zero point) knob.
9. Rinse electrode with distilled water before immersing electrode in second buffer solution (e. g. pH = 4.01 or pH = 4.00).
10. Set digital display to pH value of second buffer solution by turning "STEILHEIT" (slope) knob. Unit is thus adapted ("calibrated") to the electrode characteristics. Rinse electrode.

pH Measurement with manual temperature adjustment

11. Keep switch "FUNKTION" in position "pH man".
12. Set "TEMPERATUR" knob to the temperature of the solution to be measured.
13. For measuring the pH value, immerse electrode in solution to be measured. After adequate tempering time, the pH value is read from the display.

* "Electrode" in this context stands for pH combination electrode or for electrode combination comprising separate indicator electrode and reference electrode.

pH Measurement with automatic temperature compensation (s. remarque)

14. Set "FUNKTION" knob in position "pH auto".
15. Immerse reference thermometer Pt 100 according to DIN 43 760 into the solution to be measured. Connect the thermometer with pH meter (Sockets on the rear panel of the pH meter, see illustration 2). (The reference thermometer should immerse near by the electrode).
16. For measuring the pH value, immerse electrode in solution to be measured. After adequate tempering time, the pH value is read from the display.
Remarque: we recommend resistance thermometer type W 103 with cable-combination type KR 1 NN.

Note: The method described is the one practiced most frequently. It is realized that reproducible measuring values are always obtained, which need not necessarily be identical with the absolute pH value. In the case of higher requirements in regard to pH measurements in solutions with higher temperatures, it is necessary, in order to adjust the pH meter to the electrode (calibration), to temper the buffer solutions also to the measuring temperature and observe the pH values in accordance with DIN 19 266 applicable to this temperature.

mV Measurements

17. Commence operation of unit as described under paras 1 and 2.
18. Set "TEMPERATUR" knob in position "mV".
19. Immerse electrode combination in solution to be measured.
20. The measuring values are read from the digital display directly.

Instructions for connecting the Adapter Cable for DEAD STOP TITRATION

1. Insert the 2 mm plug into the 2 mm socket on the rear panel of the pH meter (Fig. 2, position 3).
2. Insert the 4 mm plug into the center contact jack of the standard socket according to DIN 19 262 (Fig. 2, position 1).
3. Of the two 4 mm banana plugs of the double platinum electrode,
3.1 insert one onto the 4 mm socket for the reference electrode (Fig. 2, position 2) and
3.2 insert the other one into the socket on the rear of the banana plug which is already plugged into the DIN 19 262 standard socket (the center contact jack of this socket is therewith occupied by 2 plugs, namely the adapter cable plug and one of the two electrode plugs).
4. The switch positions of the pH-meter is: Function to "mV".

The pH meter is therewith ready for making DEAD STOP TITRATIONS.

We recommend: double-platinum electrode type Pt 1200 or Pt 1400 (with ground joint socket 14) and cable-combination type L 1 NN.

Notice

For reasons of function and safety, the unit shall not, as a matter of principle, be opened.

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SCHOTT-GERÄTE GmbH