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**Mikroprozessor-
pH-Meter
CG 836 und CG 837**

Operating Instructions

**Microprocessor
pH Meters
CG 836 and CG 837**

Mode d'emploi

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
SCHOTT

Abb. 1: Vorderansicht

Fig. 1: Front view

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1 Anzeige
Display
Affichage

2 Tastenfeld
Keyboard
Touches

3 Elektrodenanschluß
Electrode Connector
Electrode de mesure ou combinée

4 Bezugselektrode und zweite Buchse für 5 und/oder 6
Reference electrode and second socket for 5 and/or 6
Electrode de référence et 2^{me} fiche pour 5 et/ou 6

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Operating Instructions for Microprocessor pH Meter CG 836 and CG 837

The microprocessor pH meters CG 836 and CG 837 can be used to make pH-value, redox voltage and temperature measurements. The measured temperature is shown in the display parallel to the pH or mV value.

The calibration data for pH measurements are preserved when the instrument is switched off.

For adjustment to match the electrode function, the data for three buffer solutions, including their temperature characteristics, are stored in the pH-meters CG 836 and CG 837.

The instrument can be operated as hand-held device with closed or removed folding cover, as table unit with folded-over cover or as pendant unit with carrying strap.

Commissioning

1. The instrument is delivered ready for operation. A battery life of about 6000 operating hours can be expected using alkaline batteries. A warning appears in good time in the display when the batteries are nearly exhausted.
2. Connect the electrode, e.g. a type **N 2042 A** combination electrode with integrated temperature sensor (resistance thermometer Pt 1000, see Fig. 1, Pos. 3 and 5).
3. Switch-on the instrument: Holds the "on/off"-key, the "pH"-key or the "mV"-key depressed. A visual display device test will run. The unit is ready for measurement thereafter.

Calibration of the pH-meter

4. The simplest procedure is to use a combination pH-electrode with integrated temperature sensor, e.g. type **N 2042 A**. The buffer solutions whose data are stored in the instrument have the following pH-values at 25 °C: CG 836: pH = 4.01, pH = 6.87 and pH = 9.18; CG 837: pH = 4.00, pH = 7.00, and pH = 10.00. The call-up sequence of the buffer solutions is arbitrary.
5. Dip the electrode into the first buffer solution used and press the corresponding "cal"-key. The messages "CAL", "WAIT", "P 1" and the corresponding pH-value appear in the display. "WAIT" flashes when the buffer solution has been identified and disappears when the reading has become stable. The message "P 2" flashes as request to continue the calibration procedure.
6. Rinse the electrode with distilled water and then repeat the procedure with the second buffer solution. The corresponding indication is "P2". Actual temperature differences of the buffer solutions are automatically normalized by the CG 836 and CG 837 to 25 °C.
7. When the pH-meter has automatically taken-over the calibration, the messages "CAL" and "OK" appear in the display, and then the electrode zero-point pH-value and the slope in %. The pH-meter with electrode is now ready for making measurements. To display the slope and zero-point again, press one of the "cal"-keys and thereafter the "pH"-key.

Single-point calibration

8. A simplified calibration can be made as described under item 5 with only one buffer solution. When the flashing message "P2" appears in the display, terminate the calibration procedure by pressing the "pH"-key. This single-point calibration corrects the electrode zero-point and keeps the slope value from the previous full calibration.

pH-measurement with combination pH-electrode N 2042 A (with integrated temperature sensor Pt 1000)

9. Press the "pH"-key to switch to the "pH"-function. The function displays "pH" and "°C auto" light up (see Fig. 1, Pos. 1).
10. Dip the electrode into the sample solution. The temperature-compensated pH-value reading appears in the display after elapse of an adequate temperature equilibration time.

pH-measurement with combination pH-electrode and separate temperature sensor Pt 1000

11. pH-measurement with automatic temperature-compensation is possible too using a combination pH-electrode and a separate Pt 1000 temperature sensor. In this case, connect one plug of the temperature sensor to the socket marked Pt 1000 and connect the other plug to the reference electrode socket (see Fig. 1, Pos. 4 and 5).

pH-measurement with manual temperature compensation

12. When using an electrode without integrated temperature sensor (e.g. combination pH-electrode **N 62**), the manually set temperature is shown in the display together with the message "°C man".
13. This temperature value can be adjusted with the keys "Δ" and "∇" to correspond to the actual measuring temperature. Manual temperature adjustment is not possible in the "mV"-function and no temperature is indicated in the display.

mV-measurement

14. Press the "mV"-key to switch the instrument to the "mV"-function. The "mV"-function display lights up.
15. Dip the electrode combination (e.g. the combination Pt-electrode, type **Pt 62**, for redox voltage measurements) into the sample solution. The digital display shows the resulting reading directly in mV.

Temperature measurement

16. Dip the Pt 1000 temperature sensor or the combination electrode **N 2042 A** (with integrated temperature sensor) into the sample medium.
17. The display shows the temperature reading. The simultaneously displayed pH- or mV-value is of no significance if only the temperature sensor is connected.

Chart recorder output

18. The instrument has a chart recorder output (Fig. 1, Pos. 4 and 6) with source resistance $R_i = 2 \text{ k}\Omega$ for connecting a servo chart recorder with floating input (pH: 118 mV/Δ pH at 25 °C, 100% slope); U: 1 mV/mV).

Hold function

19. Press the "hold"-key to hold the readings. The message "HOLD" appears in the display. Press the "on/off"-key to switch the pH-meter back to the normal measuring mode.

Automatic switch-off

20. The pH-meter switches itself off automatically 60 minutes after the last keypress. To override this automatic switch-off function, switch-on the instrument with the "on/off"-key whilst simultaneously holding the keys "cal 6.87" (for the CG 836) or "cal 7.00" (for the CG 837) and "hold" depressed. Only the message "HOLD" appears in the display as acknowledge. The instrument is ready for making measurements after pressing the "on/off"-key again. The automatic switch-off function is always disabled in "hold" mode of the instrument.

Reset function

21. If the pH-meter gives incorrect readings or displays for no apparent reason, the instrument must be switched-off. Then press the "hold"-key and the "on/off"-key at the same time. Only the message "HOLD" appears in the display. The instrument is ready for making measurements after pressing the "on/off"-key again. After switching-on in this way, the instrument is set to 98% slope and zero-point pH = 7.00. A new calibration with the electrode is necessary if the actual electrode characteristic differs from these values.

Error messages

22. Err1: Overshoot of the measuring range. This error message disappears again automatically when the reading returns to within the measuring range.
23. Err2: The response time of the electrode is too long (greater than 1 minute). This error message appears during the calibration procedure. Usually the cause is that the electrode has not been watered properly or is too old. "Err2" also appears when no electrode is connected.
24. Err3: A buffer solution which does not correspond to the pressed "cal"-key has been used. To cure this error, press the correct key.
25. Err4: The temperature during calibration lies below 0 °C or above + 100 °C. Correct the temperature to lie within this range.
26. Err%: The slope of the electrode is less than 85% or more than 105% of the theoretical value.
27. ErrΔ pH: The zero-point of the electrode lies below pH = 5.5 or above pH = 8.0.

Attention: Error messages Err1 and Err3 may also be caused by:

- a) drifting of the electrode tension too strong (ageing of electrode)
- b) exchange of an electrode for a new one with different slope and different zeropoint
- c) wrong sequence of single operating steps

In all these cases the **RESET FUNCTION** (see point 21) helps to make the micro-processor pH meters CG 836 and CG 837 function again.

Battery replacement

28. To replace the batteries, remove the screws on the rear side of the instrument. The instrument can then be opened without difficulty to replace the batteries. This procedure is necessary only at intervals of about two years, because of the very long battery life. The calibration data are preserved for about 2 minutes without batteries. If the battery replacement procedure takes longer than this or if the old batteries were completely exhausted, a reset (see item 21) is necessary after fitting the new batteries. The used-up batteries have to be disposed off in accordance with environmental protection regulations.

Safety warnings and instructions

29. The instrument may be used only for the purposes described in these operating instructions.

Only leakproof alkaline batteries (4 × Mignon, IEC-LR 6) may be fitted. Always replace the complete set of batteries at one time. Do not mix different types or makes of batteries, or batteries which have been discharged to differing extents. Do not use rechargeable batteries, because deep discharge of accumulators may destroy the instrument.

Repairs may be carried out only by properly trained persons. The warranty becomes void in the case of unauthorised tampering inside or damage to the instrument.

The nominal pH-values of the buffer solutions used for calibration must correspond exactly to the pH-values marked on the "cal"-keys.

Please also observe the operating instructions for the electrodes.

pH calibration: Dip the electrode into the 1st buffer solution, press the corresponding **[cal]** key, wait until **P2** flashes. Rinse the electrode, dip it into the 2nd buffer solution, press the corresponding **[cal]** key and wait until **0.0** appears in the display *

Temp. compensation: With Pt 1000 automatic, manual with keys **[Δ]** / **[▽]** *

pH measurement: Press the **[pH]** or **[on/off]** key * Dip the electrode into the sample solution * **mV measurement:** Press the **[mV]** key *

To hold the measured value: Press the **[hold]** key *

To cancel hold: Press the **[on/off]** key * **Reset:** Switch-off unit, press **[hold]** key and **[on/off]** key, when "Hold" appears in display press **[on/off]** key again *

Err 1: Measuring range overshoot

Err 4: Cal temperatur

Err 2: Response time too long

Err %: Slope not in range 85% to 105%

Err 3: Wrong buffer solution

Err Δ pH: pH zero point not in range 5.5 to 8.0

This brief instruction is provided as self-adhesive label. It can be mounted over the German instruction on the underside of the pH-meter.

All specifications in this instruction manual are guidance values which are valid at the time of printing. However, SCHOTT-GERÄTE may make changes at any time for technical or commercial reasons or in the case of necessity to comply with the statutory stipulations of various countries.

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SCHOTT-GERÄTE GmbH