$\Delta T \leq 0,2 \text{ K } (-50 \dots + 200 \text{ }^\circ\text{C}) \pm 1 \text{ Digit}$
 $\leq 1 \text{ K } (-150 \dots - 50 \text{ }^\circ\text{C})$
 $\leq 2 \text{ K } (-200 \dots - 150 \text{ }^\circ\text{C})$

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Digital Pocket Conductivity Meter CG 858

Commencing Operation

1. Before commencing operation for the first time, open the battery compartment on the underside of the instrument and insert the provided battery (9 V, type IEC 6 F 22). The fixing screw can be released with a coin. Use only leakproof batteries as replacement.
2. Connect the conductivity measuring cell to the sockets provided on the face end of the instrument (4 mm sockets). Measuring cells with cell constant in the range $K = 0.75 \dots 1.25 \text{ cm}^{-1}$ can be used. Connect a temperature sensor (Pt 1000 resistance thermometer) for measurements with automatic temperature compensation or for making temperature measurements.
When using a combination cell for conductivity and temperature measurement, e. g. Type No. LF 3100 T, make connection optionally to the moisture-proof 4-pole SCHOTT-GERÄTE socket connector or to the 4 mm sockets, using 4 banana plugs in the latter case.
3. Move the "LF" switch from position "0" to position "man", "abs" or "auto", according to the desired manner of temperature compensation.
4. The instrument incorporates an automatic battery checking function. A condition of low battery voltage is indicated in the liquid crystal display (LCD).

Measuring the electric conductivity and the temperature

5. The Conductivity Meter CG 858 has three conductivity measuring ranges and one temperature measuring range. (The conductivity measuring ranges are relating to the measuring cell constant $K = 1,0 \text{ cm}^{-1}$).

Setting of the measuring	
range selector switch	measuring range
200 $\mu\text{S}/\text{cm}$	0.1 ... 199.9 $\mu\text{S} \cdot \text{cm}^{-1}$
2000 $\mu\text{S}/\text{cm}$	1 ... 1999 $\mu\text{S} \cdot \text{cm}^{-1}$
20 mS/cm	0.01 ... 19.99 $\text{mS} \cdot \text{cm}^{-1}$
$^{\circ}\text{C}$	-199.9 ... + 199.9 $^{\circ}\text{C}$

6. Set the manual control knob "K [cm^{-1}]" to the measuring cell constant (cell factor). The value of the cell constant is speci-

fied on the conductivity measuring cell. The Conductivity Meter CG 858 can be set continuously to match any cell constant in the range from 0.75 to 1.25 cm^{-1} .

If the measuring cell constant is unknown, dip the measuring cell with integrated temperature sensor (combination cell with Pt 1000 resistance thermometer) into 0.1 mol/l KCl solution and set the "LF" switch to position "auto". When the display reading has become steady, adjust the control knob "K [cm^{-1}]" such that the value 12.88 mS/cm is standing in the display. The instrument is therewith matched to the measuring cell. The actual cell constant value of the measuring cell is indicated on the scale of the control knob "K [cm^{-1}]".

7. Move the measuring range selector switch to the position corresponding to the range in which the conductivity reading is expected, e.g. "2000 $\mu\text{S}/\text{cm}$ ".
8. For manual temperature compensation, set the "LF" switch to position "man" and set the manual control knob "T-Comp" to the actual temperature of the measured solution.
9. For automatic temperature compensation, set the "LF" switch to position "auto" (with a Pt 1000 temperature sensor connected, as separate sensor or in the combination cell). The manual control knob "T-Comp" has no function in this operating mode. The manual as well as the automatic temperature compensation function is matched to the temperature characteristic of natural water.
10. To measure absolute conductivity values, set the "LF" switch to position "abs". In this operating mode, the connected temperature sensor and the manual temperature compensation control knob both have no function. But the measuring cell constant must be set (see instruction step 6).
11. Dip the conductivity measuring cell and the temperature sensor into the measuring solution and then take the conductivity reading on the liquid crystal display (LCD). Always select the optimum measuring range. The temperature sensor is not required (need not be connected) in the switch setting "man".

Measuring ordinary resistors

12. Connect the resistor which is to be tested to the 4 mm sockets of the instrument which are marked "cell".
13. Set the measuring cell constant on the instrument to " $K = 1,0 \text{ cm}^{-1}$ " and move the "LF" switch to position "abs".

14. After selecting the optimum measuring range, take the conductivity reading on the liquid crystal display (LCD).

With the manual control settings "abs" and constant "K = 1,0 cm⁻¹", the numerical reading in the display directly corresponds to the conductance value of the resistor on test. The resistance value is then obtained according to the relationship:

$$\text{Resistance} = \frac{1}{\text{Conductance}} \quad R [\Omega] = \frac{1}{G \left[\frac{1}{S} \right]}$$

Measuring the temperature

15. Connect the Pt 1000 temperature sensor (as separate sensor or included in the combination cell).
16. Move the measuring range selector switch to position "°C".
17. Dip the temperature sensor into the measuring solution and then read the temperature value on the liquid crystal display (LCD).

Tolerance of the temperature measurement:

$$\begin{aligned} \Delta T &\leq 0.2 \text{ K } (-50 \text{ to } +200 \text{ °C}) \pm 1 \text{ digit} \\ &\leq 1 \text{ K } (-150 \text{ to } -50 \text{ °C}) \\ &\leq 2 \text{ K } (-200 \text{ to } -150 \text{ °C}) \end{aligned}$$

(in addition to the tolerance of the Pt 1000 resistance thermometer).

Important notes and explanatory information with regard to making electrical conductivity measurements

18. The entire metal surfaces of the measuring cell must be immersed in the sample solution. Thus make sure that the insertion depth is adequate. Air bubbles will lead to falsified readings, particularly when making continuous flow measurements.
19. In all cases when making measurements with dip-in measuring cells, wait until the reading is constant (about 3 seconds) before taking the reading.

20. Electrodes coated with platinum black should never be cleaned by mechanical means. If simple rinsing with water is inadequate, clean the electrode according to the kind of contamination, by dipping it into 1% hydrochloric acid, 1% sodium hydroxide solution or a suitable organic solvent. For plastic measuring cells with platinum black coated nickel electrodes, use acid cleaning very briefly or not at all, on account of the danger of dissolving the nickel metal.

21. Utilized abbreviations and formulae in conductometry:

x (kappa):	electric conductivity of an electrolyte solution specified in [S·m ⁻¹] or [μS·cm ⁻¹]
K:	measuring cell constant (cell factor) specified in [cm ⁻¹]
G:	conductivity specified in $S = \frac{1}{R} \left[\frac{1}{\Omega} \right]$

Please note

22. In order to avoid malfunctioning, this instrument should be opened only by competent trained persons.
23. Warranty becomes void in the case of unauthorized tampering inside this instrument or in the case of careless or deliberate damage.

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© 1985