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ORIGINALVERSION

OPERATING MANUAL

MODE D'EMPLOI

MANUAL DE INSTRUCCIONES

HandyLab 700

SI Analytics

a xylem brand

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Operating Instructions Page 23 ... 37

Important notes: Before initial operation of the HandyLab 700, please read and observe carefully the operating instructions. For safety reasons the HandyLab 700 may only be used for the purposes described in these present operating instructions. Please also observe the operating instructions for the units to be connected. All specifications in this instruction manual are guidance values which are valid at the time of printing. However, for technical or commercial reasons or in the necessity to comply with the statutory stipulations of various countries, SI Analytics may perform additions to the HandyLab 700 without changing the described properties.

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HandyLab 700 is a portable pH meter.

Guarantee

We provide guarantee for the device described for two years from the date of purchase. This guarantee covers manufacturing faults being discovered within the mentioned period of three years. Claim under guarantee covers only the restoration of functionality, not any further claim for damages or financial loss. Improper handling/use or illegitimate opening of the device results in loss of the guarantee rights.

To ascertain the guarantee liability, please return the instrument and proof of purchase together with the date of purchase freight paid or prepaid.

Registered trademarks

The following names are registered trademarks. For practical reasons they are shown without trademark symbol in this manual.

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Disposal

Please observe the applicable local or national regulations concerning the disposal of "waste electrical and electronic equipment".

Caution!

Do not operate the device when one of the following conditions applies:

- the device shows visible damage
- the device fails to perform the intended function
- prolonged storage at temperatures above 70 °C
- severe transport stresses

In this case, a professional routine test must be performed. This test should be carried out by the manufacturer.

Inserting the Batteries

With four AA batteries, the HandyLab has an operating time of over 1000 h. Open the battery compartment on the rear of the device. Be sure to observe the correct polarity when inserting the batteries (see markings in the battery chamber). Close the battery compartment cover and screw it handtight.

A battery icon in the display indicates the battery power level:

	Icon fully filled	Batteries at full capacity
	Icon partially filled	Battery capacity is sufficient
	Icon empty	Battery capacity not sufficient; calibration is possible
	Icon blinks	Max. 10 operating hours remaining, measurement is still possible Caution! It is absolutely necessary to replace the batteries.

Connecting a Sensor

The HandyLab 700 provides several connections so that many types of sensors can be used for measurement.

Note that only **one** sensor may be connected to the meter at a time.

The meter automatically recognizes a connected Memosens sensor and switches accordingly. Memosens is signaled in the display.


Caution!

Always make sure that a sensor is connected to the meter before starting measurement.

Explanation: The analog pH input of the HandyLab is an electro-meter amplifier with an extremely high-impedance. When the sensor is not in contact with the medium or not connected to the meter, electric charges on the input can generate arbitrary, stable pH or mV values which will be shown in the display.

Connections

- a M8, 4 pins for Memosens sensors
- b Temperature probe GND
- c Temperature probe
- d pH socket (DIN 19 262)

Memosens sensors have a **cable coupling**, which allows convenient replacement of sensors while the cable remains connected to the meter.

The connecting cable is connected to socket **a** (M8, 4 pins for Memosens sensors).

Switching On the Meter

When you have connected the sensor, you can switch the meter on by pressing the **on/off** or **meas** key.


When the meter is switched on with the **on/off** key, first a self test is performed and then the calibration data and settings are displayed before the meter switches to measuring mode.


When the meter is switched on with the **meas** key, it immediately switches to measuring mode.

Depending on the connected sensor and the specific measuring task, several steps for configuration and calibration must be performed as described on the following pages.

Icons

Important information about the state of the device:

Keypad

The keys of the membrane keypad have a noticeable pressure point.

They have the following functions:

on/off	Switches the meter on and displays the device and calibration data (see Start-Up)
meas	Switches the meter on / Activates measuring mode
cal	Starts calibration
set	Activates configuration/ Confirms entries
clock	Displays time and date, allows setting the clock using set
▲ ▼	When this icon is displayed, you can use the arrow keys for navigation.

Configuring

pH Configuration

Prior to measurement, a configuration should be performed to match the connected sensor and the desired measurement performance. Furthermore, you can select the suitable calibration method.

The following table gives you an overview. Factory settings are shown in **bold print**.

Measurement

↓ **set**

“Setup” display

▲	Display 1
▲	Display 2
▲	CAL Timer
▲	CAL
▲	CAL POINTS
▲	BUFFER SET (AutoCal, FREE CAL)
▼	Auto OFF
▼	Temp Unit
▼	Time Format
▼	Date Format
▼	Default

set

Select using arrow keys, confirm by pressing **set**.

pH x.xx / pH x.xxx / mV / (°C for analog pH only)		
OFF / date + time / date / time		
OFF / 1 ... 99 days		
AutoCal /Manual/DATA INPUT/(ISFET-Zero)/FREE CAL		
1 / 2 / 3 / 1-2-3 (for AutoCal, Manual, FREE CAL)		
-01-		
(De- NIST standard	1.679/4.006/6.865/9.180	
fault)		
-02-	NIST technical	1.68/4.00/7.00/10.01/12.46
-03-	SI Analytics technical	2.00/4.00/7.00/10.00
-04-	Knick Calimat	2.00/4.00/7.00/9.00/12.00
-05-	Mettler Toledo	2.00/4.01/7.00/9.21
-06-	HACH	4.01/7.00/10.01/12.00
-07-	Ciba (94)	2.06/4.00/7.00/10.00
-08-	Hamilton	2.00/4.01/7.00/10.01/12.00
-09-	Reagecon	2.00/4.00/7.00/9.00/12.00
-10-	DIN 19267	1.09/4.65/6.79/9.23/12.75
OFF / 0.1h / 1h / 6h / 12h		
°C / °F		
24h / 12h		
dd.mm.yy / mm.dd.yy		
NO / YES (reset to factory settings)		


- ▲ This icon prompts you to select a menu item using the arrow keys -
▼ the selection is confirmed by pressing **set**.

AutoCal Calibration

(Calibration with automatic buffer recognition)

The calibration method is selected in the configuration menu. Calibration is required to adjust the sensor to the meter. It is indispensable for achieving comparable and reproducible measurement results.

Note: To abort calibration, you can press **meas** at any time. This will be confirmed by the display message "CAL ABORTED". Exception: When you have selected "CAL POINTS 1-2-3" and the first calibration step has been completed, the calibration process cannot be stopped any more.

DATA INPUT Calibration

(Calibration by entering known sensor values)

The calibration method is selected in the configuration menu.

Measurement

↓ **cal**

CAL
DATA INPUT

↓

ZERO POINT

Use ▲▼ to select the value for the zero point.

↓ **cal**

SLOPE

Use ▲▼ to select the value for the slope.

↓ **cal**

The calibration data will be displayed successively:

Date and time

ZERO POINT

SLOPE

Then the meter switches to measuring mode.

Note: To abort calibration, you can press **meas** at any time.

MANUAL Calibration

(Manual calibration)

The calibration method is selected in the configuration menu.

Measurement

↓ **cal**

CAL
MANUAL

The number of calibration points has been selected in the configuration menu.

↓

CAL 1/2/3
PRESS CAL

↓ **cal**

pH display blinks
PRESS CAL

Use ▲▼ to set the temperature-corrected pH value taken from the buffer table.

↓ **cal**

mV display blinks

↓

Depending on the number of calibration points, the procedure described above for CAL 1/2/3 is repeated.

mV value blinks until calibration is completed, then successive display of:

CAL DATA

1/2/3 CAL POINTS

ZERO POINT

SLOPE

Then the meter switches to measuring mode.

Note: To abort calibration, you can press **meas** at any time. This will be confirmed by the display message "CAL ABORTED". Exception: When you have selected "CAL POINTS 1-2-3" and the first calibration step has been completed, the calibration process cannot be stopped any more.

Calibrating

FREE CAL Calibration

(Free selection of calibration method)

FREE CAL calibration is selected in the configuration menu.

Measurement

↓ **cal**

CAL
AutoCal blinks


↓ **cal**

Use ▲▼ to select the required calibration method (AutoCal, DATA INPUT or MANUAL).


Perform the selected calibration (see AutoCal, DATA INPUT or MANUAL calibration).

Error Codes and Device Messages

Error messages are indicated as "ERROR ..." on the display. Information on the sensor condition is indicated by the "Sensoface" icon (friendly, neutral, sad) possibly accompanied by an info message ("INFO ...").

Example of an error message:
ERROR 8 (identical calibration media)

Example of a "Sensoface" message:
INFO 1 (cal timer expired)

Sensoface provides information on the sensor condition (maintenance request). Measurement can still be performed. After a calibration, the corresponding Sensoface icon (friendly, neutral, sad) is shown together with the calibration data. Otherwise, Sensoface is only visible in measuring mode.

The most important error messages and "Sensoface" info messages are shown on the inside of the protective cover. A complete list of messages and their meanings is provided in the following tables.

SI Analytics

Made in Germany

on/off

meas

cal

set

clock

clock

Power-on and display of device / calibration data

Power-on / Measuring mode / Logger stop

Start calibration

View / confirm device settings

Time / date

Set time / date

HandyLab 750/ 750 EX only

RCL

STO

STO

View stored values

Save measured values

Set/start logger

😊

😐

😞

Sensor okay

Calibrate the sensor soon

Calibrate or replace sensor

Error messages

1 Value out of range

2 ORP value out of range

3 Temperature error

4 Sensor zero

5 Sensor slope

8 Identical calibration media

10 Cal media interchanged

11 Sensor unstable

14 Enter time and date

21 Sensor error

22 Sensor conflict, disconnect one sensor

See manual for other error messages.

USB (opt.)

MEMO SENS

NTC 30KΩ
P11000 00KΩ
GND

Error Codes and Device Messages

"Sensoface" Messages

The "Sensoface" icon provides information on the sensor condition:

Sensoface Meaning

Sensor is okay

Calibrate the sensor soon

Calibrate or replace the sensor

The "neutral" and "sad" Sensoface icons are accompanied by an "INFO ..." message to give a hint to the cause of deterioration.

Sensoface	Message	Cause
	INFO 1	Calibration timer
	INFO 3	Sensocheck
	INFO 5	Zero / Slope
 	INFO 6	Response time
	INFO 7	ISFET: Operating point (asymmetry potential)
	INFO 8	ISFET: Leakage current
	INFO 9	ORP offset

Error Codes and Device Messages

Error Messages

The following error messages can be shown in the display.

Message	Cause	Remedy
 blinks	Battery empty	Replace batteries
ERROR 1	pH value out of range	Check whether the measurement conditions correspond to the adjusted measuring range.
ERROR 2	ORP value out of range	
ERROR 3	Temperature value out of range	
ERROR 4	Sensor zero point too high/low	Thoroughly rinse the sensor and recalibrate. If this does not help, replace the sensor.
ERROR 5	Sensor slope too high/low	
ERROR 8	Calibration error: Identical buffers	Use a buffer solution with a different nominal value before starting the next calibration step.
ERROR 9	Calibration error: Buffer unknown	Make sure that you use the same buffer set as configured.
ERROR 10	Cal media interchanged	Repeat calibration.
ERROR 11	Measured value unstable Drift too high	Leave the sensor in the liquid until the temperature is stable. If this does not help, replace the sensor.
ERROR 14	Time and date invalid	Set time and date
ERROR 18	Configuration invalid	Restart, reset to factory settings (Setup: DEFAULT YES), configure and calibrate. If this does not help, send in the device for repair.
ERROR 19	Factory settings error	Device defective, send it in.
ERROR 21	Sensor error (Memosens)	Connect operational Memosens sensor.

Specifications

Memosens pH input (ISFET)	M8 socket, 4 pins, for Memosens lab cable	
	Display ranges ¹⁾	-2,000 ... +16,000 pH -2000 ... +2000 mV -50 ... +250 °C
Memosens Redox	M8 socket, 4 pins, for Memosens lab cable	
	Display ranges ¹⁾	-2000 ... +2000 mV -50 ... +250 °C
pH/mV (analog)	Sensor standardization ¹⁾	ORP calibration (zero adjustment)
	Permissible calibration range	ΔmV (offset) -700 ... +700 mV
	Anschluss	pH socket, DIN 19 262 (13/4 mm)
	pH range	-2 ... 16
	Decimal places ^{*)}	2 or 3
	Input resistance	1 x 10 ¹² Ω (0 ... 35 °C)
	Input current	1 x 10 ⁻¹² A (bei RT, Verdopplung alle 10 K)
	Measuring cycle	approx 1s
	Measurement error ^{2,3,4)}	< 0.01 pH, TC < 0.001 pH/K
	mV range	-1300 ... +1300
Temperature	Measuring cycle	approx 1s
	Measurement error ^{2,3,4)}	< 0,1 % meas. val. + 0,3 mV, TC < 0,03 mV/K
	2 x 4 mm dia. for integrated or separate temperature detector	
	Measuring ranges	NTC 30 kΩ -20 ... +120 °C Pt 1000 -40 ... +250 °C
	Measuring cycle	approx 1s
pH-calibration	Measurement error ^{2,3,4)}	< 0.2 K (Tamb = 23 °C); TC < 25 ppm/K
Operating modes ^{*)}	AutoCal	Calibration with automatic buffer recognition
	MANUAL	Manual calibration with entry of individual buffer values
	DATA INPUT	Data entry of zero and slope
AutoCal buffer sets^{*)}	-01- (Default)	NIST standard 1.679/4.006/6.865/9.180
	-02-	NIST technical 1.68/4.00/7.00/10.01/12.46
	-03-	SI Analytics technical 2.00/4.00/7.00/10.00
	-04-	Knick Calimat 2.00/4.00/7.00/9.00/12.00
	-05-	Mettler Toledo 2.00/4.01/7.00/9.21
	-06-	HACH 4.01/7.00/10.01/12.00
	-07-	Ciba (94) 2.06/4.00/7.00/10.00
	-08-	Hamilton 2.00/4.01/7.00/10.01/12.00
	-09-	Reagecon 2.00/4.00/7.00/9.00/12.00
	-10-	DIN 19267 1.09/4.65/6.79/9.23/12.75
Permissible calibration range	Zero point	pH 6 ... 8
	With ISFET:	-750 ... +750 mV Operating point (asymmetry)
	Slope	approx. 74 ... 104 %
Calibration timer ^{*)}	Interval 1 ... 99 days, can be switched off	

Specifications

Display	Sensoface	Gives you indications about the status of the sensor, analysis of zero-point/slope, setting time, frequency of calibration. Status indication (friendly, neutral, sad)
	Display	LCD STN 7-segment display with 3 lines and icons
	Status indicators	Battery power level
	Notices	Hourglass
	Keypad	[on/off], [cal], [meas], [set], [▲], [▼], [clock]
Diagnostics functions	Sensor data (Memosens only)	Manufacturer, sensor type, serial number, operating time
	Calibration data	Calibration date, zero, slope
	Device self-test	Automatic memory test (FLASH, EEPROM, RAM)
	Device data	Device type, software version, hardware version
Data retention	Parameters, calibration data > 10 years	
Nominal operating conditions	Ambient temperature	-10 ... +55 °C
	Transport/ Storage temperature	-25 ... +70 °C
	Relative humidity	0 ... 95 %, short-term condensing allowed
Power	Power supply	Batterien 4x AA (Mignon) Alkaline oder Lithium
	Operating time	approx 1000 h (Alkaline)
Housing	Material	PA12 GF30 + TPE
	Protection	IP66/67 with pressure compensation
	Dimensions	approx (132 x 156 x 30) mm
	Weight	approx 500 g
Licensing - Compliance - Equipment safety		
EMC	EN 61326-1 (General Requirements)	
	Emitted interference	Class B (residential area)
	Immunity to interference	Industry
	EN 61326-2-3 (Particular Requirements for Transmitters)	
RoHS conformity	According to directive 2011/65/EC	

*) User-defined

1) Ranges depending on Memosens sensor

2) According to EN 60746-1, at nominal operating conditions

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© 2014 Xylem, Inc. Version 140602
⁽²⁾ is a registered word-/trademark of Endress + Hauser Conducta GmbH & Co. KG, Gerlingen and Knick Elektronische Messgeräte GmbH & Co. KG, Berlin.
SI Analytics is a trademark of Xylem Inc. or one of its subsidiaries.
SI Analytics GmbH