LF413T LF413T-ID

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

Copyright

© 2009, SI Analytics GmbH

Reprinting - even as excerpts - is only allowed with the explicit written authorization of SI Analytics GmbH, Mainz.

Printed in Germany.

LF413T - Contents

1 Overview 4

 1.1 Structure and function 4

 1.2 Recommended fields of application 4

2 Cleaning 5

3 What to do if... 5

4 Technical data 6

1 Overview


1.1 Structure and function

Structure

1	Voltage electrode (inside, 2x)
2	Current electrode (ring, 2x)
3	Temperature sensor in a graphite enclosure
4	Shaft
5	Connection head

Automatic sensor recognition with the LF413T-ID

In the plug connector of the LF413T-ID conductivity measuring cell, data is stored for the automatic sensor recognition function. Among other, this data includes the sensor type and series number. Besides, the calibration data is written in the sensor with every calibration procedure. When the sensor is connected, the data is called up by the meter and used for measurement and measured value documentation. The correct cell constant is always used automatically if the sensor is operated with several meters because the calibration data is stored in the sensor.

To be able to use the automatic sensor recognition function, a meter is required that supports the function. More information on the automatic sensor recognition function is given in the operating manual of the meter.

1.2 Recommended fields of application

- Applications in the laboratory (e. g. research, development, quality control)
- Control measurements in food industry
- On site measurements in rivers, lakes and wastewater
- Ground water measurements

2 Cleaning

Caution

Before cleaning, disconnect the sensor from the instrument.

Cleaning

A thorough cleaning is particularly recommended for measurements of low conductivities.

Contamination	Cleaning procedure
Lime sediments	Immerse in acetic acid for 5 minutes (volume share = 10 %)
Fat/oil	Clean with warm water that contains washing-up liquid

After cleaning, thoroughly rinse with deionized water and recalibrate if necessary.

Aging of the conductivity measuring cell

Normally, the conductivity measuring cell does not age. Special measuring media (e.g. strong acids and lyes, organic solvents) or too high temperatures shorten the operational lifetime considerably or damage the measuring cell. The warranty does not cover cases where such conditions cause failure or mechanical damage.

Disposal

We recommend to dispose of the conductivity cell as electronic waste.

3 What to do if...


Error symptom	Cause	Remedy
No temperature or conductivity display	– No connection between measuring instrument and sensor – Cable defective	– Check connection between measuring instrument and sensor
Measurement delivers implausible conductivity values	– Incorrect cell constant adjusted at the measuring instrument – Measuring range exceeded – Contamination in the area of the electrodes – Electrodes damaged	– Check / correct the cell constant – Observe the application range – Clean the sensor (see section 2). – Return the sensor
Incorrect temperature display	– The temperature sensor was not immersed deep enough in the measuring solution – Temperature sensor defective	– Observe the minimum immersion depth – Return the sensor

4 Technical data

General features

Measuring principle	Four-electrode measurement
Cell constant	$0.475 \text{ cm}^{-1} \pm 1.5 \%$
Temperature sensor	integrated NTC 30 (30 k Ω / 25 °C)

Dimensions (in mm)

Weight

approx. 135 g

Materials

Shaft	Epoxy
Connection head	POM
Conductivity electrodes	Graphite
Thermistor enclosure	Graphite

Connection cable


Length	3 or 1,5 m
Diameter	6 mm
Smallest allowed bend radius	Permanent bend: 80 mm Single time or short time bend: 50 mm
Plug type	Socket, 8 pins

Pressure resistance


Sensor with closed plug connection	IP 68 (2 x 10 ⁵ Pa or 2 bar)
Cable plug	IP 67 (when plugged in)

The LF413T meets the requirements according to article 3(3) of the 97/23/EC directive ("Pressure equipment directive").

Measurement conditions	Conductivity measuring range	1 μ S/cm ... 2 S/cm
	Temperature range	-5 ... 80 °C (100 °C)
	Max. allowed overpressure	2×10^5 Pa (2 bar)
	Minimum depth of immersion	36 mm
	Maximum depth of immersion	Entire sensor+cable (up to 80 °C) Sensor shaft only (=120 mm / up to 100 °C)
	Operating position	any
Storage conditions	Recommended storing method	in air
	Storage temperature	0 ... 50 °C
Characteristic data on delivery	Temperature responding behavior	t_{99} (99 % of the final value after) < 20 s
	Precision of the temperature sensor	± 0.2 K

Pin assignment

Plug from the front:

What can Xylem do for you?

We're 12,500 people unified in a common purpose: creating innovative solutions to meet our world's water needs. Developing new technologies that will improve the way water is used, conserved, and re-used in the future is central to our work. We move, treat, analyze, and return water to the environment, and we help people use water efficiently, in their homes, buildings, factories and farms. In more than 150 countries, we have strong, long-standing relationships with customers who know us for our powerful combination of leading product brands and applications expertise, backed by a legacy of innovation.

For more information on how Xylem can help you, go to www.xyleminc.com

SI Analytics

a **xylem** brand

SI Analytics GmbH

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Tel: +49 (0)6131 / 66 5111

Fax: +49 (0)6131 / 66 5001

D-55122 Mainz

Deutschland, Germany, Allemagne, Alemania

E-Mail: si-analytics@xyleminc.com

www.si-analytics.com

SI Analytics is a trademark of Xylem Inc. or one of its subsidiaries.

© 2015 Xylem, Inc.