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KF-Soft 1.0

Operating Instructions

KF-Soft 1.0

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Operating Instructions Page 17 31

Important notes: Before initial operation of the Software KF-Soft 1.0 please read and observe carefully the operating instructions and the enclosed licence conditions. For safety reasons the Software KF-Soft 1.0 may only be used for the purposes described in these present operating instructions.

Please also observe the operating instructions for the units to be connected.

All specifications in this instruction manual are guidance values which are valid at the time of printing. However, for technical or commercial reasons or in the necessity to comply with the statutory stipulations of various countries, SCHOTT Instruments may perform additions to the Software KF-Soft 1.0 without changing the described properties.

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KF-Soft 1.0

Contents

Contents	17
Requirements	18
Hardware	18
Software	18
Introduction.....	19
Properties	19
Installation	19
Settings of the TitroLine KF Titrator	19
The titration in progress	20
The sample list	23
The sample list	23
Selection of the sample list	25
The qualifications.....	27
Configuration	29
References	30

Requirements

Hardware

The KF-Soft 1.0 requires a 133 MHz Pentium processor with 32 MB RAM. A hard disk with a minimum of additional 5 MB RAM must also be available (proper running of WINDOWS requires additional free memory space). The graphic mode has to be 800 x 600 with 16 colours or better. Each connected TitroLine *KF* Titrator requires a free serial interface COM 1 to COM xx. This COM interface must be directly addressable by WINDOWS. The software **MUST** be run using a mouse. Operation without the use of a mouse is impossible, since none of the usual menus are used.

If printing is desired, a properly connected printer has to be present.

The installation of the software requires a CD-ROM drive.

Software

The software can be run under WINDOWS 95/98/Me/NT/2000. If used with WINDOWS NT, the min. configuration is 4.0 with Service Pack 5 installed.

The use of the qualifications requires WINWORD to be installed properly. The qualifications are based on DOT templates from WINWORD 97.

Titration curves are exported in the form of "csv" files. These can be read directly into EXCEL (all versions > 4.0) or any other spreadsheet program.

For the taking over of the data to the database of KF-Soft the country-specific settings **MUST** be implemented properly. This means for a **German WINDOWS** version that the decimal **comma** has to be implemented as the decimal sign instead of the decimal **point**.

A properly installed printer driver is required for output.

Running the software under WINDOWS NT requires access to the hard disk and the WINDOWS system directories. As a rule, this requires access as a system administrator.

Introduction

Properties

The KF-Soft software is designed for use with the TitroLine *KF* Titrator. It enables titration curves to be displayed during titration and titration results to be stored in a database in GLP standard. By way of various qualifications, it also allows the TitroLine *KF* Titrator to be monitored and escorted within the QA system. Specification of the software's capabilities:

Representation of the on-line titration curve:

- Representation of the on-line titration curve
- Representation of the current measurement values during titration
- Display of the result output of the TitroLine *KF* Titrator
- Input and monitoring of limit values
- Printout of the titration in progress incl. graph
- Export of the titration in progress in the form of a csv file (spreadsheet calculation)
- Storage of all finished titrations in a database
- Selection of data on the basis of titration types of sample titration, titer, blank value
- Selection based on user, sample description, date
- Monitoring of all titer determinations
- Export of selected or individual titrations with titration curve in form of a csv file
- Printout of individually selected titrations with graph
- Printout of a table of selected titrations
- Preparing qualifications: design, installation, operation and performance qualifications
- Monitoring of the qualifications in tables
- Write-protection of completed qualifications
- The status of the qualifications is documented in a table
- Specific configuration of each interface

Installation

The CD contains the German and English versions of the software:

- Setup_dt
- Setup_engl

Start the WINDOWS Explorer, then select the German installation. Installation is triggered by double clicking on the "Setup.exe" file. Installation will now proceed automatically after the entire input has been confirmed.

The preset installation directory is "c:\programs\T-Soft". It may be customized to suit your particular needs.

Once the installation completed, the program is started by clicking the "Start" button followed by "Programs".

You can also set up the software on the desktop itself. To do so, click on the right mouse key on the desktop. On the appearing menu, click on "new" "link" and "search". Select the "Programs\KF-Soft\KF_Soft.exe" file, then click on "Open". The system will automatically add an icon to the desktop which can be used to start the software.

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The serial parameter settings on the interface between the TitroLine *KF* Titrator and the PC have to be identical to those of the computer.

In each of the methods of the TitroLine *KF* Titrator the documentation is to be set to GLP (please refer to the operating instructions of the TitroLine *KF* Titrator). If the TitroLine *KF* Titrator is set to "Display only", no results or other data can be taken over to the software.

The titration in progress

At startup, the software will display the "Titration in progress" window. The following windows are offered on the file tabs:

- Titration in progress
- Sample list
- Qualifications
- Configuration


Figure 1: „Titration in progress“ window

The "Titration in progress" window shows the graph of titration in progress including the related measurement values. If the titration was completed successfully, the results of the GLP documentation will be displayed in a window (Figure 1).

Please make sure to select the correct titrator on top left (Figure 2). The titrators including all parameters (interface, baud rate, data bits, stop bits, handshake) are defined in the "Configuration" window. Pressing "Switch-over" will cause a switch-over to the current titrator.

Figure 2: Titrator selection

The software may be started with more than one titrator connected to different COM interfaces. All dates are written in a common database. One of these titrators is selected here as the active one.

The on-line graph shows the titration curve both during conditioning and titration. The x-axis will display the time in terms of seconds, whereas the y-axis will show the titration-reagent consumption in ml. A flat curve means a titration approaching its end. The graph will be scaled automatically.

The measurement values of the KF titration will be shown directly in the form of numeric values (Figure 3).

Current Measurements:

Reagent [ml] **0.000**

Indication current [µA]

Titration time [sec] **1**

Drift [µl/min]

Limit values for result

Lower Limit Upper limit

Figure 3: Display of the current measurement values

Consumption is displayed in terms of ml of dosed reagent, the indicator current μA , the titration duration in seconds. The indicator current will approach the set limit value in the course of titration, and towards the titration end it will attempt to keep this value as constant as possible. The drift will be calculated on a basis of approximately 10 and displayed as a numeric value. The value displayed on the TitroLine *KF* Titrator will not be the same, since the numeric value is referring to a more extensive time basis with a slower course of titration.

One upper and one lower limit value can be input for each titration. If the result in its calculated form, e.g. in terms of % of water, is within the defined limit, "OK" will be displayed for 'sample' in the result table. If the result is outside these limits, the numeric value will not be evaluated as "OK".

Each result of a successfully completed titration will be represented in a window. The display will show the exact replication of a printout of the TitroLine *KF* Titrator as it would be output on a printer. The TitroLine *KF* Titrator is to be set to GLP documentation (Figure 4).

Titration Results

Figure 4: Window for the display of the original GLP printout

The titration in progress may be printed and contains the graph in the form shown on the screen, the printout of the GLP documentation, and a box for entering additional information. It is also possible to export the titration in progress in the form of a csv file. In this process the various pieces of information, separated by a comma, are written into a text file (Figure 5). The data of the titration curve is written into a separate file in order to facilitate the creation of a graph using a spreadsheet program.


Figure 5: Printing and exporting the titration in progress

You exit the software by a mouse click on the cross on the top right of the window or on the “End” button.

The sample list

The sample list

The “sample list” window can be used to view all titrations or only selected ones. The highlighted titration will be shown together with its graph (Figure 6).

The figure of the TitroLine KF Titrator has no function but to remember the origin of the analysis results.


Figure 6: Sample list window

The sample list is used to store all titrations. In this process, the data being shown may be sorted on the basis of various criteria. The display will only show such titrations as were defined in the selection window. In the example of Figure 6, only the titration performed by user “JP” was selected.

The table containing the titrations may be customized to suit individual needs. To do so, place the mouse pointer directly on the border line of the column. The mouse will take a different shape. With the left mouse key depressed, you can set the desired width of the column. Figure 7 shows the former border line of the column next to the new, slightly one. The mouse pointer in its modified form cannot be depicted here.

User	Sample ID
▶ Adolfs	09
Adolfs	07

Figure 7: Setting the column with of the sample list

Figure 8 shows a customized sample list. In this example, the width of some of the columns was reduced so far that they are not visible any more.

	User	Sample ID	Consumption	Water	Mean	Date	Time	Remarks
▶	Adolfs	09	0.995 ml	10.080 mg/ml	10.080 mg/ml	21.12.00	16:06:37	
	Adolfs	07	1.347 ml	7.443 mg/ml	7.443 mg/ml	21.12.00	15:39:51	

Figure 8: Customized column width of a sample list

These column-width settings are stored in an 'ini' file. This means it will be replicated at the next program start. All the other functions such as printing and exporting always refer to the columns in the form they are being displayed. This means that any columns which were 'hidden' will not be shown, printed, or exported!

The default setting may be reset immediately by clicking on the "Restore columns" button (Figure 9).


Figure 9: Restoring the default column width

Figure 10 shows the initial sample list in the form it is implemented as a default in the software: The width of all columns is identical and in the form as transmitted from the titrator. The spaces:

- User
- Sample ID
- Sample OK
- Released

may be used for subsequent input.

	User	Sample ID	Consumption	Water	Mean	SD	RSD
▶	Adolfs	09	0.995 ml	10.080 mg/ml	10.080 mg/ml	0.00000 mg/ml	0.00 %
	Adolfs	07	1.347 ml	7.443 mg/ml	7.443 mg/ml	0.00000 mg/ml	0.00 %

Figure 10: Restored initial column width

Figure 11 shows the functions of the sample list

Figure 11: Selection of the functions of the sample list

The table below gives a description of the function of the various buttons (Table 1):

Button	Description
Restore columns	The customized column widths are reset to the initial column widths.
Export to csv file	The selected samples are exported to a csv file. The name of the file is "KFDoku.csv". A serial number is inserted in the first column of each titration (i.e. in each line). Each line contains the results of the GLP documentation. Below the number in the first column, one file named "curvxxxx.csv" will be created for each titration. In this file name, xxxx represents the serial numbers of the individual titrations. This enables a titration curve to be prepared simply and selectively. For each export operation, the files are newly created, with the existing files being overwritten. With a large number of files, the export operation may take quite a time!
Print table	The selected table with its visible columns will be printed with identical column widths and no further formatting. The printing format is DIN A 4 broadside.
Print selected sample	The printout corresponds to that of "Titration in progress". However, any eligible comment will also be printed.

Selection of the sample list

The titration for the sample list are selected using the 'sample list' window (Figure 12).

The screenshot shows a window titled "Selection of samples". It contains the following fields and controls:

- Method:** A text input field.
- User:** A text input field.
- Sample-ID:** A text input field.
- Release:** A dropdown menu.
- Limit:** A dropdown menu.
- Date:** A date range selector with "from" and "until" labels and two corresponding text input fields.

Figure 12: Selection for the sample list

As soon as one of the boxes contains an entry, a selection for this term or partial term will be carried out. Subsequently, only the information containing the term or partial term in the corresponding box will be displayed.

If, for instance, all titer determinations are to be displayed, it will suffice to enter "Titer" (without the quotation marks) in the "Method" box. All titer determinations will be displayed, irrespective of the standard used to obtain them.

If a user's name is "JP", his samples will be displayed upon entering "J". However, the samples of all other users whose name contains a "J" will be shown as well.

With regard to release, a difference can be made as to whether a sample was already released or not. Likewise, it is possible to specify whether the samples did comply with a certain limit value or not. If the limit value is specified as a criterion, only the samples which complied with the appropriate limit value will be shown.

The date in the from/to box has to be input explicitly in the following form: „day.month.year“!

The input of values of values selected on the sample list to the selection box using the copy/paste function is a particularly convenient feature. This feature will be explained for the "date" example: On the sample list, use the mouse to select the date the samples of which are to be shown on the list. The box containing the date will be darkened. If you press <CTRL> and <c> simultaneously, the date will be taken over to the clipboard of WINDOWS. Now click on the "from" box in the "Selection" window, then press <CTRL> and <v> simultaneously. The date will be copied into this box. To continue, click on the "to" box of the from/to selection box, then press <CTRL> together with <v> again. The date will be copied into this box, too. Both boxes, i.e. "from" and "to", have to contain a date!

The titrations or samples shown will be those meeting all the selected criteria.

PLEASE NOTE: If the sample list does not contain any titrations at all, this may indicate that you have activated certain criteria which are not met for any titration at all.

The qualifications

The qualifications are intended to ensure that the TitroLine *KF* Titrator is functioning properly. The application manual contains descriptions of the individual qualifications. The signification of the qualifications is as follows:

Abbr.	Description	Meaning
DQ	Design Qualification	The DQ specifies all requirements applicable to a Karl Fischer Titration Unit for one or more sample types. As a rule, one of the DQs is required at the time of new acquisition. A re-qualification may become necessary, for instance, if a new sample type or reagent are to be used with properties which differ from those which were defined with the initial DQ.
IQ	Installation Qualification	The IQ is intended to ensure that the titrator is set up and installed properly. The IQ includes both the external conditions and the correctness of the internal settings. An IQ of this type may quite probably be necessary more than once, e.g. at the time of first acquisition, re-organisation within the laboratory, or following service or repair work.
OQ	Operational Qualification	The OQ is supposed to ensure the operation of the titrator all possible laboratory conditions. This OQ should be performed at regular intervals. It will ensure the titrator's proper functioning.
PQ	Performance Qualification	The PQ is supposed to check the operation under routine conditions. This check is intended to prevent any gradual change of the results caused by a change in the analysis system.

To support a gap-less documentation of the titrator operation, it is possible to create, alter, or complete new qualifications on a permanent basis. They will be documented together with their respective status on a list.

The qualifications are created using the "Create qualifications" window (Figure 15).

In a first step, the type of qualification to be created is selected (Figure 13):

Figure 13: Selection of qualifications

In this process, WINWORD will automatically load a format template for the qualification to be created. Please make sure that WINWORD including its path is correctly configured in the "Configuration" window! The format template will be used for the automatic creation of a WINWORD document bearing the name of the qualification and the current date. This document can be used to make any required input. The document has to be stored in WINWORD, then exit WINWORD with "File" "End". Such a qualification will now have the status of "Complete" "No" (Figure 14). This means that the document may be further processed at any time by a double mouse click on the appropriate line of the list. Also, the document may be called up or modified within the operating system.

After a document has been finally completed, the status on the line will be changed to "Yes"! At the same time, the document will be write-protected on the file level so that inadvertent subsequent modifications are excluded.

If a qualification is to be further processed, it can be loaded automatically with a double click on the selected file.

The “Author” and “Complete” fields can be edited as long as a qualification had not yet been completed.

	ID	Type	Date	Author	Finished	Filename
▶	1	DQ	23.11.00 10:50:04	ad	Ja	dq_231100_100236.doc

Figure 14: List of qualifications


Figure 15: Window for creating and maintaining qualifications

You can exit the software like any other window or file card by a click on “End”. A confirmation prompt will follow, asking you whether you really wish to exit.

Configuration

The configuration of the serial interface and the path specification of WINWORD is done in the Configuration window (Figure 16).

Figure 16: Configuration of the serial interface and the WINWORD path specification in the Configuration Window.

Considering that within the software more than one titrator may be connected to the various serial interfaces, each of the titrators has to be assigned to an appropriate interface. This means that it is quite possible that titrator no. 2 is connected to the COM 1 serial interface. Each titrator should have a unique identification. The address of a titrator is an essential characteristic. However, since each titrator is connected to a different serial interface, identical addresses may exist. But to the user the titrators have to be uniquely identifiable. This means that specific numbers or descriptions are assigned to the titrators. The internal assignment is then done on the basis of the list in Figure 17 and the COM port.

All the titrators write their sample data to one and the same database.

In this process, all the information except for parity and handshake have to be entered manually (Figure 17).

RS232-Parameters							
	Device No.	COMPort	Baudrate	Parity	Stopbits	Databits	Handshake
▶	1	1	4800	N-No	2	7	0-No
	2	2	4800	N-No	2	7	0-No
*							

Figure 17: Input to be made for the serial parameters of the individual titrations

The input for parity and handshake have to be selected from a list box. Since the qualifications are using WINWORD format templates, WINWORD including its complete path has to be entered here (Figure 18).

 The image shows a configuration form with three main input areas. The first area has 'Parity' and 'Handshake' labels, each followed by a dropdown menu showing 'N-No' and '0-No' respectively. The second area has a 'Word-Path:' label followed by a text input field containing 'd:\programme\microsoft office\office'. To the right of the text field is a 'Save' button.

Figure 18: Input to be made for parity, handshake, and WINWORD path

You can use the WINDOWS Explorer to determine the path of WINWORD by having WINDOWS search for "winword.exe". Entering the entire "winword.exe" character string manually is not necessary.

References

A

Author28

C

COM interface18
 comma18
 configuration19, 29
 Configuration20
 confirmation prompt28
 Consumption21
 country-specific settings18
 course of titration21
 criteria23
 csv file19, 22
 customized sample list24

D

database29
 date19
 decimal sign18
 default setting24
 Design Qualification27
 desktop19
 different COM interfaces20
 dosed reagent21
 DQ27
 drift21
 edited28

E

End28
 export22
 Export19
 Export to csv file25
 exported18

F

format template27

G

gap-less documentation27
 GLP documentation20
 graph of titration in progress20
 graphic mode18

H

handshake29
 hard disk18
 Hardware18
 highlighted titration23

I

indicator current21
 individual needs23
 installation18, 19
 installation directory19
 Installation Qualification27
 interface29
 IQ27

L

limit value21, 26
 limit values19
 loaded automatically27

M

measurement values19, 21
 Monitoring19
 mouse18

N

new qualifications27

O

on-line graph20
 on-line titration curve19
 Operational Qualification27
 OQ27

P

parity29
 path27
 path specification29
 Performance Qualification27
 PQ27
 Print selected sample25
 Print table25
 printer21
 printing18
 printing and exporting24
 printout21
 Printout19
 Printout of a table19
 Properties19

Q

qualifications18, 27
 Qualifications20

R

RAM18
 release26

Released.....	24
Requirements	18
Restore columns.....	24, 25
result	21
result output.....	19

S

sample description.....	19
Sample ID	24
sample list.....	23
Sample list	20
Sample OK	24
selected titrations.....	23
selection.....	26
Selection	19
selection window	23
serial interface	29
spreadsheet program	18, 22
status	19
Storage	19
switch-over	20

T

terms of seconds.....	20
The sample list.....	23
time basis	21
titer determinations	26
titration curve	20, 22
Titration curves	18
Titration in progress	20
titration-reagent consumption	20

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SCHOTT Instruments GmbH