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
**Kolbenbürette
TITRONIC® 110 *plus***

Operating Instructions

**Piston Burette
TITRONIC® 110 *plus***

Mode d'emploi

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SCHOTT

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Operating Instructions Page 26 50

Important notes: Before initial operation of the Piston Burette TITRONIC® 110 *plus* please read and observe carefully the operating instructions. For safety reasons the Piston Burette TITRONIC® 110 *plus* may only be used for the purposes described in these present operating instructions.

Please also observe the operating instructions for the units to be connected.

All specifications in this instruction manual are guidance values which are valid at the time of printing. However, for technical or commercial reasons or in the necessity to comply with the statutory stipulations of various countries, SCHOTT may perform additions to the Piston Burette TITRONIC® 110 *plus* without changing the described properties.

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TABLE OF CONTENTS

PAGE

1	Technical properties of the TITRONIC® 110 <i>plus</i> Piston Burette	26
1.1	Summary, Declaration of Conformity	26
1.2	Warning and safety information	30
2	Set-up and first operation	31
2.1	Unpacking and setting up the TITRONIC T 110 <i>plus</i> Piston Burette	31
2.2	Installing the stand	31
2.3	TA <i>plus</i> Interchangeable Unit:	32
2.4	Fitting and replacing an interchangeable unit	33
2.5	Replacing the glass cylinder	34
2.6	Initial filling and rinsing of the complete interchangeable unit	35
2.7	Combining the TITRONIC® 110 <i>plus</i> Piston Burette with other devices	36
3	Display	36
3.1	LC Display	36
3.2	Dialogue example with the TITRONIC® 110 <i>plus</i> Piston Burette	37
3.3	Status messages of the TITRONIC® 110 <i>plus</i> Piston Burette	37
4	Working with the TITRONIC® 110 <i>plus</i> Piston Burette	37
4.1	Switching the piston burette on	37
4.2	Possible applications of the TITRONIC® 110 <i>plus</i> Piston Burette	37
4.3	Titration using the TitroLine <i>alpha plus</i> Titrator	37
4.4	Titration using the TR 160 Manual Titration Controller	38
4.5	Parameterising the TITRONIC® 110 <i>plus</i> Piston Burette using a keyboard	39
4.6	Dosing using the TITRONIC® 110 <i>plus</i> Piston Burette	41
4.7	Filling and rinsing the TITRONIC® 110 <i>plus</i> Piston Burette	41
4.8	System settings of the TITRONIC® 110 <i>plus</i> Piston Burette: Rinsing	42
4.9	System settings of the TITRONIC® 110 <i>plus</i> Piston Burette: Filling speed	42
4.10	System settings of the TITRONIC® 110 <i>plus</i> Piston Burette: Dosing speed	42
5	Connections and data transfer	43
5.1	RS-232-C interfaces, inputs and outputs	44
5.2	Parameterisation of the RS interface:	44
5.3	Command list for RS communication	46
5.4	Keyboard port	47
5.5	In/Outputs, input/output socket	47
6	Maintenance and care of the TITRONIC® 110 <i>plus</i> Piston Burette	47
7	Storage and transportation	48
8	Index	49

1 Technical properties of the TITRONIC® 110 *plus* Piston Burette

1.1 Summary, Declaration of Conformity

Using the versatile TITRONIC® 110 *plus* Piston Burette, virtually all titration and dosage tasks occurring in a laboratory can be completed. The substances or mixtures of substances to be dosed must be liquid.

The examples of possible uses for the piston burette include:

- Dosage and titration within the TitriSoft Computer Titration System
- Titrating in combination with automatic titration units from SCHOTT-GERÄTE such as the TitroLine *alpha plus*
- Titrating using the TR 160 Manual Titration Controller
- Dosing with a pre-selected volume, time, or speed

The following may be used:


- All common titration solutions, with the solutions specified in the "Küster-Thiel" calculation tables (published by Walter de Gruyter) being taken as standard.
- Water and all non-aggressive, inorganic and organic fluids may be used as solvents. When handling inflammable substances the explosion protection directives of the Professional Association of the Chemical Industry are to be observed.
 - With liquids with a higher viscosity ($\geq 5 \text{ mm}^2/\text{s}$), a low boiling point, or a tendency to exhalation, the filling speed can be adapted.
 - Liquids with a viscosity above $50 \text{ mm}^2/\text{s}$ cannot be dosed.

The two RC-232-C-interfaces allow the chained arrangement of TITRONIC® 110 *plus* Piston Burettes among each other as well as the connection of other devices from SCHOTT-GERÄTE which suitable for chaining, e.g. the TitroLine *alpha plus* Titrator and the TitriSoft software in combination with a PC. A standardised connector enables a standard PC keyboard (PS2) (fig. 15, item 6) to be connected for a convenient parameterisation of various functions.

A stand can be connected directly to the TITRONIC® 110 *plus* Piston Burette for the attachment of a stirrer or an electrode clamp (Fig. 2).

Interchangeable Units are available for volumes of 50 ml, 20 ml, 10 ml, 5 ml and 1 ml (Fig. 5). The centre-piece of the interchangeable units is formed by the high-precision borosilicate glass cylinder (Fig. 5) the internal diameter of which is made with a manufacturing tolerance of $\pm 0.003 \text{ mm}$. This diameter tolerance and the low-wear high-precision spindle inside the piston burette reduce the maximum fault deviation caused by the piston burette to $\pm 0.1\%$ for the TITRONIC® 110 *plus* Piston Burette e.g. with the interchangeable units of the TA 10, TA 20, and TA 50 types.

Fig. 1 TITRONIC® 110 *plus* Piston Burette

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and the following normative
document

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Kolbenbürette

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30. Januar 2003

SCHOTT Instruments GmbH

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D-55122 Mainz

Deutschland, Germany, Allemagne

SCHOTT

Technical properties of the TITRONIC® 110 *plus* Piston Burette

(Release 30. 01.03)

Multi-functional display:	4-digit display, approx. 25 x 65 mm, 4-digit volume display with floating point, height approx. 18 mm,
Piston drive:	Pulse-controlled step motor, maximum frequency = 330 Hz,
Interchangeable Units:	1 ml, 5 ml, 10 ml, 20 ml, and 50 ml, UV protection filter, erasable front label;
Resolution:	10,000 pulses per 100 % of volume, smallest increment 0.1 µl with 1 ml interchangeable unit
Accuracy:	0.2 % to 0.1 % relative to 100 % of volume, depending on interchangeable units
Reproducibility:	0.1 % to 0.05 %, depending on interchangeable unit
Display accuracy:	0.000...9.999 ml, 10.00...99.99 ml, 100.0...999.9 ml, 1000...9999 ml
Filling speed:	100 % of volume in 20 s to 999 seconds, freely selectable
Discharge speeds:	0.01 ml/h to 100 ml/min, depending on interchangeable unit
Valve switch-over:	automatic
Inputs and outputs:	2 RS-232-C serial interfaces, 9-channel sub-miniature D jacks, keyboard port: 5-channel mini DIN socket PS2 I/O port: 15-channel sub-miniature D jack, e.g. for connecting a TR 160 corresponds to Protective Class I according to DIN 57 411 / VDE 0411, Part 1
Power supply:	
Mains connection:	230 V~ / +6 %, 50...60 Hz or 115 V~ / ± 10%, 50...60 Hz, internal conversion
Mains fuse:	Fine wire fuse 5 x 20 mm: with 230 V ~, 0.25 A, slow-blow type with 115 V ~, 0.5 A, slow-blow type
Power consumption:	35 VA
CE-Sign: 	EMC - Compatibility according to the directive 89/336/EEC of the Council; Interference emission according to the EN 61 326 standard, Part 1 Jamming resistance according to the EN 61 326 standard, Part 2 Low-voltage directive according to the Directive 73/23/EEC of the Council, as amended by the Directive 93/68/EEC of the Council Testing basis: EN 61 010, Part 1
Country of origin:	Made in Germany
Dimensions:	150 x 360 x 290 (W x H x D) with interchangeable unit
Weight:	approx. 5.1 kg with interchangeable unit approx. 4.3 kg without interchangeable unit

® Registered trademark of SCHOTT Instruments GmbH, Mainz

SCHOTT Instruments GmbH

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D-55122 Mainz

Deutschland, Germany, Allemagne

SCHOTT

Interchangeable units:

Interchangeable units

Compatibility:	The units are interchangeable with those of the T 100, T 110/220 Piston Burettes and the titrator of the TitroLine <i>alpha plus</i> type
Coding of the cylinder size:	Automatically, the mounting for the glass cylinder bears an unmistakable code
Valve:	Volume-independent cone valve made of fluorocarbon polymers (PTFE)
Cylinder:	Made of DURAN®
Hoses:	FEP hose set with UV-protection
Bracket for storage bottle:	To fit glass bottle to DIN 12 036 and DIN 12 039, rectangular glass or plastic bottle
Fastening for interchangeable unit:	With automatic checking feature
Materials:	Borosilicate glass, fluorocarbon polymers, PPS, high-grade steel,, polypropylene, acrylic glass
Dimensions:	145 x 260 x 270 mm (W x H x D) without reagents bottle
Weight:	Approx. 0.80 kg for TA 20 interchangeable unit without reagents bottle
Dosing accuracy:	According to DIN EN ISO 8655, Part 3 Correctness: 0.1 - 0.3% Precision: 0.05 - 0.1% (depending on the interchangeable unit being used)

Dosing accuracy of the TITRONIC® 110 *plus* Piston Burette when used with TA and TA *plus* Inter-changeable Units:

Inter-changeable unit type no.	Volumen [ml]	Tolerances of the Ø _i of the glass cylinders [mm]	Correctness referred to 100 % of volume [%]	Precision [%]
TA 01	1,00	± 0,003	± 0,3	0,10
TA 05 <i>plus</i>	5,00	± 0,003	± 0,15	0,07
TA 10 <i>plus</i>	10,00	± 0,003	± 0,1	0,05
TA 20 <i>plus</i>	20,00	± 0,003	± 0,1	0,05
TA 50 <i>plus</i>	50,00	± 0,003	± 0,1	0,05

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1.2 Warning and safety information

For safety reasons the TITRONIC® 110 *plus* Piston Burette may only be opened by authorised persons; for example, work on the electrical equipment, for example, may only be performed by trained specialists. If the unit is opened without authorisation, or if it is damaged maliciously or through negligence, the warranty will become void.

The TITRONIC® 110 *plus* Piston Burette conforms to protection class I. It has been built and tested in accordance with EN 57 411, Part 1, Protection Measures for Electronic Measuring Instruments and has left the factory in perfect technical condition. In order to maintain this condition and ensure safe operation, the user must observe all the instructions and warnings contained in the present operating instructions.

Before switching the unit on, please make sure that its operating voltage matches that of the mains voltage. The operating voltage is specified on the type plate. The mains plug may only be inserted into a socket with an earthed contact. The earth must not be eliminated by using an extension cable without an earthed wire. Any discontinuity in the earth wire inside or outside the device or disconnection of the earth wire may lead to the device becoming hazardous. The earth wire must not be disconnected intentionally.

When replacing any fuses, they have to be of the rated type and nominal voltage. The relevant information can be found on the back panel of the device. The use of patched up fuses or short-circuiting the fuse holders is not allowed.

If it has to be assumed that it is not possible to operate the TITRONIC® 110 *plus* Piston Burette safely, it is to be taken out of operation and secured against unintentional operation: Please switch the TITRONIC® 110 *plus* Piston Burette off, remove the mains cord, and call the service department of SCHOTT-GERÄTE.

Examples for the assumption that a safe operation is no longer possible:

- the TITRONIC® 110 *plus* Piston Burette shows visible damages,
- the TITRONIC® 110 *plus* Piston Burette does not function properly,
- serious strain is present during transportation.

The TITRONIC® 110 *plus* Piston Burette must not be stored or operated in humid rooms.

For reasons of safety, the TITRONIC® 110 *plus* Piston Burette must only be used for the range of application described in the present operating instructions. The relevant regulations applicable to the liquids to be dosed must be adhered to. The user has to ensure that the persons charged with the operation of this piston burette are experts in handling the matters occurring in context with the piston burette.

When working with fluids which are not conventional titration media, the chemical resistance of the materials of the TITRONIC T 110 *plus* Piston Burette requires particular attention (please refer to chapter 1.1, "Summary, Declaration of Conformity").

When using liquids with a high vapour pressure and/or matters or mixed matters not qualified for use under chapter 1.1, "Summary, Declaration of Conformity", the user will be responsible for the safe and proper operation of the piston burette.

When the piston rises in the cylinder a microfilm of the dosed fluid remains on the internal wall of the cylinder, but this has no influence on the dosing accuracy. However, this minimal residue of fluid may evaporate and thus penetrate into the zone beneath the piston and corrode or dissolve the materials which are in use there (also cp. chapter 6, "Maintenance and care of the TITRONIC T 110 *plus* Piston Burette").

The TITRONIC T 110 *plus* Piston Burette is equipped with integrated circuits (EPROMs). These circuits can be deleted by X-rays, radioactive rays, or UV rays. Whilst UV rays are retained completely by the housing, X-rays or other high-energy rays may penetrate the device housing and delete the program.

Please observe also the relevant operating instructions for any units to be connected to the TITRONIC T 110 *plus* Piston Burette, for example, titrators, keyboards, sample changers, manual titration controllers, personal computers.

2 Set-up and first operation

2.1 Unpacking and setting up the TITRONIC T 110 *plus* Piston Burette

The piston burette, all its accessories, and the periphery units have been carefully checked at the factory to ensure their correct function and size.

Please ensure that the small accessories must also be removed in full from the packaging.

For the scope of delivery, please refer to the enclosed parts list.

The bottom section of the piston burette may be set up on any flat surface. Please plug the mains connection cable to the mains connection jack at the back panel, and the mains plug to mains. Before switching on, please ensure that the device's operating voltage and the mains voltage are the same. The possible operating voltages are shown on the type plate (please refer to the base plate of the TITRONIC® 110 *plus* Piston Burette). Take care that the voltage selector on the base plate is set to the proper operating voltage.

2.2 Installing the stand

The stand-rod bracket (type no. TZ 1516) is contained in the scope of delivery and is to be screwed to the base plate of the piston burette using the screw. To do so, please turn the bottom section of the Piston Burette upside down, then attach the stand-rod bracket to the base. Subsequently, screw the stand rod into position. The TM 135 Magnetic Mixer is attached using the holder which comes with the device.


Fig. 2

2.3 TA *plus* Interchangeable Unit:


Fig. 3

- 1 Built-in TZ 3000 valve with automatic control
- 2 Hose between valve and titration tip
- 3 Titration tip
- 4 Inscription field
- 5 Locking feature with lower burette section
- 6 Cylinder pressure screw
- 7 UV protection
- 8 Receptacle for reagents bottle according to DIN 12 036 and DIN 12 039 or rectangular glass or plastic bottle
- 9 Adjustment of reagents bottle
- 10 Suction hose to valve
- 11 Clamped screwing connection
- 12 FEP hose for valve <-> cylinder connection

Items no. 2, 3, 10, 11, and 12 are only available with the complete TZ 3270 kit.

2.4 Fitting and replacing an interchangeable unit


An interchangeable unit is installed by simply slipping it on over the guide rail from front to back. The installation is completed when you can clearly see the mechanical retention device snap in.

If you should encounter any problems during installation, the position of the PTFE piston inside the glass cylinder is to be adjusted. To adjust this piston, please use the piston drawer provided to pull it downwards until the top edge of the PTFE piston is approx. 8 mm below the upper edge of the black pressure screw. The actual adjustment can now be done by simply pressing the interchangeable unit onto a flat surface, for example a table top.

To replace an interchangeable unit, the piston in the cylinder has to be set to the lower starting position. To do so, please trigger the filling operation (also refer to chapter 4.7, "Initial filling and rinsing of the TITRONIC® 110 *plus* Piston Burette"). Subsequently, the interchangeable unit can be easily removed towards the front by pressing the retention key on the left side of the unit (Fig. 4).

After removing the interchangeable unit, the signal tone will sound, and the display will show <<hd-->>. After installing an interchangeable unit the burette size will be recognised and displayed automatically, e.g. <<hd 10>> in the case of a 10 ml interchangeable unit.

When working with a square glass bottle, please use the clamping ring provided to prevent the reagents bottle from falling out of the holder.

Caution: After each replacement, an automatic filling incl. float balancing will take place. For this reason, please **DO NOT** remove the interchangeable unit **immediately**, but only as soon as the display shows 0.0 (cp. Fig. 3).

If a new interchangeable unit is installed without reagent, you have to initiate an "Initial filling" or "Rinsing" of the piston burette (please refer to chapter 2.6).

The front section of the interchangeable units contain a coding which communicates the size of the installed cylinder to the bottom section of the burette. This ensures accuracy of the dosed volume at any time.

Fig. 4

2.5 Replacing the glass cylinder

The hose between the cylinder and the valve is to be unscrewed from the cylinder (see Fig. 3, item 6). Pull the UV protection down, then loosen the black press screw (cp. Fig. 3, item 6) for the glass cylinder by turning it anti-clockwise, and pull it off. The cylinder with the piston inside can now be removed. A new cylinder with a piston is to be fitted in reverse order and then adjusted as described under chapter 2.4 "Fitting and replacing an Interchangeable Unit". Basically, it is to be ensured that only the intended cylinder size is fitted within an interchangeable unit, since otherwise the code in the interchangeable unit will not match the cylinder size. The result would be an incorrect dosage. For reasons of dosage and analysis accuracy, it is always recommended to replace the PTFE piston at the same time when a defect glass cylinder is replaced. This applies in particular if the glass is broken, since the sealing rings of the PTFE piston may be damaged by splinters of glass.

Please note: As a rule, the hoses and cylinder contain chemicals which may escape or be sprayed around when the unit is dismantled. The relevant safety precautions for dealing with these chemicals have to be observed.


Abb. 5 Installation of the glass cylinder

2.6 Initial filling and rinsing of the complete interchangeable unit

Initial filling of the burette unit is carried out by a rinsing program. You initiate the rinsing program using the <Rinse> key on the TZ 2835 keyboard, or by pressing the unlocking key of the TA interchangeable unit during power-up.

When starting this operation using the <Rinse> key, the display will first show the selected number of rinsing cycles. 1 to 3 cycles can be input. You start the rinsing cycle by pressing the <START> key. The number of number of rinsing cycles will be displayed during the rinsing process.

Fig. 6 Display with two rinsing cycles selected

After the start of the rinsing process, one upward and one downward move of the piston will be performed per rinsing stroke. The number of strokes can be selected from 1 to 3.

Initial filling and rinsing **cannot** be initiated using the software.

During the execution of the initial-filling or rinsing program you have to place a waste container under the titration tip.

After the start of the rinsing program, the display will additionally show a colon, and the piston will move to its bottom position.

Subsequently, the contents of the cylinder (air or old reagent) will be ejected. As soon as the piston top position is reached, fresh reagent will be drawn in. When the lowermost piston position is reached, the complete contents of the cylinder will be ejected once again. Another filling cycle will take place, and after float balancing the message <<0.000 ml>> will appear in the case of interchangeable units from 1 to 10 ml, or <<0.00 ml>> with interchangeable units of 20 and 50 ml, respectively.

Please note: The initial-filling / rinsing program is not suitable for changing titration media, solutions, and liquids which widely differ in composition and/or concentration. Such cases require a careful cleaning of the parts coming in contact with the titration agents, such as intake tube, valve, glass cylinder, and titration tip (for their disassembly procedure, please refer to chapter 2.5, "Replacing a glass cylinder").

After the execution of the program, the TITRONIC® 110 *plus* Piston Burette is ready for the next titration or dosing processes.

2.7 Combining the TITRONIC® 110 *plus* Piston Burette with other devices

The following devices can be connected to the TITRONIC® 110 *plus* Piston Burette:

- The TitroLine *alpha plus* Titrator
- The TR 160 manual titration controller
- A multi-functional keyboard of the TZ 2835 type, PC keyboards in AT mode
- A TITRONIC® *universal* Piston Burette
- TW alpha for automatic sample preparation

For a description of the connections to be made and the functioning of the unit, please refer to chapter 4, "Working with the TITRONIC® 110 *plus* Piston Burette". In addition, please note the operating instructions provided with the devices being connected to the TITRONIC® 110 *plus* Piston Burette. If you connect a device with a power supply of its own to the TITRONIC® 110 *plus* Piston Burette, please make sure both devices are switched off when making the connection.

Connected device	Connection to the 110 plus (Fig. 15)	Cable used Type no.	Explanations
TitroLine alpha T110, T200, TW alpha	RS-232-C (2)	TZ 3086	With varying addresses in daisy-chain operation
TitroLine alpha <i>plus</i> TITRONIC 110 <i>plus</i> <i>plus</i>	RS-232-C (2)	TZ 3084	With varying addresses in daisy-chain operation
TITRONIC universal	RS-232-C (2)	TZ 3087	With varying addresses in daisy-chain operation
Weighing-balance: Sartorius: A 200 S Mettler: AB.S, PG Mettler: AT, PM	RS-232-C (2)	TZ 3082 TZ 3081 TZ 3083	For other weighing-balance use special cable
TR 160	I/O digital port	Firmly connected cable	
Keyboard: TZ 2835	Keyboard connector, PS 2	Firmly connected cable	For use in stand-alone operation

3 Display

3.1 LC Display

The LCD display consists of a 4-digit, 7-segment LCD display

Fig. 7 LCD Display showing all segments.

In the process of parameterisation, the segments of the 7-segment display are used to form meaningful abbreviations.

3.2 Dialogue example with the TITRONIC® 110 *plus* Piston Burette

Display after dosing using the PC keyboard. 15.00 ml were dosed.

The image shows a digital LCD display. On the left, there is a small icon of a piston or burette tip, followed by the letter 'd'. To the right of this, the number '15.0' is displayed in a large, clear, segmented font.

Fig. 8 LCD display operating as a volume display

3.3 Status messages of the TITRONIC® 110 *plus* Piston Burette

The TITRONIC® 110 *plus* Piston Burette has no additional symbols for status reporting. While filling, the display shows a colon.

4 Working with the TITRONIC® 110 *plus* Piston Burette

4.1 Switching the piston burette on

Please do not switch the TITRONIC® 110 *plus* Piston Burette on before all devices required for work are connected.

The display sequence after power-up is as follows:

- Display of the software version number, e.g. 03 03
- Display of the interchangeable unit being used, e.g. <<hd 50>>
- Display of the set RS address, e.g. <<Ad 00>>
- Display of the interchangeable unit being used, e.g. <<hd 50>>
- Display 'Burette is filling' <<00:00>>
- - Display ready <<0.000>>

4.2 Possible applications of the TITRONIC® 110 *plus* Piston Burette

Possible applications of the TITRONIC® 110 *plus* Piston Burette

Depending on the connected devices, the TITRONIC® 110 *plus* Piston Burette is very versatile to use:

- The TITRONIC® 110 *plus* Piston Burette is controlled by other devices through the RS-232-C interface. In this case the RS-232-C interface takes priority over all other functional units (please refer to chapter 4.3, "Titrating using the TitroLine *alpha plus* Titrator").
- Titrations are performed under manual control using the TR 160 Manual Titration Controller (please refer to chapter 4.4, "Titrating using the TR 160 Manual Titration Controller").
- When connecting a keyboard in addition to the TR 160 Manual Titration Controller, it is possible to parameterise the TITRONIC® 110 *plus* Piston Burette.
- It can be used to dose a pre-defined volume.

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Titrating using titration devices equipped with a RS-232-C interface

The connected titration device (e.g. a TitroLine *alpha plus*) takes total control of the TITRONIC® 110 *plus* Piston Burette. Control by a titration device via the RS-232-C interface takes priority over all other units which are connected, e.g. over the TR 160 Manual Titration Controller.

As a rule, the data link between the TITRONIC® 110 *plus* Piston Burette and a titration device equipped with a RS-232-C interface is made through a data cable of the TZ 3084 type, with the data input 1 of the TITRONIC® 110 *plus* Piston Burette being used.


4.4 Titration using the TR 160 Manual Titration Controller

The TR 160 Manual Titration Controller connects directly to the I/O port (Fig. 15, item 3) of the TITRONIC® 110 *plus* Piston Burette. If no other devices are connected, e.g. a keyboard, the TITRONIC® 110 *plus* Piston Burette will recognise the TR 160 Manual Titration Controller as the only controlling device. With the TR 160 Manual Titration Controller connected as sole device, manual titration can begin right away..

Functioning

With the TR 160 Manual Titration Controller, titrating can be controlled manually. This involves a 10-level switch and 4 press-buttons.

All the keys are logically interlocked. If any of the key is depressed, the piston burette will not react to any other command until the operation in progress is completed.

Item 1: The level-switch

10 different titration speeds can be selected with this level switch. The titration speed may also be changed during a running titration by turning this switch. Setting <1> gives the slowest titration speed (100 % of burette volume/300 sec). Setting <10> gives the highest titration speed (100 % of burette volume/30 sec). A useful feature for manual titration is the possibility of switching directly from the highest to the lowest titration speed (setting 10 to setting 1) or vice versa. The current titration speeds are indicated on the display of the TITRONIC® 110 *plus* Piston Burette.

Item 2: <F> button

The <F> (Filling) press-button initiates the filling process of the burette from any piston position. This requires just a brief depression of the press-button. The ex-factory setting for the filling time for 100 % of volume is 30 s for all interchangeable units.

Item 3: <0> button

Press the <0> button to reset the display of the piston burette to <<0.000>> or <<0.00>>, respectively. When the titer is small, this feature enables a number of titrations to be carried out in fast succession without any refilling being required between them.

However, before changing the interchangeable unit, initiate the filling cycle (press the <F> button), to move the piston to the lower end position.

Fig. 9 TR 160 Manual Titration Controller

Item 4: < --- > button

For as long as this button is held depressed, reagent is ejected at the rate selected with the rotary switch. Thereby the ejected volume is shown continuously by the display on the piston burette. Briefly pressing the pushbutton (< 0.3 sec.) allows accurate dosing of down to 0.001 ml with the TA 01 and TA 05 interchangeable units, or down to 0.01 ml with the TA 10, TA 20, and TA 50 interchangeable units. It is also possible to titrate more than the volume of the interchangeable unit.

Item 5 : < - - -> button

If this button is pressed, volumes are ejected in steps of 0.010 ml for burette volumes 1 ml or 5 ml, or in steps of 0.10 ml for burette volumes 10 ml, 20 ml or 50 ml. The step size of the titration steps is firmly programmed and cannot be changed externally. A brief actuation of the button will suffice; prolonged pressing has no further effect.

4.5 Parameterising the TITRONIC® 110 *plus* Piston Burette using a keyboard

The TITRONIC® 110 *plus* Piston Burette can be parameterised for various tasks. This requires a multi-functional keyboard with a PS2 port.

It is also possible to use a mini keyboard. In cases of doubt, please consult your keyboard dealer whether the keyboard in question is suitable.

The PC keyboard of the TZ 2835 type from SCHOTT-GERÄTE will in any case fit for the system. This keyboard comes with an adhesive label which can be attached above the function-key bar to identify the functions of the keys clearly.


Fig. 10 PC keyboard, type no. TZ 2835

Action of the keyboards

In combination with the TITRONIC® 110 *plus* Piston Burette, the keys of the keyboards being used have specific functions assigned to which are listed in the table below (also cp. Fig. 10)

Mini PC-keyboard, type no. TZ 2835 from SCHOTT	Standard MF2-keyboard or standard mini-keyboard	A c t i o n
<START>	<F1>	Start of the dosing operation.
<STOP>	<F2>	Stop of dosing / Filling/Initial filling.
<Rinse>	<F3>	Start of initial filling, number of repetitions, 2 times min.
<FILL>	<F4>	Fill command, executed in the set time, in terms of ml/ min.
<Fill Speed>	<F5>	Call of the menu for setting the filling speed.
<Param 1>	<F6>	Call of the menu for setting the RS parameter set.
<Param 2>	<F7>	Call of the menu for setting the baud rate.
<Param 3>	<F8>	Call of the menu for setting the device address.
<Dos Speed>	<F9>	Call of the menu for entering the dosing speed.
<DOS>	<F10>	Dosing mode; if the Start key is pressed in this mode, the set volume will be dosed.
<ESC>	<ESC>	<ESC> will take you to the previous level within the menu. A max. of 2 steps is required.
<↑> <↓>	<↑> <↓>	Scrolling in overview mode.
<ROLL UP> <ROLL DOWN>	<ROLL UP> <ROLL DOWN>	Setting of the contrast of the alphanumeric display.
0...9 Decimal point	0...9 Decimal point	Parameter input.
<ENTER>	<ENTER> or <↵>	Acknowledgement of the input parameters. The present operating instructions always quote the <ENTER> key in reference to both keys.
	 or <←>	Deletion of the input figure/character above the cursor.
Letters, ASCII characters	Letters, ASCII characters	Without any function.
All other keys	_____	are without any function.

4.6 Dosing using the TITRONIC® 110 *plus* Piston Burette

Using the TITRONIC® 110 *plus* Piston Burette, it is possible to perform dosing tasks with variable parameters. The variable parameters include: the volume, the dosing speed, and the filling time.

Dosing requires a keyboard as an additional device (also refer to chapter 4.5, "Parameterisation of the TITRONIC® 110 *plus* Piston Burette using a keyboard"). When connecting the keyboard, the TITRONIC® 110 *plus* Piston Burette has to be switched off.

After power-up, the device reports with the display set to 0. Pressing the <DOS>key will take the device to the input mode for the dosing volume. The desired volume can be set using the numeric keys. The decimal point is recognised as decimal character. The numeric input can be deleted using the backspace key.

Fig. 11 Volume display after power-up.

The actual dosing procedure

Parameterisation of a dosing volume:

Input of the volume to be dosed in ml, dosing is possible from 0.000 to 9999 ml, the input has to be acknowledged with the <ENTER> key, the selected volume will be dosed immediately.

If the input volume was acknowledged using the <ENTER> key, the input volume can also be dosed using the <START> key. In principle, any input down to ± 0.001 ml is possible, since after parameterisation and storage of volume any interchangeable unit can be used. If an interchangeable unit with a volume of 20 or 50 ml is used, the last digit which was input will be cut off on the display.

Titration with manual control and predosing

If a TR 160 Manual Titration Controller is connected as an additional device for the dosing function, manually controlled titrations with predosing can be performed. To do so, the TR 160 Manual Titration Controller has to be activated after the connection and power-up of the TITRONIC® 110 *plus* Piston Burette by performing the first titration. If you move to the <DOS>function subsequently, the TR 160 Manual Titration Controller can be used to continued titration after predosing.

4.7 Filling and rinsing the TITRONIC® 110 *plus* Piston Burette

Depending on the application, the TITRONIC® 110 *plus* Piston Burette has to be filled at various times or occasions. In this process, the ejected volume from the storage bottle is replaced. This function is not suitable for "filling" an empty cylinder. This rather requires the rinsing function with a minimum of 2 rinsing cycles.

"Filling" of the TITRONIC® 110 *plus* Piston Burette is initiated:

- by pressing and then "releasing" the locking device on the interchangeable unit (Fig. 3, item 5, and Fig. 4),
- by pressing the <F4> or <FILL> key if a keyboard is connected,
- automatically after pushing on a interchangeable unit,
- by a command from the PC-controlled TitriSoft Titration System or a Titration Controller

"Rinsing" of the TITRONIC® 110 *plus* Piston Burette is initiated:

- by holding down the locking button on the TA sample changer during power-up (Fig. 3, item 5) (please refer to chapter 2.6, "Initial filling or rinsing of the complete interchangeable unit");
- by pressing the <F3> or <RINSE> key if a keyboard is connected.

In the course of the filling or rinsing process, the display of the colon indicates that these functions are being performed. Upon completion of the filling process, the display will show <<0.000 ml>> or <<0.00 ml>>.

During the rinsing process, the volume indicator will show the volume count for supervision purposes. The display shows <<r n 1>>.

After two cycles, rinsing is finished, and the display will show <<0.000 ml>> or <<0.00 ml>>.

4.8 System settings of the TITRONIC® 110 *plus* Piston Burette: Rinsing

After pressing the <Rinse> key, the selected number of rinsing cycles will be displayed. Using the numeric keys, 1 to 3 cycles can be input.

A digital display showing the text 'r 1 n 2' in a segmented font. The 'r' is on the left, followed by a space, then '1', another space, 'n', and finally '2'.

Fig. 12 Display if 2 rinsing cycles were specified

4.9 System settings of the TITRONIC® 110 *plus* Piston Burette: Filling speed

After pressing the <Fill Speed> key, the set fill speed will be displayed. An input can be made using the numeric keys. Input of the time: <numeric keys> and <ENTER>. The filling time has to be modified if the liquid to be dosed has a high viscosity or a high and vapour pressure. Filling time between the 20 s per 100% of volume and 300 s per 100% of volume are possible, as a rule, a reduction down to 300 s per 100 % of volume will be sufficient.

A digital display showing the text 'F 20' in a segmented font. The 'F' is on the left, followed by a space, then '20'.

Figure 13 Display for filling speed = 20 s

4.10 System settings of the TITRONIC® 110 *plus* Piston Burette: Dosing speed

After pressing in the <DOS SPEED>/F9 key, the set dosing speed will be displayed in terms of ml/min. Using the numeric keys, an input from 0.01 to 100 is possible. When entering the dosing speed, the max. possible speed is limited by double the cylinder volume of the interchangeable unit used per minute, for instance, 100 ml/min for a TA 50 interchangeable unit with a volume of 50 ml.

The smallest admissible speed is 0.01 ml/min. Please acknowledge your input using the <ENTER>key.

A digital display showing the text 'c 100' in a segmented font. The 'c' is on the left, followed by a space, then '100'.

Fig. 14 Display of dosing speed = 100 ml/min.

This indication is only valid for a 50 ml TA; in this case the theoretical dosing speed of 100 ml/min is reached.

5 Connections and data transfer


Fig. 15 Back panel of the TITRONIC® 110 *plus* Piston Burette

- 1 Mains connector with mains filter and fuse

2 Mains switch

3 Multi-function I/Os
- 4 Bi-directional RS-232-C-1 data ports

5 Bi-directional RS-232-C-2 data ports

6 Keyboard port

Cable	
What is being connected	
RS 1	PC with cable, length 5 m
	or with cable, length 19 m
For daisy chain, please refer to RS 2	
RS 2	TitroLine <i>alpha plus</i> , T 110 <i>plus</i>
	TW <i>alpha plus</i>
	T 110, T 200, TW <i>alpha</i>
	T <i>universal</i>
	Sartorius weighing-balance
	Mettler AT weighing-balance
	Mettler AB-S, PG etc. weighing-balance
Keyboard	Keyboard

5.1 RS-232-C interfaces, inputs and outputs

The TITRONIC® 110 *plus* Piston Burette is equipped with two serial interfaces (RS-232-C) (Fig. 15, item 4 and item 5) for data communication with other devices.

You can use these two interfaces to connect several devices to one interface of a PC. Using PC programs such as TitriSoft, a large number of titration, dosing, and sample-preparation tasks can be realised.

Interconnection is always made according to the same principle. Make a connection from the PC to interface 1 of the first device, then make a connection from interface 2 of the first device to interface 1 of the second device; subsequently, use interface 2 of the second device to connect the next device etc. In this way up to 16 devices can be connected to, and addressed by, one PC interface. The devices are addressed through one individual address which has to be set uniquely for each device. All the devices and the PC have to be set to the same transfer parameters.

The plug connectors of these interfaces are located on the back panel of the device. The “RS 1” socket forms the connection to a connected PC or titration device equipped with a serial interface. The “RS 2” socket can be used to connect further devices within the titration system, also in chained mode.

PIN assignment of the RS-232-C interfaces:

Pin no.	Signification	Pin no.	Signification
2	R x D data input	8	CTS
3	T x D data output	5	Digital mass

5.2 Parameterisation of the RS interface:

Please note: Any modification of the data parameters requires a connected keyboard. The data transfer parameters are set using parameter keys on the keyboard.

Setting the data transfer parameters

The ex-factory setting of the data transfer parameters is 4800 baud, 7 bits, 1 stop bit, and No Parity (this corresponds to mode: r 7 n 1). This basic setting is valid for operation using SCHOTT titration devices.

If you are using earlier SCHOTT titration devices, you may also use the 7 2 n STANDARD mode of these devices.

The possible parameter sets - to be set through the <Param 1> (F6) menu - are:

Parameter mode	Data bits	Stop	Parity
r 7 n 1	7	1	no (Standard setting)
r 8 n 1	8	1	no
r 7 E 1	7	1	even
r 8 E 1	8	1	even
r 7 o 1	7	1	odd
r 8 o 1	8	1	odd

Fig. 16 Example of the standard RS parameter display

Setting the baud rate :

The baud rate is set using the <Param 2> key (F7). After pressing the key, the currently set baud rate will appear. Using the 'arrow up' or 'arrow down' key, you can select the desired baud rate. You can leave the selection menu using the <ESC> key. All devices communicating with each other must have the same transfer parameters. You may choose among the following baud rates: 2400, 4800, 9600, and 19200.

Fig. 17 Example of the baud rate parameter 4800 baud.

Setting the device address

The use of a device address enables several devices to be operated within a chain. For this purpose, a RS-232-C cable, for instance of the TZ 3088 type, is used to make a connection from the computer to socket 1 of the first device of the chain. Using another RS-232-C cable, for instance of the TZ 3084 type, socket 2 of the first device is connected to socket 1 of the second device. Socket 2 of the second device can be used to connect additional devices. Please make sure that the devices in the chain have different addresses. The device address of the TITRONIC® 110 *plus* Piston Burette is set using the address parameter menu. To access this menu, please press the <<Param 3>> key (F8). Now you can use the numeric keys to enter an address from 00 to 15.

Fig. 18 Example of address parameter address 00.

5.3 Command list for RS communication

The commands have the following form:

Adresse	2-digits aa,	e.g.: 01
Command		e.g.: BP
Parameter		e.g.: 14 <CR> <LF>

Each command must be entered separately and ended by <CR> <LF>. All replies will be returned to the computer only upon completion of the respective action.

Command	Description	Reply
aaBF	Fill burette	aaY
aaBH14	Number of steps up 1 to 10000	aaY
aaBI3456	Step increment 1 to 10000	aaY
aaBN	Reset of ml display to zero	aaY
aaBR145	Number of steps down 1 to 10000	aaY
aaBP	Inquire piston position in scale units	aaBP=1234
aaBV	Inquire burette volume	aaV=12.34
aaDA13.56	Carry out dosing with addition of the volume	aaY
aaDO123.4	Dose the specified quantity in ml	aaY
aaDB123.4	Dose quantity in ml without filling and zero	aaY
aaDM34.56	Carry out dosing (ml) with sensor supervision	aaY
aaGDH11,3	Speed for dosing in ml/h	aaY
aaGDM20,3	Speed for dosing in ml/min	aaY
aaGF20	Filling time in seconds (min 20s)	aaY
aaOA1	Switch off specified output 1	aaY
aaOE1	Switch on output 1	aaY
aaOJ	Switch on output	aaY
aaON	Switch off output	aaY
aaRH	Identification request	aaIdent T200
aaRA	Inquire unit size	aaAUFSATZ:20
aaRS	Report status	aaStatus:"text"
aaRC	Inquire "last command"	aa"last comm."
aaVD123.4	Dose the specified quantity in ml	aaY
aaVE	Inquire software version no.	aaVersion:9049
aaVF	Valve in 'fill' position	aaY
aaVN	Valve in 'dose' position	aaY
aaVO	Volume inquiry from display	aaV=12.345ml
aaSEEPROM	Reset EEPROM to ex-factory values	aaY
aaSR	Stop current function	aaY

Example: Dose 12.58 ml using the TITRONIC® 110 *plus* Piston Burette having the address "03":
 "03DO12.58 <CR> <LF>". Reply from the piston burette upon completion of the dosing operation: "03Y
 <CR> <LF>".

5.4 Keyboard port

The keyboard port (ref. to Fig. 15, item 6) is used to connect the TZ 2835 PC keyboard or a PC multi-functional keyboard. The port is capable of handling PC AT data formats. When using a third manufacturer's keyboard, please make sure that the keyboard is working in the "AT format" (please refer to chapter 4.5, "Parameterisation of the TITRONIC® 110 *plus* Piston Burette using a keyboard").

5.5 In/Outputs, input/output socket

The back panel of the TITRONIC® 110 *plus* Piston Burette features a 16-channel socket serving as in- and output (please refer to Fig. 15, item 3). This I/O port is controlled either by the internal program of the TITRONIC® 110 *plus* Piston Burette, or by computer commands. Among other things, it is used to connect the TR 160 Manual Titration Controller.

The multi-functional I/O consists of 8 inputs and 1 output:

Input:	Pins no. 1; 2; 3; 4; 5; 6; 7; 8
Output:	Pin no. 9
Free:	Pins no. 10; 11; 12; 13
Ground:	Pins no. 14; 15

The inputs correspond to TTL signal and are circuited over a pull-up resistor (10 kΩ) to +5 V.

The outputs correspond to those of an open-collector component (LS 05). Therefore, the matching pin of a connected device has to be circuited via a 1 ... 10 kΩ resistor to the supply voltage of the connected device.

6 Maintenance and care of the TITRONIC® 110 *plus* Piston Burette

The preservation of the proper functioning of the piston burette requires certain inspection and maintenance work.

Maintenance intervals

Normal operation:

- In a larger sense, the valve and the glass cylinder are parts which are subject to wear and tear; this means, that the glass surface may be damaged by the use of chemicals. This wear and tear may lead to inaccurately dosed volumes. For this reason, the dosing accuracy of the device is to be verified.
- As a rule, the max. intervals for carrying out all maintenance work is 6 months.

Under particular strain:

- Once per quarter, the glass cylinder has to be inspected visually in order to check for damages, and in addition a test according to DIN EN ISO 8655, Part 3, has to be carried out. These inspections are particularly necessary if solutions were used which are likely to attack glass, e. g. alkaline, fluoride, or phosphatic solutions.
- If there is any suspicion that a solution is attacking the glass excessively, the maintenance intervals are to be reduced accordingly.
- Once per quarter the electrical contacts have to be checked for corrosion if the piston burette is used in rooms the atmosphere of which occasionally contains corrosive matters.

In case of disturbances:

- The above work has to be carried out immediately if a disturbance, an error, or another defect becomes obvious.

Maintenance work to be performed

- Check the piston and the cylinder for visible signs of mechanical damage, such as cracks, fractures and damage of the surface.
- Check the piston rod for corrosion.
- Check the hoses, screw unions, and seals for signs of visible damage, contamination, and leaks.
- Check the electrical plug contacts for corrosion and mechanical damage (on piston burette and cables).
- Defective parts must be repaired or replaced by new ones. Defect glass parts must always be replaced.
- After each maintenance operation, a test for metrological reliability according to DIN 12 650, part 6, or ISO DIN 8655, Part 3, has to be carried out.

Periods without operation

- If the interchangeable units used with the TITRONIC® 110 *plus* Piston Burette are not used over an extended period of time, the liquids contained in the system, in particular aggressive solutions, have to be drained. If the liquid is left in the system, one has to reckon with corrosion and an alteration of these solutions used over time. Since the state of the art does not offer any plastic hoses which would be absolutely free of diffusion occurrences, this precaution measures apply in particular to the hose-line section.

Cleaning

- The interchangeable units may be cleaned using a piece of cloth soaked in normal water or a household cleansing agent containing ethanol. In general, short-term immersion in, or rinsing with, water or ethanol is not harmful, but this method should only be applied in extreme cases.
- The casing of the lower section of the piston burette may also be cleaned using a piece of cloth and a household cleansing agent. The lower and rear side have to be dry-cleaned. In no case must water penetrate into the interior of the lower section.

Reproducibility of results

- The measurement results or analysis results depend on a variety of factors. Please check the plausibility of the measurement results or analysis results at regular intervals, and carry out the required reliability tests. In this regard, please adhere to the usual validation procedures.

7 Storage and transportation

If the TITRONIC® 110 *plus* Piston Burette or the interchangeable units have to be stored over some time, or to be dislocated, the use of the original packing will be the best protection of the devices. However, in many cases this packing will not be available any more, so that one will have to compose an equivalent packaging system. Sealing the lower section in a foil is hereby recommended.

The devices should be stored in a room with a temperature between +10 and +40°C, and the (relative) humidity of the air should not exceed 70 %.

8 Index

- "Parameterisation 51
- back panel 31, 32
- Back panel** 47
- baud rate** 49
- button** 42
- buttons 41
- care 31, 51
- chained mode 48
- changeable unit 8
- cleaning 37
- Cleaning** 52
- Command list 50
- cone valve 8
- Connections and data transfer 47
- Correctness 8
- cylinder 31, 34, 36, 44, 52
- Cylinder 8
- daisy-chain operation 38
- data link 40
- Declaration of Conformity 31
- device address** 49
- digit 44
- display 37, 38, 39, 43, 44, 46, 48, 50
- Display** 37, 38, 39, 46
- dosage 30
- Dosage 30
- dose 39
- dosing 50
- Dosing 30, 44
- Dosing accuracy** 8
- dosing speed 44, 46
- Dosing speed 46
- dosing tasks 44
- dosing volume** 44
- electrode clamp 30
- filling 34, 50
- Filling 43, 44
- filling speed 43, 46
- Filling speed 7, 46
- In/Outputs 51
- Installation** 36
- Installing 32
- interchangeable 37
- Interchangeable 33
- interchangeable unit 7, 8, 34, 36, 37, 39, 42, 44, 45, 46
- Interchangeable unit 8
- Interchangeable Unit 36
- interchangeable units 30, 7, 37, 52
- Interchangeable units** 8
- Interchangeable Units 30, 7, 8
- Interconnection 48
- key 37, 43, 44, 45, 46, 49
- keyboard 38, 39, 41, 42, 43, 44, 48, 51
- Keyboard 37
- Keyboard 38, 47, 51
- keyboard port 7
- Keyboardconnector 38
- keyboards 31, 43
- Magnetic Mixer 32
- Maintenance 31, 51
- Maintenance work to be performed** 52
- Manual Titration Controller 39, 40, 41, 44
- Multi-function I/Os 47
- numeric keys 46, 49
- operating voltage 31
- operating voltages 32
- parameterisation 30, 38, 44
- Parameterisation 44
- Parameterisation of the RS interface 48
- PC keyboard 30
- Periods without operation** 52
- Possible applications 39
- possible uses 30
- Precision 8
- predosing** 44
- PTFE piston 34, 36
- RC-232-C-interfaces 30
- rear side 52
- Replacing the glass cylinder 36
- results** 52
- rinsing 34, 37, 44, 46
- Rinsing 45, 46
- rinsing cycles 37, 46
- rinsing function 44
- rinsing program 37
- RS-232-C interfaces 48
- solvents 30
- speed 30
- square glass bottle 34
- Status messages 39
- stirrer 30
- Storage 52
- Summary, Declaration of Conformity 30
- System settings 46
- the filling 30
- Titrating 30, 39
- titration 30, 31, 40
- Titration 41, 51
- titration solutions 30
- titration tip 37
- TitriSoft 30, 45, 48
- TitroLine 30
- TitroLine *alpha plus* 38, 39, 40, 47
- TitroLine *alpha plus* Titrator 30
- TitroLine *alpha plus* type 8
- TR 30, 39, 41, 44
- TR 160 38
- transportation 52
- UV protection filter, 7
- valve 36, 37, 51
- Valve 50
- viscosity 30, 46
- voltage selector 32

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