**GEBRAUCHSANLEITUNG
OPERATING INSTRUCTIONS
MODE D'EMPLOI
MANUAL DE INSTRUCCIONES**

ViscoClock

VISKOSITÄTSMESSGERÄT
VISCOSITY MEASURING UNIT
APPAREIL DE MESURE DE LA VISCOSITÉ
APARATO DE MEDICIÓN DE LA VISCOSIDAD

SI Analytics
a **xylem** brand

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Important notes: Read and observe these operating instructions carefully before using the Viscosity Measuring Unit ViscoClock. For safety reasons the Viscosity Measuring Unit ViscoClock may only be used for the purpose described in these operating instructions. Please also observe the operating instructions for the units to be connected. All specifications in this instruction manual are guidance values which are valid at the time of printing. However, for technical or commercial reasons or in the necessity to comply with the statutory stipulations of various countries, SI Analytics GmbH may perform additions to the Viscosity Measuring Unit ViscoClock without changing the described properties.

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VZ 6552 Membran (5 Stück)

1	Properties and range of use of the Viscosity Measuring Unit ViscoClock....	18
1.1	Description	19
1.2	Technical data of the Viscosity Measuring Unit ViscoClock Version April 01, 2009	21
1.3	Warning and safety information	22
2	Unpacking the Viscosity Measuring Unit ViscoClock	23
3	Initial operation	24
3.1	Selection of the viscometer	24
3.2	Preparation of the sample	24
3.3	Placing the viscometer into the Viscosity Measuring Unit ViscoClock	25
3.4	Insertion of the Viscosity Measuring Unit ViscoClock into the thermostat bath.....	26
3.5	Connecting the Viscosity Measuring Unit ViscoClock	26
3.6	Connecting the Viscosity Measuring Unit ViscoClock to a computer (PC) or printer	26
3.7	Preparing the Viscosity Measuring Unit ViscoClock for measurement	27
3.8	Starting the Viscosity Measuring Unit ViscoClock	28
3.9	Reading the flow time	28
3.10	Turning the Viscosity Measuring Unit ViscoClock off	28
3.11	Viscosity calculation	29
3.12	Evaluation example.....	29
4	Maintenance	30
4.1	Purification	30
4.2	Inspection.....	30
4.3	Repair.....	30
4.4	Purging and first-time filling.....	30
5	Elimination of failures	31
6	Accessories and spare parts	32
	Declaration of Conformity	last page of the document

1 Properties and range of use of the Viscosity Measuring Unit ViscoClock

The Viscosity Measuring Unit ViscoClock, an electronic, automatic flow-time meter, can be used in combination with a SI Analytics GmbH Ubbelohde-, a Micro Ubbelohde and a Micro Ostwald viscometer for the automatic measurement of the flow time of a liquid at temperatures from - 40 °C up to a maximum of + 150 °C.

The Viscosity Measuring Unit ViscoClock automatically measures the time which a liquid requires to flow from the upper to the lower measurement mark of the measurement bulb through the capillaries of a SI Analytics GmbH viscometer.

The device can be used for any transparent liquid for which the viscometer is suitable according to its operating instructions.

The relevant regulations regarding the handling of the substances used have to be observed: The decree on hazardous matters, the Chemicals Act and the rules and information of the chemicals trade.

For the measurement of the **relative** viscosity, all calibrated or uncalibrated DIN, ISO, ASTM, Micro Ubbelohde and Micro Ostwald viscometers can be used.

For the determination of the **kinematic** (absolute) viscosity, only SI Analytics GmbH viscometers, calibrated for automatic measurements, can be used. As a principle, the measurement of the flow time should be performed in a transparent SI Analytics GmbH thermostat bath in which the Viscosity Measuring Unit ViscoClock is fixed by the manual measurement insert of the bath and which is suitable for the intended temperature range (product information: SI Analytics GmbH thermostats).

The following tempering liquids (tempering range: - 40 to + 150 °C) are permitted (→ operating instructions for the SI Analytics GmbH transparent thermostats):

- Water
- Alcohol water (e.g. ethanol, methanol)
- Glycerine
- Paraffin oil
- Silicone oil
- Glycol

SI Analytics recommend the use of silicone oils as tempering liquids for temperature ranges above 90 °C. There is no experience on hand as to the use of the Viscosity Measuring Unit ViscoClock for other purposes or for its use with other tempering liquids!

Misuse prohibited

The Viscosity Measuring Unit ViscoClock must not be used for other work than the above-mentioned. The manufacturer will not be liable for damaged resulting from such use.

Requirements for the user's staff

The Viscosity Measuring Unit ViscoClock must only be used by authorized persons whose training or knowledge and experience ensure proper operation of the device and who are informed about the dangers. On the part of the user it has to be ensured that the persons entrusted with the use of the Viscosity Measuring Unit ViscoClock are experts in the handling of substances used in the environment and in the ViscoClock or that they are supervised by specialized persons, respectively.

1.1 Description

The Viscosity Measuring Unit ViscoClock automatically measures the flow time of a liquid in a SI Analytics GmbH viscometer to determine the relative and kinematic (absolute) viscosity. The time measured is the time which a liquid requires to flow from the upper to the lower measurement mark of the measurement bulb through the capillaries of the viscometer. Measurement is done in a SI Analytics GmbH thermostat bath (SI Analytics GmbH product information), e.g. transparent thermostats of the CT 72 type series. The thermostat change is done outside the thermostat bath.

The Viscosity Measuring Unit ViscoClock consists of a stand for accommodating the viscometer and the electronic measuring unit. The two measurement marks are integrated in the plastic stand (PPA). The electronic measuring unit is built into a polypropylene casing. Owing to the use of heat-carrier liquids as bath media, colour changes may occur on the stand which, however, do not affect the measurement accuracy.

Fig.1



Fig.2

- 5 Centering prism for positioning and holding the viscometer
- 6 Stand
- 7 Lower measurement mark
- 8 Upper measurement mark
- 9 Guidance plate for viscometer centering
- 10 Fixing plate for positioning and holding the viscometer
- 11D Silicone plug for filling tube for pushing operation
- 11S Silicone plug for capillary tube in sucking operation
- 12 Handle plate
- 13 Electronic measuring unit!
- 14 5-digit LCD display in seconds
- 15 Push button <START>
- 16 Membrane rod for closing and opening of the venting tube
- 17 Silicone membrane for closing the venting tube
- 18 RS-232-C socket for connecting a printer or computer (PC)
- 19 Socket for mains jack plug
- 20 Power supply TZ 1858/100 - 240 V~, 50/60Hz

Guarantee

We provide guarantee for the unit described for two years from the date of purchase. This guarantee covers manufacturing faults being discovered within the mentioned period of two years. Claim under guarantee covers only the restoration of functionality, not any further claim for damages or financial loss.

Improper handling/use or illegitimate opening of the meter results in loss of the guarantee rights. The guarantee does not cover wear parts, as valves and pipes including the thread connections as well as the breach of glass parts. To ascertain the guarantee liability, please return the instrument and proof of purchase together with the date of purchase freight paid or prepaid.

Translation of the legally binding German version.

1.2 Technical data of the Viscosity Measuring Unit ViscoClock Version April 01, 2009

CE-Sign:	CE	EMC compatibility according to the Council Directive: 2004/108/EG; applied harmonized standards: EN 61326-1:2006 Low-voltage directive according to the Council Directive 2006/95/EG Testing basis EN 61 010, Part 1
Country of Origin:		Made in Germany
Measurable time::		up to 999.99 s; resolution: 0.01 s
Accuracy of time measurement:		$\pm 0,01 \text{ s} / \pm 1 \text{ digit}$; but not more accurate than 0.1%, indicated as a measurement uncertainty with a reliability level of 95%.
Viscosity measurement range:		0,35 ... 10 000 mm ² /s (cSt) Absolute, kinematic viscosity furthermore depends on the uncertainty of the numerical value of the viscometer constants and the measurement conditions, especially on the measurement temperature.
Viskosimetertypen:		DIN, ISO, ASTM Ubbelohde viscometers; Micro Ubbelohde viscometers and Micro Oswald viscometers
Anzeige:		5--digit LCD display 20 x 48 mm, digit height 12.7 mm, seconds display with two decimal digits behind the decimal point, resolution 0.01 s
Connections:		
Voltage supply:		Low voltage: U 9 V~/DC Plug connection: Socket for low voltage connection – jack plug-, Positive pole on the pin contact, internal contacts $\varnothing = 2.1 \text{ mm}$, USA/Japan for connecting a 100 - 240 V~, 50/60Hz TZ 1858
RS-232-C Interface:		For connecting a printer with serial RS-232-C interface or a computer (PC) for documenting the data Plug connection: 4-Channel! DIN mini round plug
Configuration of RS-232-C Interface, firmly set:		4800 baud, 7 bit word length, 1 stop bit, no parity. The measurement value is automatically transferred after each measurement. The character string consists of 4 characters ahead of the decimal point, the decimal point itself, 2 characters behind the decimal point and the final characters CR and LF.
Power supply		Corresponds to protective system III, protection class for dust and humidity IP 50. TZ 1858 power supply for 100 - 240 V~, 50/60Hz
Material:	Casing:	Polypropylene
	Stand:	Polyphalamide plastic; sealing membrane and plugs made of silicone
Dimensions :		ca. 490 x 95 x 50 mm (H x W x D)
Weight:		Approx.450 g; power supply approx. 220 g
Environment:		
	Ambient temperature:	+ 10 °C ... + 40 °C for storage and transportation
	Operating temperature:	Stand:-40 °C ... max. + 150 °C electronic measuring unit: + 10 °C ... max. + 40 °C
	Air humidity:	acc. to EN 61 010, part 1 Max. rel. humidity 80% for temperatures up to + 31 °C linear decrease down to 50% rel. hum. at a temperature of + 40 °C

Not suitable for use in hazardous environment!

The Viscosity Measuring Unit ViscoClock is suitable for all transparent thermostats by SI Analytics GmbH.

1.3 Warning and safety information

The Viscosity Measuring Unit ViscoClock corresponds to protection class III, was manufactured and tested according to DIN EN 61010, part 1, protective measures for electronic measurement devices, and has left the company in an impeccable condition as concerns safety technology. In order to maintain this condition, the user should observe the notes and warnings contained in the present operating instructions.

Please also note the corresponding Operating instructions for the devices to be connected.

The relevant regulations regarding the handling of the substances used have to be observed: The decree on hazardous matters, the Chemicals ACt and the rules and information of the chemicals trade. On the part of the user it has to be ensured that the persons entrusted with the use of the Viscosity Measuring Unit ViscoClock are experts in the handling of substances used in the environment and in the ViscoClock or that they are supervised by specialized persons, respectively.

For reasons of safety, the Viscosity Measuring Unit ViscoClock must only be opened by authorized persons; this means, e.g., that work on electrical features must only be performed by qualified electricians.

In the case of nonobservance of these provisions the Viscosity Measuring Unit ViscoClock may constitute a danger.

- Electrical accidents of persons or fire hazard.
- Moreover, in the case of unauthorized intervention in the Viscosity Measuring Unit ViscoClock as well as in the case of negligent or deliberate damage, the warranty will become void.

Prior to powering the device on it has to be ensured that the Operating voltage of the TZ 1858 power supply of the Viscosity Measuring Unit ViscoClock matches the mains voltage. The operating voltage is indicated on the power supply. Nonobservance of this provision may result in damage to the Viscosity Measuring Unit ViscoClock.

If it is assumed that safe operation is impossible *), the Viscosity Measuring Unit ViscoClock has to be put out of operation and secured against inadvertent putting to operation.

- Turn the Viscosity Measuring Unit ViscoClock off and pull the power supply out of the mains socket.
- Remove the Viscosity Measuring Unit ViscoClock from the piece of work.

*) Examples for the assumption that a safe operation is no longer possible:

- the package is damaged,
- the Viscosity Measuring Unit ViscoClock shows visible damages,
- the Viscosity Measuring Unit ViscoClock does not function properly,
- liquid has penetrated into the casing of the electronic measuring unit.

- Use original SI Analytics GmbH viscometers and thermostat baths only.
- Operate the Viscosity Measuring Unit ViscoClock only with the original SI Analytics GmbH TZ 1858 power supply for 100 - 240 V~, 50/60Hz.
- Store and operate the Viscosity Measuring Unit ViscoClock only in dry rooms (→ chapter technical data: environment).
- Protect the power supply from humidity. A moist power supply or jack plug/mains plug must not be used again prior to being inspected by a qualified electrician.
- Do not use the Viscosity Measuring Unit ViscoClock in explosion-proof rooms.
- When working with aggressive or health-endangering measurement liquids, personal protective equipment such as eye protection, protective gloves, laboratory coat and breath protection has to be worn.

Development and production is done within a system corresponding to the DIN EN ISO 9001 pp. standard.

Observe the safety data sheets of the manufacturers!

- When handling solvents please observe the pertinent memorandums of the employer's liability insurance associations (VBG 1 and VBG 113).
- When using alkaline laboratory cleaners there is a danger to health (caustification, damage to skin and eyes)!
Please observe the memorandums of the employer's liability insurance associations and the safety data sheets of the manufacturers.
- When using tampering liquids with a temperature above 50 °C there is a danger of getting burnt when touching the thermostat bath cover!
- Please observe the standard regulations and Standards for capillary-viscosity measurements.
- You can obtain a list of the national and international standards from SI Analytics GmbH.

2 Unpacking the Viscosity Measuring Unit ViscoClock

- Check packing for damage. In case of damage contact your supplier.
- Remove the Viscosity Measuring Unit ViscoClock and the TZ 1858 universal voltage power supply from the packing.

Scope of delivery: Software CD VZ 6555
Universal voltage power supply TZ 1858
Primary adapter for Europe
Primary adapter for USA
Manual pump 'press' VZ 6550

NOTE: The scope of delivery does not include viscometers.

- Place the Viscosity Measuring Unit ViscoClock in a dry environment on a plane surface. Ambient temperature: + 10 °C ... + 40 °C (→ chapter 1.2, technical data).

Owing to the small stand foot and the relatively high center of gravity we recommend storing the Viscosity Measuring Unit ViscoClock always in a horizontal position outside of the bath.

3 Initial operation

3.1 Selection of the viscometer

- Select a suitable viscometer. The flow time of the liquid to be measured depends on the size of the capillaries. The capillary size of the viscometer should be selected in such a manner that flow times greater than 100 s result. You will find assistance in the viscometer catalogue of SI Analytics GmbH.

NOTE: Please use original SI Analytics GmbH viscometers only.

List of viscometers for selecting viscometers which are suitable for measuring the flow times for the application and viscosity determination using the Viscosity Measuring Unit ViscoClock:

Viscosimeter	Type-No.	uncalibrated	Calibrated for measurements		Viscosity determination	
			manual	automatic	relative	absolute
DIN-Ubbelohde	532...	-	-	X	-	X
DIN-Ubbelohde	501...	-	X	-	X	-
DIN-Ubbelohde	530...	X	-	-	X	-
ASTM-Ubbelohde	525...	-	X	-	X	-
ASTM-Ubbelohde	526...	X	-	-	X	-
Micro-Ubbelohde	537...	-	-	X	-	X
Micro-Ubbelohde	538...	X	-	-	X	-
Micro-Ostwald	517...	-	-	X	-	X
Micro-Ostwald	516...	-	-	X	X	-

Owing to their constructional shape the type 502 41, type 502 43 and type 502 50 viscometers are not suitable.

Prior to the first operation we recommend a cleaning process according to DIN 51 562, part 1. Subsequently, the viscometer has to be dried and freed of dust!

3.2 Preparation of the sample

Prior to measuring, low-viscous samples should be filtered using an SI Analytics GmbH glass filter porosity 2 to 4 (10 ...100 µm), whereas high-viscous samples should first pass through a sieve with a mesh width of 0.3 mm (test sieve 0.2 according to DIN 4188). Samples with a setting value or pour-point which is not at least 30 °C lower than the test temperature according to DIN 51 583 or DIN 51 597, respectively, have to be heated up to 50 °C prior to the measurement

- Fill the measurement liquid through the filling tube (3) into the storage container.
→ Viscometer operating instructions.

3.3 Placing the viscometer into the Viscosity Measuring Unit ViscoClock

- Press the fixing plate (10) upwards against the handle plate (12), → Exhibit 3.

The following applies to Ubbelohde viscometers:

- Insert the viscometer from below in an angle through the guidance plate (3), and then insert it vertically into the centering prism (5).
- Lower the fixing plate (10).
- Check proper seating of the silicone plug on filling tube (3) (silicone plug 11D) and capillary tube (2) (silicone plug 11S).
- The viscometer is fastened, → Fig. 4.

The following applies to Micro Oswald viscometers:

- Insert the Micro Oswald viscometer into the stand (6) of the Viscosity Measuring Unit ViscoClock in such a manner that the capillary tube is located as is shown in Exhibits 3 + 4, but the filling tube is inserted into the oblong hole of the guidance plate (9). In this position the filling tube of the viscometer is located under the silicone membrane (17). The membrane rod (16) is located in its upper position - filling tube is open.
- Lower the fixing plate (10).
- The hand pump "sucking" VZ 6554 is a necessary prerequisite for Micro Oswald viscometers.

Fig. 3



Fig. 4


3.4 Insertion of the Viscosity Measuring Unit ViscoClock into the thermostat bath

The transparent thermostat has to be equipped with a VZ 5402 manual measuring insert in order to accommodate the Viscosity Measuring Unit ViscoClock.

- The Viscosity Measuring Unit ViscoClock, with filled viscometer, is inserted into the manual measuring insert {22} of the thermostat bath.

Fig. 5

3.5 Connecting the Viscosity Measuring Unit ViscoClock

CAUTION

Prior to plugging the power supply (29) into the socket, make sure that the mains-unit voltage matches the mains voltage.

- Plug the 9 V/ DC jack plug (19) into the back panel of the electronic measuring unit
- Plug the TZ 1858 universal voltage power supply for 100 - 240 V~, 50/60Hz, into the mains socket.

3.6 Connecting the Viscosity Measuring Unit ViscoClock to a computer (PC) or printer

There is an RS-232-C Interface (18) on the rear side of the electronic measuring unit (13).


- If desired, connect a computer (PC) or printer for documentation purposes.

Plug connection: 4-channel mini DIN round plug.

Required data cable, which is not included in the scope of delivery, → accessories, chapter 5.

Fig. 6 The transfer parameters of the RS-232-C Interface are firmly set:


- - 4800 baud
- - 7 bit word length
- - 2 stop bits
- - no parity

3.7 Preparing the Viscosity Measuring Unit ViscoClock for measurement


- Allow for the tempering time, depending on the application.
- Push membrane rod (16) downwards.
The silicone membrane (17) shuts the venting tube (1) (applies only to Ubbelohde viscometers).

Fig. 7

- Manually place the adapter (21) onto the silicone plug (11S), either
 - suck at the capillary tube (2) through the silicone plug 11D) (→ exhibit 8) or
 - push at the filling tube (3) through the silicone plug (11D) until the upper reservoir of the viscometer is at least half filled (→ viscometer operating instructions).

Fig. 8

3.8 Starting the Viscosity Measuring Unit ViscoClock


Remove the adapter (21) from the manual pump; considering that the measurement liquid will immediately start flowing down, the <START> key (15) has to be pressed immediately. Thus

the display will be set to zero.

the membrane rod (16) will be unlocked and the filling tube released (1) (applies only to Ubbelohde viscometers).

the measurement is triggered.

Fig. 9

NOTE

Do not touch the Viscosity Measuring Unit ViscoClock during measurement since the shock may cause erroneous measurements!

When the meniscus of the measurement liquid flows

- through the upper measurement mark (8), the time measurement starts, and when it flows
- through the lower measurement mark (7), the time measurement ends.

3.9 Reading the flow time

The Viscosity Measuring Unit ViscoClock takes the time as soon as the meniscus of the measurement liquid flows through the lower measurement mark (7).

- Read the flow time in terms of seconds on the display (14), if a computer (PC) or printer is connected, you can furthermore document the measurement value,

The displayed measurement value remains on the display until

- the start key (15) is pressed again or
- the power supply is interrupted,

Fig. 10

NOTE

The flow time is evaluated according to → chapter 3.22, 3.12, operating instructions for SI Analytics GmbH viscometers, and if required taking into account the tables of the Hagenbach corrections.

CAUTION

According to DIN, the next measurement may not be started until the capillaries are drained off. Otherwise there is a danger of formation of bubbles, hence of erroneous measurements!

3.10 Turning the Viscosity Measuring Unit ViscoClock off

- You turn the Viscosity Measuring Unit ViscoClock off by unplugging the power supply (20) from the mains socket.

3.11 Viscosity calculation

Concerning the calibration of a viscometer, please read the sheet with the heading "Viscometer within quality assurance systems" which is enclosed with each SI Analytics GmbH viscometer.

Viscosity is a "classical" parameter of liquids. The result obtained using capillary viscometry is kinematic viscosity, measured in the unit of mm²/s (formerly Centistokes, cSt). Considering that viscosity is highly depending on temperature, the measurement temperature has to be maintained accurately and included in the result (± 0.02 K).

The number of seconds identified for the various capillaries in the tables for Hagenbach corrections have to be subtracted from the determined flow time. Intermediate values may be interpolated.

In the case of absolute measurements the corrected flow time will directly result in kinematic viscosity in (mm²/s) by multiplying the corrected flow time by the constants identified on the viscometer (determined at 23 °C).

$$\nu = K (t - \Delta t_H)$$

Relative viscosity is the quotient of the viscosities of a solution and the pure solvent. It is especially of major importance in the food industry and the evaluation of large polymers. After Staudinger detected the contexts between the average chain length and viscosity, polymers were increasingly evaluated using viscosity determination. In practice, relative viscosity is used for the calculatory derivation of further parameters, such as specific and inherent viscosity and the viscosity J number and the k-value.

3.12 Evaluation example

Ubbelohde viscometer type no. 532 10, calibrated for automatic measurements.

Capillaries no. I

Constant = 0.009987

Flow time (averaged) = 180,00 s

Hagenbach correction for 180.00 Δt_H = 0.30 s ($\rightarrow$ correction - table in viscometer operating instructions)

Kinematic viscosity ν = $K (t - \Delta t_H)$
 = $0.009987 \cdot (180.00 - 0.30)$
 1.795 mm²/s

4 Maintenance

When using alkaline laboratory cleaners there is a danger to health (caustification, damage to skin and eyes)! Please observe the memorandums of the employer's liability insurance associations and the safety data sheets of the manufacturers.

NOTE

The cleaners may affect the glass of the viscometer.

- Please use personal protective equipment such as eye protection, protective gloves, laboratory coat and breath protection.

4.1 Purification

- After each use you should clean the stand (6) using a moist cloth and a commercially available cleaning agent.
- Use a fine, dry brush to clean the upper (8) and lower (7) measurement plane (transmitter and receiver are opposed to each other). Accurate measurement is only possible with clean measurement planes.
- The viscometer should be cleaned according to DIN 51 562 part 1. Subsequently, the viscometer must be stored dry and free of dust!

4.2 Inspection

- Check cable and plug contacts for mechanical damage and corrosion.

4.3 Repair

- For repair, please send the Viscosity Measuring Unit ViscoClock to SI Analytics GmbH.

4.4 Purging and first-time filling

At the end of its useful life the Viscosity Measuring Unit ViscoClock has to be disposed of according to the national regulations.

- When disposing components of the Viscosity Measuring Unit ViscoClock please note:
 - Metal parts should be channelled towards recycling
 - Plastic parts should be channelled towards recycling
 - Electric/electronic components should be channelled to be reused in electronic devices.

5 Elimination of failures

Please note: To ensure reliable functionality, the ViscoClock should be connected to the power supply for at least 3 minutes prior to taking test measurements.

The display is dark.

<i>Power supply has been disrupted.</i>	Check connection/function of plug.
---	------------------------------------

The specimen cannot be pushed or sucked up.

<i>The silicon stopper is not seated entirely in the filling tube (compression operation) or the capillary tube (suction operation).</i>	Check the silicon stopper for tight seating.
<i>The hand pump is not seated completely on the silicon stopper.</i>	Set the glass tube from the hand pump tightly on the silicon stopper.
<i>For suction operating mode: The silicon membrane is not on the ventilation tube.</i>	Push the membrane rod with the silicon membrane down while checking the seating of the silicon membrane on the ventilation tube.

The specimen comes up or out when pushed down in the ventilation tube.

<i>The membrane rod (release lever) with the silicon membrane is not pushed down.</i>	Push the membrane rod with the silicon membrane down.
<i>Membrane rod is pushed down, but the silicon membrane is not uniform or does not seal the ventilation tube.</i>	Check the seating of the silicon membrane on the membrane rod and pull down some if needed.

The measurement was not started or not stopped when the meniscus passed through the light barrier. The timer is not running properly or not running at all.

<i>Start button was not pressed hard enough: The membrane rod triggers and opens the ventilation tube, but the timer was not set to zero and thus was not started for the measurement.</i>	Start button must be pressed harder.
<i>The measuring time is too short.</i>	Measuring time: The measuring time should be more than 30 seconds. This is the requirement for a reliable measurement recording. To prevent increased measurement uncertainty, longer measurement times (> 100 seconds) are recommended for narrow capillaries. Measuring time must be extended by choosing another viscometer.
<i>The sample liquid is not sufficiently transparent. Check the ViscoClock with a transparent specimen.</i>	Check the ViscoClock with a transparent specimen.
<i>Check the functionality of the measurement time recorder.</i>	- Take the viscosimeter from the device. - Push down the ventilation membrane rod (release lever). - Press start button. - Run a thin object (e.g., pencil, pen) through the upper light barrier. Then start the time measurement. - The time measurement must stop again when the lower light barrier is passed through.

6 Accessories and spare parts

Accessories: TZ 3095 1.5 m data cable RS-232-C,

with plug for connecting a printer:
1st end: 4-channel mini DIN plug
2nd end: 25-channel sub-Miniature-D socket

TZ 3096 1.5 m data cable RS-232-C,

for computer (PC) with 25-channel plug:
1st end: 4-channel mini DIN plug
2nd end: 25-channel sub-Miniature-D socket

TZ 3097 1.5 m data cable RS-232-C,

for computer (PC) with 9-channel plug:
1st end: 4-channel mini DIN plug
2nd end: 9-channel sub-Miniature-D socket

TZ 3098 data cable set:
1.5 m data cable RS-232-C,

1st end: 4-channel mini DIN plug
2nd end: 9-channel sub-Miniature-D socket
9-channel plug → 25-channel socket

VZ 6550

manual pump „press“

VZ 6554

manual pump „suck“

TZ 1858

universal voltage power supply 100 – 240 V ~
9 V / DC, 550 mA

TZ 1852

primary adapter for TZ 1858, for UK

TZ 1854

primary adapter for TZ 1858, for Australia

CT 72/P or CT 72/ 2 or CT 72/4

transparent thermostat of the CT 72 product range

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