SCHOTT[®] Instruments

**Halbzelle
Pb 1100 A**

**Einstabmesskette
PB 60**

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PB 60 Lead Combination Electrode

**Accuracy when
going to press**

The use of advanced technology and the high quality standard of our products are the result of a continuous development. This may result in differences between this operating manual and your electrode. Also, we cannot guarantee that there are absolutely no errors in this manual. Therefore, we are sure you will understand that we cannot accept any legal claims resulting from the data, figures or descriptions.

**Warranty
declaration**

In case of proper handling, SI ANALYTICS guarantees the quality of the electrode for one year.

Copyright

© Mainz 2012, SI Analytics GmbH
Reprinting - even as excerpts - is only allowed with the explicit written authorization of SI Analytics GmbH, Mainz.
Printed in Germany

View Halfcell SCHOTT® Instruments Pb 1100 A

View Combination electrode SCHOTT® Instruments PB 60

Commissioning

Combination Electrode PB 60

1	Remove the protection cap.
2	Pull the closing ring downwards so that the refilling opening of the reference electrode is free.
3	Fill ELY/BR/503 electrolyte solution into the refilling opening.
4	Wipe the shaft using a clean paper towel.
5	Place the combination electrode into deionized water for approx. 10 minutes.
6	Before use, place the combination electrode into diluted standard solution for approx. 5 minutes. In order to achieve a short response time, the concentration of the standard solution should be slightly below the measuring range expected.

The combination electrode is now ready to measure. Always leave the refilling opening of the reference electrode open while measuring.

Halfcell Pb 1100 A

For measurements using the Pb 1100 A lead electrode, an additional reference electrode is required. The two electrodes together form a combination electrode.

1	Remove the protection cap.
2	Place the electrode into deionized water for approx. 10 minutes.
3	Before use, place the electrode into diluted standard solution for approx. 5 minutes. In order to achieve short response time, the concentration of the standard solution should be slightly below the measuring range expected.
4	Fill the reference electrode with ELY/BR/503 electrolyte solution.

The combination electrode is now ready to measure.

Note

The reference electrode must contain enough electrolyte solution, so that the level of the electrolyte solution is clearly higher than that of the test sample during the measurement.

Calibration

Calibrate according to the data given in the operating manual of the measuring instrument and the analysis specification.

When operating the electrode ensure that:

- the filling opening for the bridge electrolyte is open
- the inner junctions are covered with bridge electrolyte
- no air bubbles are in the bridge electrolyte
- the depth of immersion is within the optimum range:

Minimum depth of immersion:

The diaphragm must be covered

Maximum depth of immersion:

Approx. 1 cm below the fluid level of the bridge electrolyte

Before Measuring

- Before use, condition the combination electrode or electrodes respectively for approx. 5 min. in standard solution.
- Remove any air bubbles in the bridge electrolyte by slightly knocking against the shaft.
- Calibrate according to the operating manual of the meter and the analysis specification.

Sample preparation

Add 2 % ISA/FK solution. This sample conditioning solution creates optimum conditions for measuring. It provides a constant ionic strength and similar diffusion potentials at the reference electrode in standard and test sample.

Note

If you would like to have more detailed information concerning sample preparation and measuring procedures, SI Analytics provides a large number of application reports for various applications.

Response times

The response time depends on the concentration range.

It is:

- several seconds with high concentrations,
- several minutes near the detection limit.

The measured value is stable if the value does not change by more than 0.1 mV within 30 seconds.

Interferences

- Cupric, mercury and silver ions poison the membrane.
- Iron and cadmium ions poison the membrane if their concentrations are higher than that of the lead concentration.

Aging

Please note that every electrode undergoes a natural aging process. The response time increases and the slope decreases with the age of the electrode. The following factors shorten the lifetime considerably:

- Incorrect storing
- Special measuring conditions (e.g. organic solutions, frequent measuring with high concentrations of interfering ions)
- High temperatures
- High changes in temperature

The warranty does not cover damage caused by measuring conditions and mechanical damage.

Maintenance

- For combination electrodes: replenish used reference electrolyte solution.
- In case of increased response times of the electrode perform the following activities:
 - Polish the membrane surface using an SF/K polishing strip. To do so, moisten a polishing trip with deionized water, and with circular movements polish the membrane surface for approx. 30 seconds.
 - Subsequently place the electrode into diluted standard solution for approx. 5 minutes.

Storage


Combination electrode PB 60

- Between two measurements, place the combination electrode into diluted standard solution.
- Overnight up to one week: Push the closing ring upwards to close the refilling opening. Then rinse the combination electrode with deionized water, blot it dry using a clean paper towel and store it in a dry place.
- For more than a week: Remove the electrolyte solution from the reference electrode and rinse the combination electrode with deionized water, blot it dry using a clean paper towel and put on the protection cap. Store the combination electrode in a dry place.

Halfcell Pb 1100 A

- Between two measurements, place the electrode into diluted standard solution.
- For storing overnight or longer, rinse the electrode with deionized water, blot it dry using a clean paper towel and put on the protection cap (closing cap of the glass tube). Store the electrode in a dry place.

Typical calibration line of a lead electrode

Technical data
SCHOTT® Instruments ISE – Halfcell Pb 1100 A /
Combination Electrode PB 60

Measuring range 0.2 ... 20000 mg/l Pb²⁺ (10⁻⁶ ... 10⁻¹ mol/l Pb²⁺)

Reproducibility ± 4 %

pH range 4 ... 7 (see section INTERFERENCES)

Temperature range 0 ... 80 °C / up to 100 °C for a short time

**Membrane
resistance** < 1 MΩ

Length Pb 1100 A: 170 mm (including 50 mm connection head)
PB 60: 153 mm (including 33 mm connection head)

Diameter Shaft: 12 mm
Connection head: 16 mm

Error diagnostics

	Cause	Remedy
Measured value unstable	Inner junctions not sufficiently wetted with bridge electrolyte	Fill up bridge electrolyte until the inner junctions are covered with bridge electrolyte
	Inner junctions encrusted	Leave the bridge electrolyte to react on the inner junctions for some hours until the crusts have dissolved
	Cable broken	Exchange cable
Slope too low	Membrane surface contaminated	Polish the membrane surface with polishing strip S/FK (see MAINTENANCE)
	Conditioning time too short	Extend conditioning time
	Standard solutions too old	Use new standard solutions
	Inner junctions encrusted	Leave the bridge electrolyte to react on the inner junctions for some hours until the crusts have dissolved.
	Electrode defective	Exchange electrode

Accessories

Description	Type No.	Order no.
Standard solution 10 g/L Br ⁻	ES/Br	120120
Bridge electrolyte	ELY/BR/503	106575
ISA sample conditioning solution	ISA/FK	140110
Polishing strips (24 pieces)	S/FK	180130

SI Analytics

a **xylem** brand

SI Analytics GmbH

Hattenbergstrasse 10

55122 Mainz

Germany

Phone: +49 (0) 6131/66 5111

Telefax: +49 (0) 6131/66 5001

E-mail: support.si-analytics@xyleminc.com

Internet: www.si-analytics.com

T814847