


SCHOTT®
Instruments

ProLab 2000



**pH / ISE / conductivity / D.O. meter
with automatic sensor recognition
and electronic access control**

**Accuracy when
going to press**

The use of advanced technology and the high quality standard of our instruments are the result of a continuous development. This may result in differences between this operating manual and your meter. Also, we cannot guarantee that there are absolutely no errors in this manual. Therefore, we are sure you will understand that we cannot accept any legal claims resulting from the data, figures or descriptions.

Warranty

We guarantee the meter described for 3 years from the date of purchase.

The meter warranty covers manufacturing faults that are discovered within the warranty period.

The warranty does not cover components that are replaced during maintenance work, e. g. batteries.

The warranty claim extends to restoring the meter to readiness for use but not, however, to any further claim for damages. Improper handling or unauthorized opening of the meter invalidates any warranty claim.

To ascertain the warranty liability, return the meter and proof of purchase together with the date of purchase freight paid or prepaid.

**CE conformity
Radio data transmission**

SI Analytics GmbH hereby declares that the ProLab 2000 meter is in compliance with the basic requirements and the other relevant regulations of the directive 1999/5/EC.

The EC declaration of conformity can be requested from SI Analytics GmbH.

Copyright

© 2009, SI Analytics GmbH

Reprinting - even in the form of excerpts - is only allowed with the explicit written authorization of SI Analytics GmbH.

Printed in Germany.

**KONFORMITÄTSERKLÄRUNG
DECLARATION OF CONFORMITY
DÉCLARATION DE CONFORMITÉ**

Wir erklären in alleiniger
Verantwortung, dass das
Produkt

We declare under our sole
responsibility that the
product

Nous déclarons sous notre
seule responsabilité que le
produit

**pH-/ISE-/
Leitfähigkeits-/
Sauerstoff-
Messgerät
ProLab 2000**

**pH / ISE /
conductivity / DO
measuring
instrument
ProLab 2000**

**Appareil de mesure
pour pH/ISE/
conductivité/
oxygène
ProLab 2000**

auf das sich diese Erklärung
bezieht, übereinstimmt mit
den Angaben im Kapitel

to which this declaration
relates is in conformity with
the specifications in the
chapter

auquel se réfère cette
déclaration est conforme aux
indications du chapitre

Technische Daten

**pH-/ISE-/Leitfähigkeits-/Sauerstoff-Messgerät
ProLab 2000**

30. Oktober 2009

**SI Analytics GmbH
Hattenbergstr. 10
D-55122 Mainz
Deutschland, Germany, Allemagne**

30. Oktober, October 30, 30 octobre 2005
AGQSF 0000-A106-01/091030

1	Overview	9
1.1	General features	9
1.2	Keypad	10
1.3	Display	11
1.4	Socket field	12
1.5	Automatic sensor recognition	13
1.6	Electronic access control	15
2	Safety	17
2.1	Authorized use	17
2.2	General safety instructions	18
3	Commissioning	19
3.1	Scope of delivery	19
3.2	Power supply	19
3.3	Initial commissioning	20
4	Operation	21
4.1	Switch the meter on and off	21
4.2	General operating principles	22
4.2.1	Operating modes	22
4.2.2	Navigation	23
4.2.3	Navigation example 1: Setting the language	25
4.2.4	Navigation example 2: Setting the date and time	27
4.3	Access control	29
4.3.1	Administrating access authorizations	29
4.3.2	Lost your electronic key?	31
4.3.3	Lock	31
4.4	System settings (system menu)	32
4.4.1	<i>Data storage</i>	32
4.4.2	<i>Display</i>	32
4.4.3	<i>System</i>	34
4.4.4	<i>Automatic Stability control</i>	35
4.5	pH value / ORP voltage	36
4.5.1	General information	36
4.5.2	Measuring the pH value	37
4.5.3	Measuring the ORP	39
4.5.4	Settings for pH and ORP measurements	41
4.5.5	pH calibration	43
4.5.6	Calibration interval	47
4.5.7	Calibrating	48
4.6	Ion concentration	52
4.6.1	General information	52
4.6.2	Measuring the ion concentration	53
4.6.3	Settings for ISE measurements	54

4.6.4	Calibrating for ISE measurements	55
4.7	Conductivity	61
4.7.1	General information	61
4.7.2	Measuring	62
4.7.3	Temperature compensation	64
4.7.4	Settings for conductivity measuring cells	65
4.7.5	Determining the cell constant (calibration in control standard)	67
4.8	Dissolved oxygen	71
4.8.1	General information	71
4.8.2	Measuring	72
4.8.3	Settings for D.O. sensors	74
4.8.4	D.O. calibration	75
4.9	Storage	80
4.9.1	Manual storage	81
4.9.2	Automatic storage at intervals	82
4.9.3	Reading the measurement data storage	84
4.9.4	Erasing the data storage	87
4.9.5	Displaying and downloading calibration records	87
4.10	Transmitting data (to a PC or printer)	89
4.10.1	RS232 interface	89
4.10.2	USB interface (device)	90
4.10.3	Options for data transmission	91
4.10.4	Operation with MultiLab pilot	92
4.11	Reset	92
4.11.1	Resetting the sensor settings	92
4.11.2	Resetting the system settings	94
5	Maintenance, cleaning, disposal	95
5.1	Maintenance	95
5.2	Cleaning	96
5.3	Disposal	96
6	What to do if...	97
6.1	pH and ORP measurement	97
6.2	ISE measurement	99
6.3	Conductivity measurement	100
6.4	D.O. measurement	100
6.5	General errors	102
7	Technical data	103
7.1	General data	103
7.2	Measuring ranges, resolution, accuracy	104
7.2.1	pH/ORP	104
7.2.2	ISE	105
7.2.3	Conductivity	105

7.2.4	Dissolved oxygen	107
8	Lists	109
9	Index	115
	Appendix	117
A.1	Firmware update	117
A.2	Menus	119
A.2.1	<i>Data storage</i>	119
A.2.2	<i>Display</i>	120
A.2.3	<i>System</i>	120
A.2.4	<STO __>	121
A.2.5	pH/U	121
A.2.6	ISE	122
A.2.7	Conductivity	122
A.2.8	Oxygen	122

1 Overview


1.1 General features

The ProLab 2000 precision meter enables you to perform pH, ORP, conductivity, dissolved oxygen (D.O.) and ion-selective measurements rapidly and reliably.

The ProLab 2000 provides the maximum degree of ease of use, reliability and, above all, measuring certainty for all applications.


The proven calibration procedures, and stability control function (SC) and the sensor recognition function support your work with the meter.

In addition, the ProLab 2000 provides an electronic access control. Documented measurement data is thus automatically assigned to a user.






1	Keypad
2	Display
3	Reader field for electronic access control
4	Socket field

1.2 Keypad




In this operating manual, keys are indicated by brackets <.> .
 The key symbol (e.g. <MENU/OK>) generally indicates a short keystroke (under 2 sec) in this operating manual. A long keystroke (approx. 2 sec) is indicated by the underscore behind the key symbol (e.g. <MENU/OK__>).

Key functions	Key Symbol	Function
		<On/Off> Switch the meter on/off
		<MODE> Select measured parameter <MODE__> Activate operation lock
		<CAL> Call up calibration procedure <CAL__> Open menu for calibration data storage
		<SC> Switch on or off the stability control function manually.
		<STO> Store measured value <STO__> Open menu for automatic storing function
		<RCL> Open menu for manually stored measured values <RCL__> Open menu for automatically stored measured values
		<PRINT> Print <PRT__> Open menu for automatic printing
		<▲> Increment values, scroll


	<▼>	Decrement values, scroll
	<MENU/OK> <MENU/OK_>	Confirm entries Open setting menu for system settings
	<ESC>	Return to higher menu level / Cancel inputs

1.3 Display

The graphic display can display three measured parameters and the current temperature value for each measured parameter at the same time. The illumination enables to read the display even in the darkness.



1.4 Socket field



Connections:

1	Conductivity measuring cell / D.O. electrode
2	pH ions electrode
3	Reference electrode
4	Temperature sensor
5	USB interface
6	Power pack
7	RS232 interface/analog output



CAUTION


Only connect sensors to the meter that cannot return any voltages or currents that are not allowed (> SELV and > current circuit with current limiting).

Nearly all commercial sensors - especially SI Analytics sensors - fulfill these requirements.

1.5 Automatic sensor recognition

The automatic sensor recognition function enables

- to operate a sensor with different meters without recalibrating
- to operate different sensors with one meter without recalibrating
- to assign measurement data to a sensor
 - measurement datasets are always downloaded to the interface with the sensor type and sensor series number.
 - measurement datasets are always stored together with the sensor type and sensor series number.
- to assign calibration data to a sensor
 - calibration data is always downloaded to the interface with the sensor type and sensor series number.
- the automatic activation of the cell constants for conductivity sensors



To be able to use the automatic sensor recognition function a meter is required that supports the automatic sensor recognition (e.g. ProLab 2000), and a sensor (ID sensor) that is suitable for sensor recognition.

In every ID sensor, sensor data is stored that clearly identifies the sensor.

The sensor data is sent to the meter automatically via radio communication and used for sensor identification there.




Note

With the ProLab 2000 meter, you can also operate non-ID sensors. Then, however, you cannot use the advantages of the sensor recognition function.

ID sensors


SI Analytics GmbH ID sensors support the automatic sensor recognition function. Their sensor designation has the addition "ID", e.g. A 161 1M DIN-ID electrode.



Note

Information on available ID sensors is given on the Internet or directly by SI Analytics.

ID sensors connected to the ProLab 2000 meter can be recognized by the ID sensor symbol on the display of the meter.



Sensor data from ID sensors

ID sensors transmit the following sensor data to the meter:

- SENSOR ID
 - Sensor type
 - Sensor series number
- Calibration data
 - Calibration date
 - Calibration characteristics
 - Calibration interval
 - Selected buffer set (pH electrodes only)

The calibration data is updated in the ID sensor after each calibration procedure. The ID sensor symbol flashes while this is done.



Note

While the ID sensor symbol is flashing, the sensor must not be disconnected because otherwise the calibration data will not be completely transmitted. The sensor will then have no valid calibration.



Note

If non-ID sensors are used, the calibration data from the meter is used and also stored in the meter.

1.6 Electronic access control


The ProLab 2000 always documents measurement data and calibration data with an ID number (*USER ID*).

If the *USER ID* is assigned to a certain user, all measurement data is also assigned to the user in compliance with GLP.

The *USER ID* is easily and safely transmitted to the ProLab 2000 via an electronic key. Each electronic key in the form of a keyring pendant contains a key number. The key number is read by the meter via a contactless radio data connection and used as the *USER ID*.

The meter checks the access authorization for the *USER ID*. Measurements with an electronic key are only possible if the *USER ID* is registered in the meter.

The measuring data is then documented along with the registered *USER ID*.



If the *USER ID* of the electronic key is not registered in the meter, access to the meter with this electronic key is denied.

Anonymous access is always possible. For anonymous access, the *USER ID 0* is used automatically.

By labeling measurement data without electronic key with *USER ID 0*, this data can be separated from the GLP-compliant documentation of measurements.

2 Safety

This operating manual contains basic instructions that you must follow during the commissioning, operation and maintenance of the meter. Consequently, all responsible personnel must read this operating manual before working with the measuring system. The operating manual must always be available within the vicinity of the meter.

Target group

The meter was developed for work in the laboratory. Thus, we assume that, as a result of their professional training and experience, the operators will know the necessary safety precautions to take when handling chemicals.

Safety instructions

The individual chapters of this operating manual use the following safety instruction to indicate various types of danger:



CAUTION

indicates instructions that must be followed precisely in order to avoid the possibility of slight injuries or damage to the meter or the environment.

Further notes



Note

indicates notes that draw your attention to special features.



Note

indicates cross-references to other documents, e.g. operating manuals.

2.1 Authorized use

This meter is authorized exclusively for measurements of the pH, ORP, conductivity, D.O. and ion-selective measurements in the laboratory. The technical specifications as given in chapter 7 TECHNICAL DATA (page 103) must be observed. Only the operation and running of the meter according to the instructions given in this operating manual is authorized. Any other use is considered **unauthorized**.

2.2 General safety instructions

This instrument is built and inspected according to the relevant guidelines and norms for electronic measuring instruments (see page 103).

It left the factory in a safe and secure technical condition.

Function and operational safety

The smooth functioning and operational safety of the meter can only be guaranteed if the generally applicable safety measures and the specific safety instructions in this operating manual are followed during operation.

The smooth functioning and operational safety of the meter can only be guaranteed under the environmental conditions that are specified in chapter 7 TECHNICAL DATA (page 103).

If the meter was transported from a cold environment to a warm environment, the formation of condensate can lead to the faulty functioning of the meter. In this event, wait until the temperature of the meter reaches room temperature before putting the meter back into operation.

Safe operation

If safe operation is no longer possible, the meter must be taken out of service and secured against inadvertent operation!

Safe operation is no longer possible if the meter:

- has been damaged in transport
- has been stored under adverse conditions for a lengthy period of time
- is visibly damaged
- no longer operates as described in this manual.

If you are in any doubt, please contact the supplier of the meter.

Obligations of the purchaser

The purchaser of this meter must ensure that the following laws and guidelines are observed when using dangerous substances:

- EEC directives for protective labor legislation
- National protective labor legislation
- Safety regulations
- Safety datasheets of the chemical manufacturers.

3 Commissioning


3.1 Scope of delivery

- ProLab 2000 laboratory meter
- Power pack
- 4 batteries 1.5 V Micro type AAA
- 1 electronic administrator key (as keyring pendant)
- 1 electronic key (as keyring pendant)
- Cover
- USB cable (Z875)
- CD-ROM with USB driver
- Operating manual

3.2 Power supply

The power pack supplies the meter with low voltage (9 V DC).

The batteries are only used to buffer the system time if the power supply is interrupted.



CAUTION

The line voltage at the operating site must lie within the input voltage range of the original power pack (see page 103).

CAUTION

Use original power packs only (see page 103).

1	Insert the plug into the socket of the meter.
2	Connect the original power pack to an easily accessible power outlet.

3.3 Initial commissioning

Perform the following activities:


- Insert the batteries (see page 95)
- Connect the power pack (see page 19).
- Switch on the meter (see page 21)
- Set the language (see page 25)
- Set the date and time (see page 27)
- Set up the access authorization for electronic keys (see page 29)

4 Operation


4.1 Switch the meter on and off

Switching on

- 1 Press the **<On/Off>** key.
The display shows the meter designation and software version. Subsequently, the access control display appears.



- 2 Start anonymous access with **<MENU/OK>**.
The meter switches to the measuring mode.
or
Place the electronic key on the read field and leave it there for approx. 2 seconds until the key has been recognized.
The display shows the *USER ID* that was read.



- 3 Confirm the displayed data with **<MENU/OK>**.
The meter switches to the measuring mode.

Switching off

- 1 Press the **<On/Off>** key.
The meter is switched off.

4.2 General operating principles

This section contains basic information on the operation of the ProLab 2000.

Operating elements, display

An overview of the operating elements and the display is given on page 10 and page 11.

Operating modes, navigation

An overview of the operating modes and navigation of the ProLab 2000 is given on page 22 and page 23.

4.2.1 Operating modes

The meter has the following operating modes:

- Measuring
Measurement data of one to three sensors appears in the measured value display
- Calibration
The course of a calibration with calibration information, functions and settings is displayed
- Storage in memory
The meter stores measuring data automatically or manually
- Transmitting data
The meter transmits measuring data and calibration records to the serial interface automatically or manually.
- Setting
The system menu or a sensor menu with submenus, settings and functions is displayed

4.2.2 Navigation

Measured value display

In the measured value display, you can

- select a measuring window with <▲><▼> and open the relevant measuring menu with <MENU/OK>.
- open the system menu with the sensor-independent settings with <MENU/OK__>.
- change the display in the selected measuring window (e. g. pH <-> mV) by pressing <MODE>.

Menus and dialogs

The menus for settings and dialogs contain further submenus. The selection is made with the <▲><▼> keys. The current selection is displayed in a frame.

- Submenus

The name of the submenu is displayed at the upper edge of the frame. Submenus are opened by confirming with <MENU/OK>.

Example:

System	
General	
Measurement	
Interface	
Clock function	
Service information	
Reset	

- Settings

Settings are indicated by a colon. The current setting is displayed on the right-hand side. The setting mode is opened with <MENU/OK>. Subsequently, the setting can be changed with <▲><▼> and <MENU/OK>.

Example:

System	
Language:	English
Beep:	On
Illumination:	On
Contrast:	48 %

- Functions

Functions are indicated by the name of the function. They are immediately carried out by confirming with **<MENU/OK>**.

Example: Display the *Calibration record* function.

pH

<i>Calibration record</i>	
<i>Buffer:</i>	<i>NIST/DIN</i>
<i>Calibration interval:</i>	<i>7 d</i>
<i>Unit for zero point</i>	<i>pH</i>
<i>Unit for slope:</i>	<i>mV/pH</i>
■ 2.00 4.01 7.00 10.01	

Messages

Information or operating instructions are designated by the ■ symbol. They cannot be selected.

Example:

pH

■ <i>Buffer series NIST/DIN</i>	
■ <i>Immerse sensor in buffer 1</i>	
<i>Set temperature: 25 °C</i>	
<i>Continue</i>	


**Note**

The principles of navigation are explained in the two following sections by reference of examples:


- Setting the language (see page 25)
- Setting the date and time (see page 27)

4.2.3 Navigation example 1: Setting the language


- 1 Switch on the meter (see page 21)



- 2 In the measured value display:
Open the system menu with **<MENU/OK __>**
(press for approx. 2 s).
The meter is in the setting mode.



- 3 Select the *System* submenu with **<▲><▼>**.
The current selection is displayed in a frame.
- 4 Open the *System* submenu with **<MENU/OK>**.



5 | Select the *General* submenu with <▲><▼>. The current selection is displayed in a frame.

6 | Open the *General* submenu with <MENU/OK>.

System	
Language:	English
Beep:	Off
Illumination:	On
Contrast:	48 %

7 | Open the setting mode for the *Language* with <MENU/OK>.

System	
Language:	English
Beep:	Off
Illumination:	On
Contrast:	48 %

8 | Select the required language with <▲><▼>.

9 | Confirm the setting with <MENU/OK>. The setting becomes active the next time the system menu is called up.

4.2.4 Navigation example 2: Setting the date and time

The meter has a clock with a date function. The date and time are indicated in the status line of the measured value display. This indication can be switched off. When storing measured values and calibrating, the current date and time are automatically stored as well.

The correct setting of the date, time and date format is important for the following functions and displays:

- Current date and time
- Calibration date
- Identification of stored measured values.

Therefore, check the time at regular intervals.



Note

After a drop of the voltage of the buffer batteries (empty batteries), the date and time are reset and have to be adjusted again.

Setting the date, time and date format

The data format can be switched from the display of day, month, year (*dd.mm.yy*) to the display of month, day, year (*mm/dd/yy* or *mm.dd.yy*).

1	In the measured value display: Open the system menu with <MENU/OK __> (press for approx. 2 s). The meter is in the setting mode.
2	Using <▲><▼> and <MENU/OK> , select and confirm the <i>System / Clock function</i> menu. The setting menu for the date and time appears.
3	Using <▲><▼> and <MENU/OK> , select and confirm the <i>Time</i> menu. The seconds are highlighted.

System	
Time:	14:53:40
Date:	26.04.06
Date format:	dd.mm.yy

4	Change and confirm the setting with <▲><▼> and <MENU/OK> . The minutes are highlighted.
---	--

5	Change and confirm the setting with <▲><▼> and <MENU/OK>. The hours are highlighted.
6	Change and confirm the setting with <▲><▼> and <MENU/OK>. The time is set.
7	If necessary, set the <i>Date</i> and <i>Date format</i> . The setting is made similarly to that of the time.
8	If necessary, select and set the <i>Date</i> with <▲><▼> and <MENU/OK>.
9	To make further settings, switch to the next higher menu level with <ESC>. or Switch to the measured value display with <MODE>. The meter is in the measuring mode.

4.3 Access control

When the meter is delivered, only the enclosed electronic administrator key has an access authorization for the meter.


The administrator can set up access authorizations for electronic keys (see below).

Anonymous access (*USER ID*: 0) is always possible.


4.3.1 Administrating access authorizations

The ProLab 2000 provides basic functions to administrate access authorizations. The administrator can set up, erase and display access authorizations. The administrating functions are available to the administrator immediately after registration only.

1	When the meter is switched on: Switch off the meter with <On/Off> .
2	Press the <On/Off> key. The display shows the meter designation and software version. Subsequently, the access control display appears.



3	Place the electronic administrator key on the read field and leave it there for approx. 2 seconds until the key has been recognized.. The <i>USER ID</i> and the functions for access control are displayed.
---	---



Adding new USER IDs

- 4 Using <▲><▼> and <MENU/OK>, select and confirm the menu item, *Register USER ID*. An info text for registration is displayed.
- 5 Place the electronic key on the reader field. The *USER ID* is displayed and registered.
- 6 Use <MENU/OK> to complete the registration.

Displaying registered USER IDs

- 7 Using <▲><▼> and <MENU/OK>, select and confirm the *USER ID* menu item.
- 8 Using <▲><▼>, display the registered *USER IDs* individually.
- 9 If necessary (e.g. for deleting), confirm a *USER ID* with <MENU/OK>.

Erasing a registered USER ID

- 10 Using <▲><▼> and <MENU/OK>, select and confirm the *Erase* menu item. A safety query appears. After confirming the erasure, the access authorization for the currently displayed key number is erased.

Outputting a list of the registered USER IDs

- 11 Using <▲><▼> and <MENU/OK>, select and confirm the menu item, *Output to RS232/USB*. The list of all key numbers with access authorization is downloaded to the interface.

Switching to the measuring mode

- 12 Confirm the data with <MENU/OK>. The meter switches to the measuring mode.

4.3.2 Lost your electronic key?

Without an electronic key with access authorization, anonymous access is possible only: (*USER ID* 0).

Electronic key

Keys for additional users are available from SI Analytics GmbH. The access authorization for new keys is set up by the administrator (see page 29).

Electronic administrator key

If the electronic administrator key is lost, the access authorization for a new electronic administrator key can be set up in the factory only.

The addresses and telephone numbers of SI Analytics GmbH are given on the cover of this operating manual.

4.3.3 Lock

The activated lock prevents the inadvertent usage of the meter or registered *USER ID* during running operation. The lock can only be released with the currently registered electronic key or the administrator key.



Note

The lock can only be activated in the measuring mode of operation. With anonymous access (*USER ID*: 0), the lock function is not available.

Activating the lock

1	If necessary, switch to the measuring mode with <MODE> .
2	Use <MODE__> to activate the lock. The registered <i>USER ID</i> is locked against inadvertent use.

Releasing the lock

1	Press any key.
2	Place the registered electronic key or administrator key on the reader field. Place the registered electronic user key or administrator key on the read field and leave it there for approx. 2 seconds until the key has been recognized. The lock is released.

4.4 System settings (system menu)

The system menu comprises the following settings:

- *Data storage* (see page 32)
- *Display* (see page 32)
- *System* (see page 34).

4.4.1 Data storage

This menu contains all functions to display, edit and erase stored measured values and calibration records.



Note

Detailed information on the data storage functions of the ProLab 2000 is given on page 80.

4.4.2 Display

With the aid of the *Display* submenu, you can modify the measured value display to meet your requirements. When doing so, you can display or hide the following elements:

- Individual measured parameters
- Date indication in the status line
- Time indication in the status line

Settings

In the measured value display, open the system menu with **<MENU/OK__>**. After completing the settings, switch to the measured value display with **<MODE>**.

Menu item	Setting	Description
<i>Display / pH:</i>	<i>On</i> <i>Off</i>	Display of the pH/ISE measuring window. This menu item is only visible if the combination electrode and an additional sensor are connected.
<i>Display / Cond:</i>	<i>On</i> <i>Off</i>	Display of the conductivity measuring window. This menu item is only visible if a conductivity measuring cell is connected.

Menu item	Setting	Description
<i>Display / Oxi:</i>	<i>On</i> <i>Off</i>	Display of the D.O. measuring window. This menu item is only visible if a D.O. sensor is connected.
<i>Display / Time:</i>	<i>On</i> <i>Off</i>	Display of the time in the system status line
<i>Display / Date:</i>	<i>On</i> <i>Off</i>	Display of the date in the system status line

4.4.3 System

Overview

The following sensor-independent meter features can be adjusted in the system menu/*System* and its submenus:

- Menu language
- Beep on keystroke
- Display illumination
- Display contrast
- Unit of the temperature display
- Data interface
- Clock and date function
- Function to reset all sensor-independent system settings to the default condition

Settings

In the measured value display, open the system menu with **<MENU/OK__>**. After completing the settings, switch to the measured value display with **<MODE>**.

Menu item	Setting	Description
<i>System / General / Language</i>	<i>Deutsch English (further)</i>	Select the menu language
<i>System / General / Beep</i>	<i>On Off</i>	Switch on/off the beep on keystroke
<i>System / General / Illumination</i>	<i>On Off</i>	Switch the display illumination on/off
<i>System / General / Contrast</i>	<i>0 ... 100 %</i>	Change the display contrast
<i>System / Measurement / Temperature unit</i>	<i>°C °F</i>	Temperature unit, degrees Celsius or degrees Fahrenheit. All temperatures are displayed with the selected unit.
<i>System / Measurement / Stability control</i>	<i>Auto Off</i>	Switch on or off the automatic <i>Stability control</i> (see page 35).
<i>System / Interface / Baud rate</i>	<i>1200, 2400, 4800, 9600, 19200</i>	Baud rate of the data interface

Menu item	Setting	Description
<i>System / Interface / Output format</i>	<i>ASCII CSV</i>	Output format for data transmission For details, see page 89
<i>System / Interface / Output header</i>		The header is output to the interface in csv format.
<i>System / Clock function</i>	<i>Time Date Date format</i>	Settings of time and date. For details, see page 27
<i>System / Service information</i>		Meter information: Series number and software version
<i>System / Reset</i>	-	Resets the system settings to the default values. For details, see page 94

4.4.4 Automatic *Stability control*

The automatic *Stability control* function continuously checks the stability of the measurement signal. The stability has a considerable impact on the reproducibility of measured values.

You can activate or switch off the automatic *Stability control* function (see page 34).

The automatic *Stability control* function is carried out:

- as soon as the measured value is outside the allowed stability range
- when you switch over between the measured parameters with **<MODE>**.

In this case, the measured parameter flashes on the display.

4.5 pH value / ORP voltage

4.5.1 General information

You can measure the following variables:

- pH value []
- ORP [mV]



Temperature measurement

CAUTION

When connecting an earthed PC/printer, measurements cannot be performed in earthed media as incorrect values would result. The RS232 and USB interfaces are not galvanically isolated.

For reproducible pH measurements, it is essential to measure the temperature of the test sample.

You have the following options to measure the temperature:

- Automatic measurement of the temperature by the temperature sensor (NTC30 or Pt1000) integrated in electrode.
- Measurement by an external temperature sensor.
- Manual determination and input of the temperature.

The display of the temperature indicates the active temperature measuring mode:

Temperature sensor	Resolution of the temp. display	Mode
yes	0.1 °C	Automatic with pH temperature sensor
-	1 °C	Manual


Preparatory activities

Perform the following preparatory activities when you want to measure:

1	Connect a pH or ORP electrode to the meter.
2	If necessary, select the pH or mV display with <MODE> .
3	Adjust the temperature of the solutions and measure the current temperature if the measurement is made without a temperature sensor.
4	Calibrate or check the meter with the electrode.

4.5.2 Measuring the pH value

1	Perform the preparatory activities (see page 36).
2	Immerse the pH electrode in the test sample.



3	Select the pH or mV display with <MODE> .
---	--

Stability control

The Stability control function checks the stability of the measurement signal. The stability has a considerable impact on the reproducibility of measured values. The display of the measured parameter flashes until a stable measured value is available.

Independent of the setting for automatic *Stability control* (see page 35) in the *System* menu, you can start the *Stability control* function manually at any time.

1	If necessary, select a channel (measured value) with <▲><▼> .
2	Freeze the measured value with <SC> . The display switches to the display of the selected channel as necessary. In the status line, [SC] is displayed.
3	With <MENU/OK> activate the <i>Stability control</i> function. [SC] flashes while the stability control is active. As soon as a stable measured value is recognized, the current measurement data is downloaded to the interface. The measured values of all channels meeting the criterion for stability control are marked by SC.



Note

You can terminate prematurely the *Stability control* function with **<MENU/OK>** manually at any time. If the *Stability control* function is terminated prematurely, the current measurement data is not downloaded to the interface.

- | | |
|---|---|
| 4 | Using <SC> or <MODE>, release the frozen measured value. The [SC] status display disappears. The display switches to the previous indication. |
|---|---|

Criteria With identical measurement conditions, the following applies:

Measured parameter	Reproducibility	Response time
pH value	Better than 0.01	> 15 seconds


4.5.3 Measuring the ORP



Note

ORP electrodes are not calibrated. However, you can check ORP electrodes using a test solution.

- | | |
|---|---|
| 1 | Perform the preparatory activities (see page 36). |
| 2 | Immerse the ORP electrode in the test sample. |



- | | |
|---|--|
| 3 | Select the mV display with <MODE> . |
|---|--|

Stability control

The Stability control function checks the stability of the measurement signal. The stability has a considerable impact on the reproducibility of measured values. The display of the measured parameter flashes until a stable measured value is available.

Independent of the setting for automatic *Stability control* (see page 35) in the *System* menu, you can start the *Stability control* function manually at any time.

- | | |
|---|--|
| 1 | If necessary, select a channel (measured value) with <▲><▼> . |
| 2 | Freeze the measured value with <SC> .
The display switches to the display of the selected channel as necessary.
In the status line, [SC] is displayed. |
| 3 | With <MENU/OK> activate the <i>Stability control</i> function.
[SC] flashes while the stability control is active.
As soon as a stable measured value is recognized, the current measurement data is downloaded to the interface. The measured values of all channels meeting the criterion for stability control are marked by SC. |



Note

You can terminate prematurely the *Stability control* function with **<MENU/OK>** manually at any time. If the *Stability control* function is

terminated prematurely, the current measurement data is not downloaded to the interface.

- 4 Using **<SC>** or **<MODE>**, release the frozen measured value. The [SC] status display disappears. The display switches to the previous indication.

Criteria With identical measurement conditions, the following applies:

Measured parameter	Reproducibility	Response time
ORP voltage	better than 0.3 mV	> 15 seconds

4.5.4 Settings for pH and ORP measurements

Overview The following settings are possible for pH and ORP measurements:

- Resolution
- *Calibration interval*
- Buffers for calibration
- *Unit for zero point*
- *Unit for slope*
- *Calibration record* (display)

Settings The settings are made in the measuring menu of the pH/ORP measurement. To open it, activate the relevant measuring window in the measured value display and press <MENU/OK>. After completing the settings, switch to the measured value display with <MODE>.

Menu item	Possible setting	Description
<i>Calibration / Calibration record</i>	-	Displays the calibration record of the last calibration.
<i>Calibration / Buffer</i>	<i>NIST/DIN TEC MERCK 1 MERCK 2 MERCK 3 DIN 19267 Mettler US Mettler EU Fisher 1 Fluka BS Radiometer Baker Metrohm Beckmann Hamilton DC Precisa</i>	Buffer sets to be used for pH calibration. For details, see page 43.
<i>Calibration / Calibration interval</i>	<i>1 ... 999 d</i>	<i>Calibration interval</i> for the pH electrode (in days). The meter reminds you to calibrate regularly by the flashing CalClock in the measuring window.
<i>Calibration / Unit for zero point</i>	<i>mV pH</i>	Unit for the zero point.

Menu item	Possible setting	Description
<i>Calibration /Unit for slope</i>	<i>mV/pH</i> <i>%</i>	Unit of the slope. The % display refers to the Nernst slope of -59.16 mV/pH ([determined slope/Nernst slope] x 100).
<i>Man. temperature</i>	<i>-20 ... +130 °C</i>	Entry of the manually determined temperature. For measurements without temperature sensor only.
<i>Reset</i>	-	Resets all sensor settings to the delivery condition (see page 92).
<i>High resolution pH</i>	<i>On</i> <i>Off</i>	Resolution of the pH display: <i>On</i> = 0.001 <i>Off</i> = 0.01
<i>High resolution mV</i>	<i>On</i> <i>Off</i>	Resolution of the mV display: <i>On</i> = 0.1 mV <i>Off</i> = 1 mV

4.5.5 pH calibration

Why calibrate?

pH electrodes age. This changes the zero point (asymmetry) and slope of the pH electrode. As a result, an inexact measured value is displayed. Calibration determines the current values of the zero point and slope of the electrode and stores them.

Thus, you should calibrate at regular intervals.

When do you have to calibrate?

- After connecting a non-ID electrode
- If the CalClock has expired and flashes

Buffer sets for calibration

You can use the buffer sets quoted in the table for an automatic calibration. The pH values are valid for the specified temperature values. The temperature dependence of the pH values is taken into account during the calibration.

No.	Buffer set *	pH values	at
1	DIN buffers according to DIN 19266 and NIST Traceable Buffers	1.679 4.006 6.865 9.180 12.454	25 °C
2	Technical buffers	2.000 4.010 7.000 10.011	25 °C
3	Merck1*	4.000 7.000 9.000	20°C
4	Merck2 *	1.000 6.000 8.000 13.000	20°C
5	Merck3 *	4.660 6.880 9.220	20°C
6	DIN 19267	1.090 4.650 6.790 9.230	25 °C
7	Mettler Toledo USA *	1.679 4.003 7.002 10.013	25 °C

No.	Buffer set *	pH values	at
8	Mettler Toledo TEC *	1.995 4.005 7.002 9.208	25 °C
9	Fisher *	2.007 4.002 7.004 10.002	25 °C
10	Fluka BS *	4.006 6.984 8.957	25 °C
11	Radiometer *	1.678 4.005 7.000 9.180	25 °C
12	Baker *	4.006 6.991 10.008	25 °C
13	Metrohm *	3.996 7.003 8.999	25 °C
14	Beckman *	4.005 7.005 10.013	25 °C
15	Hamilton Duracal *	4.005 7.002 10.013	25 °C
16	Precisa *	3.996 7.003 8.999	25 °C

* Brand names or trade names are trademarks of their respective owners protected by law (see page 114).



Note

The buffers are selected in the sensor menu (*Buffer* setting, see page 41).

Calibration points

Calibration can be performed using one to five buffer solutions in any order (single-point to five-point calibration). The meter determines the following values and calculates the calibration line as follows:

	Determined values	Displayed calibration data
1-point	<i>ASY</i>	<ul style="list-style-type: none"> ● Zero point = <i>ASY</i> ● Slope = Nernst slope (-59.16 mV/pH at 25 °C)
2-point	<i>ASY</i> <i>SLO</i>	<ul style="list-style-type: none"> ● Zero point = <i>ASY</i> ● Slope = <i>SLO</i>
3-point to 5-point	<i>ASY</i> <i>SLO</i>	<ul style="list-style-type: none"> ● Zero point = <i>ASY</i> ● Slope = <i>SLO</i> <p>The calibration line is calculated by linear regression.</p>

**Note**

You can display the slope in the units, mV/pH or % (see page 34).

Stability control

The calibration procedure automatically activates the stability control function.

The current measurement with stability control can be terminated at any time (accepting the current value).

Calibration record

When finishing a calibration, the new calibration values are displayed as an informative message (■ symbol) first. Then you can decide whether you want to take over these values of the new calibration or whether you want to continue measuring with the old calibration data. After accepting the new calibration values the calibration record is displayed.

Display of calibration data and download to interface

You can have the data of the last calibration displayed (see page 87). Subsequently, you can transmit the displayed calibration data to the interface, e. g. to a printer or PC, with the <PRINT> key.

**Note**

The calibration record is automatically transmitted to the interface after calibrating.

Sample record:

```

Printing date 26.04.06 16:13
ProLab 2000
Ser. no. 06249876
USER ID: 1234567890

CALIBRATION pH
Calibration date 20.04.06 10:14:03
A 161 1M DIN-ID
Ser. no. A062498765
USER ID: 1234567890

Cal. interval          7 d
NIST/DIN
Buffer 1                6.86
Buffer 2                9.18
Voltage 1               7.2 mV      26.3 °C
Voltage 2              -171 mV      26.3 °C
Slope                  -59.2 mV/pH
Zero point              -0.5 mV
Sensor                  +++
    
```



Note

If non-ID sensors are used, the sensor type and series number are not output in the calibration record.

Calibration evaluation

After calibrating, the meter automatically evaluates the calibration. The zero point and slope are evaluated separately. The worse evaluation of both is taken into account. The evaluation appears on the display and in the calibration record.

Display	Calibration record	Zero point [mV]	Slope [mV/pH]
	+++	-18 ... +18	-60.5 ... -57.5
	++	-22 ... +22	-57.5 ... -56
	+	-26 ... +26	-61 ... -60.5 or -56 ... -55
	-	-30 ... +30	-62 ... -61 or -55 ... -50
Clean the electrode according to the sensor operating manual			
CalError	CalError	< -30 or > 30	... -62 or ... -50
Eliminate the error according to chapter 6 WHAT TO DO IF... (page 97)			

Preparatory activities

Perform the following preparatory activities when you want to calibrate:

1	Connect a pH electrode to the meter. The pH measuring window is displayed.
2	Keep the buffer solutions ready. Adjust the temperature of the buffer solutions, or measure the current temperature, if you measure without a temperature sensor.

4.5.6 Calibration interval

The calibration interval and calibration evaluation are indicated on the display as the CalClock.

CalClock

The remaining time of the calibration interval is indicated by the segmented ring around the calibration evaluation. This segmented ring reminds you to calibrate regularly.

After the specified calibration interval has expired the outside circle of the CalClock flashes. It is still possible to measure.

**Note**

To ensure the high measuring accuracy of the measuring system, calibrate after the calibration interval has expired.

Setting the calibration interval

The calibration interval is set to 7 (d7) days in the factory. You can change the interval (1 ... 999 days):

1	Open the menu for measurement settings with <MENU/OK> .
2	In the <i>Calibration / Calibration interval</i> menu, set the calibration interval with <▲><▼> .
3	Confirm the setting with <MENU/OK> .
4	Exit the menu with <MODE> .


4.5.7 Calibrating

Make sure the *Buffer* setting is correctly set in the sensor menu (see page 41).


Use any one to five buffer solutions of the selected buffer set in ascending or descending order.

The *NIST/DIN* calibration is described below. With other buffer sets, other nominal buffer values are displayed. Apart from that, the procedure is identical.

- | | |
|---|---|
| 1 | Press <▲><▼><MODE> to select the measured parameter pH or mV in the measured value display. |
| 2 | Start the calibration with <CAL>. The calibration display appears. |



- | | |
|---|---|
| 3 | Immerse the electrode in buffer solution 1. |
| 4 | If the <i>Set temperature</i> menu item appears, measure and enter the temperature of the buffer manually (measurement without temperature sensor). |
| 5 | Select <i>Continue</i> with <▲><▼> and press <MENU/OK>. The buffer is measured. The measured value is checked for stability (Stability control). |




Recognized nominal buffer value (referring to 25 °C)

Current electrode voltage

Current temperature value

6 Wait for the end of the measurement with Stability control or accept the calibration value with **<MENU/OK>**. The calibration display for the next buffer appears.



7 For single-point calibration, select *Exit with one point calibration* with **<▲><▼>** and confirm with **<MENU/OK>**. The calibration is completed as a single-point calibration. The new calibration values are displayed as a message (■). You have the following options:

- Accept the new calibration values with **<MENU/OK>**. Subsequently, the calibration record is displayed and downloaded to the interface at the same time.
- To switch to the measured value display without accepting the new calibration values, press **<MODE>** or **<ESC>**.




Note


For **single-point calibration**, the meter uses the Nernst slope (-59.16 mV/pH at 25 °C) and determines the zero point of the electrode.

Continuing with two-point calibration
(*Buffer: NIST/DIN*)

- | | |
|----|---|
| 8 | Thoroughly rinse the electrode with distilled water. |
| 9 | Immerse the electrode in buffer solution 2. |
| 10 | If the <i>Set temperature</i> menu item appears, measure and enter the temperature of the buffer manually (measurement without temperature sensor). |
| 11 | Select <i>Continue</i> with <▲><▼> and press <MENU/OK> . The buffer is measured. The measured value is checked for stability (Stability control). |



12 Wait for the measurement with stability control to be completed or *Terminate stability control* with <MENU/OK> and accept the calibration value.
The calibration display for the next buffer appears.



13 For two-point calibration, select *Exit with 2 point calibration* with <▲><▼> and confirm with <MENU/OK>. The calibration is completed as a two-point calibration. The new calibration values are displayed as a message (■). You have the following options:

- Accept the new calibration values with <MENU/OK>. Subsequently, the calibration record is displayed and downloaded to the interface at the same time.
- To switch to the measured value display without accepting the new calibration values, press <MODE> or <ESC>.

Continuing with three-point to five-point calibration (Buffer NIST/DIN)

14 Thoroughly rinse the electrode with distilled water.

15 Immerse the electrode in buffer solution 3.

16 If necessary, measure the temperature of buffer 3 manually, then enter and confirm it with <▲><▼> and <MENU/OK> in the *Set temperature* setting.

17 Select *Continue* with <▲><▼> and press <MENU/OK>. The buffer is measured. The measured value is checked for stability (Stability control).

pH <CAL>

- Buffer value = 4.010
- U = 184.0 mV
- Temperature = 24.8 °C

Terminate with ...

Terminate stability control

- 18 Switch to measurement of the next buffer with <MENU/OK>. or
When all buffers have been measured, use <▲><▼> to select the *Terminate with ...* menu item and confirm with <MENU/OK>.
- The calibration procedure is terminated.
The new calibration values are displayed as a message (■).
You have the following options:
- Accept the new calibration values with <MENU/OK>. Subsequently, the calibration record is displayed and downloaded to the interface at the same time.
 - To switch to the measured value display without accepting the new calibration values, press <MODE> or <ESC>.

**Note**

After a five-point calibration the calibration procedure is automatically terminated. The menu item *Terminate with ...* is not displayed.

**Note**

The calibration line is calculated by linear regression.

4.6 Ion concentration

4.6.1 General information



Note

Incorrect calibration of ion sensitive electrodes will result in incorrect measured values. Calibrate regularly before measuring.



CAUTION

When connecting an earthed PC/printer, measurements cannot be performed in earthed media as incorrect values would result. The RS232 and USB interfaces are not galvanically isolated.

Temperature measurement with ISE measurements

For reproducible measurements of the ion concentration, it is essential to measure the temperature of the test sample.

You have the following options to measure the temperature:

- Measurement of the temperature by the integrated temperature sensor of a simultaneously connected D.O. sensor or conductivity measuring cell in the test sample.
- Manual determination and input of the temperature.

The meter recognizes whether a suitable sensor is connected and automatically switches on the temperature measurement.

The display of the temperature indicates the active temperature measuring mode:

Temperature sensor, Cond or Oxi	Resolution of the temp. display	Mode
-	1 °C	Manual
yes	0.1 °C, measured temperature value flashes	The temperature value of the second sensor (Cond or Oxi) in the same test sample is taken over for measurement*

* If this is not required you can unplug the second sensor and enter the temperature manually.

Preparatory activities

Perform the following preparatory activities when you want to measure:


1	Connect the ISE combination electrode to the meter. The pH/ISE measuring window is displayed.
2	If necessary, select the ISE display (unit, mg/l) with <MODE> .
3	Measure the temperature of the test sample using a thermometer.
4	Calibrate or check the meter with the electrode.

**Note**

While no valid calibration is available, e.g. in the delivery condition, "Error" appears in the measured value display.

4.6.2 Measuring the ion concentration

1	Perform the preparatory activities according to page 52.
2	Immerse the electrode in the test sample.

**Stability control**

The Stability control function checks the stability of the measurement signal. The stability has a considerable impact on the reproducibility of measured values. The display of the measured parameter flashes until a stable measured value is available.

Independent of the setting for automatic *Stability control* (see page 35) in the *System* menu, you can start the *Stability control* function manually at any time.

1	If necessary, select a channel (measured value) with <▲><▼> .
---	--

- 2 | Freeze the measured value with **<SC>**.
The display switches to the display of the selected channel as necessary.
In the status line, [SC] is displayed.
- 3 | With **<MENU/OK>** activate the *Stability control* function.
[SC] flashes while the stability control is active.
As soon as a stable measured value is recognized, the current measurement data is downloaded to the interface. The measured values of all channels meeting the criterion for stability control are marked by SC.



Note

You can terminate prematurely the manual *Stability control* function with **<MENU/OK>** manually at any time. If the manual *Stability control* function is terminated prematurely, the current measurement data is not downloaded to the interface.

- 4 | Using **<SC>** or **<MODE>**, release the frozen measured value.
The [SC] status display disappears. The display switches to the previous indication.

Criteria

With identical measurement conditions, the following applies:

Measuring signal	Reproducibility	Response time
Electrode voltage	better than 0.1 mV	> 30 seconds

Temperature while calibrating and measuring

For precise ISE measurements the temperature difference between measurement and calibration should not be greater than 2 K. Therefore, adjust the temperature of the standard and measuring solutions accordingly. If the temperature difference is greater the *[TempErr]* warning appears in the measured value display.

4.6.3 Settings for ISE measurements

Overview

The following settings are possible for ISE measurements:

- *Calibration record* (display)

Settings

The settings are made in the measuring menu of the ISE measurement. To open it, activate the relevant measuring window in the measured value display and press the **<MENU/OK>** key. After completing the settings, switch to the measured value display with **<MODE>**.

Menu item	Possible setting	Description
<i>Calibration / Calibration record</i>	-	Displays the calibration record of the last calibration.
<i>Man. temperature</i>	-20 ... +130 °C	Entry of the manually determined temperature. For measurements without temperature sensor only.

4.6.4 Calibrating for ISE measurements

Why calibrate?

Ion-selective electrodes age and are temperature-dependent. This changes the slope. As a result, an inexact measured value is displayed. Calibration determines the current value of the slope of the electrode and stores it in the instrument.

Thus, you should calibrate before each measurement and at regular intervals.

When to calibrate?

- Before any measurement if possible
- After connecting another ISE electrode
- When the CalClock flashes, e. g. after a voltage interruption

Standard solutions

Use two or three different standard solutions. For three-point calibration, the standard solutions have to be selected in either increasing or decreasing order.

Standard solution	Values [mg/l]
Std 1	0.01; 0.02; 0.05; 0.1; 0.2; 0.5; 1; 2; 5; 10; 20; 50; 100; 200; 500; 1000
Std 2	If Std 2 > Std 1, Std 3 must be > Std 2
Std 3	If Std 2 < Std 1, Std 3 must be < Std 2



Note

The measurement precision is also dependent on the selected standard solutions. Therefore, the selected standard solutions should cover the expected value range of the subsequent concentration measurement.

Temperature while calibrating and measuring

For precise ISE measurements the temperature difference between measurement and calibration should not be greater than 2 K. Therefore, adjust the temperature of the standard and measuring solutions accordingly. If the temperature difference is greater the *[TempErr]* warning appears in the measured value display.

ISE Cal

This is the conventional **two-point** or **three-point calibration procedure** that uses 2 or 3 freely selectable standard solutions. The concentration expected in the measurement determines the concentration of the calibration standards.

Stability control

During calibration, the stability control is automatically activated. The current measurement with stability control can be terminated at any time (accepting the current value).

Calibration record

When finishing a calibration, the new calibration values are displayed as an informative message (■ symbol) first. Then you can decide whether you want to take over these values of the new calibration or whether you want to continue measuring with the old calibration data. After accepting the new calibration values the calibration record is displayed.

Display of calibration data and download to interface

You can have the data of the last calibration displayed (see page 87). Subsequently, you can transmit the displayed calibration data to the interface, e. g. to a printer or PC, with the **<PRINT>** key.

**Note**

The calibration record is automatically transmitted to the interface after calibrating.

Sample record:


```
Printing date 26.04.06 16:13
ProLab 2000
Ser. no. 06249876
USER ID: 1234567890

Calibration ISE
Calibration date 20.04.06 10:14:03
USER ID: 1234567890

Standard 1          0.010 mg/l
Standard 2          0.020 mg/l
Voltage 1           0.0 mV 24.0 °C
Voltage 2           9.0 mV 24.0 °C
Slope               29.9 mV
Sensor              +++
```

Calibration evaluation

After calibrating, the meter automatically evaluates the calibration.

Display	Calibration record	Magnitude of the slope [mV]
	+++	50.0 ... 70.0 or 25.0 ... 35.0
<i>Error</i> Eliminate the error according to chapter 6 WHAT TO DO IF... (page 97)	<i>Error</i>	< 50 or > 70 or < 25 or > 35

Preparatory activities

Perform the following preparatory activities when you want to calibrate:

1	Connect the ISE combination electrode to the meter. The pH/mV/ISE measuring window is displayed.
2	Keep the standard solutions ready.
3	Measure the temperature of the standard solutions using a thermometer.

Carrying out an ISE calibration

Proceed as follows to calibrate the meter:

1	In the measured value display, select the ISE measuring window with <▲><▼> and <MODE>.
2	Start the calibration with <CAL>. The calibration display appears.

ISE <CAL>

■ Immerse sensor in std. 1


Set temperature: 24 °C

Continue


Set standard

3	Thoroughly rinse the electrode with distilled water.
4	Immerse the electrode in standard solution 1.
5	Using <▲><▼>, select the <i>Set standard</i> setting and press <MENU/OK> .
6	Set the concentration of the standard solution with <▲><▼> and press <MENU/OK>.

7	Measure the temperature of the standard solution using a thermometer.
8	Select the <i>Set temperature</i> setting with <▲><▼> and press <MENU/OK>.
9	Set the temperature with <▲><▼> and press <MENU/OK>.
10	Select <i>Continue</i> with <▲><▼> and press <MENU/OK>. The standard solution is measured. The measured value is checked for stability (Stability control).




11	Wait for the end of the measurement with Stability control or accept the calibration value with <MENU/OK>. The calibration display for the next standard solution appears.
----	--




Continuing with two-point calibration

12	Thoroughly rinse the electrode with distilled water.
13	Immerse the electrode in standard solution 2.
14	Using <▲><▼>, select the <i>Set standard</i> setting and press <MENU/OK> .
15	Set the concentration of the standard solution with <▲><▼> and press <MENU/OK>.
16	Measure the temperature of the standard solution using a thermometer.

- | | |
|----|---|
| 17 | Select the <i>Set temperature</i> setting with <▲><▼> and press <MENU/OK>. |
| 18 | Set the temperature with <▲><▼> and press <MENU/OK>. |
| 19 | Select <i>Continue</i> with <▲><▼> and press <MENU/OK>. The standard solution is measured. The measured value is checked for stability (Stability control). |



- | | |
|----|--|
| 20 | Wait for the end of the measurement with Stability control or accept the calibration value with <MENU/OK>. The calibration display for the next standard solution appears. |
|----|--|



- | | |
|----|---|
| 21 | For two-point calibration, select <i>Exit with 2 point calibration</i> with <▲><▼> and confirm with <MENU/OK>. The calibration is completed as a two-point calibration. The new calibration values are displayed as a message (■). You have the following options: <ul style="list-style-type: none"> ● Accept the new calibration values with <MENU/OK>. Subsequently, the calibration record is displayed and downloaded to the interface at the same time. ● To switch to the measured value display <u>without</u> accepting the new calibration values, press <MODE> or <ESC>. |
|----|---|

Continuing with three-point calibration

Repeat the steps 12 to 20 in the same way for the third standard solution. After finishing the last calibration step, the new calibration values are displayed as a message (■).

You have the following options:

- Accept the new calibration values with **<MENU/OK>**. Subsequently, the calibration record is displayed and downloaded to the interface at the same time.
- To switch to the measured value display without accepting the new calibration values, press **<MODE>** or **<ESC>**.

**Note**

The calibration line is calculated by linear regression.

4.7 Conductivity

4.7.1 General information

You can measure the following variables:

- Conductivity
- Resistivity
- Salinity
- Total dissolved solids (TDS)

The meter provides the following functions:

- AutoRange (automatic switchover of the measurement range). If a measuring range is exceeded, AutoRange causes the meter to automatically change to the next higher measuring range and back again. Therefore, the instrument always measures in the measuring range with the highest possible resolution.
- Stability control to check the stability of the measurement signal. This ensures the reproducibility of the measuring signal. The display of the measured parameter flashes until a stable measured value is available.

Temperature measurement



The LF413T(-ID), LFOx1400(-ID) and LF313T(-ID) conductivity measuring cells have an integrated temperature sensor.

CAUTION

When connecting an earthed PC/printer, measurements cannot be performed in earthed media as incorrect values would result. The RS232 and USB interfaces are not galvanically isolated.

Preparatory activities

Perform the following preparatory activities when you want to measure:

1	Connect a conductivity measuring cell to the meter. The conductivity measuring window is displayed.
2	Check whether the <i>Measuring cell</i> and cell constant settings are suitable for the connected conductivity measuring cell. If necessary, correct the settings.




Note

The selection of the measuring cell and setting of the cell constant is made in the conductivity measuring menu (see page 65). The cell constant to be set must either be taken from the operating manual of the measuring cell or is printed on the measuring cell.

4.7.2 Measuring

You can carry out conductivity measurements as follows:

1	Perform the preparatory activities according to page 61.
2	Immerse the conductivity measuring cell in the test sample.



Selecting the displayed measured variable

You can switch between the following displays with **<MODE>**:

- Conductivity [$\mu\text{S}/\text{cm}$] / [mS/cm]
- Resistivity [$\Omega\cdot\text{cm}$] / [$\text{M}\Omega\cdot\text{cm}$]
- Salinity SaL []
- Total dissolved solids TDS [mg/l]

The factor to calculate the total dissolved solids is set to 1.00 in the factory. You can adjust this factor to meet your requirements in the range of 0.40 to 1.00. The factor can be set in the TDS measuring menu.

Stability control

The Stability control function checks the stability of the measurement signal. The stability has a considerable impact on the reproducibility of measured values. The display of the measured parameter flashes until a stable measured value is available.

Independent of the setting for automatic *Stability control* (see page 35) in the *System* menu, you can start the *Stability control* function manually at any time.

1	If necessary, select a channel (measured value) with <▲><▼> .
2	Freeze the measured value with <SC> . The display switches to the display of the selected channel as necessary. In the status line, [SC] is displayed.

- 3 With **<MENU/OK>** activate the *Stability control* function. [SC] flashes while the stability control is active. As soon as a stable measured value is recognized, the current measurement data is downloaded to the interface. The measured values of all channels meeting the criterion for stability control are marked by SC.



Note

You can terminate prematurely the manual *Stability control* function with **<MENU/OK>** manually at any time. If the manual *Stability control* function is terminated prematurely, the current measurement data is not downloaded to the interface.

- 4 Using **<SC>** or **<MODE>**, release the frozen measured value. The [SC] status display disappears. The display switches to the previous indication.

Criterion

The meter checks the stability of the measured value on the basis of the temperature measurement. For identical measurement conditions, the following applies:

Measured parameter	Reproducibility	Response time
Temperature	better than 0.02 °C	> 10 seconds

4.7.3 Temperature compensation

The calculation of the temperature compensation is based on the preset reference temperature, 20 °C or 25 °C. It appears on the display as *Tr20* or *Tr25*.

You can select one of the following temperature compensation methods:

- **Nonlinear temperature compensation (*nLF*)** according to EN 27 888
- **Linear temperature compensation (*lin*)** with selectable coefficients of 0.000 ... 3.000 %/K
- No temperature compensation (off)



Note

The reference temperature and temperature compensation are set in the conductivity measuring menu (see page 65).

Application tips

Select the following temperature compensations given in the table according to the respective test sample:

Test sample	Temperature compensation	Display indicator
Natural water (ground water, surface water, drinking water)	<i>nLF</i> according to EN 27 888	<i>nLF</i>
Ultrapure water	<i>nLF</i> according to EN 27 888	<i>nLF</i>
Other aqueous solutions	<i>lin</i> Set linear temperature coefficient 0.001 ... 3.000 %/K	<i>lin</i>
Salinity (seawater)	Automatic <i>nLF</i> according to IOT	<i>Sal, nLF</i>

4.7.4 Settings for conductivity measuring cells

Overview

The following settings are possible for conductivity measuring cells:

- Measured parameter
- Reference temperature
- Temperature compensation
- TDS factor
- Calibration interval
- Measuring cell / cell constant

Settings

The settings are made in the measuring menu of the conductivity measurement. To open the settings, activate the relevant measuring window in the measured value display and press the **<MENU/OK>** key. After completing the settings, switch to the measured value display with **<MODE>**.

Menu item	Possible setting	Description
<i>Calibration / Calibration interval</i>	<i>1 ... 999 d</i>	<i>Calibration interval</i> for the measuring cell (in days). The meter reminds you to calibrate regularly by the flashing CalClock in the measuring window.
<i>Calibration / Calibration record</i>	-	Displays the calibration record of the last calibration.
<i>Reset</i>	-	Resets all sensor settings to the delivery condition (see page 92).
<i>Temp. comp. (TC) / Reference temp.</i>	<i>20°C 25 °C</i>	Reference temperature This setting is only available when the Conductivity or Resistivity display is set.
<i>Temp. comp. (TC) / Compensation</i>	<i>nLF lin Off</i>	Procedure for temperature compensation (see page 64). This setting is only available when the Conductivity or Resistivity display is set.

Menu item	Possible setting	Description
<i>Temp. comp. (TC) / Linear coeff.</i>	<i>0.000 ... 3.000 %/K</i>	Coefficient of the linear temperature compensation. This setting is only available when the linear temperature compensation is set.
<i>TDS factor</i>	<i>0.40 ... 1.00</i>	TDS factor This setting is only available when the TDS display is set.
<i>Measuring cell</i>	<i>Cal</i>	<i>Measuring cell used</i> Measuring cells the cell constant of which is determined by calibration in the KCL control standard solution (Calibration ranges see page 67 and page 106). The currently valid cell constant is displayed in the status line.
	<i>C0.01</i>	Measuring cell with a nominal cell constant of 0.010 cm^{-1} . The cell constant is permanently set.
	<i>LF313T</i>	Measuring cell LF313T(-ID), nominal cell constant 0.100 cm^{-1} . The cell constant can be adjusted in the range 0.090 to 0.110 cm^{-1} .
	<i>man</i>	Any measuring cells with freely adjustable cell constants in the range 0.250 to 25.000 cm^{-1} .
<i>Cell constant</i>	0.090 to 0.110 cm^{-1}	Display and setting option for the cell constant of the measuring cell <i>LF313T(-ID)</i> .
<i>Man. cell const.</i>	0.250 to 25.000 cm^{-1}	Display and setting option of the cell constant of any measuring cells (<i>man</i>).

4.7.5 Determining the cell constant (calibration in control standard)

Why determine the cell constant?

Due to aging, the cell constant slightly changes. As a result, an inexact measured value is displayed. Calibration determines the current value of the cell constant and stores this value in the instrument. Therefore, calibrate at regular intervals. The calibration data is stored in the meter. The calibration data of non-ID sensors is stored in the meter. For ID sensors, the calibration data is stored in the sensor.

You can determine the cell constant of the conductivity measuring cell in the ranges $0.450 \dots 0.500 \text{ cm}^{-1}$, $0.585 \dots 0.715 \text{ cm}^{-1}$ or $0.800 \dots 1.200 \text{ cm}^{-1}$ by calibrating in the 0.01 mol/l KCl control standard or set it manually in the range $0.250 \dots 2.500 \text{ cm}^{-1}$ or $0.090 \dots 0.110 \text{ cm}^{-1}$. Besides, the fixed cell constant 0.010 cm^{-1} can be selected.

In the delivery condition, the calibrated cell constant of the meters is set to 0.475 cm^{-1} (conductivity measuring cells LF413T(-ID) and LFOx1400(-ID)).

CalClock

The expiry of the CalClock reminds you to clean the measuring cell. It is still possible to measure.

The CalClock is set to 180 days in the factory. You can change the interval (see page 65).

Stability control

During calibration, the stability control is automatically activated.

Display of calibration data and download to interface

You can have the data of the last calibration displayed (see page 87). Subsequently, you can transmit the displayed calibration data to the interface, e. g. to a PC, with the **<PRINT>** key.



Note

The calibration record is automatically transmitted to the interface after calibrating.

Sample record:

```

Printing date 26.04.06 16:13
ProLab 2000
Ser. no. 06249876
USER ID: 1234567890

Calibration Cond
Calibration date 20.04.06 10:14:03
LFOx1400-ID
Ser. no. A062498765
USER ID: 1234567890

Cal interval          150 d
Cell constant.       0.475 1/cm      25.0 °C
Sensor               +++
    
```



Note

If non-ID sensors are used, the sensor type and series number are not output in the calibration record.

Calibration evaluation

After the calibration, the meter automatically evaluates the current status of the calibration. The evaluation appears on the display and in the calibration record.

Display	Calibration record	Cell constant [cm ⁻¹]
	+++	within the ranges 0.450 ... 0.500 cm ⁻¹ 0.585 ... 0.715 cm ⁻¹ or 0.800 ... 1.200 cm ⁻¹
<i>Error</i> Eliminate the error according to chapter 6 WHAT TO DO IF... (page 97)	<i>Error</i>	outside the ranges 0.450 ... 0.500 cm ⁻¹ or 0.800 ... 1.200 cm ⁻¹

Determining the cell constant

Proceed as follows to determine the cell constant:

- 1 | Connect a conductivity measuring cell to the meter.
- 2 | In the measured value display, select the conductivity measuring window with <▲><▼> and <MODE>.
- 3 | Start the calibration with <CAL>. The calibration display appears.

⌘ <CAL>

- Cell constant = 0.475 1/cm
- Immerse sensor in solution

Continue

- 4 Immerse the conductivity measuring cell in the control standard solution, 0.01 mol/l KCl.
- 5 Press <MENU/OK>. The measurement with stability control to determine the cell constant starts.

⌘ <CAL>

- Cell constant = 0.481 1/cm
- Temperature =24.8 °C

Terminate stability control

- 6 Wait for the end of the measurement with Stability control or accept the calibration value with <MENU/OK>. The determined cell constant is displayed as a message (■).

⌘ <CAL>

- Cell constant = 0.481 1/cm
- Sensor+++

Accept

- 7 You have the following options:
- Accept the new calibration values with **<MENU/OK>**. Subsequently, the calibration record is displayed and downloaded to the interface at the same time.
 - To switch to the measured value display without accepting the new calibration values, press **<MODE>** or **<ESC>**.

4.8 Dissolved oxygen

4.8.1 General information

You can measure the following variables:

- D.O. concentration
- D.O. saturation index ("D.O. saturation")
- D.O. partial pressure

D.O. measurements with the ProLab 2000 can be carried out using the D.O. sensors, LFOx1400(-ID), 9009/6.(-ID) or Ox1400Stir(-ID). The stirrer of the Ox1400Stir(-ID) D.O. sensor has to be supplied with voltage separately using the NT/pH Mix 540 power supply. The meter automatically recognizes the type of the connected D.O. sensor.

The meter provides the following functions:

- AutoRange (automatic switchover of the measurement range). If a measuring range is exceeded, AutoRange causes the meter to automatically change to the next higher measuring range and back again. Therefore, the instrument always measures in the measuring range with the highest possible resolution.
- Stability control to check the stability of the measurement signal. This ensures the reproducibility of the measuring signal. The display of the measured parameter flashes until a stable measured value is available.



Note

Incorrect calibration of D.O. sensors will result in incorrect measured values. Calibrate at regular intervals.



CAUTION

When connecting an earthed PC/printer, measurements cannot be performed in earthed media as incorrect values would result. The RS232 and USB interfaces are not galvanically isolated.

Temperature sensor

The D.O. sensor has an integrated temperature sensor that always measures the current temperature of the test sample.


Preparatory activities

Perform the following preparatory activities when you want to measure:

1	Connect the D.O. sensor to the meter. The D.O. measuring window is displayed.
2	Calibrate or check the meter with the sensor.

4.8.2 Measuring

1	Perform the preparatory activities according to page 71.
2	Immerse the D.O. sensor in the test sample.



Selecting the displayed measured variable

You can switch between the following displays with **<MODE>**:

- D.O. concentration [mg/l]
- D.O. saturation [%]
- D.O. partial pressure [mbar].

Salinity correction

When measuring the concentration of solutions with a salt content of more than 1 g/l, a salinity correction is required. For this, you have to measure and input the salinity of the measured medium first. When the salinity correction is switched on, the *[SAL]* indicator is displayed in the measuring window.



Note

The salinity correction is switched on or off and the salinity is entered in the measuring menu of the D.O. measurement (see page 74).



Note

The LFOx1400(-ID) double sensor can perform the salinity correction automatically. The conductivity module of the sensor measures the salt content of the test sample simultaneously with the D.O. The meter takes the measured value into account.

Stability control

The Stability control function checks the stability of the measurement signal. The stability has a considerable impact on the reproducibility of measured values. The display of the measured parameter flashes until a stable measured value is available.

Independent of the setting for automatic *Stability control* (see page 35) in the *System* menu, you can start the *Stability control* function manually at any time.

1	If necessary, select a channel (measured value) with <▲><▼> .
---	--

- 2 | Freeze the measured value with **<SC>**.
The display switches to the display of the selected channel as necessary.
In the status line, [SC] is displayed.
- 3 | With **<MENU/OK>** activate the *Stability control* function.
[SC] flashes while the stability control is active.
As soon as a stable measured value is recognized, the current measurement data is downloaded to the interface. The measured values of all channels meeting the criterion for stability control are marked by SC.



Note

You can terminate prematurely the manual *Stability control* function with **<MENU/OK>** manually at any time. If the manual *Stability control* function is terminated prematurely, the current measurement data is not downloaded to the interface.

- 4 | Using **<SC>** or **<MODE>**, release the frozen measured value.
The [SC] status display disappears. The display switches to the previous indication.

Criteria

With identical measurement conditions, the following applies:

Measuring mode	Reproducibility	Response time
D.O. concentration	Better than 0.05 mg/l	> 10 seconds
D.O. saturation	Better than 0.6 %	> 10 seconds
D.O. partial pressure	Better than 1.2 mbar	> 10 seconds

4.8.3 Settings for D.O. sensors

Overview

The following settings are possible for D.O. sensors:

- Salinity correction
- Salinity (salinity equivalent)
- Calibration interval
- Comparison measurement

Settings

The settings are made in the measuring menu of the D.O. measurement. To open the settings, activate the relevant measuring window in the measured value display and press the <MENU/OK> key. After completing the settings, switch to the measured value display with <MODE>.

Menu item	Possible setting	Description
<i>Calibration / Calibration interval</i>	1 ... 999 d	<i>Calibration interval</i> for the D.O. sensor (in days). The meter reminds you to calibrate regularly by the flashing CalClock in the measuring window.
<i>Calibration / Comparison meas.</i>	<i>On</i> <i>Off</i>	Enables to adjust the measured value with the aid of a comparison measurement, e.g. Winkler titration. For details, see page 75.
<i>Calibration / Calibration record</i>	-	Displays the calibration record of the last calibration.
<i>Reset</i>	-	Resets all sensor settings to the delivery condition (see page 92).
<i>Sal automatic</i>	<i>On</i> <i>Off</i>	Automatic salt content correction for concentration measurements. Note: This function is available with the LFOx1400(-ID) double sensor <u>only</u> .
<i>Sal correction</i>	<i>On</i> <i>Off</i>	Manual salt content correction for concentration measurements.
<i>Salinity</i>	0.0 ... 70.0	Salinity or salinity equivalent for the salt content correction. This function is only available if the manual salt content correction is switched on.

4.8.4 D.O. calibration

Why calibrate?

D.O. sensors age. This changes the slope of the D.O. sensor. Calibration determines the current slope of the sensor and stores this value in the instrument.

When to calibrate?

- After connecting another non-ID D.O. sensor
- When the CalClock flashes (after the calibration interval has expired).

Calibration datasets

The ProLab 2000 administrates two sets of calibration data:

- Set 1 for the type, "Ox":
 - 9009/6. or
 - Ox1400Stir
- Set 2 for the type, "FOx1400":
 - LFOx1400

Sensors of different types can be calibrated separately from each other. When a sensor of one type is calibrated, the calibration data of the other types remain stored. The ProLab 2000 recognizes the type of the connected sensor and automatically uses the correct calibration data.

Calibration procedures

The ProLab 2000 provides 2 calibration procedures:

- Calibration in water vapor-saturated air.
Use an air calibration vessel for calibration.
- Calibration via a comparison measurement (e.g. Winkler titration according to DIN EN 25813 or ISO 5813). At the same time, the relative slope is adapted to the comparison measurement by a correction factor. When the correction factor is active, the *[Factor]* indicator appears in the measuring window.

Stability control

During calibration, the stability control is automatically activated.

Display of calibration data and download to interface

You can have the data of the last calibration displayed (see page 87). Subsequently, you can transmit the displayed calibration data to the interface, e. g. to a printer or PC, with the **<PRINT>** key.



Note

The calibration record is automatically transmitted to the interface after calibrating.

Sample record:

```

Printing date 26.04.06 16:13
ProLab 2000
Ser. no. 06249876
USER ID: 1234567890

CALIBRATION Oxi
Calibration date 20.04.06 10:14:03
LFOx1400-ID
Ser. no. A062498765
USER ID: 1234567890

Cal. interval      14d
Relative slope     0.97
Sensor            +++
    
```



Note

If non-ID sensors are used, the sensor type and series number are not output in the calibration record.

Calibration evaluation

After the calibration, the meter evaluates the current status of the sensor against the relative slope. The evaluation appears on the display and in the calibration record. The relative slope has no effect on the measuring accuracy. Low values indicate that the electrolyte will soon be depleted and the sensor will have to be regenerated.

Display	Calibration record	Relative slope
	+++	S = 0.8 ... 1.25
	++	S = 0.7 ... 0.8
	+	S = 0.6 ... 0.7
<i>Error</i>	<i>Error</i>	S < 0.6 or S > 1.25

Eliminate the error according to chapter 6 WHAT TO DO IF... (page 97)

Calibration in water vapor saturated air (air calibration vessel)

For this calibration procedure, the *Comparison meas.* setting must be set to *Off* in the measuring menu.


Proceed as follows to calibrate the meter:




Note

The sponge in the air calibration vessel must be moist (not wet). Leave the sensor in the air calibration vessel for a time long enough to adjust.


1	Put the D.O. sensor into the air calibration vessel.
2	Connect the D.O. sensor to the meter.
3	In the measured value display, select the D.O. measuring window with <▲> <▼> and <MODE> .
4	Start the calibration with <CAL> . The calibration display for the relevant sensor type appears.



5	Put the D.O. sensor into the air calibration vessel.
6	Press <MENU/OK> . The measurement with stability control to determine the relative slope starts.



7	Wait for the end of the measurement with Stability control or accept the calibration value with <MENU/OK> . The determined relative slope is displayed as a message (■).
---	--



- 8 You have the following options:
- Accept the new calibration values with **<MENU/OK>**. Subsequently, the calibration record is displayed and downloaded to the interface at the same time.
 - To switch to the measured value display without accepting the new calibration values, press **<MODE>** or **<ESC>**.

Calibration via a Comparison meas.

For this calibration procedure, the *Comparison meas.* setting must be set to *On* in the measuring menu.




Note

Before calibrating via a comparison measurement, the sensor should be calibrated in the air calibration vessel.

Proceed as follows to calibrate the meter:


- 1 Connect the D.O. sensor to the meter.
- 2 Immerse the D.O. sensor in the reference solution.
- 3 In the measured value display, select the D.O. measuring window with **<▲><▼>** and **<MODE>**.
- 4 Start the calibration with **<CAL>**. The calibration display appears.



- 5 Press **<MENU/OK>**.
The measurement with stability control to determine the D.O. concentration starts.

LFOx1400-ID <CAL>

■ Concentration = 8.5 mg/l
■ Temperature =24.8 °C



Terminate stability control

- 6 Wait for the end of the measurement with Stability control or accept the measured value with **<MENU/OK>**.
The determined D.O. content is displayed as a message (■).

LFOx1400-ID <CAL>

■ Concentration = 8.49 mg/l
■ Temperature =24.8 °C

1.000

Accept

- 7 Open the setting of the correction factor with **<MENU/OK>**.
- 8 Using **<▲><▼>**, set the correction factor to adjust the displayed concentration value to the nominal value (value of the comparison measurement). Subsequently, accept the correction factor with **<MENU/OK>**.
- 9 Select *Accept* with **<▲><▼>** and press **<MENU/OK>**.
The meter switches to the measured value display. The *[Factor]* indicator appears in the measuring window.

4.9 Storage

You can transmit measured values (datasets) to the data storage in two ways:

- Manual storage (see page 81)
- Automatic storage at intervals (see page 82)

Each storage process transmits the current dataset to the interface at the same time.

Measurement dataset

A complete dataset consists of:

- Date/time
- *USER ID*
- Meter designation with series number
- SENSOR ID
 - Sensor designation
 - Sensor series number
- *SAMPLE ID*
- Measured value of the connected sensor
- Measured temperature value of the connected sensor
- Info on stability control: *SC* appears with the measured value if the Stability control criterion was met while storing (stable measured value). Otherwise, there is no *SC* display.
- Measurement conditions

Storage locations


The meter has separate data storages for manually stored measured values and automatically stored measured values.

Storage	Maximum number of datasets
Manual storage	308
Automatic storage	2244

4.9.1 Manual storage

Proceed as follows to transmit to the data storage and simultaneously download to the interface a measurement dataset:

- 1 Open the menu for manual storage with **<STO>**.



<STO> 1 of 646
- 26.04.2006 11:24:16
- pH 7.000 24.8 °C SC


SAMPLE ID: 1

Continue

- 2 If necessary, use **<▲>****<▼>** and **<MENU/OK>** to change and confirm the *SAMPLE ID* (1 ... 999).
The dataset is stored. The meter switches to the measured value display.

If the data storage is full

The following window appears if all storage locations are occupied:



Warning

Data storage full. Erase?

Yes

No

You have the following options:


- To erase the entire data storage, confirm *Yes*.
- To cancel the storage process and switch to the measured value display, confirm *No*. Then you can e.g. transmit the stored data to a PC (see page 84) and subsequently erase the storage (see page 87).

4.9.2 Automatic storage at intervals

The storage interval (*Interval*) determines the chronological interval between automatic storage processes. Each storage process transmits the current dataset to the interface at the same time.

Configuring the automatic storage function

- 1 Open the menu for automatic storage with <STO __>.




Settings

You can configure the automatic storage function with the following settings:

Menu item	Possible setting	Description
SAMPLE ID	1 ... 999	Determine the sample designation for the series of datasets.
Interval	1 s, 5 s, 10 s, 30 s, 1 min, 5 min, 10 min, 15 min, 30 min, 60 min	Storage interval. The lower limit of the storage interval can be restricted by the number of free storage locations. The upper limit is restricted by the storage duration.
Duration	1 min ... x min	Storage duration. Specifies after which time the automatic storage should be terminated. The lower limit of the storage duration is restricted by the storage interval. The upper limit is restricted by the number of free storage locations.


Data storage administration

The two lower display lines indicate the use of the data storage calculated in advance for the selected settings:



Starting the automatic storage function

To start the automatic storage function, select *Continue* with \blacktriangle and \blacktriangledown and confirm with \blacktriangle . The meter switches to the measured value display.



The active automatic storage function can be recognized from the progress bar in the status line. The progress bar indicates how much of the adjusted storage duration has already expired.

Terminating the automatic storage function prematurely

Proceed as follows to switch off the automatic storage function before the adjusted storage duration has expired:

- 1 With **<STO__>** open the following window.

Warning

Stop automatic storage?

Yes

No

- 2 Select and confirm *Yes* with **<▲><▼>** and **<MENU/OK>**. The meter switches to the measured value display. The automatic storage function is terminated.

4.9.3 Reading the measurement data storage

You can select the contents of the manual or automatic measurement data storage by means of different filter criterions and

- read them out on the display, and
- download them to the interface.

Each measurement data storage has a separate erasure function for the entire contents, independent of the filter settings.

Editing the data storage

The data storage is edited in the *System / Data storage* menu. In the measured value display, open the system menu with **<MENU/OK__>**. With the **<RCL>** or **<RCL__>** key you can move directly to the menu for the manual or automatic data storage.



Note

The settings are explained here using the manual storage as an example. The same settings and functions are available for the automatic storage.

Settings	Menu item	Setting/ function	Description
	<i>Data storage / Manual data storage / Display</i>	-	<p>Displays in pages all measuring datasets that correspond to the filter settings (see menu <i>Data filter</i>).</p> <p>Further options:</p> <ul style="list-style-type: none"> ● Scroll through the datasets with <▲><▼>. ● Output the displayed dataset to the interface with <PRINT>. ● Quit the display with <ESC>.
	<i>Data storage / Manual data storage / Output to RS232/USB</i>	-	<p>Outputs to the interface all measuring datasets that correspond to the filter settings. The download takes place in ascending order of the storage location number.</p> <p>The process can take several minutes. To terminate the process prematurely, press <ESC>.</p>
	<i>Data storage / Manual data storage / Data filter</i>	-> see <i>explanations below this table, page 86</i>	Allows to set certain filter criteria in order to display datasets and download them to the interface.
	<i>Data storage / Manual data storage / Erase</i>	-	<p>Erases the entire contents of the selected measuring data storage, independent of the filter settings.</p> <p>Note: All calibration data remains stored when performing this action.</p>

Data filter	Menu item	Setting/function	Description
	<i>Filter</i>	<i>No filter</i>	Filter criteria: Data filter switched off
		<i>Date & SAMPLE ID</i>	Selection according to period and <i>SAMPLE ID</i> .
		<i>SAMPLE ID</i>	Selection according to <i>SAMPLE ID</i>
		<i>Date</i>	Selection according to period
	<i>From</i>	<i>TT.MM.JJ</i>	Selects all datasets within the specified period.
	<i>Until</i>	<i>TT.MM.JJ</i>	
	<i>SAMPLE ID</i>	0 ... 999	Selects all datasets with the specified <i>SAMPLE ID</i> .

Display presentation of a dataset

```

30.04.2006 11:24:16 (1)
ProLab 2000: 06249876
USER ID: 1234567890
A 161 1M-BNC-ID: A062498765
SAMPLE ID: 1
- pH 7.000 24.8 °C SC
    
```

Sample printout

```

26.04.06 16:13
ProLab 2000
Ser. no. 06249876
USER ID: 1234567890

SAMPLE ID: 1
A 161 1M-BNC-ID
Ser. no. A062498765
LFOx1400-ID
Ser. no. A062498765
pH 6.87 24.0 °C, SC +++
Oxi 7.76 mg/l 25.0 °C, SC, ++
Sal = 0.7
Cond 1389 uS/cm 25.0 °C, SC, +++
C = 0.958 l/cm, Tref25, nLF

-----

26.04.06 16:03
ProLab 2000
Ser. no. 06249876
USER ID: 1234567890
    
```

Quitting the display

To quit the display of stored measuring datasets, you have the following options:

- Switch directly to the measured value display with **<MODE>**.
- Quit the display and move to the next higher menu level with **<ESC>** or **<MENU/OK>**.

4.9.4 Erasing the data storage

How to erase the measured value storage is described on page 84.

4.9.5 Displaying and downloading calibration records

You can display the calibration data

- of one selected sensor, or
- of all sensors connected

and then download them to the interface.

Displaying the calibration record of a selected sensor

The calibration record of the last calibration can be found under the *Calibration / Calibration record* menu item in the respective measuring menu. To open it, activate the relevant measuring window in the measured value display and press **<MENU/OK>**.

Displaying the calibration records of all sensors connected

The calibration records of the last calibration can be found under the *Data storage / Calibration data storage* menu item in the system menu. In the measured value display, open the system menu with **<MENU/OK__>**.

With the **<CAL__>** key you can go directly to the *Calibration data storage* menu.

Menu item	Setting/ function	Description
<i>Data storage / Calibration data storage /Display</i>	-	Displays in pages all calibration records of all sensors. Further options: <ul style="list-style-type: none"> ● Scroll through the calibration records with <▲><▼>. ● Output the displayed calibration record to the interface with <PRINT>. ● Quit the display with <ESC> or <MENU/OK>. ● Switch directly to the measured value display with <MODE>.
<i>Data storage / Calibration data storage / Output to RS232/USB</i>	-	Outputs to the interface the calibration records of all sensors.

Sample printout

```

Printing date 26.04.06 16:13
ProLab 2000
Ser. no. 06249876
USER ID: 1234567890

CALIBRATION pH
Calibration date 20.04.06 10:14:03
A 161 1M DIN-ID
Ser. no. A062498765
USER ID: 1234567890

Cal. interval      7 d
DIN
Buffer 1           6.86
Buffer 2           9.18
Voltage 1          7.2 mV      26.3 °C
Voltage 2          -171 mV     26.3 °C
Slope              -59.2 mV/pH
Asymmetry          -0.5 mV
Sensor             +++

CALIBRATION ISE
Calibration date 20.04.06 10:14:03
Standard 1         0.010 mg/l
Standard 2         0.020 mg/l
Voltage 1          38.5 mV      25.0 °C
Voltage 2          58.0 mV      23.0 °C
    
```

4.10 Transmitting data (to a PC or printer)

The meter has two interfaces:

- RS232 interface (serial port)
- USB interface (device)

Via both interfaces, you can transmit data to a PC and update the meter software.

The RS232 interface enables to transmit data to an external serial printer.



Note

The relevant interface cable has to be connected if you want to download data to an interface (USB or RS232).

It is not possible to download data to both interfaces (USB and RS232) at the same time. After connecting a meter to the USB socket the RS232 interface is inactive. The RS232 interface is active if no meter is connected to the USB interface.



CAUTION


The interfaces are not galvanically separated.

When connecting an earthed PC/printer, measurements cannot be performed in earthed media as incorrect values would result.

4.10.1 RS232 interface

1	Connect the interface to the PC or printer via the cable Z390 (PC) or Z393 (ext. printer).
2	If necessary, disconnect a connected USB cable from the meter.
3	Set up the following transmission data on the PC/printer: <ul style="list-style-type: none"> – Baud rate: selectable from 1200, 2400, 4800, 9600, – Handshake: RTS/CTS + Xon/Xoff – PC only: – Parity: none – Data bits: 8 – Stop bits: 1

**Socket assignment
(RS232)**



RS 232

- 1 *
- 2 RxD
- 3 TxD
- 4 *
- 5 SGnd
- 6 CTS

* not used

4.10.2 USB interface (device)

Connect the interface to the PC via the supplied Z875 USB cable. The data output automatically switches to *USB*. The RS232 interface is deactivated.

**Installation of the USB
driver on the PC**

System requirements of the PC for installation of the USB driver:

- PC with Pentium processor or higher with at least one free USB connection and CD-ROM drive
- Windows 2000, XP, Vista.

1	Insert the supplied installation CD in the CD drive of your PC.
2	Install the USB driver on the PC. Follow the Windows installation instructions as necessary.
3	Connect the ProLab 2000 to the PC via the USB interface. The meter is listed as a virtual COM interface among the connections in the Windows instrument manager.

4.10.3 Options for data transmission

Via the RS 232 interface, you can transmit data to a PC or external printer. The following table shows which data is transmitted to the interface in which way:

Data	Control	Operation / description
Current measured values of the connected sensor	Manual	<ul style="list-style-type: none"> ● With <PRINT>. ● Simultaneously with every manual storage process (see page 81).
	automatic, at intervals	<ul style="list-style-type: none"> ● With <PRT__>. Then you can set the transmission interval. ● Simultaneously with every automatic storage process (see page 82).
Stored measured values	Manual	<ul style="list-style-type: none"> ● Displayed dataset: with <PRINT> after recall from data storage. ● All datasets according to the filter criteria via the <i>Output to RS232/USB</i> function. <p>For details, see page 84.</p>
Calibration records	Manual	<ul style="list-style-type: none"> ● Calibration record of a sensor: Calibration record of a sensor with <PRINT> (after calling up from the data storage or at the end of a calibration procedure). ● All calibration records after calling up from the data storage via the <i>Output to RS232/USB</i> function. <p>For details, see page 87.</p>
	automatic	<ul style="list-style-type: none"> ● For the respective sensor at the end of a calibration procedure.



Note

The following rule applies: Except for the menus, the display content is generally downloaded to the interface with **<PRINT>** (displayed measured values, measurement datasets, calibration records).

4.10.4 Operation with MultiLab pilot

With the aid of the MultiLab pilot software, you can record and evaluate measuring data with a PC. The data is transmitted after the meter was connected to the RS232 serial interface (COM port) or the USB interface of the PC.



Note

More detailed information can be found in the MultiLab pilot operating manual.

4.11 Reset

You can reset (initialize) all sensor settings and sensor-independent settings separately from each other.

4.11.1 Resetting the sensor settings



Note

The calibration data is reset to the default settings together with the measuring parameters. Recalibrate after performing a reset.

pH

The following settings for pH measurements are reset to the default settings with the *Reset* function:

Setting	Default settings
<i>Cal. interval</i>	7 d
<i>Unit for slope</i>	mV/pH
Measured parameter	pH
<i>High resolution pH</i>	<i>On</i>
<i>High resolution mV</i>	<i>On</i>
Zero point	0 mV
Slope	-59.16 mV
Temperature, manual	25 °C

The sensor settings are reset under the *Reset* menu item in the respective measuring menu. To open it, activate the relevant measuring window in the measured value display and press the **<MENU/OK>** key.

Cond The following settings for conductivity measurements are reset to the default settings with the *Reset* function:

Setting	Default settings
<i>Cal. interval</i>	150 d
Measured parameter	∅
Cell constant (c)	0.475 cm ⁻¹ (calibrated) 0.475 cm ⁻¹ (set up)
Temperature compensation	nLF
Reference temperature	25 °C
Temperature coefficient (TC) of the linear temperature compensation	2.000 %/K
TDS factor	1.00

The sensor settings are reset under the *Reset* menu item in the respective measuring menu. To open it, activate the relevant measuring window in the measured value display and press the <MENU/OK> key.

Oxi The following settings for D.O. measurements are reset to the default settings with the *Reset* function:

Setting	Default settings
<i>Cal. interval</i>	14d
<i>Comparison meas.</i>	Off
Measured parameter	D.O. concentration
Relative slope (S_{Rel})	1.00
Salinity (value)	0.0
Salinity (function)	Off

4.11.2 Resetting the system settings

The following system settings can be reset to the delivery status:

Setting	Default settings
<i>Language</i>	<i>Deutsch</i>
<i>Temperature unit</i>	°C
<i>Beep</i>	<i>On</i>
<i>Baud rate</i>	4800 Baud
<i>Output format</i>	ASCII
<i>Contrast</i>	40 %
<i>Illumination</i>	<i>On</i>

The system settings are reset under the *System / Reset* menu item in the system menu. In the measured value display, open the system menu with <MENU/OK__>.

5 Maintenance, cleaning, disposal

5.1 Maintenance

The only maintenance activity required is replacing the batteries.




Note

For maintenance of the measuring cells and electrodes refer to the relevant operating manuals.

Replacing the batteries

- | | |
|---|---|
| 1 | Using a screw driver, unscrew the screws (2) of the lid of the battery compartment. |
|---|---|



- | | |
|---|---|
| 2 | Open the battery compartment (1) on the underside of the meter. |
| 3 | Remove the four batteries from the battery compartment. |
| 4 | Place four new batteries (type Micro AAA) in the battery compartment. |
| 5 | Close the battery compartment (1).
The date (day) flashes on the display. |
| 6 | Using a screw driver, tighten the screws (2) of the lid of the battery compartment. |
| 7 | Set the date and time (see page 27) |

**CAUTION**

Make sure that the poles of the batteries are the right way round. The \pm signs on the batteries must correspond to the \pm signs in the battery compartment.

Only use leakproof alkaline manganese batteries.

5.2 Cleaning

Occasionally wipe the outside of the meter with a damp, lint-free cloth. Disinfect the housing with isopropanol as required.

**CAUTION**

The housing components are made out of synthetic materials (polyurethane, ABS, PMMA), metal and glass. Thus, avoid contact of the synthetic material with acetone and similar detergents that contain solvents. Remove any splashes immediately.

5.3 Disposal**Note**

This instrument contains batteries. Batteries that have been removed must only be disposed of at the recycling facility set up for this purpose or via the retail outlet.

It is illegal to dispose of the instrument in household refuse.

6 What to do if...

6.1 pH and ORP measurement

Error message <i>Err1</i>	Cause	Remedy
	pH electrode:	
	– Measured value outside the measuring range	– Use a suitable electrode
	– Air bubble in front of the junction	– Remove air bubble
	– Air in the junction	– Extract air or moisten junction
	– Cable broken	– Replace electrode
	– Gel electrolyte dried out	– Replace electrode

Error message <i>Err2</i>	Cause	Remedy
	– No electrode connected	– Connect electrode
	– Setting time during calibration too long	– Adjust temperature if necessary – Recalibrate

Error message <i>Err4</i>	Cause	Remedy
	– Temperature not stable during calibration.	– Adjust temperature if necessary – Recalibrate

Error message <i>CalError</i>	Cause	Remedy
	pH electrode:	
	– The values determined for zero point and slope of the electrode are outside the allowed limits.	– Recalibrate
	– Junction contaminated	– Clean junction
	– Electrode broken	– Replace electrode

Cause	Remedy
Buffer solutions:	
– Incorrect buffer solutions	– – Change buffer solutions
– Buffer solutions too old	– Use only once. Note the shelf life
– Buffer solutions depleted	– Change buffer solutions

No stable measured value

Cause	Remedy
<i>pH electrode:</i>	
– Junction contaminated	– Clean junction
– Membrane contaminated	– Clean membrane
<i>Test sample:</i>	
– pH value not stable	– Measure with air excluded if necessary
– Temperature not stable	– Adjust temperature if necessary
<i>Electrode + test sample:</i>	
– Conductivity too low	– Use a suitable electrode
– Temperature too high	– Use a suitable electrode
– Organic liquids	– Use a suitable electrode

Obviously incorrect measured values

Cause	Remedy
<i>pH electrode:</i>	
– Not connected	– Connect electrode
– Cable broken	– Replace electrode
– pH electrode unsuitable	– Use a suitable electrode
– Temperature difference between buffer and test sample too high	– Adjust temperature of buffer or sample solutions

	Cause	Remedy
	– Measurement procedure not suitable	– Follow special procedure
CalClock flashes	Cause	Remedy
	– Calibration interval expired	– Recalibrate the measuring system

6.2 ISE measurement

Error message <i>Err1</i>	Cause	Remedy
	– Measuring range exceeded	– Dilute test sample
Obviously incorrect measured values	Cause	Remedy
	– Electrode not connected	– Connect electrode
	– Cable broken	– Replace electrode
Error message, Error (invalid calibration)	Cause	Remedy
	<i>ISE electrode:</i>	
	– Moisture in the plug	– Dry plug
	– Electrode obsolete	– Replace electrode
	– Electrode unsuitable for the range to be measured	– Use suitable electrode
	– Socket damp	– Dry socket
	<i>Calibration procedure:</i>	
	– Wrong sequence of standards for three point calibration	– Select correct sequence
	– Calibration standards do not have the correct temperature (max. ± 2 °C temperature difference)	– Adjust the temperature of the calibration standards

Warning [TempErr]

Cause	Remedy
– Temperature difference between measurement and calibration greater than 2 K.	– Adjust the temperature of the test sample

Warning [ISEErr]

Cause	Remedy
– Electrode voltage outside calibrated range	– Recalibrate

6.3 Conductivity measurement**Conductivity measuring cell was not recognized**

Cause	Remedy
– Measuring cell not connected	– Connect measuring cell
– Cable broken	– Replace measuring cell

Error message, Error (invalid calibration)

Cause	Remedy
– Measuring cell contaminated	– Clean cell and replace it if necessary
– Unsuitable control standard	– Use control standard, 0.01 mol/l KCl

CalClock flashes

Cause	Remedy
– Calibration interval expired	– Recalibrate the measuring system

6.4 D.O. measurement**D.O. sensor was not recognized**

Cause	Remedy
– Sensor not connected	– Connect probe
– Cable broken	– Replace probe

Error message, <i>OFL</i>	Cause	Remedy
	– Measuring range exceeded	– not possible
Error message, <i>Error</i> (invalid calibration)	Cause	Remedy
	<i>D.O. sensor:</i>	
	– Electrolyte solution depleted	– Regenerate probe
	– Membrane contaminated	– Clean membrane
	– Electrode system poisoned	– Regenerate probe
Error message, <i>LEAK</i>	Cause	Remedy
	<i>D.O. sensor:</i>	
	– Membrane damaged	– Exchange the membrane cap
	– Membrane head screwed on not tight enough	– Screw membrane head tight
Measured value flashes for a long time (No stable measured value)	Cause	Remedy
	<i>D.O. sensor:</i>	
	– Membrane contaminated	– Clean membrane
Measured value too low	Cause	Remedy
	<i>D.O. sensor:</i>	
	– Insufficient flow	– Provide flow to the probe
Measured value too high	Cause	Remedy
	– High amount of dissolved substances	– Correct solubility function using the salinity equivalent
	– Air bubbles bump on the membrane with high velocity	– Avoid direct flow to the membrane

Cause	Remedy
– The carbon dioxide pressure is too high (> 1 bar)	– Measuring not possible

CalClock flashes

Cause	Remedy
– Calibration interval expired	– Recalibrate the measuring system

6.5 General errors**Meter does not react to keystroke**

Cause	Remedy
– Operating condition undefined or EMC load unallowed	– Processor reset: Press the <On/Off> and <SC> key simultaneously.

Time is lost

Cause	Remedy
– Buffer batteries empty	– Exchange buffer batteries (see page 95)

Translation of the legally binding German version

Stand October 30, 2009

7 Technical data

7.1 General data

Dimensions	approx. 240 x 280 x 70 mm	
Weight	approx. 2.5 kg (without power pack)	
Mechanical structure	Type of protection	IP 54
Electrical safety	Protective class	III
Test certificates	cETLus, CE	
Ambient conditions	Storage	- 25 °C ... + 65 °C
	Operation	+ 5 °C ... + 45 °C
	Climatic class	2
Power supply	Power pack	FRIWO FW7555M/09, 15.1432.500-00 Friwo Part. No. 1883259 Input: 100 ... 240 V ~ / 50 ... 60 Hz / 400 mA Output: 9 V = / 1.5 A Connection max. overvoltage category II Primary plugs contained in the scope of delivery: Euro, US, UK and Australian.
	Batteries (to buffer the system clock in the case of a power failure)	4 x 1.5 V alkali-manganese batteries, Type AAA
Serial interface	Connection of the PC cable (Z390) or printer cable (Z393)	
	Baud rate	adjustable: 1200, 2400, 4800, 9600, 19200 Baud
	Type	RS232, data output
	Data bits	8
	Stop bits	2
	Parity	None
	Handshake	RTS/CTS
	Cable length	Max. 15m
	USB interface	Automatic switch-over when a USB cable is connected.
Type		USB 1.1 (device)
Cable length		max. 3 m
Guidelines and norms used	EMC	EC guideline 2004/108/EC EN 61326-1 Class B FCC Class A

Translation of the legally binding German version

Stand October 30, 2009

Instrument safety	EC guideline 2006/95/EC EN 61010-1 ANSI/UL 61010-1 CAN/CSA-C22.2 No. 61010-1
Radio data transmission	EC guideline 1999/5/EC EN 300 330-2 EN 50364 EN 60950-1
Climatic class	VDI/VDE 3540
IP protection class	EN 60529

7.2 Measuring ranges, resolution, accuracy

7.2.1 pH/ORP

Measuring ranges, resolution

Variable	Measuring range	Resolution
pH	- 2.000 ... + 20.000 - 2.00 ... + 20.00	0.001 0.01
U [mV]	- 1999.9 ... + 1999.9 - 1999 ... + 1999	0.1 1
T [°C]	- 10.0 ... + 120.0	0.1

Manual temperature input

Variable	Range	Increment
T _{manual} [°C]	- 10 ... + 130	1

Accuracy (± 1 digit)

Variable	Accuracy	Temperature of the test sample
<i>pH</i>	± 0.003 ± 0.01	+ 15 °C ... + 35 °C
<i>U [mV] / range</i>		
- 999.9 ... + 999.9	± 0.2	+ 15 °C ... + 35 °C
- 1999 ... + 1999	± 1	+ 15 °C ... + 35 °C
<i>T [°C] / temperature sensor</i>		
	± 0.1	0 °C ... + 55 °C

Translation of the legally binding German version

Stand October 30, 2009

7.2.2 ISE

Measuring ranges, resolution	Variable	Measuring range	Resolution
	ISE [mg/l]	0.000 ... 10.000 0.00 ... 100.00 0.0 ... 1000.0 0 ... 2000	0.001 0.01 0.1 1

Manual temperature input	Variable	Range	Increment
	T_{manual} [°C]	- 20 ... + 130	1

7.2.3 Conductivity

Measuring ranges, resolution	Variable	Measuring range	Resolution
	κ [$\mu\text{S}/\text{cm}$]	0.000 ... 2.000*	0.001
0.00 ... 20.00**		0.01	
0.0 ... 200.0		0.1	
0 ... 2000		1	
κ [mS/cm]	0.00 ... 20.00	0.01	
	0.0 ... 200.0	0.1	
	0 ... 2000	1	
Specific resistance [kOhm*cm]	0.000 ... 2.000	0.001	
	0.00 ... 20.00	0.01	
	0.0 ... 200.0	0.1	
	0 ... 2000	1	
Resistivity [MOhm*cm]	0.00 ... 20.00	0.01	
	0.0 ... 200.0	0.1	
	0 ... 2000	1	
SAL	0.0 ... 70.0 according to the IOT table	0.1	
TDS [mg/l]	0 ... 2000 Factor can be set between 0.40 and 1.00	1	
T [°C]	- 5.0 ... + 105.0	0.1	

* only possible with cells of the cell constant, 0.010 cm^{-1} ** only possible with cells of the cell constant, 0.010 cm^{-1} or $0.090 \dots 0.110 \text{ cm}^{-1}$

Translation of the legally binding German version

Stand October 30, 2009

Cell constants	Cell constant C	Values
	Can be calibrated in the ranges	0.450 ... 0.500 cm ⁻¹ 0.585 ... 0.715 0.800 ... 1.200 cm ⁻¹
	Adjustable	0.010 cm ⁻¹ (fixed) 0.090 ... 0.110 cm ⁻¹ 0.250 ... 25.000 cm ⁻¹

Reference temperature	Reference temperature	Values
	Adjustable	20 °C (Tr20) 25 °C (Tr25)

Accuracy (± 1 digit)	Variable	Accuracy	Temperature of the test sample
<i>κ</i> / Temperature compensation			
	None (Off)	± 0.5 %	
	Nonlinear (nLF)	± 0.5 %	0 °C ... + 35 °C according to EN 27 888
		± 0.5 %	+ 35 °C ... + 50 °C Extended nLF function according to WTW measurements
	Linear (lin)	± 0.5 %	+ 10 °C ... + 75 °C
SAL / range			
	0.0 ... 42.0	± 0.1	+ 5 °C ... + 25 °C
		± 0.2	+ 25 °C ... + 30 °C
TDS [mg/l]			
		± 1	
T [°C] / temperature sensor			
	NTC 30	± 0.2	0 °C ... + 55 °C
	PT 1000	± 0.3	0 °C ... + 55 °C

Translation of the legally binding German version

Stand October 30, 2009

7.2.4 Dissolved oxygen**Measuring ranges,
resolution**

Variable	Measuring range	Resolution
D.O. concentration [mg/l]	0 ... 20.00 0 ... 90.0	0.01 0.1
Saturation [%]	0 ... 200.0 0 ... 600	0.1 1
D.O. partial pressure [mbar]	0 ... 200.0 0 ... 1250	0.1 1
T [°C]	0 ... 50.0	0.1

Accuracy (± 1 digit)

Variable	Accuracy
D.O. concentration [mg/l]	± 0.5 % of measured value at ambient temperature + 5 °C ... + 30 °C
Saturation [%]	± 0.5 % of measured value when measuring in the range of ± 10 °C around the calibration temperature
D.O. partial pressure [mbar]	± 0.5 % of measured value at ambient temperature + 5 °C ... + 30 °C

Correction functions

Temperature compensation	Accuracy better than 2 % at 0 ... + 40 °C
Salinity correction	0 ... 70.0 SAL
Air pressure correction	Automatic through integrated pressure sensor in the range of 500 ... 1100 mbar

FCC Class A Equipment Statement

Note: This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.

Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

8 Lists

This chapter provides additional information and orientation aids.

Abbreviations The list of abbreviations explains the indicators and abbreviations that appear on the display and in the manual.

Specialist terms The glossary briefly explains the meaning of the specialist terms. However, terms that should already be familiar to the target group are not described here.

Trademarks used The list comprises the trademarks used in the present document and their owners.

Index The index will help you to find the topics that you are looking for.

Abbreviations

κ	Conductivity value (international γ)
°C	Temperature unit, degrees Celsius
°F	Temperature unit, degrees Fahrenheit
C	Cell constant [cm ⁻¹] (internat. k)
Cal	Automatic calibration using a selected buffer set
d	Day
h	Hour
j	Year
Lin	Linear temperature compensation
LoBat	Batteries almost empty (Low battery)
m	Month
mV	Voltage unit
mV/pH	Unit of the electrode slope (internat. mV)
nLF	Nonlinear temperature compensation
OxiCal	Automatic calibration for D.O. measurements
pH	pH value
S	Slope (internat. k)
SAL	Salinity
SELV	Safety Extra Low Voltage
TDS	Total Dissolved Solids
Tr20	Reference temperature of 20 °C
Tr25	Reference temperature of 25 °C
U	Voltage

Glossary

Adjusting	To manipulate a measuring system so that the relevant value (e. g. the displayed value) differs as little as possible from the correct value or a value that is regarded as correct, or that the difference remains within the tolerance.
AutoRange	Name of the automatic selection of the measuring range.
Calibration	Comparing the value from a measuring system (e. g. the displayed value) to the correct value or a value that is regarded as correct. Often, this expression is also used when the measuring system is adjusted at the same time (see adjusting).
Cell constant, k	Characteristic quantity of a conductivity measuring cell, depending on the geometry.
Conductivity	Short form of the expression, specific electrical conductivity. It is a measured value of the ability of a substance to conduct an electric current. In water analysis, the electrical conductivity is a dimension for the ionized substances in a solution.
D.O. partial pressure	Pressure caused by the oxygen in a gas mixture or liquid.
D.O. saturation	Short name for the relative D.O. saturation. Note: The D.O. saturation value of air-saturated water and the D.O. saturation value of oxygen-saturated water are different.
Electrode zero point	The zero point of a pH electrode is the pH value at which the electromotive force of the pH electrode at a specified temperature is zero. Normally, this is at 25 °C.
Electromotive force of an electrode	The electromotive force (voltage) U of the electrode is the measurable electromotive force of an electrode in a solution. It equals the sum of all the galvanic voltages of the electrode. Its dependency on the pH results in the electrode function, which is characterized by the parameters, slope and zero point.
Junction	The junction is a porous body in the housing wall of reference electrodes or electrolyte bridges. It forms the electrical contact between two solutions and makes electrolyte exchange more difficult. The expression, junction, is also used for ground or junction-less transitions.
Measured parameter	The measured parameter is the physical dimension determined by measuring, e. g. pH, conductivity or D.O. concentration.
Measured value	The measured value is the special value of a measured parameter to be determined. It is given as a combination of the numerical value and unit (e. g. 3 m; 0.5 s; 5.2 A; 373.15 K).

Molality	Molality is the quantity (in Mol) of a dissolved substance in 1000 g solvent.
Offset potential	The measurable potential of a symmetrical electrode, the membrane of which is immersed in a solution with the pH of the nominal electrode zero point. The asymmetry is part of the offset potential.
ORP voltage	The ORP is caused by oxidizing or reducing substances dissolved in water if these substances become effective on an electrode surface (e. g. a gold or platinum surface).
pH value	The pH is a measure of the acidic or basic effect of an aqueous solution. It corresponds to the negative decadic logarithm of the molal hydrogen ions activity divided by the unit of the molality. The practical pH value is the value of a pH measurement.
Potentiometry	Name of a measuring technique. The signal (depending on the measured parameter) of the electrode is the electrical potential. The electrical current remains constant.
Reference temperature	Fixed temperature value to compare temperature-dependent measured values. For conductivity measurements, the measured value is converted to a conductivity value at a reference temperature of 20 °C or 25 °C.
Reset	Restoring the original condition of all settings of a measuring system.
Resistance	Short name for the electrolytic resistivity. It corresponds to the reciprocal value of the electrical conductivity.
Resolution	Smallest difference between two measured values that can be displayed by a meter.
Salinity	The absolute salinity S_A of seawater corresponds to the relationship of the mass of dissolved salts to the mass of the solution (in g/Kg). In practice, this dimension cannot be measured directly. Therefore, the practical salinity is used for oceanographic monitoring. It is determined by measuring the electrical conductivity.
Salt content	General designation for the quantity of salt dissolved in water.
Slope	The slope of a linear calibration function.
Slope (relative)	Designation used by WTW in the D.O. measuring technique. It expresses the relation of the slope value to the value of a theoretical reference probe of the same construction type.
Stability control	Name for a function to check the stability of the measured value.
Standard solution	The standard solution is a solution where the measured value is known by definition. It is used to calibrate a measuring system.
TDS	Total dissolved solids

Temperature coefficient	Value of the slope of a linear temperature function.
Temperature compensation	Name of a function that considers the temperature influence on the measurement and converts it accordingly. Depending on the measured parameter to be determined, the temperature compensation functions in different ways. For conductimetric measurements, the measured value is converted to a defined reference temperature. For potentiometric measurements, the slope value is adjusted to the temperature of the test sample but the measured value is not converted.
Temperature function	Name of a mathematical function expressing the temperature behavior of a test sample, a sensor or part of a sensor.
Test sample	Designation of the test sample ready to be measured. Normally, a test sample is made by processing the original sample. The test sample and original sample are identical if the test sample was not processed.
Zero point	Designation for the offset potential of a pH electrode. It is the measurable potential of a symmetrical electrode, the membrane of which is immersed in a solution with the pH of the nominal electrode zero point (pH = 7).

Trademarks used

Trademark	Owner of the trademark
Merck	Merck KGaA
Mettler Toledo	Mettler Toledo
Fisher	Fisher Scientific Company
Fluka	Fluka AG
Radiometer	Radiometer
Baker	Mallinckrodt Baker, Inc.
Metrohm	Metrohm AG
Beckman	Beckman Instruments, Inc.
Hamilton	Hamilton Company Corporation
Precisa	Precisa Instruments AG
Windows	Microsoft Corporation

9 Index

A

Air calibration vessel	75
Authorized use	17

B

Battery compartment	95
---------------------------	----

C

Calibration	
ISE	55
Conductivity	67
Oxygen	75
pH	43
Calibration evaluation	
ISE	56
Conductivity	68
Oxygen	76
pH	46
Calibration interval	47
Calibration points	
ISE	55
pH	45
calibration records	87
Cell constant	67
Cleaning	96
Comparison measurement (D.O.)	75
Connecting sensors	12

D

Data filter	86
Dataset	80
Date and time	27
Default settings	
Measured parameter	92
System settings	94
Display illumination	11

F

Filter	86
Firmware update	117

I

Initial commissioning	20
Initialize	92
ISE standard solutions	55

K

Key functions	10
Keys	10

M

Measured value display	23
Measured value memory	
Edit	84
Measurement accuracy	47
Measurement data storage	
Delete	84
Storage locations	80
Measurement dataset	80
Measurement window	11
Measuring	
ISE	53
Conductivity	62
Dissolved oxygen	72
ORP voltage	39
pH	37
Measuring menu	
ISE	54
Conductivity	65
Oxygen	74
pH/ORP	41
Measuring ranges	
Conductivity	105
Memory administration	83
Menus (navigation)	23
Messages	24
MultiLab pilot	92

O

Obligations of the purchaser	18
Operational safety	18

P

pH buffer sets	43
Precautions	17
Print	91

R

Reset	92
RS232 interface	90

S

Safe operation	18
Safety	17
Sample display	11
Saving	80
at intervals	82
automatic	82
manual	81
Scope of delivery	19
Single-point calibration	
pH	49
Slope	
ISE	55
pH	43, 67
Relative (O ₂ sensor)	75
Socket field	12
Stability control	
automatic	35
manual	
Conductivity	62
Dissolved oxygen	72
ISE	53
ORP	39
pH	37
Status line	11
Storage interval	82
System menu	32
General information	34
Sample display	32

T

Temperature compensation	64
Temperature measurement	
ISE	52, 71
Conductivity	61
pH	36
Three-point calibration	
ISE	59
pH	50
Transmitting data	89
Transmitting measured values	89
Two-point calibration	
ISE	58
pH	49

Z

Zero point of the pH electrode	43, 67
--------------------------------------	--------

Appendix

A.1 Firmware update

General information

With the Update_ProLab1000_2000 program and a PC you can update the firmware of the ProLab 2000 to the newest version. Available firmware updates can be found on the Internet.

Connect the meter to a PC for the update.

For the update via the USB interface, the following is required:

- a free USB interface (virtual COM port) on your PC
- the driver for the USB interface (installation of the USB driver from the installation CD-ROM, see page 89)
- the USB cable (included in the scope of delivery of the ProLab 2000).

For the update via the RS232 interface, the following is required:

- a free RS232 interface on your PC
- the RS232 cable, Z390.

Program installation

- | | |
|---|---|
| 1 | With the installation program, "Install_Update_ProLab1000_2000_Vx_yy_English.exe", install the firmware update program on a PC. |
|---|---|

Program start

- | | |
|---|---|
| 2 | Start the "Update_ProLab1000_2000" program from the Windows start menu. |
| 3 | If necessary, change the language in the language menu. |

Firmware update

- | | |
|---|---|
| 4 | Connect the ProLab 2000 to a USB interface (virtual COM port) of the PC with the aid of the USB interface cable.
or
Connect the ProLab 2000 to a serial interface (COM port) of the PC with the aid of the Z390 interface cable. |
| 5 | Switch on the ProLab 2000. |
| 6 | Start the updating procedure with OK in the firmware update program. |
| 7 | Follow the instructions of the firmware update program. During the programming process, a corresponding message and a progress bar (in %) is displayed. The programming process takes approx. two minutes. A terminatory message is displayed after a successful programming process. The firmware update is completed. |


- 8 | Disconnect the ProLab 2000 from the PC.
The ProLab 2000 is ready for operation.

After switching the meter off and on you can check whether the meter has taken over the new software version.

A.2 Menus

A.2.1 Data storage

Storage & config (<MENU/OK__>)



A.2.2 Display

Storage & config (<MENU/OK__>)

- | Data storage (see page 119)
- | Display
 - | pH : On / Off
 - | Cond : On / Off
 - | Oxi : On / Off
 - | Time : On / Off
 - | Date : On / Off
- | System (see page 120)

A.2.3 System

Storage & config (<MENU/OK__>)

- | Data storage (see page 119)
- | Display (see page 120)
- | System
 - | General
 - | Language
 - | Deutsch
 - | English
 - | (etc.)
 - | Beep : On / Off
 - | Illumination : On / Off
 - | Contrast : 0 ... 100 %
 - | Measurement
 - | Temperature unit : ° C / ° F
 - | Stability control : Auto / Off
 - | Interface
 - | Baud rate : 1200, 2400, 4800, 9600, 19200
 - | Output format : ASCII / CSV
 - | Output header
 - | Clock function
 - | Time
 - | Date
 - | Date format
 - | Service information
 - | Reset

A.2.4 <STO __>

<STO __> (<STO __>)

- | SAMPLE ID : 1 ... 999
- | Interval : 1 s, 5 s, 10 s, 30 s,
 1 min, 5 min, 10 min, 15 min, 30 min, 60 min
- | Duration : 1 min ... x min

A.2.5 pH/U

pH/U (<MENU/OK>)

- | Calibration
 - | Calibration record
 - | Buffer
 - | NIST/DIN
 - | TEC
 - | MERCK 1
 - | MERCK 2
 - | MERCK 3
 - | DIN 19267
 - | Mettler US
 - | Mettler EU
 - | Fisher 1
 - | Fluka BS
 - | Radiometer
 - | Baker
 - | Metrohm
 - | Beckmann
 - | Hamilton DC
 - | Precisa
 - | Calibration interval : 1 ... 999 d
 - | Unit for zero point : mV / pH
 - | Unit for slope : mV/pH / %
 - | Man. temperature : -20 ... +130 °C
- | Reset
- | High resolution pH : On / Off
- | High resolution mV : On / Off

A.2.6 ISE

ISE (<MENU/OK>)

- | Calibration
 - | Calibration record
- | Man. temperature : -20 ... +130 °C

A.2.7 Conductivity

Conductivity, Resistance, TDS, SAL (<MENU/OK>)

- | Calibration
 - | Calibration record
 - | Calibration interval : 1 ... 999 d
- | Reset
- | Temp. comp. (TC)
 - | Reference temp. : 20 °C / 25 °C
 - | Compensation
 - | nLF
 - | lin
 - | Off
- | TDS factor : 0.40 ... 1.00
- | Measuring cell : On /Off
 - | cal
 - | C0.01
 - | LF313T
 - | man
- | Cell constant : 0,090 ... 0,110 cm⁻¹
- | Man. cell const. : 0,250 ... 25,000 cm⁻¹

A.2.8 OxygenO₂ (<MENU/OK>)

- | Calibration
 - | Calibration record
 - | Comparison meas. : On /Off
 - | Calibration interval : 1 ... 999 d
- | Reset
- | Sal automatic : On /Off
- | Sal correction : On /Off

| Salinity: 0.0 ... 70.0

SI Analytics GmbH
Postfach 2443
D-55014 Mainz
Hattenbergstr. 10
D-55122 Mainz

Telefon +49 (0) 61 31/66 5111
Telefax +49 (0) 61 31/66 5001
Email: support@si-analytics.com
Internet: www.si-analytics.com

