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
LabPump

Operating Instructions

LabPump

SCHOTT

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Important notes: Before initial operation of the LabPump please read and observe carefully the operating instructions. For safety reasons the pump LabPump may only be used for the purposes described in these present operating instructions.

Please also observe the operating instructions for the units to be connected.

All specifications in this instruction manual are guidance values which are valid at the time of printing. However, for technical or commercial reasons or in the necessity to comply with the statutory stipulations of various countries, SCHOTT may perform additions to the LabPump without changing the described properties.

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
TABLE OF CONTENTS

PAGE

1	Front and rear panel Pump LabPump	6
2	Mission	7
3	Commissioning	7
4	Technical properties	7
5	Error Corrective Action.....	7
6	Warning and safety instructions	8

Front view Pump LabPump:

2 Mission

The range of application for the pump is the delivery of air, gases and vapours. The temperature range of the media to be delivered lies between + 5°C and + 40°C.

The Pump LabPump will be used with hand and half automatic measurements for sucking and pumping up the solutions in viscometers:

- Filling viscometers.
- Rinsing with next sample.
- Sucking up between the manual measurements.
- Discharging viscometers without demounting out of the thermostat.

Because all used materials and connections of the Pump LabPump are built in PPS (Pumphead)), PTFE or stainless steel, the pump is suitable for applications with aggressive mediums.

The LabPump pump is not permitted to be used in an explosion-endangered environment.

The mission for use at half automatic processing of samples in use with a viscosity measuring unit, from SCHOTT-GERÄTE is possible up to a viscosity range of 30 000 mm²/s at 25°C. For half automatic working will be used a PTFE tubing set with waste bottle type no. VZ 5624. When working with a multiple suction, the hose kit VZ 5614 is employed.

3 Commissioning

With the power switch (green) on the front panel the Pump LabPump will be switched on and will be connected with the main power.

Also on the front panel operates the switch for starting the Pump LabPump.

With the aid of the included foot switch it is possible to switch the pump on and off, as and when required. In order to use the foot switch, it is connected to the pump instead of the mains cable. In order to operate the pump, switch on the pump on the mains switch and start the pumping action by using the foot switch when required.

With the turning valve knob on the front the suction rate will be setted in accordance to the symbol (turning right is corresponding to bigger suction rate, turning left is corresponding to lower suction rate).

4 Technical properties

Main power:	Power consumption 100 VA / 50 /60 Hz LabPump 230V: 230 V~ / 0,65 A LabPump 115V: 115 V~ / 1,2 A
Power connection:	Euro type main power plug with fuse (according to IEC 320-C14)
Protection system:	IP 20
Overload protection:	via thermal switch
pump capacity:	5,5 L/ min
final vacuum:	160 mbar absolute
max. pressure:	2,5 bar ü
Materials:	Pumphead: PPS membran: EPDM PTFE coated valve: FFPM
Weight:	approx. 1,9 kg

5 Error Corrective Action

ATTENTION: before working on the pump, disconnect the power supply by removing the mains switch from the mains supply.

Error: Pump does not deliver:

- The thermal switch has acted due to overheating of the pump.
- > Remove the pump from the mains supply, let it cool down, find the cause for the overheating and remove cause.
- Connections or pipelines are blocked.
- Liquid (condensation) has collected in the pump head.
- > Operate the pump with the media air for a few minutes.
- > Mount pump at the highest point of the system.
- Membrane or valve plates have been worn out.
- > Send in pump for repairs.

Error: Delivery volume, pressure or vacuum are too low:

- Leakages at the connections, the pipelines or at the head of the pump
- Membrane or valve plates are defective or soiled.

NOTE:

- In the event that you need to send in your pump for repairs, we would kindly ask you to tell us the media being delivered. Our Service engineer needs to be informed about delivery media, in particular in the case of aggressive.
- If you have delivered hazardous or aggressive gases, we would ask you to clean the pump before sending it off.

6 Warning and safety instructions

For safety and functional reasons, the Pump LabPump must basically be opened by authorized persons only; e. g. work at the electric appliance must only be carried out by trained specialists. With unauthorized operation of the device as well as careless or deliberate damaging, the warranty is no longer valid.

The Pump LabPump has protection class I. It has been built and tested according to protective arrangements for electronic measuring instruments DIN 57 411, part 1 / VDE 0411, part 1, and has left our plant in perfect condition concerning safety regulations. To maintain this condition and to guarantee safe operation the user has to observe the instructions and warning marks which are included in these operating instructions.

Before switching on, it has to be guaranteed that the operating voltage and supply voltage correspond. Operating voltage is indicated on the type plate. The mains plug must only be inserted in a socket outlet with earthing contact. The protective effect must not be invalidated by an extension line without earthed lead. Any interruption of the earthed lead inside or outside the device or detaching the earthed lead connection could make the device dangerous. Deliberate interruption is not permissible.

It has to be guaranteed that fuses and rated current strength indicated are used as substitute. The use of patched fuses or short-circuiting of fuse holder is not permissible.

The built-in safety devices should never be made inoperative.

If there is reason to believe that non-dangerous operation is not possible, the device should be stopped and secured against unintentional operation. Turn off the Pump LabPump, remove the line cord and call the SCHOTT-GERÄTE service department.

It is reasonable to assume that safe operation of the Pump LabPump is no longer possible,

- if the Pump LabPump is visibly damaged,
- if the Pump LabPump does not operate properly,
- if a major damage during transportation has happened.

The Pump LabPump must not be stored or operated in damp rooms.

For security reasons, the device must only be used for purposes described in these instructions.

The LabPump pump is not permitted to be used in an explosion-endangered environment.

Please note the corresponding operating instructions for to be connected devices.

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Supplier's Certificate

We certify that the equipment was verified according EN ISO 9001, part 8.2.4 "Monitoring and measurement of product" and that the specified requirements for the product are met.

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SCHOTT-GERÄTE GmbH