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
Originalversion

OPERATING MANUAL

LabPump

Pumpe | Pump

SI Analytics

a **xylem** brand

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Important notes:

The operating manual is part of the product. Before initial operation of the unit, please carefully read and observe the operating instructions and keep it. For safety reasons the product may only be used for the purposes described in these present operating instructions. Please also observe the operating instructions for the units to be connected

All specifications in this operating manual are guidance values which are valid at the time of printing. However, for technical or commercial reasons or in the necessity to comply with the statutory stipulations of various countries, the manufacturer may perform additions to the product without changing the described properties.

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TABLE OF CONTENTS

1	Warning and safety information	11
2	Working with the LabPump.....	12
2.1	Setup of the device.....	12
2.2	Application scope	12
2.3	Commissioning.....	13
2.4	Specifications	13
2.5	Error Corrective	13

Version 161028 US

Notes to the manual

The provided manual will allow you the proper and safe handling of the product.
For maximum security, observe the safety and warning instructions in the operating manual!

The pictogram has the following meaning:

 Warning of a general danger to personnel and equipment.
Non-compliance results (can result) in injury or material damage.

 Important information for device use.

 Refers to another part of the operating manual.

Status at time of printing

Advanced technology and the high quality of our products are guaranteed by a continuous development. This may result in differences between this operating manual and your product. We cannot exclude mistakes. We are sure you understand that no legal claims can be derived from the information, illustrations and descriptions.

A potentially more recent version of this manual is available on our internet website at www.si-analytics.com. The German version is the original version and binding in all specifications.

Garantieerklärung

We provide guarantee for our device of one year from the date of purchase. This guarantee covers manufacturing faults being discovered within the mentioned period of one year. Claim under guarantee covers only the sensor itself, not any further claim for damages or financial loss. Warranty claims shall not include minor deviation from the agreed quality, of only minor impairment of usefulness, of usual wear and any damage that occurs after the transfer of risk from faulty handling, excessive strain, unsuitable equipment or due to special external influences.

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1 Warning and safety information

The LabPump corresponds to protection class I.

It was manufactured and tested according to DIN EN 61 010, Part 1, "**Protective Measures for electronic measurement devices**" and control devices and has left the factory in an impeccable condition as concerns safety technology. In order to maintain this condition and to ensure safe operation, the user should observe the notes and warning information contained in the present operating instructions. Development and production is done within a system which meets the requirements laid down in the DIN EN ISO 9001 standard.

 For reasons of safety, the device and the power supply must be opened by authorised persons only; this means, for instance, that work on electrical equipment must only be performed by qualified specialists. **In the case of nonobservance of these provisions the device and the power supply may constitute a danger: electrical accidents of persons or fire hazard.** Moreover, in the case of unauthorised intervention in the device or the power supply, as well as in the case of negligently or deliberately caused damage, the warranty will become void.

 Prior to switching the device on it has to be ensured that the operating voltage matches the mains voltage. The operating voltage is indicated on the specification plate (Underside of the device and the power supply). The mains plug must only be inserted in a socket outlet with earthing contact. The protective effect must not be invalidated by an extension line without earthed lead. Any interruption of the earthed lead inside or outside the device or detaching the earthed lead connection could make the device dangerous. Deliberate interruption is not permissible. **Nonobservance of this provision may result in damage to the device and the power supply, or in personal injury or damage to property!**

 **If it has to be assumed that safe operation is impossible, the device has to be put out of operation and secured against inadvertent putting to operation.** In this case please switch the device off, pull plug of the mains cable out of the power supply, and remove the device from the place of work.

Examples for the assumption that a safe operation is no longer possible,

- if the package is damaged,
- if the LabPump shows visible damages,
- if the power supply shows visible damages,
- if the LabPump does not function properly,
- if the LabPump has been altered technologically or if unauthorized personnel tried or succeeded to open the instrument as attempt to repair it.

In case that the user operates such a device, all thereof resulting risks are on the user.


 The LabPump is not permitted to be used in an explosion-endangered environment.

 **The relevant regulations regarding the handling of the substances used have to be observed:** The Decree on Hazardous Matters, the Chemicals Act, and the rules and information of the chemicals trade. It has to be ensured on the side of the user that the persons entrusted with the use of the device are experts in the handling of substances used in the environment and in the device or that they are supervised by specialised persons, respectively.

 During all work with titration solutions: **Please wear protective glasses!**

2 Working with the LabPump

2.1 Setup of the device

2.2 Application scope

The range of application for the pump is the delivery of air, gases and vapours. The temperature range of the media to be delivered lies between + 5°C and + 40°C.

The Pump LabPump will be used with hand and half automatic measurements for sucking and pumping up the solutions in viscometers:

- Filling viscometers
- Rinsing with next sample
- Sucking up between the manual measurements
- Discharging viscometers without demounting out of the thermostat.

Because all used materials and connections of the Pump LabPump are built in PPS (Pump head), PTFE or stainless steel, the pump is suitable for applications with aggressive mediums.

⚠ The LabPump pump is not permitted to be used in an explosion-endangered environment.

The mission for use at half automatic processing of samples in use with a viscosity measuring unit, from SI Analytics®, is possible up to a viscosity range of 30 000 mm²/s at 25°C. For half automatic working will be used a PTFE tubing set with waste bottle type no. VZ 5624. When working with a multiple suction, the hose kit VZ 5614 is employed.

2.3 Commissioning

With the power switch on the front panel the LabPump will be switched on and will be connected with the main power.

Also on the front panel operates the switch for starting the LabPump.

With the aid of the included foot switch it is possible to switch the pump on and off, as and when required. In order to use the foot switch, it is connected to the pump instead of the mains cable. In order to operate the pump, switch on the pump on the mains switch and start the pumping action by using the foot switch when required.

i With the turning valve knob on the front the suction rate will be setted in accordance to the symbol:

- turning right is corresponding to bigger suction rate,
- turning left is corresponding to lower suction rate.

2.4 Specifications

Main power:	Power consumption 100 VA / 50/ 60 Hz LabPump 230V: 230 V~ / 0,65 A LabPump 115V: 115 V~ / 1,2 A
Power connection:	Euro type main power plug with fuse (according to IEC 320-C14)
Protection system:	IP 20
Overload protection:	via thermal switch
Capacity:	5.5 L/ min
final vacuum:	160 mbar absolut
max. pressure:	2.5 bar ü
Materials:	Pump head: PPS Membrane: EPDM PTFE coated Valve: FFPM
Weight:	approx. 1.9 kg

2.5 Error Corrective

 before working on the pump, disconnect the power supply by removing the mains switch from the mains supply!

Error: Pump does not deliver:

- The thermal switch has acted due to overheating of the pump.
 - Remove the pump from the mains supply, let it cool down, find the cause for the overheating and remove cause.
- Connections or pipelines are blocked.
- Liquid (condensation) has collected in the pump head.
 - Operate the pump with the media air for a few minutes.
 - Mount pump at the highest point of the system.
- Membrane or valve plates have been worn out.
 - Send in pump for repairs.

Error: Delivery volume, pressure or vacuum are too low:

- Leakages at the connections, the pipelines or at the head of the pump.
- Membrane or valve plates are defective or soiled.

i

In the event that you need to send in your pump for repairs (Service see backside of this manual), we would kindly ask you to tell us the media being delivered. Our Service engineer needs to be informed about delivery media, in particular in the case of aggressive. If you have delivered hazardous or aggressive gases, we would ask you to clean the pump before sending it off.

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Supplier's Certificate

We certify that the above equipment has been tested in accordance with DIN EN ISO 9001, Part 8.2.4 "Monitoring and measurement of product" and that the specified quality requirements for the product have been met.

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Copyright
Xylem Analytics Germany GmbH
Xylem Analytics Germany Sales GmbH & Co.KG