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**Mikro - Ostwald
Viskosimeter**

Operating instructions

**Micro Ostwald
viscometer**

Mode d'emploi

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Manual de instrucciones

**Micro viscosímetro
de Ostwald**

SI Analytics

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Operating Instructions Page 6 8

Important notes: Before initial operation of the micro Ostwald viscometer please read and observe carefully the operating instructions. For safety reasons the micro Ostwald viscometer may only be used for the purposes described in these present operating instructions.

Please also observe the operating instructions for the units to be connected.

All specifications in this instruction manual are guidance values which are valid at the time of printing. However, for technical or commercial reasons or in the necessity to comply with the statutory stipulations of various countries, SI Analytics GmbH may perform additions to the micro Ostwald viscometer without changing the described properties.

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
CONTENTS

PAGE

1	Preparation of the sample	7
2	Selection of the capillary	7
3	Cleaning the viscometer	7
4	Measuring operation	7
5	Calculation of viscosity	8
6	Example of calculation	8
7	Table of the Hagenbach Correction (HC)	8

Operating Instructions

Ostwald Micro Viscometer

1 Preparation of sample

Low-viscosity samples must be filtered through a SCHOTT glass filter porosity 2 to 4 (10 ... 100 μm) before the measurement; high-viscosity samples, through a sieve with 0.3 mm mesh width (test sieve cloth 0.2, DIN 4188). Samples, whose stock value in accordance with DIN 51 583 or pour point in accordance with DIN 51 597 is not at least 30 °C lower than the test temperature, must be heated up to 50 °C before the measurement.

2 Selection of capillary

The size of the capillary is to be selected in such a way that the uncertainty inherent in the Hagenbach Correction does not exceed the allowable error for the time measurement (see table). For precision measurements should, therefore, be chosen no flow times below 30 seconds. If necessary, a viscometer with a smaller capillary is to be used.

3 Cleaning of viscometer

Before the first use, a cleaning with 15 % H_2O_2 and 15 % HCl is recommended. The viscometer should then be rinsed with a suitable solvent. It must be completely dry and dust-free.

4 Measuring operation

Exactly 2 ml of the sample liquid are filled with a pipette into the wider one of the viscometer tubes.

Hang viscometer with its stand Type No. 053 97 into a Glass-Panelled Thermostatic Bath from SI Analytics GmbH.

To obtain the measuring accuracy of the viscometer, the constant temperature bath should maintain the set temperature at a constant ± 0.01 °C (**glass panelled thermostatic baths from SI Analytics GmbH**). Differences in temperature of only 0.1 °C may cause an error of as much as 0.6 % in mineral oils. Measuring should take place only after an equilibration time of approx. 5 minutes.

The liquid is siphoned above the upper measuring mark M1. Then the flow time between the two timing marks M1 and M2 is measured. The measuring operation can be repeated as often as necessary.

When using the viscometers in automated viscosity-measuring units (AVS®) by SI Analytics GmbH, the viscosity is measured automatically.

Subjective measuring errors are eliminated, and the efflux times measured are available as a printed-out document. According to the type and number of samples to be measured, an optimum measuring device may be assembled which can be expanded to an automatic sampler. Rinsing and filling of viscometers and measuring of sample with subsequent changing of sample is performed automatically. The timing marks required for manual measuring are replaced by light barriers.

The accuracy obtained with the AVS® automated measuring system is greater since certain parameters such as errors in reading, clock errors, etc. are eliminated.

5 Calculation of viscosity

The seconds contribution given in the table for the Hagenbach Corrections is to be subtracted from the measured flow time for the various capillaries. Intermediate values can be interpolated.

With absolute measurements, the corrected flow time multiplied by the viscometer constant K, produces the kinematic viscosity [mm²/s]* directly.

$$\nu = K (t - \vartheta)$$

The viscometer constant K is mentioned in the enclosed production certificate.

6 Example of calculation

Micro Ostwald Viscometer

Type No. 516 10

Capillary I

Constant = 0.0100

Flow time (averaged) = 40.00 s

Hagenbach Correction for 40,00 s $\vartheta = 0.18$ s

Kinematic viscosity $\nu = K (t - \vartheta)$
 = 0.0100 (40,00 – 0,18)
 = 0.3982 [mm²/s]*

* previously centistokes [cSt]; 1 cSt = 1 mm²/s

7 Table of Hagenbach Corrections (HC) for:

Micro Ostwald Viscometers

Type No. 516 .., 517 .., 518 ..

Correction seconds¹⁾:

Flow time [s]	Capillary no. I	Capillary no. Ic
30	0.32	0.05
35	0.23	0.03
40	0.18	0.03
45	0.14	0.02
50	0.11	—
55	0.09	—
60	0.08	—
65	0.06	—
70	0.06	—
75	0.05	—
80	0.04	—

¹⁾ The correction seconds stated are related to the respective theoretical constant.

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