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Ostwald-Viskosimeter

Operating Instructions

Viscometer Ostwald

Mode d'emploi

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SI Analytics

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Operating Instructions Page 5

Important notes: Before initial operation of the Ostwald viscometer please read and observe carefully the operating instructions. For safety reasons the Ostwald viscometer may only be used for the purposes described in these present operating instructions.

Please also observe the operating instructions for the units to be connected.

All specifications in this instruction manual are guidance values which are valid at the time of printing. However, for technical or commercial reasons or in the necessity to comply with the statutory stipulations of various countries, SI Analytics GmbH may perform additions to the Ostwald viscometer without changing the described properties.

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1	Description	5
2	Preparation of the sample	5
3	Cleaning of viscometer	5
4	Measuring operation	5
5	Table	5

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Drückhöhe: 115 mm
Füllmenge: 3 ml

Operating Instructions Viscometer Ostwald

1 Description

The Viscometer Ostwald is not only used for viscosity measurement, but also very frequently for determination of the molecular weight of high density polymers by means of viscometry. Its dimensions are adjusted to the capillary diameters now commonly used for this application, and identified by different colours.

2 Preparation of the sample

Low-viscosity samples shall be passed through a glass filter porosity 2 ... 4 (10 ... 100 µm) before measuring, and high-viscosity samples through a strainer of 0.3 mm mesh (test strainer 0.2, DIN 4 188). Samples whose setting point according to DIN 51 583 or pourpoint according to DIN 51 597 is not at least 30 °C below, the test temperature must be heated up to 50 °C before measuring.

3 Cleaning of viscometer

Before first use, a cleaning with 15 % H₂O₂ and 15 % HCl is recommended. The viscometer should then be rinsed with a suitable solvent. The viscometer must be completely dry and dust-free and can thereby be used for automatic measurements.

4 Measuring operation

3 ml of the sample liquid are filled with a pipette into the wider one of the viscometer tubes. Thereafter, the viscometer with mounting ref. no. 064 99 is suspended in a constant temperature bath.

About 10 minutes are required for equalisation of temperature. The liquid is siphoned above the upper measuring mark. Then the flow time between the two timing marks is measured. The measuring operation can be repeated as often as necessary.

With regard to the preferred application of Ostwald Viscometers for the determination of the molecular weight, the determination of the viscometer constant by precision calibration can be dispensed with.

As various solvents are used for the determination of the molecular weight, the corresponding capillaries are to be selected.

The Hagenbach correction is $\vartheta = \frac{0.12}{K \cdot t}$

The calculated correction seconds are to be subtracted from the flow time t.

Kinematic viscosity $\nu = K (t - \vartheta)$

K is either taken from the following table, or, for accurate measurements, determined by means of calibrated reference viscometers such as Ubbelohde Viscometers of SI Analytics GmbH (ref. no. 501 ..), or by means of liquids of known viscosity.

5 Table

Ref. no.	Capillary Ø i.d. (mm)	Flow time of water (s) approx.	Constant K (approx.)	Measuring above mm ² /s (cSt) approx.
509 03	0.3	250	0.004	0.3
509 04	0.4	75	0.01	1
509 05	0.5	30	0.03	2.5
509 06	0.6	15	0.07	5.5
509 07	0.7	10	0.1	10

Common features for all sizes: Capillary length: 100 mm Hydrostatic head: 115 mm
 Flow rate: 0.5 ml Filling capacity: 3 ml

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