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Laborkocher

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Operating Instructions

Laboratory hot plates

with glass-ceramic heating surface

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Operating Instructions Page 12 22

Important notes: Before of the laboratory hot plate is first put into operation, please carefully read and observe the operating instructions. For reasons of safety, the laboratory hot plate with glass-ceramic material may only be used for the purposes described in these operating instructions.

All the information provided in this instruction manual is applicable at the time of printing. However, for technical or commercial reasons or in order to comply with national legislation in the various countries, SI Analytics may need to make additions and/or changes to the laboratory hot plate with glass-ceramic material, without this having any impact on the properties described.

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CONTENTS**PAGE**

1 Properties	13
General	13
Warranty.....	13
Physical properties of the glass-ceramic material	13
Chemical properties of the glass-ceramic material.....	13
Technical characteristics.....	14
2 Warning and safety instructions	15
3 Set-up and putting into operation	17
General information, correct use.....	17
Initial cleaning and heating	18
4 Working with the laboratory hot plates	18
Control panel of the SLK 1 and SLK 2 laboratory hot plates.....	18
Detailed explanations on the Chapter: "Working with the laboratory hot plates".....	19
5 Care, cleaning and maintenance.....	21
Care	21
Cleaning	21
Maintenance.....	21
6 Malfunction, possible defects	22
7 Accessories	22
Declaration of conformity.....	at the documents end

Notes on the operating instructions

The present manual is designed to enable you to use the laboratory hot plates safely for their intended purpose.

To ensure maximum possible safety, you should always observe all safety and warning instructions. The pictogram used has the following meaning:

Warning of a general hazard for people and materials.

Failure to follow these instructions may result in physical injury or material damage.

1 Properties

General

The laboratory hot plates with glass-ceramic heating surface are used primarily to heat up aqueous solutions not containing any inflammable components in glass vessels. As the interior of the laboratory hot plates and the glass-ceramic heating surface become hot, all other applications should be checked carefully to establish whether they could give rise to hazardous situations.

The glass-ceramic material of the SI Analytics laboratory hot plates has special thermal properties. The glass-ceramic heating surface is virtually free of thermal expansion and can withstand temperature fluctuations from - 200 °C to + 700 °C. Ice-cold water can be poured over a hot glass-ceramic heating surface without causing any damage.

Besides its thermal properties, the glass-ceramic material also has the chemical stability and surface quality of glass. Owing to its high infrared permeability, heat energy is transmitted quickly with only little loss. The residual-heat indicator of the laboratory hot plates warns of burning hazards. The residual heat can be used in an energy-saving manner to pre-heat another vessel.

The glass-ceramic surfaces remain level, are pore-free and do not have any restricting frames or corners where dirt could accumulate. The laboratory hot plates of SI Analytics are therefore easy to clean.

The back panel of the laboratory hot plates contains a borehole with an M8 thread to take a 10 mm Ø tripod rod (not included in the scope of supply, Order No. Z 601).

The functions of the laboratory hot plates are controlled by sensor keys.

Limitations on use are described in the chapter "Warning and safety instructions".

Warranty

We hereby guarantee the designated piece of equipment against any manufacturing defects arising within two years from the date of purchase. The warranty covers restoration of functionality, but no further claims to compensation.

In case of improper use or unauthorised opening of the device, the warranty will be rendered invalid. Parts subject to wear, such as for example the radiant heater elements, are excluded from the warranty, as are breakage of the glass-ceramic plate and corrosion damage resulting from improper use in an aggressive atmosphere.

To establish obligations under the guarantee, we kindly request you to send the equipment to us carriage or postage-free together with the sales receipt showing the date of purchase.

Physical properties of the glass-ceramic material

Coefficient of thermal expansion α :	+ 20 ... + 300 °C: $- 0.25 \cdot 10^{-6} \text{ K}^{-1}$ + 20 ... + 700 °C: $+ 0.10 \cdot 10^{-6} \text{ K}^{-1}$
Density ρ :	2.58 g / cm ³
Module of elasticity E:	92 kN / mm ²
Knoop hardness under 1N load:	> 575 HK (test according to ISO 9385)
Maximum temperature stability: Tmax permanent:	700 °C; temporary (max. 1 min.): $\leq 850 \text{ °C}$
Static load capacity: surface load:	0.1 kp / cm ² (approx. 1 N / cm ²)

Chemical properties of the glass-ceramic material

Hydrolytic resistance	DIN ISO 719:	Hydrolytic class 1
Alkali resistance	DIN 52 322 / ISO 685:	Class 2
Acid resistance	DIN 12 116:	Class 3

The **SLK 1** and **SLK 2** laboratory hot plates including their radiant heater elements in the respective heating zone are suitable for quickly heating up liquids. The highest heat intensity, i.e. << 9 >>, corresponds to a heating output of 1.2 kW or 1.8 kW, respectively.

An over-temperature protection facility prevents the laboratory hot plates from overheating.

SI Analytics

Technical characteristics SLK 1 and SLK 2 laboratory hot plates

Status: 11 January 2011

CE Mark:		Electromagnetic compatibility according to Council Directive 2004/108/EEC; low-voltage directive according to Council Directive 2006/95/EEC.
Country of origin:		Germany
Display:		7-segment light-emitting diode display (LED), height 10 mm
Mains connection:		Integrated European-type plug: cold-device plug DIN VDE 0625 EN 60 320 , IEC 320 / C 14 / C 18 (DIN 49 457 B9),
Safety measure:		Device of Protection Class II, manufactured and tested in accordance with safety measures for electric devices; not suitable for use in explosive environments.
Power supply:		Mains: 230 V~, 50 / 60 Hz
Casing material:		Composite HUP 25/25 RN (SMC on unsaturated polyester-resin basis)
Climate:		
Ambient temperature:		+ 10 ... + 40 °C for operation and storage
Atmospheric humidity:		in accordance with EN 61 010, Part 1: maximum relative humidity 80 % for temperatures up to 31 °C, decreasing linearly to 50 % relative humidity at a temperature of 40 °C

Overview of the two SI Analytics laboratory hot plates:

Technical characteristics	SLK 1	SLK 2
Dimensions approx. w x h x d [mm]	295 x 110 x 390	295 x 110 x 390
Glass-ceramic surface approx. w x d [mm]	280 x 280	280 x 280
Glass-ceramic heating surface approx. [mm]	Ø 165	Ø 200
Type of heating element	Infrared	Infrared
Heating output [W]	1,200	1,800
Weight, approx. [kg]	3.6	3.6

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2 Warning and safety instructions

For safety and functional reasons, the laboratory hot plate may as a general rule only be opened by authorised persons; this means, for example, that only trained specialists may carry out any work on the electrical system. In case of unauthorised manipulation of the laboratory hot plate and negligent or deliberate damage, the warranty will lapse.

The laboratory hot plate complies with Protection Class II. It has been manufactured and tested in accordance with DIN VDE 0411, Part 1 (EN 61 010, Part 1) and DIN VDE 0411, Part I 2-010 (EN 61 010, Part 2-010), safety measures for electric heating devices. It has left the plant in perfect safety condition. In order to maintain this condition and ensure safe operation, the user must observe the information and warnings given in these operating instructions.

Prior to switching on, care must be taken to ensure that the voltage indicated on the rating plate of the laboratory hot plate matches the mains voltage.

Please make sure that connecting cables do not come into contact with the hot heating zone or any other hot objects. The laboratory hot plate must be kept at least 0.5 m away from inflammable substances.

The decoration marks the centre of the heating surface. However, the surface outside the decorated area also becomes hot.

Please prevent hard objects from falling on the glass-ceramic surface. Point impact loads may lead to breakage of the glass-ceramic surface.

If cracks, fissures or breaks should appear in the glass-ceramic surface and if there is reason to assume that safe operation is no longer possible, the laboratory hot plate must be put out of operation (remove mains cable) and secured against inadvertent usage.

The temperature of the heating surfaces can reach max. 555 °C. Therefore:

CAUTION: Risk of deflagration, explosion and fire when heating flammable liquids!
Use only to heat up liquids with a flash point exceeding 580 °C!

Be aware of increased risk of ignition of hot liquids!

Always cover vessel to prevent hot liquids from coming into contact with the heating surface (e.g. by liquids splashing or boiling over or through escaping gases)!

The SI Analytics laboratory hot plate must be set up and connected in such a way that maximum possible safety is guaranteed for people and materials.

If the laboratory hot plate is used outside the laboratory, e.g. in the office, workshop or home, it is important to bear in mind that the device may give rise to hazards typical of such areas. Do not use the laboratory hot plates to heat rooms (improper use!).

 The surface of the heating devices becomes hot during operation, so please be careful: such devices should **generally be kept out of reach of children**, in particular **infants**.

 Food prepared using fat or oil (e.g. French fries) should only be prepared under supervision. Food should not be prepared in aluminium foil or plastic containers on the hot heating surfaces! Such substances, especially sugar in either solid or liquid form (foods containing sugar), will melt and become stuck and could cause cracks or breaks in the glass-ceramic surface when cooling down.

Avoid repeatedly burning in residues or boiled-over liquids. Calcareous deposits will harm the glass-ceramic surface.

If the heating surface can no longer be switched off owing to a defect in the sensor control, immediately disconnect the laboratory hot plate from the mains and remove the mains cable.

The laboratory hot plate is not to be used as a storage or working surface.

The laboratory hot plate **must not** be placed in furniture niches during operation.

The laboratory hot plate **must never be immersed in water.**

Explosion hazard!

Never use the laboratory hot plate in places where there is risk of explosion!

Danger of electric shocks!

Never use the laboratory hot plate in areas with running water and high humidity, e.g. bathrooms.

Fire hazard!

For safety reasons, keep the laboratory hot plate at least 50 cm away from inflammable material!

Risk of tripping!

Never route connection cables in highly frequented areas!

Place the laboratory hot plate on a stable, horizontal surface. The area under the hot plate must be non-combustible!

Choose a flat, clean, dry and non-slip surface.

For safety reasons, the laboratory hot plate may only be used for the purposes described in these operating instructions.

⚠ When working with harmful and aggressive media:

Risk of poisoning and chemical burns! Risk of damage to the device by aggressive gases or vapours.

Only use the laboratory hot plate underneath an exhaust system!

Operation under a recirculating air conduit may accelerate damage to the device!

When using the optional compressed-air facility (Order No. 28 541 6595), the laboratory hot plate will be supplied with compressed air.

⚠ Failure to comply with the above information may cause the laboratory hot plate to give rise to hazards such as electrical accidents to human beings or fire hazards. In the case of unauthorised manipulation of the laboratory hot plate and negligent or deliberate damage, the warranty will lapse. ⚠

Laboratory hot plate SLK 1

3 Set-up and putting into operation

When unpacking the device, make sure you also remove the accessories from the packaging. Please contact the manufacturer if any parts are missing, and the forwarding agent or supplier in case of any visible damage.

Place the laboratory hot plate on a flat, horizontal and non-slip surface, which should not be surrounded by any combustible materials. A space of at least 0.5 m must be left clear on all sides of the hot plate.

Insert a 10 mm Ø tripod rod into the holder on the back (M8 female thread) and tighten the tripod rod manually.

Before connecting to the mains, check that the details on the rating plate agree with the mains voltage! Plug the mains cable into the integrated European-type socket in the back panel of the laboratory hot plate. After connecting the cable to the mains, the sensor, glass-ceramic and ambient conditions will be automatically calibrated. The “KeyLock” LED lamp will light up during this process.

During calibration, the glass-ceramic surface in the vicinity of the sensors must not be covered with any objects, nor may the area around the sensors be exposed to any strong light source such as sunlight, powerful halogen lamps, etc. The “Keylock” LED lamp will briefly light up to indicate that the laboratory hot plate is being supplied with mains voltage.

General information, correct use

The decoration marks the centre of the heating surface. However, the surface outside the decorated area also becomes hot. You can quickly interrupt the heating process by simply moving the vessel out of the hot zone to colder adjacent areas.

The laboratory hot plate is equipped with a 2-stage residual-heat indicator << H >> and << h >>. The << H >> display lights up immediately after the hot laboratory hot plate has been switched off. After a short time, the display will change to << h >> indicating that there is no longer any significant residual heat. Likewise, residual heat will be indicated after about 10 seconds at the << 0 >> heat intensity setting. As long as the << H >> symbol is illuminated, the switched off heating zone of the laboratory hot plate is still hot and can be used to save energy. This residual-heat indicator warns of the risk of burning. Once the hot plate has cooled down sufficiently, the residual-heat indicator will go off. Caution! If you pull the mains plug out of the mains socket or in case of a power failure, the residual-heat indicator will stop functioning. As soon as the heating process is resumed, the residual-heat indicator will be reactivated.

Caution: Danger of burning!

Always disconnect the laboratory hot plate after use.

After 5 hours of heating at the << 9 >> setting without the sensor surface being touched, the safety cut-out circuit will switch off the laboratory hot plate. When used at a heat intensity between << 1 >> and << 8 >>, the laboratory hot plate will be switched off by the safety cut-out circuit after 120 hours to prevent the glass-ceramic heating surface from overheating. The over-temperature protection facility protects the laboratory hot plate against overheating.

If several sensor keys are touched within 1 second, the laboratory hot plate interprets this as an operating error, which is ignored (safety measure). An interval of at least 2 seconds must be left between touching keys.

Please observe the instructions for cleaning and care. In particular, you should keep the fields of the sensor keys clean, as soiling here may be misinterpreted as finger contact.

If sugar, plastic or aluminium foil accidentally come into contact with the hot spots of the heating surfaces, **do not** under any circumstances switch the laboratory hot plate off, but **immediately** remove these substances with a razor-blade scraper.

Caution: Danger of burning!


After using the razor-blade scraper, retract the blade properly. Danger of injury! Then clean the heating surfaces when cool.

Initial cleaning and heating

When the laboratory hot plate is heated up for the first time or after a long period out of use, a temporary smell will develop, which is caused by small particles of dust burning and water inside the insulation system evaporating. To remove any accessible dust particles from the glass-ceramic heating area, use a mild detergent to clean the glass-ceramic surface. Then rub the surface dry with a soft cloth.

4 Working with the laboratory hot plates

Control panel of the SLK 1 and SLK 2 laboratory hot plates

- 1 = Heating symbol.
- 2 = Display for heat intensity settings 0 ... 9.
- 3 = LED lamp for "Programming heat intensity" activated.
- 4 = Sensor keys for heat intensity / timer operating time:
 - setting a higher heat intensity/ longer timer operating time.
 - setting a lower heat intensity/ shorter timer operating time.
- 5 = LED lamp for "Programming timer" "Heating zone" activated.
- 6 = Sensor key "Timer On / Off", Programming timer.
- 7 = Timer display indicating remaining time (in seconds), 00 ... 99.
- 8 = LED lamp for "Programming timer" "Egg timer" activated.
- 9 = Sensor key "On / Off".
- 10 = Sensor key: "Lock keys".
- 11 = LED lamp for "Lock keys" activated.

The symbols or for residual heat appear in the LED heat intensity display.


The symbol for hot surfaces is located on the upper side of the glass-ceramic heating surface.

Detailed explanations on the Chapter: “Working with the laboratory hot plates”.

“On / Off” sensor key, for switching the laboratory hot plate on and off (On / Off switch).

- To switch on the laboratory hot plate, press this key for at least 2 seconds. An acoustic signal will indicate successful operation.
- If no other sensor is touched within 10 seconds of switching on the laboratory hot plate, the laboratory hot plate will automatically switch off.
- To switch off, press this key for at least 2 seconds. An acoustic signal will indicate successful operation.
- When switching off the laboratory hot plate, all settings with the exception of those of the “Egg timer” and an activated key lock will for safety reasons be switched down to 0.

Sensor keys for increasing <↑> or decreasing <↓> the heat intensity and/or timer operating time.

- An acoustic signal confirms every successfully completed operation.
- If both keys are touched at the same time, the system will switch to << 0 >>.
- << H >> in the display indicates the presence of residual heat (see Chapter on Putting into operation).

⚠ Danger of burning ⚠

Display of the heat intensity settings 0 ... 9.

- setting << 1 >> = lowest heating output
- setting << 9 >> = highest heating output

Lock keys:

By pressing the key <Lock keys>, the control system can be put into a safe operating mode; in this mode no changes can be made to the heating levels.

In locked mode, the laboratory hot plate can only be switched off.

To lock the keys, press <Lock sensor key> for at least 2 seconds.

Locking is confirmed by an acoustic signal.

To indicate that the key is locked, a LED lamp will light up above <Lock sensor key>.

If a sensor key is pressed in this mode, there will be no reaction.

The key lock will remain activated after the laboratory hot plate has been switched off.

Unlock keys:

To unlock the control system again, press <Lock sensor key> for at least 2 seconds. Release is confirmed by an acoustic signal.

The LED lamp above <Lock sensor key> goes out again.

The laboratory hot plate can now be used again without any restriction.

Timer functions

The laboratory hot plate contains a timer with the two following functions:

- A timer function for the heating zone. The heating zone automatically switches off after the time set in the timer has expired. The timer for the heating zone can only be set, when a heat intensity between 1 and 9 has been selected.

- A freely available “egg timer” function. The “egg timer” function operates regardless of whether the heating zone has been switched on or off. An acoustic signal is emitted when the set time has expired. The time set on the “egg timer” does not in any way affect the operation of the heating zone.

In order to use either of the timer functions, the timer must first be activated by pressing the <Timer sensor key>.

When the <Timer sensor key> is pressed for the first time, the “egg timer function” can be programmed.

To indicate this mode, only the LED lamp to the right of the timer time display will light up.

When the <Timer sensor key> is pressed for a second time and the heating zone is switched on, the heating zone can be programmed. This mode is indicated by an additional flashing LED lamp to the left of the timer time display.

The egg timer:

The "egg timer" can be activated by pressing the <Timer sensor key> once, irrespective of whether the heating zone has been switched on or not.

The "00" display and the illuminated LED lamp to the right of the LED timer display indicate that the egg timer is ready for setting.

Setting and starting the egg timer:

The operating time of the egg timer can be set by pressing the  /  sensor keys.

After the sensor keys have been pressed for the last time, the egg timer begins to count down the time.

The egg timer is **not** stopped by switching off the laboratory hot plate. It continues to count down the time.

The remaining time is shown in the timer display, until it has expired, when an acoustic signal is emitted.

To switch off the acoustic signal, just press any sensor key.

The heating zone timer:

The heating zone timer can only be set when the heating zone is switched on (any heat intensity between 1 and 9 must be set). The heating zone timer is activated by pressing the timer sensor key twice. The "00" appearing in the LED timer display, the illuminated LED lamp to the right of the LED timer display and the flashing LED lamp to the left of the LED timer display indicate that the heating zone timer is ready for setting.

Setting and starting the heating zone timer:

The timer setting for the heating zone can now be set by pressing the  /  sensor keys.

The LED lamp to the left of the LED timer display indicates that the heating zone timer is ready for setting.

After the sensor keys have been pressed for the last time, the left-hand LED lamp will start to flash and the timer display will show the remaining operating time of the heating zone timer, until it has expired, when an acoustic signal is emitted. To switch off the acoustic signal, just press any sensor key.

To display and delete timer operating times

The remaining operating time of the heating zone timer is always displayed when the heating zone is switched on. If no heating zone is switched on, any remaining time set in the egg timer is displayed.

If the timer sensor key is touched, the remaining operating times for the heating zone timer or for the egg timer can be displayed alternately.

The respective modes are displayed by various combinations of the two LED lamps to the left and to the right of the timer display:

When the LED lamps on the left and right light up continuously, the remaining operating time for the egg timer is shown in the timer display.

When the LED lamp on the left flashes and the LED lamp on the right lights up continuously, the remaining operating time for the heating zone timer is shown in the timer display.

To immediately delete the heating zone timer, set the heat intensity at 0 or switch off the control system with the On / Off sensor.

The "egg timer" is not deleted by switching off, it will continue to run in the background until the time has expired or it is manually deleted.

To delete settings in the heating zone timer or egg timer in operating mode, the heating zone timer or egg timer must first be selected with the <Timer key>.

Then the respective timer can be set to 00 by pressing the  /  sensor keys.

5 Care, cleaning and maintenance

Care

If used for its intended purpose, the laboratory hot plate does not need any special care. To facilitate cleaning, it is advisable to:

- avoid liquids boiling over!
- avoid residues being burnt hard!

Cleaning

The glass-ceramic material is largely chemically resistant. The heating surface remains level at all times, is free of pores and therefore easy to clean. If the heating surface is slightly soiled, the best way to clean it is with hot water and a few drops of a mild household detergent when the heating surface is lukewarm or cold.

To remove encrusted residues, lime and water marks, as well as stains with a metallic shimmer you should use a household cleaning agent suitable for hobs with glass-ceramic heating surfaces. Use a wet cloth or sponge to thoroughly remove any remaining cleaning agent, as some of the commercially available cleaners can be corrosive at higher temperatures. If residues have already become encrusted, remove them manually with a razor-blade scraper. Any **plastic objects** and **aluminium foil** as well as **substances containing sugar** coming into contact with the hot surface **must be removed immediately** with the scraper.

Never use abrasive cleaning agents! In most cases, it is sufficient to clean the hot plate with a damp cloth and some mild household detergent after each use. Then rub dry.

Any encrusted and boiled-over food should best be soaked first with a wet cloth. Then remove the residues with a glass scraper.

Note:

Sugar and melted plastic should immediately be removed from the hot heating surface.

Maintenance

The laboratory hot plate does not require any maintenance if used for its intended purpose.

- Before each use, first check the mains cable to ensure that it is in perfect condition. Do not use the laboratory hot plate with a damaged cable!

CAUTION! Danger of electric shocks!

Do not use laboratory hot plate with damaged cables!

6 Malfunction, possible defects

Malfunction / defects; causes; remedies

Malfunction / Defect	Cause	Remedy
Laboratory hot plate does not switch on; heating zone does not heat up.	Sensor field not clean	Clean sensor field with a damp cloth, rub dry
	No mains voltage	Check mains cable including plugs; check fuse of the circuit from the mains socket, replace if necessary
	Heating zone defective; electronics defective.	Disconnect laboratory hot plate from power supply; disconnect from mains; send laboratory hot plate in for repair
No display	7-segment display or electronics defective.	Caution: Danger of burning! Disconnect laboratory hot plate from mains; send laboratory hot plate in for repair
Heating zone cannot be switched off	Electronics defective.	Disconnect laboratory hot plate from mains; send laboratory hot plate in for repair
"F" appears in the display	Sensor field covered.	Make sure that the sensor field is not covered.
"F1", "F2", "F3" or "F4" appears in the display	Error in calibration of sensors.	Disconnect laboratory hot plate from mains and plug in again. Do not cover the sensor areas or expose them to strong light.
"Ax" appears in the display	Problems in calibration with one sensor key.	Repeat calibration process by disconnecting the laboratory hot plate from the mains and plugging it in again.

7 Accessories

Tripod rod incl. M 8 mounting nut
(stainless steel; length 450 mm; Ø 10 mm):

Z 601

Order No. 28 541 6492

Set of magnetic stirring bars for standard applications
(AlNiCo5; round cross-section; PTFE-coated;
1 pce per set: 15, 20, 30, 40, 50, 60, 70, 80 mm):

Z 603

Order No. 28 541 6554

Magnetic stirring bar for medium quantities
(SmCo; round cross-section; PTFE-coated;
5 pces ea. 9 x 15 mm):

Z 604

Order No. 28 541 6562

Magnetic stirring bar for larger quantities
(SmCo; elliptical cross-section; PTFE-coated;
1 pce ea. 19 x 75 mm):

Z 605

Order No. 28 541 6579

For safety and guarantee reasons, do not use any other accessories!

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We declare under our sole responsibility that the products

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**Laboratory hot
plates
with
glass-ceramic
material**

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to which this declaration refers, is in conformity with the normative document

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Technical Characteristics
Laboratory hot plate with glass-ceramic material

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Supplier's Certificate

We certify that the above equipment has been tested in accordance with DIN EN ISO 9001, Part 8.2.4 "Monitoring and measurement of product" and that the specified quality requirements for the product have been met.

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SI Analytics GmbH