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Important notes: Please carefully read and observe the user manual before the first startup of using the Viscosity Index.

All fees included in this manual are valid data at the time of printing and refer to the respective version number. However, SI Analytics GmbH may provide additions or amendments for technical and commercial reasons to the application without affecting the described properties and functions.

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TABLE OF CONTENT

PAGE

1	General information regarding the Viscosity Index.....	19
1.1	Function.....	20
2	Initial startup	21
2.1	System requirements	21
2.2	Installation	22
2.3	First application start	24
2.4	Requirements for WinVisco	25
3	Working with Viscosity-Index.....	25
3.1	The interface	25
3.2	Functions.....	26
3.3	Preferences.....	29
3.4	Restoring deleted measurement data.....	33
3.5	Support.....	33
4	Problem solving / troubleshooting.....	34

1 General information regarding the Viscosity Index

The Viscosity Index calculates a viscosity index (hereinafter referred to as VI) based on two kinetic viscosities. The basis for calculating the VI is represented by application WinVisco. It provides measurement results (including kinematic viscosities) for individual series of measurements. However, two viscosities are required for the calculation of a VI, for which WinVisco is incapable of the calculation and storage. Instead, the measurement results can only be stored in CSV files. An external program can access these to calculate the index. Viscosity Index finds its use here.

1.1 Function

Viscosity Index is an easy to understand and use application that allows the user to automatically read the saved CSV files from a directory and combine it to a measurement series by including the sample name. The viscosity index is then calculated automatically.

All read data are stored in an access database and can be accessed anytime. Furthermore, there is the option of manually printing the measurement series, as well as to activate the automatic printing mode, where measurement series are printed on the default printer.

An additional option to determine the national language increases the user friendliness and clarity of the application.

2 Initial startup

2.1 System requirements

Viscosity Index can be used on the following systems:

Hardware:

Components	Minimum requirement
Processor	Intel Pentium 4 with 1.5 GHz
RAM	512 MB
Graphics	Resolution 1024x768 Color mode 16 Bit
Drive	CD-ROM for installation
Printer	Installed and configured default printer

Software:

Components	Prerequisite
Operating system	Microsoft Windows • 2000 SP4 • XP SP3 • 7
.NET Framework	Microsoft .NET Framework 2.0 SP1

Special software requirements (only affects Windows 2000)

Components	Version
Microsoft Windows Installer	3.1 (or higher)
Microsoft Internet Explorer	6 (or higher)
Microsoft Updaterollup	1 (for Windows 2000 SP4)
Microsoft Data Access Components (MDAC)	2.8

2.1.1 Storage


Approximately 25 MB free memory are required on the hard drive to install the software, depending on the system (only application data without a database). The respective specific memory requirements increase with the size of the database used.

2.1.2 Printer


Viscosity Index uses various printing modes whereby a trouble-free installation of a printer is a necessary condition.

2.2 Installation

The installation is performed by performing *setup.exe*. The Setup Wizard scans your computer for the required components. If you are missing ([here](#): under Windows 2000, missing .NET Framework 2.0), you can download and install these, if an Internet connection exists. Otherwise, Framework must be installed manually.

When you accept the license terms, the required components are downloaded and installed (see the image below), the actual installation wizard appears.


The installation wizard guides you through the installation of Viscosity Index. Follow the instructions on the screen and confirm each step by clicking "Next".

Please note:

The following message may appear in some Windows versions (specifically Windows 2000) when calling up the installation routine (*setup.exe*).

The reason of the message, despite an installed Service Pack 4 (Windows 2000) is a missing update (*Update Rollup 1 for Windows 2000 SP4*) that Microsoft had subsequently released in September 2005. We recommend installing the update. You will receive the update on the following page: [Microsoft Website](#).


Alternatively, you can start the installation by running the *Setup.msi* (possible approach, but **not recommended**, because any missing and required components cannot be identified and installed later).

2.3 First application start

Before Viscosity Index can enter its first service, a number of adjustments must be made, since no configurations are yet available after a successful installation. The following dialog message with a request for setting the software appears.

After confirming the information, the software settings required for the operation of Viscosity Index must be made:

Viscosity-Index uses its own included database to store the data. The specified path within the setting interface usually corresponds with the correct path to the database. If changes are made in this case, the file must be copied to the specified location!

Viscosity Index starts its operation by accepting the settings. All subsequent settings can now be executed in the Preference interface of Viscosity Index ("*Working with Viscosity Index*, [Preferences](#)"). After you complete the settings, this message is no longer displayed when the program starts.

Please note:

If the required settings are not made, the application cannot start. You can also make the necessary settings within the application while the software is running (see item 4.3, "[Preferences](#)").

Please note:

The directory specification for the CSV files should not refer to the installation path of the application (typically C:\Program Files\SI Analytics GmbH\Viscosity-Index\)! Create a sub-directory, or enter an alternate path.

2.4 Requirements for WinVisco

A prerequisite for the generation and the presence of CSV files is a correctly configured "*application.ini*" of WinVisco. According to the configuration file ("*ApplicationIni V379_8.pdf*"), line 5 of the file must be expanded by the respective CSV directory, in which the files are stored.

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Please note:

The CSV directory (line 5) may **not** correspond with the directory of line 8!!!

Background:

WinVisco automatically creates a complete collection of all measurements and also stores these in a CSV file. The storage location for this is listed in line 8. If line 8 is empty, this is not incorrect. In this case, the complete collection is stored in the installation directory of WinVisco.

However, this complete collection is **not** suitable to analyze the viscosity index. The directory information of **line 5** and **line 8** should therefore **vary**! In addition, the CSV directory of line 5 should preferably **not** be located in the program directory (**not** C:\Program Files\ etc.), since you will rarely will have the required access rights as Windows user.

3 Working with Viscosity-Index

After successfully specifying all the required settings, Viscosity Index is now ready for use. The default language is English and can also be adjusted in the settings (see item 3.3 "[Preferences](#)").

3.1 The interface

The interface consists of three groups:


- Menu bar
- Data section
- Control panel

3.1.1 Menu bar


The menu bar provides you with the following options:

- Issue print jobs
- End program

- Elementary functions
 - o Copy
 - o Cut
 - o Paste
 - o Delete

- Restore print data
- Changing application specific settings

- General information for the Viscosity Index with contact options for potential questions or problems

3.1.2 Data section

The data area contains a table to display the scanned and processed data.

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Measurement results that have been read from the CSV files and written to the database are displayed here. Incomplete measurement series (orange) and incomplete or erroneous information (red) is visually clearly highlighted by the application.

3.1.3 Control panel

The control panel includes functions for processing the measurement series.

You will find options for combining incomplete measurement series with a new sample name, a search function and the option to determine the print mode.

3.2 Functions

3.2.1 Combining measurement series

Measurement series can be combined if they are not complete.

Prerequisites:

- An 80% match of the two respective sample names
- Both measurement times are within the set tolerance (see item 3.3.3 "[Delay for last measurement](#)")

Precisely two rows must be marked to combine measurement series. Mark two measurement series for combining, by...

- 1) marking the first by a mouse click,
- 2) pressing the *CTRL* button on the keypad,
- 3) to then mark the second measurement series (mouse click).

	Test Probe	1	38,8	29.03.2008 08
	Test Probe 2	4	38,8	21.04.2008 13
▶	Test Probe 3			
	Test Probe 4	3	38,8	01.01.2008 17

After two measurement series have been marked, they can be added to a measurement series by confirming the process by clicking on "Merge".

Merge-settings:

Sensitivity: %

New Heading-ID:

* necessary

Merge

Optionally, a new sample name can be specified under which the measurement series are to be combined. If no sample name is provided, one of the other two names will be automatically selected. After combining, the new record is displayed at the end of the table, the previous incomplete measurement series are removed from the listing.

Any errors in the application of this option are displayed by an info message below the data section. Example:


Two rows are to be merged.

Merge-settings:

3.2.2 Search of measurement series

Measurement series can be searched for under the sample name, the viscosity at 40°C, and a viscosity of 100°C. Linking several search criteria (such as the sample name and viscosity) is naturally possible.

The following control element is used to search for values:

Search:

Heading samples:

Viscosity 40°C


Viscosity 100°C

Search

A successful search request is visually confirmed by a green highlighted record.

Messung_A	1	30,0	01.03.2008 01
Messung_B	1	38,8	25.03.2008 09
Messung_C	2	38,8	25.03.2008 09

Unsuccessful search requests are also displayed by an info message below the data section.

3.2.3 Selection of the print mode


Viscosity Index offers two printing modes:

- Automatic printing
- Manual printing

3.2.3.1 Manual printing

The manual print mode is automatically activated at the program start of the Viscosity Index. In order to print measurement series, you must

- 1) mark a measurement series (mouse click on the respective measurement series)
- 2) and confirm the printing process with a click on the "Print" button (alternative: Menu item "File, Print" or "CTRL-P")


Please note:

If the automatic print mode is activated, and you try to manually print a measurement series, you will see a message in which you are informed that the mode will be changed. The automatic printing operation is stopped and disabled.

Additional information can be found under item 3.2.3.2 "[Automatic printing](#)".

3.2.3.2 Automatic printing

In addition to the manual print mode, it is possible to print measurement series automatically. To use this function, it must be enabled, of course.

In this print mode, the data section is sequentially scanned every 10 seconds for data. If complete measurements are available, they will be printed on the default printer. Incomplete measurements are ignored, perhaps you can combine these manually into a record (see item 3.2.1 "[Combining measurement series](#)").

3.2.3.3 Deleting measurement series


Viscosity Index offers the ability to manually delete measurement series as an additional function. To delete a measurement series, select it and press the Delete key on your keyboard (German keyboard layout: *ENTF* button, American keyboard layout: *DEL* button). A dialog window is displayed that is waiting for the confirmation of the deletion.

You can also delete several measurement series at once. As already described in [Combining measurement series](#), select the respective lines via the *CTRL* button.

- ❗ Deleting measurement series triggers the marking of the records within the database as "already printed". These data are no longer displayed in the data section. However, they can be viewed under menu item "[Restoring](#)" of the options menu and restored on request, provided the database has not been cleaned in the meantime.


3.3 Preferences

Viscosity Index offers a series of setting options, which you can reach via the menu item "*Options, Preferences*".


You can therefore activate the following settings:

- Language of the application
- Polling time to browse the CSV directory
- Delay for last measurement
- Options for cleaning the database
- Report parameters (your company address, company logo, ...)
- etc.

To obtain more information about the respective functions, move the mouse pointer over one of the setting options, it displays a tool tip with further details.

3.3.1 Language

The national language of the application can be defined under the selection dialog "*Language*". The following languages are available:

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3.3.2 Polling time

The polling time describes the time after how many seconds the CSV directory will be searched periodicaly for new CSV files.

In the default application settings, a time of 30 seconds is specified. Maximum value of this option is a value of 9999 seconds (~ 166.65 minutes, ~ 2.77 hours). Minimum value of this option is a value of 10 seconds.

ⓘ Please note:

The less time is selected, the greater the loading and utilization of the system! Other parallel processes could be affected by the load during the operation and response. In computationally powerful systems, time can be reduced.

3.3.3 Delay for last measurement

The interval describes the interval between two records in hours. This is relevant when merging two measurement series.

Example 1:

Set interval: 5 hours

Date of the 1st measurement series: 01.01.2008 10:00:00

Date of the 2nd measurement series: 01.01.2008 17:30:00

If you now want to combine this measurement series, an appropriate information message will be displayed below the data section:

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Example 2:

Set interval: 5 hours
Date of the 1st measurement series: 01.01.2008 10:00:00
Date of the 2nd measurement series: 01.01.2008 12:00:00 PM

If you would want to combine these measurement series, no info message would appear. Both measurement series would be combined successfully. Additional information can be found under item 3.2.1 "[Combining measurement series](#)".

3.3.4 Clear database

Viscosity Index allows you to clean the database of already printed records. Especially when a large number of records (depending on the number of the read CSV files and the use period of the software), the cleaning of the database can increase the speed of the system. SI Analytics recommends to clean the database regularly.

3.3.5 Database path


Describes the directory where the database is located. Viscosity Index is delivered with a pre-created database. This is located in the program application directory.

3.3.6 Directory specification of CSV files

Specifying a directory where the CSV files of the WinVisco can be found is of the highest importance for the application. Only on the basis of this read in data can values be written to the database and added to the data section within the application. If no, or a non-existing directory is specified, Viscosity Index creates this directory automatically in response. This path must correspond with that of line 5 of the "application.ini"!


You can also determine which files should be handled after they are read. Viscosity Index shifts the files to sub-directory "*edited*" in the specified CSV directory by default (deleting disabled). If necessary, these files can be deleted when exceeding the allotted memory consumption, or immediately after the reading.

-  However, the immediate deletion of the files is not recommended. If, despite a built-in fault management, Viscosity Index should experience complications, these files can no longer be recovered.


3.3.7 Report details

You have the option to expand the report with your own company information and a logo. Additionally, the alignment of the logo can be determined. Our manufacturer data are specified in the delivery state.


3.3.8 CSV export

Data used by Viscosity Index can be exported on demand into a separate CSV file. The content of this file is the heading-ID, the viscosity at 40°C and 100°C and the viscosity index. The export data can be exported to a single file or in multiple files (depending on the respective date - one file per day), as needed. You must provide the export directory accordingly!

3.4 Restoring deleted measurement data

If measurement series have been accidentally deleted, you can restore them manually. To do this, open the recovery form under the menu item "Options".

In order to restore a measurement record, or to explicitly mark this as printed, set the respective hook inside the "print" column.

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By specifying the respective heading-ID, a measurement record can be searched easily. The application visualizes the search results accordingly:

Proben-ID	Gerät 1	Temperatur 1	Datum 1	Viskosität 1	Gerät 2	Temperatur 2	Datum 2	Viskosität 2	Gedruckt	ID
Test Probe 5	2	38,8	21.04.2008 1...	4,91		0		0	☐	11
Test Probe 4	3	38,8	01.01.2008 1...	55,02				0	☐	7
Test Probe 3				0	2	98,2	21.04.2008 14...	13,91	☐	6
Test Probe 2	4	38,8	21.04.2008 1...	5,27				0	☐	5
Test Probe	1	38,8	27.03.2008 0...	20,2	3	98,2	27.03.2008 09...	73,3	☐	4

3.5 Support

Should any questions or problems arise during the application of Viscosity Index, you can always use our Support. The necessary contact information can be found under menu item "Help, About".

4 Problem solving / troubleshooting

Even during the development of Viscosity Index, we observed the highest quality and integrity of our product. At the time of the delivery of the product, there were no known bugs.

A list of possible dialogue references and the resulting potential causes can be seen below:

Message:

Potential cause:

Error during calculation of viscosity index.

The number format of a kinetic viscosity does not correspond with any valid format. Check the viscosity value within the CSV file.

Measurements could not be deleted from the database. Repeat the process at a later date.

The database was not found under the set path. Check the information for accuracy.

Cannot read database.

The database was not found under the set path. Check settings.

Cannot write records into the database.

The database was not found under the set path. Check settings.

Cannot read path for database.

- 1) Text file with directory information for the database does not exist.
- 2) Path information for the database is incorrect.

Error while searching for new CSV files.

- 1) Directory information for CSV files does not exist. Check settings.
- 2) Setting file was not found or could not be read. Check the existence of the file (settings.ini)

Some fields have not been filled in.

One or more input fields in the settings interface were not completed. Please enter all required data.

The source file for viscosity calculation has been changed. Incorrect results are not excluded!

The source file for the viscosity calculation is subject to the DIN ISO standard 2909 and is checked for changes by a security algorithm, in order to avoid miscalculations. Changing the file contents will be detected and clarified by this error message. Obtain the original file again from the manufacturer.

The CSV file has an unknown column for the database.

Column name:

File name:

Updates to a newer or older version of WinVisco, could result in to the algorithm of the Viscosity Index unknown column names within the exported CSV file of WinVisco. In this case, a message is displayed (message will only been shown from Viscosity Index Version 1.14.6.13 ongoing).

The value of this column cannot be added to the database and will be ignored.

The values of the unknown column will be ignored and not included into the database. In general, however, these are columns that are irrelevant to the calculation of the viscosity index. An update of the Viscosity Index is still recommended.

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Supplier's Certificate

We certify that the above equipment has been tested in accordance with DIN EN ISO 9001, Part 8.2.4 "Monitoring and measurement of product" and that the specified quality requirements for the product have been met.

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SI Analytics is a trademark of Xylem Inc. or one of its subsidiaries.
© 2014 Xylem, Inc. Version 140616 M EDV xxxxxxxx