GEBRAUCHSANLEITUNG
ORIGINALVERSION
OPERATING MANUAL
MODE D'EMPLOI
MANUAL DE INSTRUCCIONES

ViscoPump II

VISKOSITÄTSMESSGERÄT
VISCOSITY MEASURING UNIT
L'APPAREIL DE MESURE DE LA VISCOSITE
UNIDAD DE MEDICION DE LA VISCOSIDAD

SI Analytics
a **xylem** brand

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Operating Instructions 21 ... 38

Important notes:

The operating instructions is an important part of the device. Please read the operating instructions before the first use and keep them carefully. For safety reasons the product may only be used for the purposes described in these present operating instructions. Please also observe the operating instructions for the units/products to be connected.

All specifications in this instruction manual are guidance values which are valid at the time of printing. However, for technical or commercial reasons or in the necessity to comply with the statutory stipulations of various countries, SI Analytics GmbH may perform additions to the product without changing the described properties.

Guarantee:

We provide guarantee for the device described for two years from the date of purchase. This guarantee covers manufacturing faults being discovered within the mentioned period of two years. Claim under guarantee covers only the restoration of functionality, not any further claim for damages or financial loss. Improper handling/use or illegitimate opening of the device results in loss of the guarantee rights. The guarantee does not cover wear parts.

To ascertain the guarantee liability, please return the instrument and proof of purchase together with the date of purchase freight paid or prepaid

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Table of Content

1	Basics.....	23
1.1	Approved/Intended Use.....	23
1.2	Safety guidelines	23
1.3	Technical data AVS® 370	24
1.4	Functioning of the device.....	26
1.5	Capillary viscometry	26
1.6	Measurements principles	27
2	Initial operation	28
2.1	Unpacking.....	28
2.2	Installation and dismantling of the module	28
2.3	Connecting the devices	29
2.4	Connecting the viscometers and other devices.....	32
2.5	Trouble shooting.....	34
3	Data transfer	35
3.1	RS-232-C Interfaces.....	35
3.2	Interface configuration.....	35
4	Maintenance and care of the device and the viscometers	36
4.1	Maintenance work to be carried out.....	36
4.2	Maintenance and care of the VZ 7215 absorbent bottle	36
4.3	Periods without operation	37
4.4	Reproducibility of results.....	37
4.5	Viscometers within quality assurance systems	37
4.6	Storage and transportation	38
4.7	Recycling und disposal	38

Copyright

© 2014, SI Analytics GmbH

Reprinting - even as excerpts - is only allowed with the explicit written authorization of SI Analytics GmbH, Mainz.

Printed in Germany.

1 Basics

1.1 Approved/Intended Use

The ViscoPump II module is a functional unit used to determine absolute and relative viscosity of liquid media. The module is controlled exclusively through the use of a special program to be run either on a computer or on a suitable control unit (such as the AVS®470 Viscosity Measuring Unit). The module forms an integral part a viscosity measuring device and is designed for being built into the respective device or controller unit.

1.2 Safety guidelines

The Development and production of the ViscoPump II module is done within a system which meets the requirements laid down in the DIN EN ISO 9001 standard. It was manufactured and tested according to DIN EN 61 010, Part 1, Protective Measures for electronic measurement devices and control devices and has left the factory in an impeccable condition as concerns safety technology. In order to maintain this condition and to ensure safe operation, the user should observe the notes and warning information which marked by "⚠" contained in the present operating instructions.

The ViscoPump II modules control the entire measurement process!

(among other things the temperature pre-adaptation process of the samples in the viscometers, the process of pumping the liquid up into the storage containers of the viscometers, the measurement of the flow times etc.)

If the module is incorporated into an AVS®370 or AVS®470 Viscosity Measuring Unit, please verify after start-up that the proper operating mode is set, i.e. "sucking" or "pressing" action; these operating modes can be seen on the corresponding LED on the front panel of the AVS®470, or in the context menu of the WinVisco software of the AVS®370 Viscosity Measuring Unit, respectively. Please refer to the operating instructions of the device concerned to make the proper setting for the application, and make sure to use the suitable hose set.

For reasons of safety, the ViscoPump II modul must be opened by authorised persons only!

This means, for instance, that work on electrical equipment must only be performed by qualified specialists. In the case of unauthorised intervention in the Visco Pump II modul, as well as in the case of negligently or deliberately caused damage, the warranty will become void.

The installed safety units must never be disabled!

If it has to be assumed that safe operation is impossible, the ViscoPump II modul has to be put out of the device and secured against inadvertent putting to operation!

(In this case please switch the basic unit off, pull plug of the mains cable out of the mains socket, remove the modul and return it to the SI Analytics service department).

Examples for the assumption that a safe operation is no longer possible,

- the package is damaged,
- the unit shows visible damages,
- the unit does not function properly
- liquid has penetrated into the casing.

The ViscoPump II modul must not be operated or stored in moist rooms and - for safety reasons - must be used solely for the applications described in this operating manual!

The relevant regulations regarding the handling of the substances used have to be observed (the Decree on Hazardous Matters, the Chemicals Act, and the rules and information of the chemicals trade)! Please wear protective glasses ever!

It has to be ensured on the side of the user that the persons entrusted with the use of the unit are experts in the handling of substances used in the environment and in the unit or that they are supervised by specialised persons, respectively.

Please read the separate operating instructions of the transparent thermostats as well.

1.3 Technical data AVS® 370

Translation of the legally binding German version

(Release 8. July 2014)

Note

The data given below is only valid if the ViscoPump II module is built into a certified device (e.g. AVS®370, AVS®470, AVS®Pro)

CE-sign:	CE	EMC compatibility according to the Council Directive: 2004/108/EG; applied harmonized standards: EN 61326-1:2006 Low-voltage directive according to the Council Directive 2006/95/EG applied harmonized standards: EN 61 010, Part 1
Country of origin:		Germany
Measurement parameters:		Flow time in seconds [s]
Capture of measurement value:		Flow time: Optoelectronic or thermo-resistive capture of the meniscus passage through the measurement planes of the viscometers
Measurement ranges:		
Viscosity:		"pressing" action 0.35 ... 1.800 mm ² /s (cSt) measuring temperature of approx. 20 ... 25 °C "sucking" action 0.35 ... 5.800 mm ² /s (cSt) measuring temperature of approx. 20 ... 25 °C
Time:		0.01 to 9999.99s, Resolution 0.01 s
Pumping pressure:		Fully automatically controlled "sucking" action to approx. -160 mbar Fully automatically controlled "pressing" action to approx. +160 mbar
Measurement ranges:		Precision (reproducibility and comparability) in accordance with DIN 51562, Part 1
Time measurement:		± 0.01 s ± 1 digit, but not more accurate than 0.01 % measurement uncertainty in the determination of absolute, kinematic viscosity furthermore depends on the uncertainty of the numerical value of the viscometer constants and the measurement conditions, especially as concerns the measurement temperature.

Data transfer parameters:

Data interface:	bi-direktional serielle interface according to EIA RS-232-C
Data format:	word length 7 bits, 2 stop bits, 4800 baud, no parity

Connections

a) Back panel	96-channel VG bar for contact making with the module plugged in a controller unit, including voltage supply and 2 serial interfaces
Gas exhaust pipe:	Stainless-steel tube used to expel or suck the gases delivered by the rotary pump

Front panel (to be connected to the front panel of ViscoPump II modules)

Pneumatic connectors:	Pressure/suction venting, to be connected to viscometer.
Overflow guard for VZ 8552 suction line	

Capacitive sensor: Round-plug connector, 4-channel, Screwed cap according to DIN 45321

Type VZ 8511 ViscoPump II for optoelectronic sampling

Measurement stand AVS®/S: Round-plug connector with bayonet catch, DIN 5 channels for AVS®/S (measurement stand), 5-channel DIN socket

Type VZ 8512 ViscoPump II for TC sampling

TC viscometer: 4-channel DIN socket

Operating indicator: Green LED (voltage indicator)

Various add-on devices: EIA-RS-232-C interface with 9-channel Submin D plug

Casing: No casing, module in plug-in design

Dimensions: approx. 50 x 173 x 140 mm (WxHxD) overall dimensions
Weight: approx. 0.8 kg
Front panel: 2 mm anodised aluminium

Ambient conditions: Ambient temperature: +10 ... +40°C for operation and storage

Air humidity according to EN 61 010, Part 1:
max. rel. hum. 80% for temperatures up to 31°C, I
nearly decreasing down to 50% rel. hum at a temperature of 40°C

1.4 Functioning of the device

The ViscoPump II Modul is used to perform flow-time measurements in capillary viscometers.

The available capillary viscometers enable viscosity measurements 0.35 of approx 5,000 mm²/s (cSt) to be carried out at a measuring temperature of approx. 20 ... 25 °C.

The ViscoPump II Modul can be fitted out in two ways to be used for meniscus sensing. Connecting TC viscometers to the ViscoPump II VZ 8512 module will also enable the measurement of black and opaque liquids. As an alternative option it is possible to use viscometers in combination with the light-optical ViscoPump II VZ 8511 module viscometers for meniscus sensing with the measurement stand, e.g. the AVS®/S.

Time recording extends up to 9999.99 seconds with a resolution of 0.01 s. Serial data transmission is used to transfer the time readings to the respective controller unit.

Prior to the measurement as such, the liquid to be measured is sucked upwards inside the capillary viscometer through two measurement planes (N2 and N1) which are designed as light barriers or thermistor sensors, depending on the viscometer type (fig. 1 and 2).

When using Ubbelohde viscometer, the design of the program ensures that the suspended spherical level will form prior to the start of the measurement.

1.5 Capillary viscometry

Capillary viscometry is the most accurate method for the determination of the viscosity of liquids with a Newtonian flowing behaviour. The measurement as such consists in a time measurement. The time measured is that which a specific quantity of liquid requires to pass through a capillary having a defined width and length. Conventionally, this process is watched with the human eye, and the flow time is measured manually using a stop watch.

In the case of the AVS® 370 Viscosity Measuring Unit, as with all viscometry measuring devices from SI Analytics GmbH, the liquid meniscus is captured on the measurement planes, either in an optoelectronic manner by means of light barriers, or else on a thermo-resistivity basis by thermistors.

1.6 Measurements principles

a) Optoelectronic sensing of the liquid meniscus

The near-infrared light which is generated in LEDs located in the upper section of the measurement stand is conducted through a glass-fibre light-conductor cable onto the measurement planes. The light shines through the viscometer before it arrives at another light-conductor cable located on the opposite side; inside this second cable, the light is conducted to a receiver in the upper section of the measurement stand.

While the liquid meniscus passes through the measurement planes, the lens-like effect of the meniscus causes a short-term darkening of the light beam, followed by a magnification. This process generates a measurement signal which can be evaluated accurately.


Fig. 1 Viscometer for optoelectronic measurements

b) Viscometer with thermistor sensors (TC viscometer)

In the case of TC viscometers, glass-coated thermistors serving as sensors are molten in on the level of the measurement planes. While the meniscus passes through the measurement planes, the differences in the thermal conductivity properties of air and liquid lead to a change in the heat balance. The thermistors of the TC viscometers are molten hermetically tightly into the glass coating of the viscometer, so that the viscometers located inside are chemically resistant to all kinds of substances other than strong leaches, fluoric solutions, or concentrated hot phosphate solutions.


Fig. 2 Viscometer for optoelectronic measurements

2 Initial operation

2.1 Unpacking

The ViscoPump II module is delivered in a folding packaging with a plastic-foam protection. Please take the module out of this packaging and make sure that all items specified on the packing list are present.

2.2 Installation and dismantling of the module

2.2.1 Installation

Warning

Please switch off the basic / controller unit first! Prior to installing or replacing this functional unit, please make sure that the mains plug is pulled off of the mains socket!

The ViscoPump II module is installed by inserting it into the slot provided. In this process, please see that both the printed board and the metal shielding plate are inserted correctly into the corresponding guidance rails, and make sure that the socket bar locks in place. Subsequently, complete the installation by tightening the front-panel screws.

Notice

The ViscoPump II modules control the entire measurement process, among other things the temperature pre-adaptation process of the samples in the viscometers, the process of pumping the liquid up into the storage containers of the viscometers, the measurement of the flow times etc!

If the module is incorporated into an AVS®370 or AVS®470 Viscosity Measuring Unit, please verify after start-up that the proper operating mode is set, i.e. "sucking" or "pressing" action; these operating modes can be seen on the corresponding LED on the front panel of the AVS®470, or in the context menu of the WinVisco software of the AVS®370 Viscosity Measuring Unit, respectively. Please refer to the operating instructions of the device concerned to make the proper setting for the application, and make sure to use the suitable hose set.

2.2.2 Dismantling

To dismantle or replace the ViscoPump II module, please proceed as follows:

- Remove the pneumatic and electrical connections from the front panel of the ViscoPump II module to be replaced.
- Loosen the screws located at the corners of the front panel.
- Use the upper and lower insertion handles to leverage the ViscoPump II module out of its rear-side plug connection.
- Pull the ViscoPump II module out of the AVS 470 Viscosity Measuring Unit.
- After inserting the new ViscoPump II module, please secure it again with the screws of the front panel. Re-establish the electrical and pneumatic connections.

2.3 Connecting the devices

2.3.1 Suitable viscometer types, racks, and measurement stands

Viscometer Type	Rack Type no.	Measurement Type no.
Ubbelohde (DIN)		
532...	053 92	AVS®/S AVS®/SK
530...		
501...		
541...		
Ubbelohde (ASTM)		
525...	053 92	AVS®/S AVS®/SK
526...		
527...		
545...		
Micro-Ubbelohde		
536...	053 92	AVS®/SAVS®/SK
537...		
538...		
Ubbelohde dilute-solution viscometers		
531...	----	AVS®/SK-V
Cannon-Fenske-routine		
513...	----	AVS®/SK-CF
520...		
Micro-Ostwald		
517...	053 97	AVS®/S AVS®/SK
Ubbelohde (TC)		
567...	053 93	----
568...		
569		
562...		
563...		
564...		
Micro-Ubbelohde (TC)		
572...	053 93	----
573...		
574...		

2.3.2 Deployable hose combinations

Hose combinations Type no.	Description	Application
VZ 5505	Silicone hose kit, oppressive, for Ubbelohde viscometers (3 legs), and Cannon-Fenske and Ostwald viscometers	Standard, but please observe: the sample can leak from the capillary tube during a malfunction
VZ 5505 + VZ 8526	Silicone hose kit, suctioning, for Ubbelohde viscometers (3 legs)	Standard, safer than oppressive operation, as the sample cannot leak from the capillary tube. Unsuitable for volatile samples.
VZ 8523	PTFE hose set, suction, for Ubbelohde viscometer (3 legs)	For aggressive samples that attack silicone, such as sulphuric acid. Hose lengths are designed to fit the AVS®370 on the support console VZ 8571. All hose sets for suctioning operation can be combined with the suction set and the sample fill set.
VZ 8524	PTFE hose kit, suctioning, with soda lime filter VZ 7215 for Ubbelohde viscometers (3 legs)	For aggressive samples whose vapors are absorbed by soda lime filters to protect the ViscoPump. Depending on the sample, an active carbon filter VZ 7216 can be used instead of the soda lime filter VZ 7215. Hose lengths are designed to fit the AVS®370 on the support console VZ 8571.
VZ 8530	PTFE hose kit, suctioning, with soda lime filter VZ 7215 for Ubbelohde viscometers (4 legs)	For Ubbelohde viscometers with additional 4th pipe for filling and cleaning. Hose lengths are designed to fit the AVS®370 on the support console VZ 8571.
VZ 5606	For TC viscometers (3 and 4 legs) with screw connections: Silicone hose fittings with connecting cable. For oppressive operation	For TC viscometers (3 and 4 legs). Typical applications are measurements of oils.

Note


The hose combinations are to be selected on the basis of the required application, e.g. VZ 8521, VZ 8523, and VZ 8524.


Note

When using Micro Ubbelohde viscometers with TV sensors, the ignition temperature of the media to be measured has to be taken into account:


- It has to be higher than 250 °C.


DIN-Ubbelohde-viscometer

Micro-Ubbelohde-viscometer

Canon-Fenske-Routine-viscometer

Micro-Ostwald-viscometer



Ubbelohde-viscometer
with TC sensorsMicro-Ubbelohde-viscometer
with TC sensors

Fig. 3 Suitable viscometer types

2.4 Connecting the viscometers and other devices

The ViscoPump II modul allows the use of most various viscometer types: DIN, ASTM, Ubbelohde and Micro Ubbelohde viscometers as well as Cannon-Fenske Routine, Micro TC and Micro Ostwald viscometers.

Owing to careful manufacture and quality-assurance procedures, all viscometers from SI Analytics meet the highest accuracy standards.

The K viscometer constant is determined individually by way of a calibration of each glass capillary viscometer. Owing to the use of high-quality measurement and testing equipment and the application of national standard gauges, SI Analytics guarantees an absolutely precisely reproducible calibration. For Ubbelohde viscometers having the same constant, the same correction seconds (Hagenbach correction) are valid. Gauging by the user is not necessary, since the corrections correspond to the theoretical values as taken from the operating instructions for the viscometers. This statement is true for both Ubbelohde viscometers of normal size as well as for micro viscometers.

Note

It is also possible to connect or control other devices (such as absorption traps, overflow guards etc). Depending on the intended use of the ViscoPump II modul, it may be highly recommendable to connect these devices, please refer to the items below.

2.4.1 TC viscometers with thermistor sensors

Fill the viscometer (approx. 18 - 20 ml), then place it in the thermostat bath.

Connect the ViscoPump II modul VZ 8512 and the TC viscometer using the hose/cable combination which comes with the device. To do so, place the device in the holders, then attach the quadruple plug of the cable to the viscometer and the ViscoPump II module (first plug, then screw); subsequently, make the screwed connections in accordance with the numbers indicated on the hoses and the rack. In the case of "pressing" operation, the capillary tube remains open, for "sucking" operation the filling tube is to be left open. Please observe the colour codes (red = suction, black = venting) when connecting the pneumatic screwed connections to the ViscoPump II module for TC viscometers.

2.4.2 Viscometers using light-barrier sensing

Use the hose/cable combination to make an electrical and pneumatic connection between the ViscoPump II modul VZ 8511 and the measurement stand. The plugs are firmly connected to the sockets by rotating the union sleeve. Please observe the colour codes (red = sucking, black = venting) when screwing the threaded pneumatic connections into the ViscoPump II module. Please insert the selected capillary viscometer into the fixating rack as is shown in fig 4., then fill it. Insert the fixating rack together with the viscometer into the measurement stand (with the cut-out at the bottom sheet pointing forwards). The cut-out will latch into the lug provided. Pressing the viscometer slightly against the fixating rack will latch it into the holding spring located on the measurement stand.


Fig. 4 Inserting or replacing a viscometer with light-barrier sensing

2.4.3 Connection of VZ 7215 absorption traps

In the “sucking” mode (vacuum), volatile components can enter the ViscoPump II module. This is particularly problematic for corrosive solvents such as formic acid or dichloroacetic acid.

 Notice

For these cases, a hose fitting “suctioning” VZ 8524 must be used, which includes the absorption traps VZ 7215 and appropriate connecting hoses.

The absorption traps which prevent contaminations from penetrating into the pneumatic system of the ViscoPump have to be inspected at regular intervals. If sodium lime is used as an absorption agent with acidic solvents, the colour condition of the indicator is to be checked on a daily basis. As soon as this condition has shifted to BLUE in the half of the absorber material, this is the very last moment to replace the material for safety reasons.

 Notice

If such a colour shift cannot be observed over an extended period of time, this may be attributable to the fact that an acidic over-saturation of the material has caused a de-colouration; this may then appear as “normal”, but it will definitely result in the destruction of the pneumatic system after some time.

This situation is explicitly excluded from the warranty coverage!

For non-corrosive solvents and oils, which contain volatile constituents, absorption traps with activated carbon filling are available. When using activated carbon as an absorption agent (e.g. with solvents or used mineral oils), a replacement should be made at intervals between 1 and 2 weeks; this depends on the load factor which, in turn, is a function of the volatility of the materials.

2.4.4 Connection of the VZ 8552 overflow guard

We urgently recommend the connection of the VZ 8552 overflow guard (available as an option) for the suction-mode operation of the ViscoPump II module. The connection of the VZ 8552 overflow guard (capacitive sensor for the safety bottle) excludes over-pumping in suction mode (contamination of the ViscoPump II module). The holder on the safety bottle accommodates the capacitive sensor.

When using the ViscoPump II module VZ 8511 (meniscus sensing by light barriers) the holder for the safety bottle is to be attached to the measurement stand, e.g. the AVS®/S.

When using the ViscoPump II module VZ 8512 (thermo-resistive measurement) the holder for the safety bottle is to be attached to the “viscometer gallows” provided for the TC viscometer 5732.

Should any liquid be over-pumped into the safety bottle, the safety sensor will trigger a stop. After emptying the safety bottle, the lateral LED on the capacitive sensor will go out. You may continue with the measurements. The electrical connection of the VZ 8552 overflow guard is made using DIN plugs on the front side of the respective module of the ViscoPump II.

 Note

The sensitivity of the capacitive sensor has to be adapted to the medium being used.

To do so, please use the enclosed screw driver to adjust the lateral set screw in such a manner that the capacitive sensor in the built-in condition (i.e. without medium) are just close from responding (i.e. the LED is on).

 Warning

Make sure that the Basic Measuring Unit is always the first device you switch off!

Prior to replacing any functional unit, please be sure to the mains plug MUST be pulled out of the mains socket.

Caution: Liquid dripping off may be hazardous to the user!

2.5 Trouble shooting

Check whether the Basic Measuring Device is switched on:

Trouble

Air bubbles in viscometer

Trouble shooting

- Is the filling quantity sufficient? Check, fill viscometer if required.
- Is the viscometer properly connected?
 - in the case of pressing operation, please check whether the filling tube is connected; if necessary connect properly.
 - for operation in suction mode, please check whether the capillary tube is connected; if necessary, connect properly.
 - please check whether the venting port is tightly connected; if necessary re-tighten its screwed connection.

Excessive pumping of measuring medium into the thermostat bath or into the safety bottle:

- Is the viscometer properly connected?
 - for pressing operation?
 - for operation in suction mode?

When using AVS measuring tripods:

- check position of the rack within the stand
- check the electrical connection from the viscometer to the ViscoPump type II module
- Is the green LED on the measuring stand illuminated?

When using TC viscometers:

- Is the viscometer properly connected?

3 Data transfer

3.1 RS-232-C Interfaces

The ViscoPump II modul has three RS-232 C interfaces.

The interfaces on the internal data bus (96-channel VG bar) are used for communication with the respective controller unit or computer. The interface on the front panel serves to connect optional devices such as thermostats etc. The connected are treated as sub-devices of the ViscoPump II module.

3.2 Interface configuration

The interface parameters cannot be changed. All transmission parameters are firmly set to the following values:

Parity:	None
Stopbits:	1
Datenbits:	7

Changes to the transmission parameters of the RS-232-C interface are only possible on a software basis. Setting is done through the program being used or through the software of the controlling device. Please refer to the operating instructions of the respective device or program being used for a description on how to make the required settings. The parameter settings of the connected peripherals must be identical, otherwise communication via this interface will be impossible.

All other devices to be connected must have the same parameter settings.

The other possible settings include:

RS parameter

baud:	bit:	stopp:	parity:
2400	7	2	No
4800	8	1	No
9600	7	1	Odd
	8	1	Odd
	7	1	Even
	8	1	Even

4 Maintenance and care of the device and the viscometers

Maintaining the proper functioning requires certain inspection and maintenance work.

Maintenance and service work includes:

- Visual check
- Interface function, ViscoPump II, and rinsing burettes
- Once per quarter, the electrical contacts have to be inspected for corrosion, if the viscosity measuring unit is used in premises with an occasional occurrence of corrosive matters in their atmosphere.

Maintenance intervals

Normal operation	intervals for carrying out all work is 6 months
Under particular strain	the max. intervals for carrying out all maintenance work are 4 weeks
In case of disturbances	If any disturbance, malfunction, or other defect becomes obvious, the work has to be carried out immediately.

4.1 Maintenance work to be carried out

- Check the hoses, screwed connections for signs of visible damage, contamination, and leaks.
- Check the electrical plug contacts for corrosion and mechanical damage (on the AVS® 370 Viscosit Measuring Unit and on the cables).
- If necessary, the exterior of the casing of the viscosity measuring unit can be cleaned with a piece of cloth soaked with a household cleaning agent. The lower and rear sections have to be dry-treated.
- In no case must liquid penetrate into the interior of the lower section.
- Defective parts must be repaired or replaced with new ones.
- Defective glass parts must always be replaced.

4.2 Maintenance and care of the VZ 7215 absorbent bottle

The VZ 7215 absorption traps which prevent contaminations from penetrating into the pneumatic system of the ViscoPump have to be inspected at regular intervals.

Notice

The measurements will not work if there is an incorrect connection. There is a risk of the sample leaking from the viscosimeter or being suctioned into the ViscoPump II.

Notice

If sodium lime is used acidic solvents as an absorption agent, the colour condition of the indicator is to be checked on a daily basis. As soon as this condition has shifted to BLUE in the half of the absorber material, this is the very last moment to replace the material for safety reasons.

If such a colour shift cannot be observed over an extended period of time, this may be attributable to the fact that an acidic over-saturation of the material has caused a de-colouration; this may then appear as "normal", but it will definitely result in the destruction of the pneumatic system after some time.

This situation is explicitly excluded from the warranty coverage!

Note

When using activated carbon as an absorption agent (e.g. with solvents or used mineral oils), a replacement should be made at intervals between 1 and 2 weeks; this depends on the load factor which, in turn, is a function of the volatility of the materials.

4.3 Periods without operation

If the capillary viscometers are not used over a long period of time, the liquids contained in the system, in particular aggressive solutions, have to be drained. If the liquid is left in the system, one has to reckon that the solutions used will alter in the course of time and attack the glass, in particular the capillaries.

 Notice

Cleaning agents should be matched to the previous samples or impurities.

In many cases, an aqueous cleaning agent (glass cleaners, detergents) or organic solvents (such as acetone or hydrocarbons) are sufficient.

 Caution

Strong oxidizing cleaning agents such as chromic acid may only be used by trained personnel and must be suitably disposed for safety and environmental reasons - the current guidelines for handling hazardous materials must be observed.

 Note

In the last rinse cycle, the viscometer should be rinsed with a suitable solvent with a low boiling point (such as acetone), and dried by an air flow, which is preferably generated by underpressure (for example, water jet pump). The viscometer is dry and dust-free by this treatment and can thus be used for manual and automatic measurements.

4.4 Reproducibility of results

The measurement or analysis results depend on a variety of factors. Please check the plausibility of the measurement results or analysis results at regular intervals, and carry out the required reliability tests. In this regard, please adhere to the usual validation procedures and especially to the "Viscometers within quality assurance systems" chapter.

4.5 Viscometers within quality assurance systems

Recommendations for companies that have introduced a quality assurance system in accordance with the DIN EN ISO 9001 standards. In this quality assurance system, an inspection of the measuring equipment is planned. The intervals and required accuracy can be defined by each company according to its own requirements. The standard DIN/ISO 10 012, Part 1 serves as a guideline in this matter. We recommend regular inspection of the viscometers in defined intervals.

Inspection of the viscometer constants:

a) Calibration using comparative measurements with reference measuring standards

Comparative measurements must be performed with a viscometer (reference measuring standard) which was tested at the PTB (Federal German Physical-Technical Institute) and provided with a constant. During this comparative measurement, the viscometer to be inspected and the PTB - tested viscometer were placed simultaneously in the same thermostat bath. The test liquid tested, the viscosity of which must not be known exactly, is filled into both viscometers, tempered and the flow-through time then measured. The constants of the viscometers to be inspected are then calculated according to the following equation:

$$K = \frac{K_{\text{PTB}} \cdot t_{\text{PTB}}}{t}$$

K constant of the tested viscometer

K_{PTB} constant of the standard reference viscometer

t low time (HC) of the tested viscometer (corrected by Hagenbach-Couette)

t_{PTB} flow time (HC) of the standard reference viscometer (corrected by Hagenbach-Couette)

Within the quality assurance system in accordance with DIN EN ISO 9001, traceability of the measuring equipment to national measuring standards is demanded. This traceability can be achieved by inspecting the comparative viscometers (reference measuring standards) at regular intervals at the PTB. The time intervals are defined according to the specifications made in the quality assurance system of the user.

b) Calibration of the capillary viscometer with normal oils of the PTB

During this calibration, a normal oil from the PTB with known viscosity is used as a reference measuring standard. The measurement is performed by means of flow-through measurement of the PTB normal oil in the viscometer to be inspected in a temperature bath, the temperature of which must correspond precisely to the test temperature of the PTB. In this case, it is extremely important to make sure that the temperature is absolutely correct. In case of temperature variation, this will always result in a constant for the viscometer that deviates from the constant applied. A temperature variation of 0.01 K, for instance, will result in a measuring error of up to 0.01 %. The calibration of the deviating temperature into the viscometer constant is not permitted.

c) Inspection by SI Analytics with quality certificate in accordance with DIN 55 350-18, 4.2.2

The inspection at SI Analytics GmbH is carried out by means of comparative measurements using viscometers as reference measuring standards that were tested at the PTB (corresponds to Item 1).

Note

General Information on the stability of viscometer constants

Each inspection (even with a certificate) can guarantee the technical measuring direction only for a limited period of time. The constants of viscometers made of the borosilicate glass DURAN®, however, can remain unchanged for long periods of time if the viscometers are kept away from altered influences. Especially extreme changes can be expected, for instance, during the use of liquids that attack glass, in particular hot caustic soda hydrated (NaOH) orduring glass-blowing repairs (even for apparently insignificant repairs).

Liquids whose components adhere to the glass wall also cause errors. In such cases, regular cleaning is required whereby the corrosive action cleaning agent on the glass must be eliminated.

For this reason, we recommend that the user should write up a special processing instructions for all important measurements and include them in his quality assurance manual in accordance with DIN EN ISO 9001. In all cases the user is responsible for the correctness of his measuring and testing equipment and is not released from his responsibility for quality (cp. DIN 55 350, Part 18).

4.6 Storage and transportation

If the AVS® 370 Viscosity Measuring Unit has to be stored over some time, or to be dislocated, using the original packing will be the best protection of the device. However, in many cases this packing will not be available any more, so that one will have to compose an equivalent packaging system. Sealing the device in a foil is recommended.

If viscometers are to be stored temporarily, or transported, all liquids contained in the system, especially aggressive solutions, have to be drained.

If the viscosimeters are to be stored intermediately or transported again, the fluids contained in the system, particularly aggressive solutions, must be removed.

4.7 Recycling und disposal

The present viscosity measuring unit and its packing material were mainly made from materials which can be disposed of and recycled in an environmentally friendly manner.

Should you have any questions regarding disposal, please contact SI Analytics GmbH.

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SI Analytics is a trademark of Xylem Inc. or one of its subsidiaries.
© 2014 Xylem, Inc. Version 141118 **M 602 827 5**