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Originalversion

## OPERATING MANUAL

## MODE D'EMPLOI

## MANUAL DE INSTRUCCIONES

# CK 300/ CK 310

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FLOW-THROUGH COOLER

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REFRIGERADOR DE CIRCULACIÓN

**SI Analytics**

a xylem brand

Gebrauchsanleitung .....	3 .....	14
--------------------------	---------	----

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**Important notes:** Before initial operation of the Flow-Through Cooler CK 300 / CK 310 please carefully read and observe the operating instructions. For safety reasons the Flow-Through Cooler CK 300 / CK 310 may only be used for the purposes described in these present operating instructions.

Please also observe the operating instructions for the units to be connected.

All specifications in this instruction manual are guidance values which are valid at the time of printing. However, for technical or commercial reasons or in the necessity to comply with the statutory stipulations of various countries, SI Analytics GmbH may perform additions to the Flow-Through Cooler CK 300 / CK 310 without changing the described properties.

Mode d'emploi .....	27 .....	38
---------------------	----------	----

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---

## Operating Instructions Flow-Through Cooler CK300 / CK310

### Table of contents

QUALITY ASSURANCE .....	16
Guarantee .....	16
YOUR CONTACTS AT SI ANALYTICS GMBH .....	16
UNPACKING AND CHECKING .....	16
INTENDED USE .....	17
Description .....	17
OPERATOR RESPONSIBILITY – SAFETY NOTES .....	17
Operator responsibility .....	17
Safety notes .....	18
Explanation of safety notes .....	18
Safety recommendations .....	20
1. TECHNICAL DATA OF FLOW-THROUGH COOLER CK300/CK310 .....	21
2. OPERATING CONTROLS AND FUNCTIONAL ELEMENTS .....	22
3. PREPARATIONS .....	23
3.1. Location .....	23
3.2. Bath liquids .....	23
3.3. Electrical connections .....	24
3.4. Connecting the hoses .....	24
4. OPERATING PROCEDURES .....	25
4.1. Power connection .....	25
4.2. Switching on .....	25
5. TROUBLESHOOTING .....	25
6. CLEANING THE UNIT / MAINTENANCE .....	26

Version 141021 US

**Copyright** © 2014, SI Analytics GmbH

Reprinting - even as excerpts - is only allowed with the explicit written authorization of SI Analytics GmbH, Mainz. Printed in Germany.

## Quality assurance

Dear customer,

SI Analytics GmbH implements a **Quality Management System** certified according to EN ISO 9001.

This guarantees the presence of organizational structures which are necessary to ensure that our products are developed, manufactured and managed according to our customers expectations. Internal and external audits are carried out on a regular basis to ensure that our QMS system is fully functional.

Please inform us if, despite our precautionary measures, you should find any product defects. You can thus help us to avoid such faults in future.

## Guarantee

We provide guarantee for the device described for two years from the date of purchase. This guarantee covers manufacturing faults being discovered within the mentioned period of two years. Claim under guarantee covers only the restoration of functionality, not any further claim for damages or financial loss.

Improper handling/use or illegitimate opening of the device results in loss of the guarantee rights. The guarantee does not cover wear parts. To ascertain the guarantee liability, please return the instrument and proof of purchase together with the date of purchase freight paid or prepaid.


## Your contacts at SI Analytics GmbH

Please get in contact with us or the authorized agent who supplied you with the unit if you have any further questions.

SI Analytics GmbH

Hattenbergstrasse 10  
D-55122 Mainz

Telefon +49 61 31 / 66 51 11  
Telefax +49 61 31 / 66 50 01  
E-Mail: [si-analytics@xyleminc.com](mailto:si-analytics@xyleminc.com)  
[www.si-analytics.com](http://www.si-analytics.com)

The following specifications should be given when product enquiries are made:

- **Unit name** printed on the front of the unit and as specified on the name plate.

## Unpacking and checking

Unpack the unit and accessories and check for damages incurred during transit. These should be reported to the responsible carrier, railway, or postal authority, and a request for a damage report should be made. These instructions must be followed fully for us to guarantee our full support of your claim for protecting against loss from concealed damage. The form required for filing such a claim will be provided by the carrier.

## Intended use

The flow cooler CK 300/CK 310 cools liquid media (water, silicone oil or mineral oil, ethanol and polyglycols).

Precise information about the bath media to be cooled and their operating temperature can be found in Section 3.2 ("Tempering liquids") in this operating manual.

The flow-through Cooler CK300/CK 310 are not suitable for direct cooling of foods, semi-luxury foods and tobacco, or pharmaceutical and medical products. Cooling means the unprotected contact of the object with the bath medium (bath fluid).

## Description

The flow-through cooler CK300/CK310 cools down the liquid which is pumped through the built-in heat exchanger. Therefore, an external pump is required like the one built-in the immersion circulators from SI Analytics GmbH.

A typical application for the flow-through cooler CK300/CK310 is the usage in the external liquid circuit of a transparent bath.

## Operator responsibility – Safety notes

The products of **SI Analytics GmbH** ensure safe operation when installed, operated, and maintained according to common safety regulations. This section explains the potential dangers that may arise when operating the circulator and also specifies the most important safety precautions to preclude these dangers as far as possible.

### Operator responsibility

Make sure all persons tasked with operating, installing, and maintaining the unit have read and understand the safety information and operating instructions.

#### **Set up, maintenance, repair**

Only qualified personnel is authorized to perform configuration, installation, maintenance and repairs of the thermostat.

#### **Handling**

Make sure the product is checked for proper condition regularly (depending on the conditions of use).

---

## Safety notes

These notes are intended to enhance your own safety consciousness.

We have set the highest quality standards for ourselves and this unit during development and production. Every unit meets relevant safety regulations.

**The correct unit usage and proper handling is however solely your responsibility.**

## Explanation of safety notes

In addition to the safety warnings listed, warnings are posted throughout the operating manual. These warnings are designated by an exclamation mark inside an equilateral triangle. "Warning of a dangerous situation (Attention! Please follow the documentation)." The danger is classified using a signal word. Read and follow these important instructions for averting dangers.

**Warning**

Describes a **possibly** highly dangerous situation. If these instructions are not followed, serious injury and danger to life could result.

**Caution**

Describes a **possibly** dangerous situation. If this is not avoided, slight or minor injuries could result. A warning of possible property damage may also be contained in the text.

**Notice**

Describes a **possibly** harmful situation. If this is not avoided, the product or anything in its surroundings can be damaged.

**Note!**

Draws attention to something special.

**Important!**

Indicates usage tips and other useful information.

This icon is used in the operating instructions to indicate flashing values or parameters which have to be set or confirmed.

### **Before initial operation**

These operating instructions must be carefully studied! It contains important information on the connection to the local mains supply, correct unit usage and safe handling.

Check for transportation damage during unpacking.

Do not try to start up a damaged unit before the damage has been repaired or you have ascertained the effect of the damage. Get in contact with supplier and/or carrier for settlement of damage claims

Repairs, alterations or modifications must only be carried out by specialist personnel. Considerable damage can be caused by improper repairs. The SI Analytics GmbH service department is at your disposal for repair work.

### **Initial operation**

This device is not designed according to the IEC 60601-1 standard and should not be operated in rooms used for medical purposes and/or in the vicinity of patients.

Only connect the unit to a power socket with an earthing contact (PE – protective earth). Electrical connections and any other work must be performed by qualified personnel only. The power supply plug serves as a safe disconnecting device from the line and must always be easily accessible. Never operate units with damaged mains power cables.

### **Operation**

Only use this unit solely for the intended application.

Never operate the unit without bath fluid in the bath! Only use the heat transfer liquids recommended by SI Analytics GmbH

Do not operate the unit with wet or oily hands.

Do not expose the unit to spray water or immerse it in water.

Employ suitable connecting tubing.

Make sure that the tubes are securely attached.

### **Cleaning**

Do not clean the unit with solvents (fire risk!) – a wet cloth soaked in household detergent is normally sufficient.

### **Transportation**

Always empty the bath before moving the unit. Transport the unit with care. Sudden jolts or drops may cause damage in the interior of the unit.

**Notice:**

The bath fluid has to be checked at regular intervals (cleanness and chemical condition). The cooling circuit can be strongly damaged, especially by chemical contamination.

**Notice:**

A distance of at least 20 cm on each side must be maintained for ventilation, allowing internal heat to be conducted away from the unit.

---

## Safety recommendations

### The following measures were taken for the protection of the operator:

- Protection Class I according to VDE 0106 T1 (IEC536) i.e. protection against electric shocks by grounding all parts which carry the risk of electric contact.
- Protection IP 21 according to EN 60529, i.e. regarding the protection against accidentally touching live parts and damage by foreign matter, it has been ensured that foreign bodies with a thickness or diameter of more than 12 mm cannot penetrate.

Notice:

### No special precautions were taken against the penetration of water and dust.

- The device should therefore not be used in a dusty atmosphere or in the neighborhood of spray water.
- Do not insert wires or tools in any of the openings.


Notice:

### Shock hazard

Complete separation from mains is required when:

- all dangers caused by this device are to be avoided,
- cleaning is carried out;
- repairs or maintenance work is about to be carried out

Complete separation means:

## Pull out the mains plug!

### Disposal

Valid in EU countries

See the current official journal of the European Union – WEEE directive.

Directive of the European Parliament and of the Council on waste electrical and electronic equipment (WEEE).

This directive requires electrical and electronic equipment marked with a crossed-out trash can to be disposed of separately in an environmentally friendly manner.

Contact an authorized waste management company in your country. Disposal with household waste (unsorted waste) or similar collections of municipal waste is not permitted!

## 1. Technical data of Flow-Through Cooler CK300/CK310

Translation of the legally binding German version

(Version as of 03.07.2014)

Country of origin: Germany / Made in Germany

Temperature ranges:

Working temperature range (20 °C ambient temperature)	[°C]:	+ 5 ... + 80
Ambient temperature range	[°C]:	+ 5 ... + 40

Cooling capacity: 360 W at 20 °C working temperature

Refrigerant: R134a

Connections:

mechanical:  
nozzles 2 x 12 mm dia. for Immersion Thermostat (CT 72)

electrical:  
6-pole connector for control cable (Immersion Thermostat CT72)

Mains power connection: 230 V~, 50/60 Hz; on request 115 V~, 60 Hz

Total power consumption: [W]: 400

Measurements: Total, incl line cord connection:  
(W x D x H) [cm]: 21x40x46

Weight: Approx. [kg]: 22,4


Environmental conditions acc. to EN61010-1

: Ambient temperature: +5 ... +40 °C  
Max. rel. humidity 80 % for temperatures up to +31 °C,  
decreasing linearly to 50 % rel. humidity at +40 °C

Protection class: IP 21 acc. EN 60 529

Power supply: acc. to class 1,  
**not for use in explosive atmosphere**

## 2. Operating controls and functional elements

Example: CK300

1  Mains power switch, illuminated

2  6-pin connecting socket for adapter cable

3  Tube connection - IN

4  Tube connection - OUT

5  Circuit breaker

6  Mains power cable with plug

### 3. Preparations

#### 3.1. Location

- The instrument should be set up at a frost-proof and dry location.
- The ambient temperature must not exceed 40 °C.
- A distance of at least 20 cm on each side must be maintained for ventilation, allowing internal heat to be conducted away from the unit.
- Wait at least one hour after setting up a unit before placing it into operation. Mishandling during transport (e.g., laying on side) may have caused oil to enter the evaporator, capillary tube, or injection valve, considerably reducing cooling performance

#### 3.2. Bath liquids

1. SI Analytics GmbH recommends the use of silicone oil or mineral oil with a maximum viscosity of  $30 \text{ mm}^2 \times \text{S}^{-1}$   
Temperature range: -20 °C to +80 °C
2. Mixture water/glycol (3:1  $\approx$  -10 °C)  
Temperature range: +5 °C bis +80 C°
3. Water (mixing ratio: 70 % deionized water and 30 % tap water)  
Temperature range: +5 °C bis +80 C
- 4.

**Notice:**

Carefully read the safety data sheet of the bath liquid used, particularly with regard to the flash point.

**Glycol: - Low toxicity-**

Strictly observe the safety data and handling instructions from the manufacturer.


**Notice:**

No liability for use with water. Danger of freezing at working temperatures  $<5 \text{ }^\circ\text{C}$ .

**Do not use distilled or de-ionized water.**

(Mixing ratio: 70% deionized water and 30% tap water).

### 3.3. Electrical connections

For the connection use the shielded control lead supplied with the Transparent Thermostat.

CK300/CK310

CT72

10 Connect the control lead between the connectors


### 3.4. Connecting the hoses

Recommended in combination with CT 53, CT54, CT72

11 Connect the tube connection OUT to pump connector Return of the immersion thermostat

12 Connect the tube connection IN to pump connector Feed of the immersion thermostat

Nozzles 2 x 12 mm dia

**Notice:**

- Make sure that the tubes are securely attached.
- Employ suitable connecting tubing (temperature range).
- Avoid mechanical work (torsion, bending) to be applied to tubing.
- Avoid sharp bends in the tubing, and maintain a sufficient distance from surrounding walls.
- Regularly check the tubing for material defects (e.g. for cracks).

## 4. Operating procedures

### 4.1. Power connection

Caution

- Only connect the unit to a power socket with earthing contact (PE – protective earth).
- The power supply plug serves as safe disconnecting device from the line and must be always easily accessible.
- Never operate equipment with damaged mains power cables.
- Regularly check the mains power cables for material defects (e.g. for cracks).
- We disclaim all liability for damage caused by incorrect line voltages!
- Check to make sure that the line voltage matches the supply voltage specified on the identification plate.  
Deviations of  $\pm 10\%$  are permissible.

### 4.2. Switching on

- Turn on the mains power switch (1) to place the unit into operation. The control light in the switch will illuminate.

The compressor starts only when the immersion thermostat mounted on the transparent bath is switched on by pressing the Start key. The flow-through cooler receives the corresponding signal via the control lead.

(see 4.3. Electrical connections)

**The flow-through cooler is switched via the control lead during normal operation and also if a failure appears.**

## 5. Troubleshooting

### Malfunction of compressor:

The cooling compressor is equipped with an overload protection device that will be triggered by overheating or excessive current consumption. Possible causes include insufficient ventilation or contamination of the condenser.

After a cool-down phase, the motor is automatically switched on again.

### Interruption of the cooling loop by a bended tube.

In such case it is essential that the bath liquid is suitable for lower temperatures than the active working temperature. There is every risk for the heat exchanger of the flow-through cooler to freeze (e.g. using a water-glycol mixture).


Both disturbances might cause the temperature in the transparent bath to rise. The immersion thermostat provides a warning function that instantly reports a temperature deviation if set correctly.

## 6. Cleaning the unit / Maintenance

Caution

- Always turn off the unit and disconnect the mains cable from the power source before cleaning the unit.
- Prevent humidity from entering into the circulator.
- Electrical connections and any other work must be performed by qualified personnel only.

The Flow-Through Cooler is designed for continuous operation under normal conditions. Periodic maintenance is not required.

Regularly check the condenser for dirt contamination. Clean the ribbed condenser, because dust and dirt will reduce cooling performance of the unit.

Cleaning the Cooling Compressor:

- Switch off the unit, disconnect mains power cable.
- Remove the hood.
- Clean the ribbed condenser with a vacuum cleaner.
- Replace the hood or the ventilation grid.
- Switch on the unit.

Clean the outside of the unit using a wet cloth and low surface tension water.

Repairs:

Before asking for a service technician or returning a circulator for repair, please contact our service station.

SI Analytics GmbH  
Hattenbergstrasse 10  
D-55122 Mainz  
Telefon +49 61 31 / 66 51 11  
Telefax +49 61 31 / 66 50 01  
E-Mail: [si-analytics@xyleminc.com](mailto:si-analytics@xyleminc.com)  
[www.si-analytics.com](http://www.si-analytics.com)

When returning a unit, take care of careful and adequate packing. SI Analytics GmbH is not responsible for damages that might occur from insufficient packin

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