Testing high polymers e.g. PETP or PBTP

Testing of polymers

Measurement problem

One of the major quality features of synthetic materials is the mean molecular weight of the polymer molecules. The molecular weight characterises the chain length of the molecules which has a decisive influence on the processing properties of a synthetic material.

The strain exerted on the plastic by the processing process may lead to changes in the polymer changes (as a rule a decay of the chains). Under certain circumstances the properties of the finished part might be changed to such an extent that it's no longer suitable for its intended purpose.

In the research and development of polymers new polymers are being developed and produced. In this process, too, the chain length of the polymer molecules is of essential importance as to the characterisation of the finished product.

This results in the following measurement tasks:

Polymer research and development

- Determination of the mean chain length or mean polymerisation degree of the polymer molecules
- Objectives:
 - Characterisation of the finished product
 - Optimisation of its chemical and physical properties
 - Rating of polymerisation installations
 - Determination of process parameters

Polymer chemistry (polymer production)

- Determination of the mean chain length or mean degree of polymerisation of the finished product (raw granules)
- Objectives:
 - Characterisation of the finished product
 - Quality assurance
 - Optimisation of the process parameters
 - Prevention of the production of spoiled batches

Polymer processing

- Characterisation of the properties and the capabilities of the starting material (raw granules)
- Objectives:
 - Rating of plants for polymer processing
 - Determination of optimum process parameters
- Determination of the chemical and physical properties of the finished part
- Objectives:
 - Quality assurance
 - Optimisation of the process parameters

Solution:

That determination of the chain length, processing properties and quality of a synthetic material is done in a form of viscosity measurement on solutions of the plastic in suitable solvents using capillary viscometers (solution viscometry). The following table will inform you about the solvents, viscometers, the application and the relevant standards.

| Polymer | Abbrev. | Solvent | Capillaries | DIN |
|--|------------|--|-------------------------------------|----------------------------|
| Polyamide | PA | Formic acid (90%) Sulphuric acid (96%) m-cresol | I, Micro Ic II, Micro IIc II | 53 727 53 727 53 727 |
| Polycarbonate | PC | Dichlormethane | Oc | 7744, Part 2 |
| Polyethylenterephthalate Polybutylenterephthalate | PET PBT | Phenol / 1,2-dichloro benzene (1:1 parts by weight) 2-chlorophenol m-Kresol Dichloro acetic acid | Ic Ic II II, Mikro IIc | 53 728, Teil 3 |
| Polyvinyl chloride | PVC | Cyclohexanone Tetrahydrofurane | I Ic | 53 726 |
| Polyethylene Polypropylene | PE PP | Decahydronaphthalene (Decalin) | I | 53 728, Sheet 4 |
| Polystyrol | PS | Toluene o-Xylol 1,2-Dichlor benene | I I I | 7741, Part 2 |
| Polymethylmethacrylate | PMMA | Chloroform Acetophenone | Oc I | 7745, Part 2 |
| Cellulose acetate | CA | Dichlormethane/ Methanole (9:1 part by volume) | Oc | 53 7228, Sheet 1 |

Further information about viscometry you'll find in our brochure "Theory and practice of capillary viscometry".

The **viscosity number** (for a definition, please refer the following table) gives information about the processibility of plastic material. It plays a decisive role within a framework of quality control of the granules. In addition it's important to verify the viscosity number of the finished plastic part.

In most cases the indication of the viscosity number or of the **relative viscosity** (please refer the following table) is sufficient as a quality criterion of established plants. This requires the determination of the solvent and of the plastic solution (concentration mostly 0,5 g/100 ml).

Instead of the viscosity number the determination frequently involves the **K value after Fikentscher**.

The determination of the mean molecular weight of the polymer molecules is done via the **limiting viscosity number** (please refer the following table). It's of particular importance in the range of research and development of polymers and of procedures and installations for their production and processing. In addition it's an important feature as regards quality assurance in case of special applications, such as plastic recycling and the processing of recycled plastics.

To determine the limiting viscosity number polymer solutions of different concentrations are produced (so called dilute series). The limiting viscosity number results from extrapolation of the viscosity numbers to a concentration = 0.

Definition and description of the terms

| Value | Description |
|---|--|
| η | dynamic viscosity |
| $\nu = \eta / \rho$ | kinematic viscosity |
| $\eta_r = \eta / \eta_S$ | relative viscosity, viscosity ratio |
| $(\eta - \eta_S) / \eta_S = \eta_r - 1$ | relative viscosity change, specific viscosity |
| $J_v = 1 / c \cdot (\eta - \eta_S) / \eta_S$ | Staudinger function, viscosity number |
| $\text{Ln} (\eta / \eta_S) / c$ | inherent viscosity |
| $J_g = \lim_{c \rightarrow 0} [1 / c \cdot (\eta - \eta_S) / \eta_S]$ | Staudinger index, limiting viscosity number, intrinsic viscosity |

Necessary basic equipment

The basic equipment should have an appropriate standard of a modern chemical laboratory.

That means, the bench top surface should be solvent resistant and should be ventilation equipped (exhaust cabine) for existing solvent evaporations, probably.

The picture shows an AVSPro automatic sampler system for the determination of the solution viscometry on high polymers of SCHOTT Instruments GmbH in a suitable laboratory environment.

Necessary measuring equipment

| | |
|---|---|
| ViscoSystem | Type: AVS 360, AVS 370 or AVS 470 |
| Glass panelled thermostat | Type: CT 53/ CT 54 |
| Flow through cooler | Type: CK 300 |
| Thermometer | Type: VZ 7101 (for 25°C) |
| Measuring stand | Type: AVS/SK (PVDF) |
| PTFE tube cable combination | Type: VZ 5622 |
| Overflow safety device | Type: VZ 6245 (AVS 360), VZ 8552 (AVS 370 or AVS 470) |
| Ubbelohde viscometer (standard version) | Type: 530 20 or suitable depending to the application |
| Micro-Ubbelohde viscometer | Type: 538 23 or suitable depending to the application |
| Software | Type: WinVisco 370 |
| Data cable | Type: TZ 1599 (AVS 360), VZ 7116 (AVS 370) |
| as option: | half automatic discharging (Please request for a detailed offer) |

Necessary laboratory equipment:

(order numbers according to SCHOTT laboratory glass catalogue)

Analytical balance (resolution min. 1 mg)

Magnetic stirrer and/or **shaker**

Water jet or another **vacuum pump** with connection to a filtering flask

Weighing glasses (no.: 24 210 ..)

Glass funnels, div. diameters (no.: 21 331 ..)

Spatula, forceps

Volumetric flask 50 ml, 100 ml (no.: 21 678 17, 21 678 24)

Filtering device (250 ml) with conical gasket (no.: 21 201 36, 29 202 27)

Filter crucibles 30 or 50 ml, e.g. porosity 2 (no.: 25 851 22, 25 851 32)

Screw cap flasks of glass 100 ml (no.: 21 801 24)

Pasteur pipets

div. **magnetic stirrer bars (PTFE coated)**

pH paper

div. **beakers** (no.: 21 106 .. or 21 116 ..)

Waste container for acids

Necessary solvents:

concentrated sulfuric acid, m-kresol, Dichlorbenzene-phenol, Dichlore acidic acid (dep. on application)

(refer the table on page 3)

distilled or diionized water

Acetone (p.a.) for drying the viscometers, quickly

The sample preparation never should be done without e.g. analytical balance, heatable magnetic stirrer or shaker.

The reason:

Only a very accurate sample concentration grants for reasonable measuring results!

Sample preparation:

In laboratory practice are two different styles of sample preparation usual, momentary.

For example – Polyester (PETP / PBTP):

1. A sample ($0,25 \pm 0,002$) g of PETP or PBTB should be weight in with an accuracy of 0.1 mg and filled on quantity in a volume flask of 50 ml. After this approx. 25 ml of solvent will be added until the polymer sample is completely diluted at approx. 80°C under warming up in a water bath or under shaking. Please remark that longer heating may cause a thermal reduction of the PETB or PBTB condition. A protection again too long dilution times it would be possible to use a thermal oil bath with higher temperature. Nevertheless, the combination between dilution time and dilution temperature is to check if there was reduction of the polymer structure, time by time.

For arbitral reasons it seems recommended to make an agreement about dilution time and temperature.

After cooling down the solution to approx. 20°C the volumetric flask will be filled up to the measuring mark with solvent which has the same temperature. The finished solution must be mixed by shaking.

The four steps of traditional sample preparation with volumetric flasks.

2. In the last time there are in use motor burettes for solvent dosing, more and more (see herefore our available T 110 *plus*). For this reason the sample will be weight into any sample bottle and the necessary volume of solvent will be dosed with the motor burette according to the weight.

The advantage of this method is that here its possible to use a heatable magnetic stirrer or a thermal block for diluting the polymer.

The advantage is clear!

The uncomfortable and in parts dangerous manual handling of aggressive and sometimes corrosive solvents is missing nearly completely!

For the further operation the validity is the same as described under 1..

The following graphic shows the most modern style of sample preparation.


With this style most of possible failure reasons will be eliminated in practice.

The weight data will be transferred from the analytical balance directly to a sample preparation system (intelligent burette or PC controlled system), directly.

The sample preparation system calculates the necessary volume of solvent for the polymer weight and eliminates failures (accurate sample concentration) for the sample preparation.

The reason:

The very time intensive and difficult weighting of the necessary polymer granules will be eliminated (it should be accurate to **1 mg [0,001 g]**).

Experiences out of the industrial practice confirm that using motor burettes minimises the possible failure of a incorrect or possible wrong solution concentration.

The following procedure description is authentic and advises as example.

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Frequently occurring errors when using glass capillary viscometers

| Error symptom | Error causes | Possible error elimination |
|---|---|---|
| Syst. measurement error: Flow time too large with short flow times | After-flow error, Hagenbach correction too small | Experimental determination of the Hagenbach correction using substances having a similar viscosity and a surface tension as the measurement product |
| Syst. measurement error: Flow time too small with short flow times | After-flow error, Hagenbach correction too large | As above, better: Choose a viscometer with a smaller capillary diameter |
| Syst. measurement error: Flow time too small with (Ostwald-, Cannon-Fenske- or BS/IP/RF tube reverse flow viscometer) | Substance quantity filled in too small | Empty, clean and refill viscometer |
| Syst. measurement error: Flow time too long (Ostwald-, Cannon-Fenske- or BS/IP/RF tube reverse flow viscometer) | Substance quantity filled in was too big | as above |
| Syst. measurement error: Flow time too small with short flow times (Ubbelohde viscometer) | Disturbance of the suspended level | Select viscometer with a smaller capillary diameter |
| Syst. measurement error: Flow time too short | Temperature of the bath liquid too high | Check temperature; if necessary, readjust the thermostat |
| Syst. measurement error: Flow time too long | Contamination in the capillary | Empty and clean viscometer; repeat the measurement |
| | Temperature of the bath liquid too low | Check temperature; if necessary, readjust the thermostat |
| Drift of the flow times | Drift of the bath temperature | Protect the thermostat from direct radiation exposure; if necessary replace the thermostat |
| | Temperature adjustment of the measuring substance not completed | Continue temperature adjustment until the time values are stabil |
| | Evaporation of a high volatile component; reaction of the product being analysed with the air | Apply pressing operation mode |

| Error symptom | Error causes | Possible error elimination |
|---|--|---|
| Increased stochastic scattering of the measurement value | Contamination in the viscometer | Empty and clean viscometer; Repeat measurement |
| | Contamination in the sample solution | Empty and clean viscometer; repeat measurement with fresh filtrated sample; use probably a filter with smaller porosity |
| | Air bubbles in the viscometer | In case calf pure matters with chemical and physical heat resistance, drive out bubbles by shorter increase of temperature |
| | Air bubbles in the viscometer | Clean and empty viscometer; during refilling, ensure absence of bubbles |
| Excessive stochastic scattering occurring during measurements using optoelectric barriers (Possibility of total malfunction) | Contamination of the optical sensors | Remove the measuring stand out of the thermostatic bath; clean optical system using nondenatured alcohol on a soft cloth |
| | Errors triggered by the optoelectric barriers as a result of the formation of bubbles, foam or liquid lamellae | Use a TC Ubbelohde-, Ostwald- or Cannon-Fenske-Routine viscometer |
| Excessive stochastic scattering occurring during automatic measurements using TC viscometers (possibility of total malfunction) | Incrustation of the sensors (in case of thermally instable media) | Transparent media: use optical detection Opaque media: Use reverse flow viscometer |
| | Wear and tear of the sensors | Replace viscometer |
| Increased stochastic scattering in the case of short flow times (Ubbelohde viscometer) | Beginning deformation of the suspended level | Select a viscometer with a smaller capillary diameter |
| Periodically fluctuating flow times | Heating up or cooling down phases of the thermostat is too long | Set the heating and cooling of the thermostat in such a manner that at least two complete cycles are completed during one viscosity measurement cycle |
| | No timely stability of the bath liquid temperature (defective thermostat) | Replace the thermostat |
| Malfunction caused by air bubbles during the sucking-in process of the liquid into the delivery vessel | Substance quantity filled in was too small | Ubbelohde viscometer: fill up the measurement substance Others: empty and clean viscometer, repeat measurement |

Source of supply for standards:

You may download the given standards in this description at Beuth publishing
Over the internet.

The internet address is as follows:

www2.beuth.de

Source material: DIN 51 562, Part 1 – 4
 DIN 53 727;
 DIN 53 728, Part 3
 DIN 53 733;
 ASTM D 445
 ASTM D 446
 ISO 3105
 ISO 3106
 div. documentation of SCHOTT Instruments GmbH
 div. documentation of O. Hofbeck

Please have appreciation if we decline requests for copies of any national or international standards or other assumed documentations.

It's not legal to duplicate and transfer copies of standards.

A transfer of this document is only allowed in its genuine version.

If you want to get more information about our laboratory program, please visit our webside or order our laboratory catalogue. You'll find detailed descriptions about our viscometry-, titration-, pH-electrodes and pH-measuring instruments and accessories.

Ottmar Hofbeck
Product management viscometry
ottmar.hofbeck@schottinstruments.com
050206

SCHOTT Instruments GmbH
Hattenbergstr. 10
D – 55122 Mainz
www.schottinstruments.com