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Upgrade 2010 for photoLab® 6000 / spectroFlex Series

What is New in Version 2.00?

New COD methods for WTW test kits COD1-TC, COD2-TC, COD3-TC

For the COD tests COD-1 TC, COD-2 TC and COD-3 TC (order nr. 251 990, 251 991, 251 992) without barcode and lot certificate, method data are available now.

Color Measurement

ADMI spectrophotometric method acc. APHA 2120 F with measurement at 31 wavelengths: 400, 410,...700nm

Determination of Ammonia NH₃

Method for determination of ammonia content via pH and temperature. Before measuring Ammonium concentration with WTW kit model 14752, pH and temperature have to be determined and entered during method procedure. Ammonia finally will be calculated and displayed. Measuring range: 0.010-3.000 mg/l NH₄-N, NH₃ depending on pH and temperature.

Reagent-free Nitrate Determination

Determination of Nitrate acc. to Standard Methods 4500-B by absorbance measurement at 220 und 275 nm.

Measuring range: 0.0 - 7.0 mg/l

Printing data to PDF

This solves potential incompatibilities of printers: Data can be printed to PDF and transferred to the standard PC environment via USB-Stick. This requires the setting *PDF file* for the PRINT key function in Menu *Setting - Data transfer / Printer*

User-defined Methods – Method Transfer

With Version 2.00, user-defined methods, which have been programmed on a model of photoLab® 6000 Serie, can be copied to USB and transferred to another instrument of this Series. Eingabe von methodenspezifische Grenzwerte für Overflow / Underflow-Anzeige

Zero adjustment for programmed methods

Due to many new options for measuring routines some changes are effective for the zero adjustment of the programmed methods. This goes along with many advantages for the user (see separate paragraph). By means of the well-proven menu navigation the changes are transparent for present users and easy to accommodate.

Multi-Wavelength Method

Enhancements for user defined programming options, incl. complex routines.

Setting of preconditions (optional) for validity of measurement results

Entering of measurement range limits

Operation Manuals:

The English and German version will be posted on WTW download page within December, the other languages following within early 2011. Analytical description might follow a little bit time shifted.

New Languages:

Swedish, Serbian, Thai, Vietnamese, Macedonian

Zero Adjustment of programmed Methods

Due to enhanced functionality of methods and special application methods, the zero adjustment for concentration measurement has been adapted. These changes are beneficial for many applications. The last valid zero adjustment is displayed on screen.

The following table shows the difference and effects of zero adjustment between Version 1.xx and Version 2.00.

	Zero Adjustment V. 1.xx	Zero Adjustment V. 2.00
	With ZERO-BLANK key the zero adjustment for all wavelengths of programmed methods is performed at once.	With ZERO-BLANK key the zero adjustment only for the wavelength of the selected methods is performed.
Factory Zero	A factory zero of rectangular cell (10, 20, 50 mm) and the 16 mm round cell for all wavelengths of programmed methods is performed and stored.	A factory zero of rectangular cell (10, 20, 50 mm) and the 16 mm round cell for the 12 Wellenlängen of programmed methods using barcodes is performed and stored. (Cell tests or reagent test kits with AutoSelector).
First use in Concentration	The user will be asked, if the factory zero should be taken over. It will be effective for the wavelengths of all programmed methods.	The user will be asked, if the factory zero should be taken over. It will be effective for the 12 wavelengths used by barcoded method data.
Zero adjustment by user in concentration mode	Pressing the ZERO-BLANK key a new zero will be performed with the actually used cell and selected cell size. This zero is effective for all wavelengths used in programmed methods at once.	Pressing the ZERO-BLANK key a new zero will be performed with the actually used cell and selected cell size. This zero is effective for the actual wavelength of selected method.
Effect	Using a special cell (Quartz, Half-micro or plastic) the zero adjustment for all programmed methods is changed. When using different cuvette for different methods, a zero must be performed each time, which then again changes the zero for all wavelength at once.	The zero adjustment takes place for the selected method / wavelength only . Thus, different cell types can be used method specifically . Also an adjustment with other solutions than water can be used method specifically, (if refraction index of the solution differs not strongly from water).
Information about zero adjustment	Selecting a method, user is asked each time, if existing zero (date xx, time xx) shall be continued to use.	Date and Time of last valid zero adjustment is displayed in the upper right corner.
Effect	Enquiry for each method selection	Die letzte Null wird immer angezeigt. No enquiry