
**Bedienungsanleitung
Operating manual
Mode d'emploi
Instrucciones de operación**

pH 1970i

Portables pH/mV-Meßgerät

Seite 3

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**Safety
instructions**

Safety

The individual chapters of this operating manual use safety instructions such as that shown below to indicate various hazards or dangers:

Warning


indicates instructions that must be followed precisely in order to avoid the possibility of slight injuries or damage to the instrument or the environment.

pH 1970i - List of contents


Safety	17
Display and socket field	18
Battery / mains operation	19
Switching on the measuring instrument	20
Measuring	20
Calibration	21
Saving data	25
Outputting the data memory	25
Erasing the data storage	25
Transmitting data	25
Recorder (Analog output)	26
Configuring	26
Resetting to default settings	27
Disposal	28
Technical data	29

Display and socket field

Display

Socket field

Sensor / Instrument	Socket / Position
pH depth armature The sockets, 1 and 2 must not be assigned at the same time.	1
pH electrode or pH combination electrode	2
pH combination electrode with temperature sensor	2 and 4
Reference electrode	3
Temperature sensor, external	3 and 4
Printer, recorder or PC (serial interface, RS232)	5
Plug-in power supply unit	6
Watertight valve for internal pressure equalization	7

**Warning**

Only connect electrodes to the measuring instrument that do not return any unallowed voltages or currents ($> \text{SELV}$ and $> \text{current circuit with current limiting}$).

Almost all electrodes - in particular WTW electrodes - fulfill these conditions.

Battery / mains operation


The plug-in power supply supplies the measuring instrument with low voltage (12 V DC). At the same time, the rechargeable battery is charged (approx. 16 hours). The *LoBat* display indicator appears when the rechargeable battery is nearly empty.

**Warning**

The line voltage at the operating site must lie within the input voltage range of the original plug-in power supply (see TECHNICAL DATA).

**Warning**

Only use original plug-in power supplies (see TECHNICAL DATA).

- Connect the plug (1) to the socket (2) of the measuring instrument.
- Connect the original WTW plug-in power supply (3) to an easily accessible mains socket.

Switching on the measuring instrument

- Press the **<ON/OFF>** key.
The *Display test* appears briefly on the display.
After this, the measuring instrument automatically switches to the measuring mode. The display shows the relevant measured value.

Measuring

Overview of the measuring modes:

Special functions:

AutoRead (drift control)

The AutoRead function checks the stability of the measurement signal (except for the measurement of the ORP voltage). Activate AutoRead with **<AR>**. Press **<RUN/ENTER>** to start the AutoRead measurement. During the AutoRead measurement, *AR* flashes on the display until a stable measured value is reached. This can be terminated at any time taking over the current measured value with **<RUN/ENTER>**.

Temperature measurement during pH measurements

You can use pH electrodes with and without temperature sensors. The measuring instrument recognizes which temperature sensors are connected and switches automatically to the correct mode for the temperature measurement (display *TP*).

Manual temperature input: Set the temperature value with **<▲>** **<▼>**.

Calibration

AutoCal TEC

This process is specially adapted to the WTW technical buffer solutions as a fully automatic **single, two-point** or **three-point calibration**. The buffer solutions are automatically recognized by the measuring instrument. Depending on the instrument setting, the instrument displays the relevant buffer nominal value or the current electrode voltage in mV.

Valid buffers (values at 25 °C): 2.00 / 4.01 / 7.00 / 10.01

Note

The calibration for pH 10.01 is optimized for the WTW technical buffer solution TEP 10 Trace or TPL 10 Trace. Other buffer solutions can lead to an erroneous calibration. The correct buffer solutions are given in the WTW catalog or in the Internet.

AutoCal DIN

as AutoCal TEC, but matched to buffer solutions according to DIN 19266.

Valid buffers (values at 25 °C): 1.679 / 4.006 / 6.865 / 9.180

ConCal

is the conventional **two-point calibration** using two buffer solutions (pH 7.0 ± 0.5 and any other buffer solution) or a **single-point calibration** using any buffer solution as a high-speed method.

Calibration evaluation

After calibrating, the measuring instrument automatically evaluates the calibration. The asymmetry and slope are evaluated separately. The worst evaluation appears on the display.

Display

					 Invalid calibration
Asymmetry [mV]	-15 ... +15	-20 ... +20	-25 ... +25	-30 ... +30	< -30 or > 30
Slope [mV/pH]	-60.5 ... -58	-58 ... -57	-61 ... -60.5 or -57 ... -56	-62 ... -61 or -56 ... -50	< -62 or > -50


Note

If a printer is connected to the interface, a calibration protocol is automatically printed after a valid calibration.

AutoCal TEC / AutoCal DIN procedure:

The following example shows the AutoCal TEC calibration. In the AutoCal DIN calibration, *Cd...* appears on the display instead of *Ct...*. Otherwise, both procedures are identical.

The selection of the calibration procedure, AutoCal TEC (display *Ct1*) or AutoCal DIN (display *Cd1*) is made with the **<CAL>** key.

- Connect the pH electrode to the measuring instrument.
- Press the **<CAL>** key repeatedly until the *Ct1* display indicator and the *Auto Cal TEC* function display appear. The sensor symbol displays the evaluation of the last calibration (or no sensor symbol in the delivery state or after the measuring parameters have been reset).
- Immerse the pH electrode in the first buffer solution.
- When measuring without a temperature sensor:
Enter the temperature of the first buffer solution with **<▲>** **<▼>**.
- Press **<RUN/ENTER>**. The AutoRead measurement begins. The electrode voltage (mV) or the buffer nominal value (see CONFIGURING) appears on the display. When the measured value is stable, *Ct2* appears.


**Note**

At this point, the AutoCal TEC calibration can be terminated with **<M>**. This corresponds to a **single-point calibration**. To do this, the instrument uses the Nernst slope (-59.2 mV/pH at 25 °C) and determines the asymmetry of the electrode.


- To continue the **two-point calibration**, immerse the pH electrode in the second buffer solution.
- When measuring without a temperature sensor:
Enter the temperature of the second buffer solution with **<▲>** **<▼>**.
- Press **<RUN/ENTER>**. The AutoRead measurement begins. The electrode voltage (mV) or the buffer nominal value appears on the display.
- When the measured value is stable, the instrument displays the value of the slope and the evaluation of the two-point calibration.
- Press **<RUN/ENTER>**. The instrument displays the value of the asymmetry.

**Note**

At this point, the AutoCal TEC calibration can be terminated with **<M>**. This corresponds to a **two-point calibration**.

- Immerse the pH electrode in the second buffer solution to continue the **two-point calibration**.
- When measuring without a temperature sensor:
Enter the temperature of the third buffer solution with <▲> <▼>.
- Press <RUN/ENTER>. The AutoRead measurement begins. The electrode voltage (mV) or the buffer nominal value appears on the display.
- When the measured value is stable, the instrument shows the value of the slope and the evaluation of the three-point calibration.
- Press <RUN/ENTER>. The instrument displays the value of the asymmetry.
- Switch to the measuring mode with <M>.

ConCal procedure:

- Connect the pH electrode to the measuring instrument.
- Press the **<CAL>** key repeatedly until the **ASY** display indicator and the **Cal** function display appear. The sensor symbol displays the evaluation of the last calibration (or no sensor symbol in the delivery state or after the measuring parameters have been reset).
- Immerse the pH electrode in the first buffer solution (pH 7.0 ± 0.5 in two-point calibration).
- When measuring without a temperature sensor:
Enter the temperature of the first buffer solution with **<▲>** **<▼>** while pressing the **<RUN/ENTER>** key.
- Press the **<RUN/ENTER>** key.
The measured pH value appears on the display.
- Set the measured value to the nominal pH value of the buffer solution (at the current temperature) with **<▲>** **<▼>**.
- When the measured value is stable, press the **<RUN/ENTER>** key. The value of the asymmetry appears.
- Press the **<RUN/ENTER>** key. **SLO** appears.

**Note**

At this point, the ConCal calibration can be terminated with **<M>**. This corresponds to a **single-point calibration**. To do this, the instrument uses the Nernst slope (-59.2 mV/pH at 25 °C) and determines the asymmetry of the electrode.

- To continue the **two-point calibration**, immerse the pH electrode in the second buffer solution.
- When measuring without a temperature sensor:
Enter the temperature of the second buffer solution with **<▲>** **<▼>** while pressing the **<RUN/ENTER>** key.
- Press **<RUN/ENTER>**.
- Set the measured value to the nominal pH value of the buffer solution (at the current temperature) with **<▲>** **<▼>**.
- When the measured value is stable, press the **<RUN/ENTER>** key. The instrument displays the value of the slope and the calibration evaluation.
- Press **<RUN/ENTER>**. The instrument displays the value of the asymmetry.
- Switch to the measuring mode with **<M>**.

Saving data

The measuring instrument has an internal data memory. It can store up to 800 datasets.

Saving data manually

Press the **<STO>** key in the measuring mode (display *No.* with the number of the next free memory location). Then, press **<RUN/ENTER>** and enter the ID number with **<▲>** **<▼>**. Finish the save with **<RUN/ENTER>**. The instrument changes to the measuring mode.

Saving data automatically *Int 1*

The save interval (*Int 1*) determines the chronological interval between automatic save processes.

To set up the save interval, press **<STO>** while pressing the **<RUN/ENTER>** key (display *Int 1*) and set the interval with **<▲>** **<▼>** key. Then press **<RUN/ENTER>**. The number of the available free memory locations is displayed. Press **<RUN/ENTER>** once more and enter the ID number with **<▲>** **<▼>**. Switch to the last active measuring mode with **<RUN/ENTER>**. The automatic save is switched on (display *Auto Store*).

Outputting the data memory

You can output the data memory with the **<RCL>** key. By repeatedly pressing the **<RCL>** key, you reach the following functions:

- StO dISP* Output stored data on the display
- StO SEr* Output stored data on the serial interface
- CAL dISP* Output calibration data on the display
- CAL SEr* Output calibration data on the serial interface

Start the output with **<RUN/ENTER>**.

Erasing the data storage

Switch off the measuring instrument. Press the **<STO>** key and keep it depressed. Press the **<ON/OFF>** key. Confirm the erasing procedure with **<RUN/ENTER>** or cancel the erasing procedure by pressing any other key.

Transmitting data

Transmitting data manually

Press the **<RUN/ENTER>** key in the measuring mode. This manually triggers a data transmission of the current measured value to the serial interface at any time - independently of the selected intervals.

Transmitting data automatically *Int 2*

The interval for the data transmission (*Int 2*) determines the chronological interval between automatic data transmissions. After the selected interval expires, the current data record is transmitted to the interface. To set up the transmission interval, press **<RCL>** while holding down the **<RUN/ENTER>** key (display *Int 2*). Then, set the interval with **<▲>** **<▼>**.

Recorder (Analog output)

Via the analog output, you can transmit data to a recorder. Connect the analog output to the recorder using the AK323/S interface cable. The data output automatically switches to *Recorder output*.

Note

The output on the analog output corresponds to the value shown on the display (pH 7.00 is output as 700 mV).

Configuring


Note

You can leave the configuration menu at any time with **<M>**. The parameters that have already been changed are stored.

Procedure for configuring (factory settings appear in bold typeface):

- Switch off the measuring instrument.
- Press **<ON/OFF>>** while pressing the **<M>** key.
- Select the required Baud rate with **<▲>** **<▼>**.
Selection: 1200, 2400, **4800**, 9600 Baud.
- Then press the **<RUN/ENTER>** key.
- Select the required display during the pH calibration with **<▲>** **<▼>**.
Selection: **Buffer nominal value** (/pH) or current electrode voltage (mV).
- Then press the **<RUN/ENTER>** key.
- Set the required calibration interval with **<▲>** **<▼>**.
Selection: 1 ... **7** ... 999 d.
- Then press the **<RUN/ENTER>** key.
- Select the date and time step-by-step with **<▲>** **<▼>**. In between, press the **<RUN/ENTER>** key each time.
- When the last **<RUN/ENTER>** key has been pressed, the instrument switches to the last active measuring mode.

Resetting to default settings

You can reset the measurement parameters and the configuration to the delivery status separately from one another (initialization).

Measurement parameters

The following measurement parameters can be reset to the delivery status:


Measuring mode	pH
Slope	-59.16 mV/pH
Asymmetry	0 mV
Manual temperature input	25 °C

Configuration parameters

The following configuration parameters (*InI*) can be reset to the delivery status:

Baud rate	4800
Display during pH calibration	Buffer nominal value
Interval 1 (automatically saved)	OFF
Interval 2 (for data transmission)	OFF

Reset procedure:

- Press **<CAL>** while pressing the **<RUN/ENTER>** key.
- To reset the measurement parameters, use **<▲>** **<▼>** to select **YES** and confirm with **<RUN/ENTER>**, or
- Continue to the configuration parameters without resetting (*InI*) with the **<RUN/ENTER>** key.
- After the *InI* configuration parameters, the instrument switches to the last active measurement mode.

Disposal

Remove the rechargeable battery from the instrument and dispose of it at a suitable facility according to local legal requirements.


Disassembling the rechargeable battery

Ni-Cd

Cd

- Screw off the fastening elements (1).
- Remove the instrument from the enclosure by vigorously pressing against the socket field. Then take out the rechargeable battery (2) and cut off the battery cable.

Technical data

Dimensions	approx. 90 x 200 x 190 mm	
Weight	approx. 1.5 kg (without plug-in power supply)	
Mechanical structure	Type of protection	IP 67
Electrical safety	Protective class	III
Test certificates	CE	
Ambient conditions	Storage	- 25 °C ... + 65 °C
	Operation	-10 °C ... + 55 °C
	Climatic class	2
Power supply	Rechargeable battery	Nickel-cadmium(NiCad) rechargeable battery
	Operational life	approx. 600 hours with one charging
	Plug-in power supply unit (charging device)	The following applies to all plug-in power supplies: Connection max. Overvoltage category II Plug-in power supply unit (Euro, US , UK, Australian plug) FRIWO FW7555M/09, 15.1432 Friwo Part. No. 1883259 Input: 100 ... 240 V ~ / 50 ... 60 Hz / 400 mA Output: 9 V = / 1,5 A
Serial interface (AK 340/B or AK 325/S cable)	Type	RS232, data output
	Baud rate	adjustable: 1200, 2400, 4800, 9600 Baud
	Data bits	8
	Stop bits	2
	Parity	None
	Handshake	RTS/CTS + Xon/Xoff
	Cable length	Max. 15m

**Analog output
(AK 323/S cable)**

Automatic switchover when the recorder is connected by the cable, AK 323/S.

pH output signal	-200 ... +1999 mV for the range - 2.00 ... + 19.99
mV output signal	-1999 ... +1999 mV for the range -1999 ... +1999 mV for the range
Accuracy	± 0.5 % of display value
Internal resistance	< 5 Ohm (current limited to max. 0.2 mA output current)

**Guidelines
and norms used**

EMC	EC guideline 89/336/EEC EN 61326-1:1998 EN 61000-3-2 A14:2000 EN 61000-3-3:1995 FCC Class A
Instrument safety	EC guideline 73/23/EEC EN 61010-1 A2:1995
Climatic class	VDI/VDE 3540
IP protection	EN 60529:1991

Measuring ranges, resolutions, accuracies

Variable	Measuring range	Resolution	Accuracy (± 1 Digit)
pH	- 2.00 ... + 19.99	0.01	± 0.01 *
U [mV]	- 199.9 ... + 199.9 - 1999 ... + 1999	0.1 1	± 0.5 ± 1
T [°C]	- 5.0 ... + 105.0	0.1	± 0.1

**Manual
temperature
input**

Variable	Range	Increment
T _{manual} [°C]	- 20 ... + 130	1

* when measuring in a range of ± 2 pH around a calibration point

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Copyright
Wissenschaftlich-Technische Werkstätten GmbH