


pH/Oxi 3400i



pH/Oxygen measuring instrument

**Accuracy when
going to press**

The use of advanced technology and the high quality standard of our instruments are the result of continuous development. This may result in differences between this operating manual and your instrument. Also, we cannot guarantee that there are absolutely no errors in this manual. Therefore, we are sure you will understand that we cannot accept any legal claims resulting from the data, figures or descriptions.

Warranty

We guarantee the instrument described for 3 years from the date of purchase.

The instrument warranty covers manufacturing faults that are discovered within the warranty period. The warranty does not cover components that are replaced during maintenance work, e.g. batteries.

The warranty claim extends to restoring the instrument to readiness for use but not, however, to any further claim for damages. Improper handling or unauthorized opening of the instrument invalidates any warranty claim.

To ascertain the warranty liability, return the instrument and proof of purchase together with the date of purchase freight paid or prepaid.

Copyright

© Weilheim 2007, WTW GmbH

Reproduction in whole - or even in part - is prohibited without the express written permission of WTW GmbH, Weilheim.

Printed in Germany.

1	Overview	5
1.1	General features	5
1.2	SETs of equipment	6
1.3	Keypad	8
1.4	Display	9
1.5	Jack field	9
1.6	Operating structure	10
2	Safety	11
2.1	Authorized use	11
2.2	General safety instructions	12
3	Commissioning	13
3.1	Scope of delivery	13
3.2	Initial commissioning	13
4	Operation	15
4.1	Switching on the measuring instrument	15
4.2	pH value / ORP voltage	16
4.2.1	General information	16
4.2.2	Measuring the pH value	18
4.2.3	Measuring the ORP voltage	20
4.3	pH calibration	21
4.3.1	AutoCal TEC	23
4.3.2	AutoCal DIN	27
4.3.3	ConCal	30
4.4	Dissolved oxygen	33
4.4.1	General information	33
4.4.2	Measuring the D. O. concentration	34
4.4.3	Measuring the D. O. saturation	35
4.4.4	AutoRead AR (Drift control)	36
4.4.5	D. O. calibration	37
4.4.6	Entering the salt content (salinity)	40
4.5	Calibration intervals (Int 3, Int 4)	41
4.6	Saving	42
4.6.1	Saving manually	42
4.6.2	Switching on AutoStore (Int 1)	44
4.6.3	Outputting the data storage	46
4.6.4	Clearing the memory	51
4.7	Transmitting data	52
4.7.1	Data transmission interval (Int 2)	52
4.7.2	PC/external printer (RS232 interface)	53
4.7.3	Remote control	54
4.8	Configuration	55
4.9	Reset	58


5	Maintenance, cleaning, disposal	61
5.1	Maintenance	61
5.2	Cleaning	62
5.3	Disposal	62
6	What to do if... ..	63
6.1	pH system messages	63
6.2	Oxi system messages	65
6.3	General errors	66
7	Technical data	67
8	Lists	71

1 Overview

1.1 General features

This compact precision handheld meter enables you to carry out pH / ORP measurements and oxygen measurements rapidly and reliably. The handheld meter provides the maximum degree of operating comfort, reliability and measuring certainty for all applications.

The proven MultiCal[®] calibration procedures and the special Auto-Read function support you in your work with the handheld meter.



1	Keypad
2	Sample display
3	Jack field



Note

If you need further information or application notes, you can obtain the following material from WTW:


- Application reports
- Primers
- Safety datasheets.

You will find information on available literature in the WTW catalog or via the Internet.

1.2 SETs of equipment

The measuring instrument is also available as part of individual SETs of equipment.


You will find additional information on this and other accessories in the WTW catalog or via the Internet.



Set (sample configuration):

1	pH/Oxi 3400i measuring instrument
2	50 ml ELY/G electrolyte solution for D.O. probes
3	Stand
4	RL/G cleaning solution for D.O. probes, 50 ml
5	50 ml pH buffer solution, STP 4
6	50 ml pH buffer solution, STP 7
7	Storing solution for pH electrodes
8	CellOx 325 D.O. probe with calibration vessel
9	SF300 abrasive film for D.O. probes
10	pH combination electrode
11	Plug-in power supply (optional)
12	3 WP 90/3 exchange membrane caps


1.3 Keypad




Key functions

	Select the measured variable <M> : – pH value / ORP voltage – D. O. concentration/saturation/partial pressure
	Save a measured value <STO>
	Switch measuring instrument on/off <ON/OFF>
	Calibrate the currently set measured variable <CAL>
	Display/transmit measured values <RCL>
	Activate/deactivate the AutoRead function <AR>
	Select the measuring mode, increase values, scroll <▲ >
	Select the measuring mode, decrease values, scroll <▼ >
	Confirm entries, start AutoRead <RUN/ENTER>

1.4 Display




1.5 Jack field



1	D. O. probe
2	pH electrode
3	pH temperature sensor
4	Plug-in power supply (optional)
5	Serial interface RS 232

1.6 Operating structure

The following overview diagram shows which keys you have to press to select between the different measuring modes:



Note

When connecting a D. O. probe with temperature sensor, the measuring instrument recognizes the probe and automatically switches to the last selected D. O. measuring mode. As soon as the probe is disconnected, the instrument switches to the pH or ORP measuring mode.

2 Safety

This operating manual contains basic instructions that you must follow during the commissioning, operation and maintenance of the measuring instrument. Consequently, all responsible personnel must read this operating manual before working with the measuring system. The operating manual must always be available within the vicinity of the measuring system.

Target group

The measuring instrument was developed for work in the field and in the laboratory.

We assume that, as a result of their professional training and experience, the operators will know the necessary safety precautions to take when handling chemicals.

Safety instructions

The individual chapters of this operating manual use the following safety labels to indicate various types of danger:



Warning

indicates instructions that must be followed precisely in order to avoid the possibility of slight injuries or damage to the instrument or the environment.

Further notes



Note

indicates notes that draw your attention to special features.



Note

indicates cross-references to other documents, e.g. operating manuals.

2.1 Authorized use

The authorized use of the measuring instrument consists exclusively of the pH and ORP measurement or the measurement of the oxygen content in the field and laboratory.

The technical specifications as given in chapter 7 TECHNICAL DATA must be observed. Only the operation and running of the measuring instrument according to the instructions given in this operating manual is authorized.

Any other use is considered to be **unauthorized**.

2.2 General safety instructions

This instrument is built and inspected according to the relevant guidelines and norms for electronic measuring instruments (see chapter 7 TECHNICAL DATA).

It left the factory in a safe and secure technical condition.

Function and operating safety

The smooth functioning and operational safety of the measuring instrument can only be guaranteed if the generally applicable safety measures and the specific safety instructions in this operating manual are followed during operation.

The smooth functioning and operational safety of the measuring instrument can only be guaranteed under the environmental conditions that are specified in chapter 6 WHAT TO DO IF....

If the instrument was transported from a cold environment to a warm environment, the formation of condensate can lead to the faulty functioning of the instrument. In this event, wait until the temperature of the instrument reaches room temperature before putting the instrument back into operation.

Safe operation

If safe operation is no longer possible, the instrument must be taken out of service and secured against inadvertent operation!

Safe operation is no longer possible if the measuring instrument:

- has been damaged in transport
- has been stored under adverse conditions for a lengthy period of time
- is visibly damaged
- no longer operates as described in this manual.

If you are in any doubt, please contact the supplier of the instrument.

Obligations of the purchaser

The purchaser of the measuring instrument must ensure that the following laws and guidelines are observed when using dangerous substances:

- EEC directives for protective labor legislation
- National protective labor legislation
- Safety regulations
- Safety datasheets of the chemical manufacturers.

3 Commissioning

3.1 Scope of delivery

- 3400i handheld meter
- Plug-in power supply (optional)
- Operating manual and short operating manual
- 4 batteries, 1.5 V Mignon type AA (in the instrument)

For details of scope of delivery of SETs, see chapter 1.2 SETs of EQUIPMENT and WTW catalog.

3.2 Initial commissioning

Perform the following activities:

- Set the date and time
- Connect the plug-in power supply (optional)

Setting the date and time

1	Press the <M> key and hold it down.
2	Press the <ON/OFF> key. The display test appears briefly on the display.
3	Press the <RUN/ENTER> key repeatedly until the date appears on the display.
4	Set the date of the current day with <▲> <▼> .
5	Confirm with <RUN/ENTER> . The date (month) flashes in the display.
6	Set the current month with <▲> <▼> .
7	Confirm with <RUN/ENTER> . The year appears on the display.
8	Set the current year with <▲> <▼> .
9	Confirm with <RUN/ENTER> . The hours flash on the display.
10	Set the current time with <▲> <▼> .
11	Confirm with <RUN/ENTER> . The minutes flash on the display.
12	Set the current time with <▲> <▼> .

- | | |
|----|---|
| 13 | Confirm with <RUN/ENTER> .
The instrument switches to a measuring mode. |
| 14 | Switch the instrument off using <ON/OFF> . |

Connecting the original plug-in power supply unit




You can either operate the measuring instrument with batteries or with the plug-in power supply. The plug-in power supply supplies the measuring instrument with low voltage (12 VDC). This saves the batteries.

Warning

The line voltage at the operating site must lie within the input voltage range of the original plug-in power supply (see chapter 7 TECHNICAL DATA).

Warning

Use original plug-in power supplies only (see chapter 7 TECHNICAL DATA).



- | | |
|---|---|
| 1 | Plug the jack (1) into the socket (2) of the measuring instrument. |
| 2 | Connect the original WTW plug-in power supply (3) to an easily accessible mains socket. |

4 Operation

4.1 Switching on the measuring instrument

- | | |
|---|---|
| 1 | Press the <ON/OFF> key.
The display test appears briefly on the display.
After this, the measuring instrument automatically switches to the measuring mode. |
|---|---|

Measuring mode when switching on

Connected sensor	Measuring mode
No sensor or pH/ORP electrode	pH or ORP measurement (depending on the last setting)
D. O. probe	Last selected measuring mode
2 sensors	Last selected measuring mode



Note

The measuring instrument has an energy saving feature to avoid unnecessary battery depletion. The energy saving feature switches the measuring instrument off if no key has been pressed for an hour.

The energy saving feature is not active

- if the power is supplied by the plug-in power supply,
- if the AutoStore function is active,
- if the communication cable and a PC with a running communication program are connected,
- if the printer cable is connected (for external printers).

4.2 pH value / ORP voltage

4.2.1 General information

Preparatory activities

Perform the following preparatory activities when you want to measure:

1	Connect the pH electrode to the measuring instrument. If necessary, press the <M> key repeatedly until the status display <i>pH</i> (pH measurement) or <i>U</i> (measurement of the ORP voltage) appears.
2	Adjust the temperature of the buffer solutions or test solutions, or measure the current temperature, if you measure without a temperature sensor.
3	Calibrate or check the measuring instrument with the electrode.
4	Using <▲ > <▼> , toggle between the <i>pH</i> or <i>mV</i> measuring modes.



Note

Incorrect calibration of pH electrodes leads to incorrect measured values. Calibrate regularly before measuring. You can only connect electrodes of the NTC30 type or without temperature sensor.



Warning

When connecting an earthed PC/printer, measurements cannot be performed in earthed media as incorrect values would result. The RS232 interface is not galvanically isolated.

Temperature measurement in pH measurements

You can perform pH measurements with or without a temperature sensor as well as with the temperature sensor of a D. O. probe. The measuring instrument recognizes which sensors are connected and automatically switches to the correct mode for the temperature measurement.

The following cases are distinguishable.

Temperature sensor		Display	Mode
<i>pH</i>	<i>Oxi</i>		
yes	-	<i>TP</i>	Automatic with pH temperature sensor
yes	yes	<i>TP</i>	
-	-		Manual
-	yes	<i>TP</i> flashes	The temperature value of the 2nd sensor (Oxi) in the same test sample is used for the pH measurement*

* If you do not wish that, you can:

- either disconnect the 2nd sensor and use the manual temperature input or
- use an electrode with a temperature sensor.

If a temperature sensor is connected, it is indicated on the display by *TP*.




Note

When calibrating without a temperature sensor (no *TP* display indicator displayed), enter the current temperature of the respective buffer solution manually using the <▲> <▼> keys while keeping the <RUN/ENTER> key depressed.

4.2.2 Measuring the pH value

1	Perform the preparatory activities according to section 4.2.1.
2	Immerse the pH electrode in the test sample.
3	Press the <▲> <▼> keys until <i>pH</i> appears in the status display. The pH value appears on the display.



4	When measuring without a connected temperature sensor: Options: ● Determine the current temperature using a thermometer and, while keeping the <RUN/ENTER> key depressed, enter this temperature value with <▲> <▼>. ● <i>TP</i> display indicator not displayed, socket for the second probe is free: Connect the second probe (Oxi) and immerse it in the same test sample. <i>TP</i> flashes, the temperature is automatically measured using the second probe.
---	--

AutoRead AR (drift control)

The AutoRead function (drift control) checks the stability of the measurement signal. The stability has a considerable impact on the reproducibility of the measured values.

For identical measurement conditions, the following criteria apply:

Reproducibility	Response time
Better than 0.02	> 30 seconds

1	Call up the pH measuring mode with <M> and <▲> <▼> .
2	Activate the AutoRead function with <AR>. The current measured value is frozen (hold function).

3	Start AutoRead with <RUN/ENTER> . <i>AR</i> flashes until a stable measured value is reached. This measured value is transmitted to the interface.
4	If necessary, start the next AutoRead measurement with <RUN/ENTER> .
5	To terminate the AutoRead function: Press the <AR> key.


**Note**

The current AutoRead measurement can be terminated at any time (accepting the current value) by pressing **<RUN/ENTER>**.

4.2.3 Measuring the ORP voltage

In conjunction with an ORP electrode, e.g. SenTix ORP, the measuring instrument can measure the ORP voltage (mV) of a solution.

1	Perform the preparatory activities according to section 4.2.1.
2	Submerge the ORP electrode in the sample.
3	Press the <▲> <▼> keys until the U status display appears. The ORP voltage (mV) of the test sample appears on the display.
4	Wait for a stable measured value.

**Note**

ORP electrodes are not calibrated. However, you can check ORP electrodes using a test solution.

4.3 pH calibration

Why calibrate?

pH electrodes age. This changes the asymmetry (zero point) and slope of the pH electrode. As a result, an inexact measured value is displayed. Calibration determines the current values of the asymmetry and slope of the electrode and stores them in the measuring instrument. Thus, you should calibrate at regular intervals.

When to calibrate?

- After connecting another electrode
- When the sensor symbol flashes
 - (after the calibration interval has expired)
 - after a voltage interruption, e. g. after changing the batteries

You can select one of 3 calibration procedures:

AutoCal TEC

is specially matched to the WTW technical buffer solutions as a fully automatic **two-point calibration**. The buffer solutions are automatically recognized by the measuring instrument. Depending on the instrument setting (see section 4.8 CONFIGURATION), the instrument displays the relevant buffer nominal value or the current electrode voltage in mV. The calibration can be terminated after the first buffer solution. This corresponds to a **single-point calibration**. When doing so, the instrument uses the Nernst slope (-59,2 mV/pH at 25 °C) and determines the asymmetry of the electrode.

AutoCal DIN

is specially adapted to the permanently programmed buffer solutions in accordance with DIN 19266 as a fully automatic **two-point calibration**. The operating sequence of the AutoCal DIN calibration corresponds to that of the AutoCal TEC calibration. The calibration can only be terminated after the first buffer solution (**single point calibration**).

ConCal

This function is a conventional **two-point calibration** using two buffer solutions ($\text{pH } 7.0 \pm 0.5$ and any other buffer solution) or a **single-point calibration** using any buffer solution, which is used as a high-speed method.

AutoRead

When calibrating with AutoCal TEC and AutoCal DIN, the AutoRead function is automatically activated.

The current AutoRead measurement can be terminated at any time (accepting the current value) by pressing **<RUN/ENTER>**.

Displaying the calibration data

You can view the data of the last calibration on the display. The proceeding is described on page 48.

Printing the calibration record

The calibration protocol contains the calibration data of the current calibration. You can transmit the calibration protocol to a printer via the serial interface (see page 50).



Note

You can automatically print a calibration protocol after the calibration. To do so, connect a printer to the interface according to section 4.7.2 before calibrating. After a valid calibration, the record is printed.

Sample printout:

```
CALIBRATION PROTOCOL
02.03.01      14:19
Device No.: 12345678
Calibration pH
Cal time: 01.03.01 / 15:20
Cal interval: 7d
AutoCal TEC   Tauto
Buffer 1      2.00
Buffer 2      4.01
Buffer 3      7.00 *
Buffer 4      10.00
C1           184.1 mV  25.0°C
C2            3.0 mV  25.0°C
S1           -59.4 mV/pH
ASY1         -    4 mV
Probe:       +++
```

Calibration evaluation

After the calibration, the measuring instrument automatically evaluates the calibration. The asymmetry and slope are evaluated separately. The worst evaluation appears on the display.

Display	Asymmetry [mV]	Slope [mV/pH]
	-15 ... +15	-60.5 ... -58
	-20 ... +20	-58 ... -57
	-25 ... +25	-61 ... -60.5 or -57 ... -56
 Clean the electrode according to the electrode operating manual	-30 ... +30	-62 ... -61 or -56 ... -50
 Eliminate the error according to chap- ter 6 WHAT TO DO IF...	< -30 or > 30	... -62 or ... -50

Preparatory activities

1	Connect the pH electrode to the measuring instrument. If necessary, press the <M> key repeatedly until the status display <i>pH</i> (pH measurement) or <i>U</i> (measurement of the ORP voltage) appears.
2	Keep the buffer solutions ready.
3	Adjust the temperature of the solution and measure the current temperature if the measurement is made without the use of a temperature sensor (the <i>TP</i> display indicator is missing from the display).

4.3.1 AutoCal TEC

For this procedure, use any two WTW technical buffer solutions (pH values at 25 °C: 2.00 / 4.01 / 7.00 / 10.01).



Note

The calibration for pH 10.01 is optimized for the WTW technical buffer solution TEP 10 Trace or TPL 10 Trace. Other buffer solutions can lead to an erroneous calibration. The correct buffer solutions are given in the WTW catalog or in the Internet.



Note


The buffer solutions are automatically recognized by the measuring instrument. Depending on the instrument setting (see section 4.8 CONFIGURATION), the instrument displays the relevant buffer nominal value or the current electrode voltage in mV.





Note

Skip the steps 2 and 7 if you use a pH electrode with a temperature sensor or the temperature sensor of a D. O. probe.


1	Press the <CAL> key repeatedly until the <i>Ct1</i> display indicator and the function display <i>AutoCal TEC</i> appears. The sensor symbol displays the evaluation of the last calibration (or no sensor symbol in the delivery state or after the measurement parameters have been reset).
---	--



- | | |
|---|--|
| 2 | If required, enter the temperature of the first buffer solution with $\blacktriangle$ $\blacktriangledown$ while keeping the <RUN/ENTER> key depressed. |
| 3 | Immerse the pH electrode in the first buffer solution. |
| 4 | Press the <RUN/ENTER> key.
The <i>AR</i> display indicator flashes.
The electrode voltage (mV) or the buffer nominal value appears on the display. Example: |




- | | |
|---|--|
| 5 | When the measured value is stable, <i>Ct2</i> appears. |
|---|--|

**Note**


At this point, the AutoCal TEC calibration can be terminated with **<M>**. This corresponds to a **single-point calibration**. When doing so, the in-

strument uses the Nernst slope ($-59,2 \text{ mV/pH}$ at 25°C) and determines the asymmetry of the electrode.


6	Thoroughly rinse the electrode with deionized water.
7	If required, enter the temperature of the second buffer solution with $\blacktriangle$ $\blacktriangledown$ while keeping the <RUN/ENTER> key depressed.
8	Immerse the pH electrode in the second buffer solution.
9	Press the <RUN/ENTER> key. The <i>AR</i> display indicator flashes. The electrode voltage (mV) or the buffer nominal value appears on the display. Example:



10	When the measured value is stable, <i>AR</i> disappears. The value of the slope (mV/pH) appears on the display. The probe symbol shows the evaluation of the current calibration.
----	--



11	Press the <RUN/ENTER> key. The value of the asymmetry (mV) appears on the display.
----	---



- | | |
|----|--|
| 12 | Switch to the measuring mode with <M> . |
|----|--|

4.3.2 AutoCal DIN


For this procedure, use two different standard buffer solutions according to DIN 19266 (type A, C, D or F with pH values at 25 °C of: 1.679 / 4.006 / 6.865 / 9.180).




Note

Skip the steps 2 and 7 if you use a pH electrode with a temperature sensor or the temperature sensor of a D. O. probe.


- 1 Press the **<CAL>** key repeatedly until the display *Cd1* and the function display *AutoCal DIN* appear. The sensor symbol displays the evaluation of the last calibration (or no sensor symbol in the delivery state or after the measurement parameter has been reset).



- 2 If required, enter the temperature of the first buffer solution with **<▲>** **<▼>** while keeping the **<RUN/ENTER>** key depressed.
- 3 Immerse the pH electrode in the first buffer solution.
- 4 Press the **<RUN/ENTER>** key.
The *AR* display indicator flashes.
The electrode voltage (mV) or the buffer nominal value appears on the display. Example:




- | | |
|---|--|
| 5 | When the measured value is stable, <i>Cd2</i> appears. |
|---|--|


**Note**

At this point, the AutoCal DIN calibration can be terminated with **<M>**. This corresponds to a **single-point calibration**. When doing so, the instrument uses the Nernst slope (-59,2 mV/pH at 25 °C) and determines the asymmetry of the electrode.


- | | |
|---|--|
| 6 | Thoroughly rinse the electrode with deionized water. |
| 7 | If required, enter the temperature of the second buffer solution with <▲> <▼> while keeping the <RUN/ENTER> key depressed. |
| 8 | Immerse the pH electrode in the second buffer solution. |
| 9 | Press the <RUN/ENTER> key.
The AR display indicator flashes.
The electrode voltage (mV) or the buffer nominal value appears on the display. Example: |



- 10 When the measured value is stable, *AR* disappears. The value of the slope (mV/pH) appears on the display. The probe symbol shows the evaluation of the current calibration.



- 11 Press the **<RUN/ENTER>** key. The value of the asymmetry (mV) appears on the display.



- 12 Switch to the measuring mode with **<M>**.

4.3.3 ConCal

Two-point calibration

For this procedure, use two buffer solutions:


- pH 7.0 ± 0.5
- any other buffer solution



Note

Skip the steps 2 and 8 if you use a pH electrode with a temperature sensor or the temperature sensor of a D. O. probe.

1	Press the <CAL> key repeatedly until the <i>ASY</i> display and the <i>Cal</i> function display appears. The sensor symbol displays the evaluation of the last calibration (or no sensor symbol in the delivery state or after the measurement parameter has been reset).
---	--




2	If required, enter the temperature of the first buffer solution with <▲> <▼> while keeping the <RUN/ENTER> key depressed.
3	Immerse the pH electrode in the first buffer solution (pH 7.0 ± 0.5 in two-point calibration).
4	Press the <RUN/ENTER> key. The measured pH value appears on the display.
5	Set the measured value to the nominal pH value of the buffer solution (at the current temperature) with <▲> <▼> .
6	When the measured value is stable, press the <RUN/ENTER> key. The value of the asymmetry appears.
7	Press the <RUN/ENTER> key. <i>SLO</i> appears.




Note


At this point, the ConCal calibration can be broken off with **<M>**. This corresponds to a **single-point calibration**. When doing so, the instrument uses the Nernst slope ($-59,2 \text{ mV/pH}$ at 25°C) and the fixed asymmetry of the electrode.



8	Thoroughly rinse the electrode with deionized water.
9	If required, enter the temperature of the second buffer solution with <▲> <▼> while keeping the <RUN/ENTER> key depressed.
10	Immerse the pH electrode in the second buffer solution.
11	Press the <RUN/ENTER> key.
12	Set the measured value to the nominal pH value of the buffer solution (at the current temperature) with <▲> <▼> .
13	When the measured value is stable, press the <RUN/ENTER> key. The value of the slope (mV/pH) appears on the display. The probe symbol shows the evaluation of the current calibration.



14	Press the <RUN/ENTER> key. The value of the asymmetry (mV) appears on the display.
----	---



- | | |
|----|--|
| 15 | Switch to the measuring mode with <M>. |
|----|--|

4.4 Dissolved oxygen

4.4.1 General information



Note

D. O. measurements with the pH/ Oxi3400i can only be carried out using a CellOx 325 or StirrOx G D. O. probe. The stirrer of the StirrOx G D. O. probe has to be supplied with voltage separately using the NT/pH Mix 540 power supply.

You can measure the following variables:

- D. O. concentration
- D. O. saturation

The pH/Oxi 3400i handheld meter is equipped with the following functions:

- AutoRange (automatic measurement range selection),
- The AutoRead function (drift control) for checking the stability of the measurement signal. This ensures the reproducibility of the measurement signal. For details of how to switch the AutoRead function on/off, see page 36.

Warning

When connecting an earthed PC/printer, measurements cannot be performed in earthed media as incorrect values would result. The RS232 interface is not galvanically isolated.

Preparatory activities

Perform the following preparatory activities when you want to measure:

1	Connect the D. O. probe to the measuring instrument. The instrument automatically switches over to D. O. measurement (O_2 status display). If the D. O. probe is already connected, press the <M> repeatedly until the O_2 status display appears.
2	Calibrate or check the measuring instrument with the probe. How to calibrate is described in section 4.4.5 from page 37.
3	Immerse the D. O. probe in the test sample.
4	Use <▲> <▼> to toggle between the measuring modes, D. O. concentration (mg/L) and D. O. saturation (%).



Note

Incorrect calibration of D. O. probes will result in incorrect measured values. Calibrate at regular intervals.

Temperature sensor

The D. O. probe has an integrated temperature sensor that always measures the current temperature of the test sample.

4.4.2 Measuring the D. O. concentration


When measuring the concentration of test samples with a salt content of more than 1 g/l, a salinity correction is required.

**Note**

How to enter the current salt content is described in section 4.4.6 ENTERING THE SALT CONTENT (SALINITY) on page 40. Switching the salt content correction on or off, see below.


To measure the D. O. concentration with and without salt content correction, proceed as follows:

1	Perform the preparatory activities according to section 4.4.1.
2	Immerse the D. O. probe in the test sample.
3	Press the <▲> <▼> key repeatedly until the D. O. concentration in <i>mg/l</i> appears on the display.

**Switching on/off salt content correction**

Proceed as follows to switch on the salt content correction:

1	Perform the preparatory activities according to section 4.4.1.
2	While pressing the <RUN/ENTER> key, switch on the salt content correction with <▲> . The <i>SAL</i> display indicator appears on the display. The specified salt content is taken into consideration during the measurement.




- | | |
|---|--|
| 3 | While pressing the <RUN/ENTER> key, switch off the salt content correction with <▼>. The SAL display indicator is no longer displayed. |
|---|--|

4.4.3 Measuring the D. O. saturation

You can measure the D. O. saturation as follows:


- | | |
|---|--|
| 1 | Perform the preparatory activities according to section 4.4.1. |
| 2 | Immerse the D. O. probe in the test sample. |
| 3 | Press the <▲> <▼> key repeatedly until the D. O. saturation in % appears on the display. |



4.4.4 AutoRead AR (Drift control)

The AutoRead function (drift control) checks the stability of the measurement signal. The stability has a considerable impact on the reproducibility of the measured values.

1	Call up the measuring mode with <M> and/or <▲> <▼> .
2	Activate the AutoRead function with <AR> . The current measured value is frozen (hold function).
3	Start AutoRead with <RUN/ENTER> . <i>AR</i> flashes until a stable measured value is reached. This measured value is transmitted to the interface.



4	If necessary, start the next AutoRead measurement with <RUN/ENTER> .
5	To terminate the AutoRead function: Press the <AR> key.

Criteria With identical measurement conditions, the following applies:

Measured parameter	Reproducibility	Response time
D. O. concentration	better than 0.05 mg/l	> 10 seconds
D. O. saturation index	better than 0.6 %	> 10 seconds

4.4.5 D. O. calibration

Why calibrate?	D. O. probes age. This changes the slope of the D. O. probe. Calibration determines the current slope of the probe and stores this value in the instrument.
When to calibrate?	● After connecting another D. O. probe ● When the probe symbol flashes (after the calibration interval has expired).
Calibration procedures	The calibration is performed in water vapor-saturated air. Use the OxiCal [®] -SL air calibration vessel (accessories) for the calibration.
AutoRead	The calibration procedure automatically activates the AutoRead function. The <i>AR</i> display indicator flashes. The calibration process is finished when <i>AR</i> stops flashing.
Displaying the calibration data	You can view the data of the last calibration on the display. The proceeding is described on page 48.
Printing the calibration protocol	The calibration record contains the calibration data of the current calibration. You can transmit the calibration protocol to a printer via the serial interface (see page 50).

**Note**

You can automatically print a calibration protocol after the calibration. To do so, connect a printer to the interface according to section 4.7.2 before calibrating. After a valid calibration, the record is printed.

Sample printout:

```
CALIBRATION PROTOCOL
02.03.01      14:19
Device No.: 12345678
CALIBRATION 02
Cal time: 02.03.01 / 14:19
Cal interval: 14d
OxiCal      Tauto AR
Relative Slope: 0,88
Probe:      +++
```

Probe evaluation


After the calibration, the measuring instrument evaluates the current status of the probe against the relative slope. The evaluation appears on the display. The relative slope has no effect on the measuring accuracy. Low values indicate that the electrolytic filling will soon be depleted and the probe has to be regenerated.

Display	Relative slope
	S = 0.8 ... 1.25
	S = 0.7 ... 0.8
	S = 0.6 ... 0.7
 Eliminate the error according to chapter 6 WHAT TO DO IF...	S < 0.6 or S > 1.25

Starting the calibration


Proceed as follows to calibrate the instrument:

- | | |
|---|--|
| 1 | Perform the preparatory activities according to section 4.4.1. |
| 2 | Keep the OxiCal®-SL air calibration vessel ready. |




Note
The sponge in the air calibration vessel must be moist (not wet). Follow the instructions in the OxiCal®-SL operating manual.


- | | |
|---|--|
| 3 | Put the D. O. probe into the air calibration vessel. |
| 4 | Press the <CAL> key repeatedly until the calibration mode appears. The sensor symbol displays the evaluation of the last calibration (or no sensor symbol in the delivery state or after the measurement parameter has been reset). |



- | | |
|---|--|
| 5 | Press the <RUN/ENTER> key. AutoRead is active, <i>AR</i> flashes. |
|---|--|



- | | |
|---|---|
| 6 | As soon as a stable value is reached, the <i>AR</i> display indicator stops flashing. The calibration is finished now. The value of the relative slope appears on the display. The probe symbol shows the probe evaluation (see page 38). |
|---|---|



- | | |
|---|--|
| 7 | Switch to the measuring mode with <M> . |
|---|--|

**Note**


In chapter 6 WHAT TO DO IF... from page 63, you will find the measures to take for error elimination.

4.4.6 Entering the salt content (salinity)

A salt content correction is required in the oxygen concentration measurement of test samples with a salt content of more than 1 g/l. To do this, you have to enter the salinity equivalent of the test sample (range 0.0 - 70.0) and to switch on the salinity correction.

Entering the salinity

- | | |
|---|--|
| 1 | Determine the salinity of the test sample. |
| 2 | Press the <CAL> key repeatedly until Sal appears on the display. |



- | | |
|---|---|
| 3 | Enter the salt content with <▲> <▼> . |
| 4 | Switch to the measuring mode with <M> . |

**Note**

How to switch on the salt content correction is described on page 34.



4.5 Calibration intervals (Int 3, Int 4)

For each measured variable, a time interval is stored. When it has expired, you will be reminded to calibrate. After a calibration interval has expired, the probe symbol of the relevant measured variable flashes. It is still possible to measure. By calibrating the relevant probe, the function is reset and the interval starts anew.


The following calibration intervals are set in the factory:

Measured parameter	Designation	Default setting
pH/ORP	Int 3	7 days (d 7)
Dissolved oxygen	Int 4	14 days (d 14)

Setting the calibration interval

You can change each of these intervals (1 ... 999 days):

1	Switch off the measuring instrument.
2	Press the <M> key and hold it down.
3	Press the <ON/OFF> key. The display test appears briefly on the display. After this, the measuring instrument automatically switches over to configuration.
4	Press the <RUN/ENTER> key repeatedly until Int 3 or Int 4 and the required measured variable pH or O ₂ appears. Example:



5	Set the required interval (in days) until the next calibration with <▲> <▼>.
6	Confirm with <RUN/ENTER>.
7	Switch to the measuring mode with <M>.

4.6 Saving

The measuring instrument has an internal data memory. It can store up to 500 data records.

A complete data record consists of:

- Number of the storage location
- Date/time
- Measured values of the probes connected
- Temperature
- Temperature measuring procedure
- ID number


You can transmit measured values (data records) to the data storage in two ways:

- Save manually
- Switch on AutoStore (Int 1), (see page 44).


4.6.1 Saving manually

You can transmit a measured value to the data storage as follows:

- 1 Press the **<STO>** key.
The current number (location number *No.*) of the next free storage location appears under the current measured value on the display.



- 2 Confirm with **<RUN/ENTER>**.
The display switches to entering the ID number.



- | | |
|---|--|
| 3 | Using <▲> <▼>, set the required ID number (1 ... 999). |
| 4 | Confirm with <RUN/ENTER>.
The measured values are stored. The instrument changes to the measuring mode. |

Message **StoFull**

This message appears when all of the 500 storage locations are occupied.

You have the following options:

Saving the current measured value. The oldest measured value (storage location 1) will be overwritten by this	Press <RUN/ENTER>
Returning to the measuring mode without saving	press any key
Outputting the data storage	see section 4.6.3
Clearing the memory	see section 4.6.4

4.6.2 Switching on AutoStore (Int 1)

The save interval (Int 1) determines the chronological interval between automatic save processes. After the fixed interval has expired, the current data record is transmitted to the internal storage and to the interface.

Setting the save interval

The default setting for the save interval (Int 1) is OFF.


By this, the *AutoStore* function is switched off.

To switch the function on, set an interval (5 s, 10 s, 30 s, 1 min, 5 min, 10 min, 15 min, 30 min, 60 min):


- 1 Press the **<RUN/ENTER>** key and hold it down.
- 2 Press the **<STO>** key. *Int 1* appears on the display.



- 3 Set the required interval between the saving procedures with **<▲>** **<▼>** (Selection: 5 s, 10 s, 30 s, 1 min, 5 min, 10 min, 15 min, 30 min, 60 min).
- 4 Confirm with **<RUN/ENTER>**.
The number of free memory locations appears on the display.



- 5 Confirm with **<RUN/ENTER>**.
The prompt for the ID number appears on the display.



6	Set the required ID number with <▲> <▼>.
7	Confirm with <RUN/ENTER>. The measuring instrument switches to the last active measuring mode and starts the measuring and saving procedure. AutoStore flashes on the display.

As soon as all of the 500 storage locations are occupied, AutoStore is terminated (Int 1 = OFF). If there are not enough storage locations available for your measurements:

- Output and backup the data storage (see page 46) and
- Clear the memory (see page 51).



Note

The AutoStore function is interrupted if you start other functions, e.g. output the data storage. After the function is finished, the AutoStore function is continued. By this, however, temporal gaps in the recording of the measured values will occur.

Switching off AutoStore

Switch AutoStore off by:

- setting the save interval (Int 1) to OFF, or
- switching the measuring instrument off and then on again.

4.6.3 Outputting the data storage


You can output the contents of the data storage:

- Stored data on the display
- Calibration data on the display
- Stored data on the serial interface
- Calibration protocol to the interface

Outputting stored data on the display


1

Press the <RCL> key repeatedly until *StO dISP* appears on the display.



2

Press the <RUN/ENTER> key.
A measured value appears on the display.
The storage location of the data record is displayed for approx. 2 s, then the respective temperature appears.



You can perform the following activities:

Display further data of the data record (ID number, date, time, storage location)	Press <RUN/ENTER>
Toggle between two saved measured variables	Press <RUN/ENTER> + <M>
Advance one data record (storage location)	Press <▲>
Go back one data record (storage location)	Press <▼>

**Note**

If you want to search for a certain element of the data record (e.g. date), proceed as follows:

- | | |
|---|--|
| 1 | Using <RUN/ENTER> , select the element (e.g. date). |
| 2 | Press <▲> or <▼> repeatedly until the required element appears on the display.
After approx. 2 s the temperature of the displayed measured value appears. |

Outputting stored data to the interface

- | | |
|---|--|
| 1 | Press the <RCL> key repeatedly until <i>Sto SEr</i> appears on the display. |
|---|--|



- | | |
|---|---|
| 2 | Press the <RUN/ENTER> key.
The complete contents of the storage are transmitted to the interface. During the data transmission the instrument increments the storage numbers. After the data transmission, the instrument automatically switches to the last active measurement mode. |
|---|---|

**Note**

You can cancel the transmission with **<M>** or **<RUN/ENTER>**.

The transmitted data contains the entire contents of the storage in incrementing order of the location numbers.

Sample printout:

```

No.      1:
  09.03.99      17:10
pH 10.01      25 °C
Tman          AR
Ident : 1

No.      2:
  09.03.99      17:11
pH 10.01      24.7 °C
Tauto          AR
Ident : 1

No.      3:
  09.03.99      17:12
    305.2 mV
Tauto
Ident : 13
...


```

Outputting the calibration data on the display


- 1 Press the **<RCL>** key repeatedly until *CAL disp* appears on the display.




- 2 Press the **<RUN/ENTER>** key.
The data of the last calibration of all measured variables appears in the following sequence:
 - pH: Slope *SLO* and asymmetry *ASY*
 - Dissolved oxygen: Relative slope *SLO*
 Information concerning the calibration procedure is output as well.



- 3 Press <RUN/ENTER> to display the value of the asymmetry (mV).




- 4 Press <RUN/ENTER> to display the relative slope of the D. O. probe.



- 5 With <M> you can switch back to the last active measuring mode.

Outputting the calibration protocol on the interface

- 1 Press the **<RCL>** key repeatedly until *CAL SEr* appears on the display.



- 2 Press the **<RUN/ENTER>** key.
The calibration protocol for all measured variables is transmitted to the interface. After the data transmission, the instrument automatically switches to the last active measurement mode.

Sample printout:

```

CALIBRATION PROTOCOL
 02.03.01      14:19
Device No.: 12345678

Calibration pH
Cal time: 01.10.01 / 15:20
Cal interval: 7d
AutoCal TEC    Tauto
Buffer 1      2.00
Buffer 2      4.01
Buffer 3      7.00 *
Buffer 4      10.00
C1            174.1 mV 25.0°C
C2            -133.3 mV 25.0°C
S1            -59.4 mV/pH
ASY1         -    4 mV
Probe:       +++

CALIBRATION 02
Cal time: 02.03.01 / 14:19
Cal interval: 14d
OxiCal        Tauto AR
Relative Slope: 0,88
Probe:       +++

CALIBRATION ISE
Cal time: 01.10.01 / 15:30
Std 1      10,00 mg/l
Std 2      20,00 mg/l
Std 3      50,00 mg/l

```

4.6.4 Clearing the memory

With this function, you can delete stored data records. 500 storage locations will then be available again.

**Note**

The *Clear memory* function only appears when there are data records stored in the memory. Otherwise, the measuring instrument automatically switches to the last active measuring mode.

Proceed as follows to clear all data records:

1	Switch off the measuring instrument.
2	Press the <STO> key and hold it down.
3	Press the <ON/OFF> key. The display test appears briefly on the display.



4	Confirm the clearing process with <RUN/ENTER> . Pressing any other key prevents the clearing, the data records will remain stored.
---	--

**Note**

The calibration data remain stored and can be called up.

4.7 Transmitting data

You have the following possibilities of transmitting data:

- One of the following options:
 - With the *AutoStore* function (page 44), measured values are periodically saved internally (save interval Int 1) and output on the interface.
 - With the *Data transmission interval* function (Int 2), measured values are periodically output on the interface (see below).
- With the *Output data storage* function (page 46), calibration data or saved measured values are output on the interface.

4.7.1 Data transmission interval (Int 2)

The interval to the data transmission (Int 2) determines the chronological interval between automatic data transmissions. After the selected interval expires, the current data record is transmitted to the interface.



Setting the Data transmission interval

Note

The setting of the interval (Int 2) is only effective if the save interval (*AutoStore* function) is switched off!

The default setting for the interval is OFF.

To start the data transmission, set an interval (5 s, 10 s, 30 s, 1 min, 5 min, 10 min, 15 min, 30 min, 60 min):

- | | |
|---|--|
| 1 | Press the <RUN/ENTER> key and hold it down. |
| 2 | Press the <RCL> key. <i>Int 2</i> appears on the display. |



- | | |
|---|--|
| 3 | Set the required interval between the saving procedures with <▲> <▼> . |
|---|--|

- 4 Confirm with **<RUN/ENTER>**.
The measuring instrument switches to the last active measuring mode.

**Note**

When the *AutoStore* function is active at the same time, the data transmission is performed according to the setting of the save interval (Int 1). Set the save interval (Int 1) to OFF to activate the Data transmission *interval* (Int 2).

4.7.2 PC/external printer (RS232 interface)


Via the RS 232 interface, you can transmit the data to a PC or an external printer. Use the AK340/B (PC) or AK325/S (ext. printer) cable to connect the interface to the instruments. The data output automatically switches to the RS232 interface.

**Warning**

The RS232 interface is not galvanically isolated. When connecting an earthed PC/printer, measurements cannot be performed in earthed media as incorrect values would result.

Set up the following transmission data on the PC/printer:

Baud rate	selectable between: 1200, 2400, 4800 , 9600
Handshake	RTS/CTS + Xon/Xoff
Parity	none
Data bits	8
Stop bits	2
Cable length	Max. 15 m

Socket assignment

RS 232
REC

- 1 CTS
- 2 RxD
- 3 Ground
- 4 TxD

4.7.3 Remote control

The measuring instrument can be remotely controlled from a PC. This requires the KOM pilot communication kit. It is available as an accessory. The instrument is then controlled via commands that simulate keystrokes and request the current display contents.

**Note**

A more detailed description is provided within the scope of delivery of the communication kit.

4.8 Configuration

You can adapt the measuring instrument to your individual requirements. To do this, the following parameters can be changed (the status on delivery is marked in bold):

Baud rate	1200, 2400, 4800 , 9600
Display during the pH calibration	Buffer nominal value , or current electrode voltage
Air pressure display	Current value in mbar (no input possible)
Intervals Calibration (Int 3) in days (d)	– pH: 1 ... 7 ... 999 d <i>Int 3</i> – O ₂ 1 ... 14 ... 999 d <i>Int 4</i>
Date/time	Any



Note


You can leave the configuration menu at any time with **<M>**. The parameters that have already been changed are stored.

1	Switch off the measuring instrument.
2	Press the <M> key and hold it down.
3	Press the <ON/OFF> key. The display test appears briefly on the display. The measuring instrument then switches automatically to the setting of the baud rate.

Baudrate




4	Select the required Baud rate with <▲> <▼> .
5	Confirm with <RUN/ENTER> . <i>CAL dISP</i> appears on the display.




- | | |
|---|---|
| 6 | Select the required display during the pH calibration with <▲> <▼> .
<i>mV</i> : Display of the current electrode voltage
<i>/pH</i> : Display of the buffer nominal value. |
| 7 | Confirm with <RUN/ENTER> . On the display, <i>P mbar</i> and the current air pressure appears. |

Displaying the air pressure




- | | |
|---|---|
| 8 | Confirm with <RUN/ENTER> . <i>pH Int 3</i> and the measured variable pH appear on the display. |
|---|---|

Calibration intervals pH




- | | |
|----|---|
| 9 | Set the required interval (in days) with <▲> <▼> . |
| 10 | Confirm with <RUN/ENTER> . <i>Int 4</i> and the measured variable <i>O₂</i> appears on the display. |



- | | |
|----|--|
| 11 | Set the required interval (in days) with <▲> <▼>. |
| 12 | Confirm with <RUN/ENTER>. The date flashes on the display. |

Date and time



- | | |
|----|---|
| 13 | Set the date of the current day with <▲> <▼>. |
| 14 | Confirm with <RUN/ENTER>.
The date (month) flashes on the display. |
| 15 | Set the current month with <▲> <▼>. |
| 16 | Confirm with <RUN/ENTER>.
The year appears on the display. |
| 17 | Set the current year with <▲> <▼>. |
| 18 | Confirm with <RUN/ENTER>.
The hours flash on the display. |
| 19 | Set the current time with <▲> <▼>. |
| 20 | Confirm with <RUN/ENTER>.
The minutes flash on the display. |
| 21 | Set the current time with <▲> <▼>. |
| 22 | Confirm with <RUN/ENTER>.
The measuring instrument switches to the last active measuring mode. |

4.9 Reset

You can reset (initialize) the measurement parameters and the configuration parameters separately from one another.

Measurement parameters

The following measuring parameters are reset to the delivery status:

pH	Measuring mode	pH
	Slope	-59.16 mV/pH
	Asymmetry	0 mV
	Manual temperature input	25 °C
Oxi	Measuring mode	D. O. concentration
	Relative slope	1.00
	Salinity (value)	0.0
	Salinity (function)	off



Note

The calibration data gets lost when the measuring parameters are reset. Recalibrate after performing a reset.


Configuration parameters

The following configuration parameters (InI) are reset to the delivery status:


Baud rate	4800
Display during the pH calibration	Buffer nominal value
Interval 1 (automatic save)	OFF
Interval 2 (for data transmission)	OFF

Resetting the measuring parameters

- | | |
|---|--|
| 1 | Press the <RUN/ENTER> key and hold it down. |
| 2 | Press the <CAL> key. The setting to reset the pH measuring parameters appears on the display. |



- | | |
|---|--|
| 3 | Using <▲> <▼> , toggle between <i>no</i> and <i>YES</i> .
<i>YES</i> : Resetting the pH measuring parameters
<i>no</i> : Retaining settings. |
| 4 | Confirm with <RUN/ENTER> .
The measuring instrument switches to the setting to reset the O ₂ measuring parameters. |



- | | |
|---|---|
| 5 | Using <▲> <▼> , toggle between <i>no</i> and <i>YES</i> .
<i>YES</i> : Resetting the O ₂ measuring parameters.
<i>no</i> : Retaining settings. |
| 6 | Confirm with <RUN/ENTER> . The measuring instrument switches to the configuration parameters. |

Resetting the configuration parameters




- | | |
|---|--|
| 7 | Using <▲> <▼>, toggle between <i>no</i> and <i>YES</i> .
<i>YES</i> : Resetting the configuration parameters
<i>no</i> : Retaining settings. |
| 8 | Confirm with <RUN/ENTER>.
The measuring instrument automatically switches to the last active measuring mode. |

5 Maintenance, cleaning, disposal

5.1 Maintenance

The measuring instrument is almost maintenance-free.
The only maintenance task is replacing the batteries.
LoBat on the display indicates that the batteries should be changed.
The batteries are then largely depleted.

Replacing the batteries



1	Open the housing after the instrument has been switched off: – Undo the four screws on the underside of the instrument – Pull down the lower cover (1).
2	If necessary, take the four depleted batteries (2) out of the battery compartment.
3	Place four new batteries (type Mignon AA) in the battery compartment.
4	Close the lower cover (1).



Warning

Make sure that the poles of the batteries are the right way round.
The $\pm$ signs on the batteries must correspond to the $\pm$ signs in the battery compartment. Only use leakproof alkaline manganese batteries.

**Note**

For maintenance of electrodes and sensors, follow the corresponding operating manual.

5.2 Cleaning

Occasionally wipe the outside of the measuring instrument with a damp, lint-free cloth. Disinfect the housing with isopropanol as required.

**Warning**

The housing is made of a synthetic material (ABS). Thus, avoid contact with acetone and similar detergents that contain solvents. Remove any splashes immediately.

5.3 Disposal

Packing

This measuring instrument is sent out in a protective transport packing. We recommend: Keep the packing material. The original packing protects the instrument against damage during transport.

Batteries

This note refers to the battery regulation that applies in the Federal Republic of Germany. We would ask end-consumers in other countries to follow their local statutory provisions.

**Note**

This instrument contains batteries. Batteries that have been removed must only be disposed of at the recycling facility set up for this purpose or via the retail outlet.

It is illegal to dispose of them in household refuse.

Measuring instrument

Dispose of the measuring instrument as electronic waste at an appropriate collection point. It is illegal to dispose of the instrument in household refuse.

6 What to do if...

6.1 pH system messages

Error message

Cause	Remedy
pH electrode:	
– Not connected	– Connect electrode
– Air bubble in front of the diaphragm	– Remove air bubble
– Air in the diaphragm	– Extract air or moisten diaphragm
– Cable broken	– Replace electrode
– Gel electrolyte dried out	– Replace electrode

Error message

Cause	Remedy
Electrode	
– Diaphragm contaminated	– Clean diaphragm
– Membrane contaminated	– Clean membrane
– Moisture in the plug	– Dry plug
– Electrolyte out of date	– Replenish electrolyte or replace electrode
– Electrode worn out	– Replace electrode
– Electrode broken	– Replace electrode
Measuring instrument:	
– Incorrect calibration procedure	– Select correct procedure
– Incorrect solution temperature (without temperature sensor)	– Set up correct temperature
– Socket damp	– Dry socket

No stable measured value	Buffer solutions	
	– Incorrect buffer solutions	– Change calibration procedure
	– Buffer solutions too old	– Use only once. Note the shelf life
	– Buffer solutions depleted	– Change solutions
	Cause	
	Remedy	
	pH electrode:	
	– Diaphragm contaminated	– Clean diaphragm
	– Membrane contaminated	– Clean membrane
	Test sample:	
Obviously incorrect measured values	– pH value not stable	– Measure with air excluded if necessary
	– Temperature not stable	– Adjust temperature if necessary
	Electrode + test sample:	
	– Conductivity too low	– Use suitable electrode
	– Temperature too high	– Use suitable electrode
	– Organic liquids	– Use suitable electrode
	Cause	
	Remedy	
	pH electrode:	
	– pH electrode unsuitable	– Use suitable electrode
	– Temperature difference between buffer and test sample too high	– Adjust temperature of buffers or test samples
	– Measurement procedure not suitable	– Follow special procedure

6.2 Oxi system messages

Error message

Cause	Remedy
Display range exceeded	
D. O. probe:	
– Not connected	– Connect the probe
– Cable broken	– Replace probe
– Depleted	– Replace probe
– Short-circuit between gold and lead electrode	– Clean probe and replace it if necessary

Error message

Cause	Remedy
Invalid calibration	
D. O. probe:	
– Electrolyte solution depleted	– Regenerate probe
– Membrane contaminated	– Clean membrane
– Electrode system poisoned	– Regenerate probe
– Worn out	– Replace probe
– Broken	– Replace probe

Error message

Cause	Remedy
Membrane damaged	
– Membrane damaged	– Regenerate probe
– Membrane head not screwed on tight enough	– Screw membrane head tight

AR flashes continuously

Cause	Remedy
No stable measured value	
– Membrane contaminated	– Clean membrane

Measured value too low**Cause**

- Insufficient flow

Remedy

- Provide flow to the probe

Measured value too high**Cause**

- High amount of dissolved substances
- Air bubbles bump on the membrane with high velocity
- The carbon dioxide pressure is too high (> 1 bar)

Remedy

- Correct solubility function using the salinity equivalent
- Avoid direct flow to the membrane
- Measuring not possible

6.3 General errors**Display LoBat****Cause**

- Batteries almost empty

Remedy

- Exchange the batteries (see section 5.1 MAINTENANCE)

Instrument does not react to keystroke**Cause**

- Operating condition undefined or EMC load unallowed

Remedy

- Processor reset: Switch the instrument on while pressing the <CAL> key

Display to**Cause**

- Time-out of the interface

Remedy

- Check the instrument that is connected

Probe symbol flashes**Cause**

- Calibration interval expired

Remedy

- Recalibrate the measuring system

Message StoFull**Cause**

- All memory locations are full

Remedy

- Output data storage and clear data storage

7 Technical data

Dimensions and weight	Length [mm]	172	
	Width [mm]	80	
	Height [mm]	37	
	Weight [kg]	Approx. 0.3	
Mechanical structure	Type of protection	IP 66	
Electrical safety	Protective class	III	
Ambient conditions	Storage	- 25 °C ... + 65 °C	
	Operation	-10 °C ... + 55 °C	
	Climatic class	2	
pH/ORP measuring ranges		Measuring range	Resolution
	pH	- 2.00 ... + 19.99	0.01
	U [mV]	- 1999 ... + 1999	1
	T [°C]	- 5.0 ... + 105.0	0.1
Precision of pH/ORP (± 1 digit)	pH (after calibration)	± 0.01	
	U [mV]	± 1	
	T [°C]	± 0.1	
pH temperature input:	- 20 °C ... + 130 °C		

Oxi measuring ranges

	mg/l	%	°C
Measuring range I Resolution	0 ... 19.99 0.01	0 ... 199.9 0.1	
Measuring range II Resolution	0 ... 90.0 0.1	0 ... 600 1	
Temperature resolution		0.1	- 5,0 ... + 105.0 0.1

**Precision of
Oxi (± 1 digit)**

Concentration	± 0.5 % of the measured value at an ambient temperature of 5 °C ... 30 °C
Saturation	± 0.5 % of the measured value when measuring in the range of ± 10 K around the calibration temperature
Temperature compensation	< 2 % at 0 ... 40 °C
T [°C]	± 0.1

**Correction functions of
Oxi**

Salinity correction	0 ... 70.0 SAL
Air pressure correction	Automatic through installation of pressure sensor in the range 500 ... 1100 mbar

Serial interface

Type	RS232, data output
Baud rate	Selectable 1200, 2400, 4800, 9600 baud
Data bits	8
Stop bit	2
Parity	None
Handshake	RTS/CTS + Xon/Xoff
Cable length	Max. 15m

Power supply

Batteries	4 x 1.5 V alkali-manganese batteries, type AA
Operational life	Approx. 3000 operating hours
Mains (optional)	The following specifications apply to all plug-in power supplies: Max. overvoltage connection category II Plug-in power supply unit (Euro, US , UK, Australian plug) FRIWO FW7555M/09, 15.1432 Friwo Part. No. 1822089 Input: 100 ... 240 V ~ / 50 ... 60 Hz / 400 mA Output: 9 V = / 1,5 A Plug-in power supply with Euro plug: FRIWO FW1199, 11.7864 Friwo Part. No. 1762613 Input: 230 V ~ / 50 Hz / 5.6 VA Output: 12 V = / 130 mA / 1.56 VA Plug-in power supply with US plug: FRIWO FW1199, 11.7880 Friwo Part. No. 1794043 Input: 120 V ~ / 60 Hz / 6 VA Output: 12 V = / 150 mA Plug-in power supply with UK plug: FRIWO FW1199, 11.7872 Friwo Part No. 1816491 Input: 230V ~ / 50 Hz / 5.6 VA Output: 12 V = / 130 mA / 1.56 VA

**Guidelines
and norms used**

EMC	E.C. guideline 89/336/EEC EN 61326-1:1997 EN 61000-3-2 A14:2000 EN 61000-3-3:1995 FCC Class A
Instrument safety	E.C. guideline 73/23/EEC EN 61010-1 A2:1995
Climatic class	VDI/VDE 3540
Type of protection	EN 60529:1991

FCC Class A Equipment Statement

Note: This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.

Test certificates

cETLus, CE

8 Lists

This chapter provides additional information and orientation aids.

Abbreviations

The list of abbreviations explains the indicators and the abbreviations that appear on the display and in the manual.

Specialist terms

The glossary briefly explains the meaning of the specialist terms. However, terms that should already be familiar to the target group are not described here.

Index

The index will help you to find the topics that you are looking for.

Abbreviations

AR	AutoRead (drift control)
ARng	Automatic range switching Measuring instrument measures with highest resolution
ASY	Asymmetry
AutoCal TEC	Automatic pH calibration with WTW technical buffer solutions according to DIN 19267
°C	Temperature unit, degrees Celsius
Cal	Calibration
Cm...	Display indicator during calibration for pH measurements. Indicates the selection of buffer data records for buffer solutions of the Merck company
ConCal	Conventional single-point or two-point calibration for pH measurements
Ct...	Display indicator during calibration for pH measurements. Indicates the selection of the buffer data records for WTW technical buffer solutions
E3	Error message see chapter 6 WHAT TO DO IF...
InI	Initialization Resets individual basic functions to the status they had on delivery
LoBat	Batteries almost empty (Low Battery)
mV	Voltage unit
mV/pH	Unit of the electrode slope (internat. mV)
OFL	Display range exceeded (Overflow)
OxiCal	Automatic calibration for D. O. measurements
pH	pH value
S	Slope (internat. k)
SAL	Salinity
SELV	Safety Extra Low Voltage
SLO	Slope setting on calibration

TP	Temperature measurement active (Temperature Probe)
U _{ASY}	Asymmetry

Glossary

Adjusting	To manipulate a measuring system so that the relevant value (e. g. the displayed value) differs as little as possible from the correct value or a value that is regarded as correct, or that the difference remains within the tolerance.
Amperometry	Name of a measuring technique. The signal (depending on the measured parameter) of the probe is the electric current. The electrical voltage remains constant.
Asymmetry	Designation for the offset potential of a pH electrode. It is the measurable potential of a symmetrical electrode, the membrane of which is immersed in a solution with the pH of the nominal electrode zero point (WTW electrodes: pH = 7).
AutoRange	Name of the automatic selection of the measuring range.
AutoRead	WTW name for a function to check the stability of the measured value.
Calibration	Comparing the value from a measuring system (e. g. the displayed value) to the correct value or a value that is regarded as correct. Often, this expression is also used when the measuring system is adjusted at the same time (see adjusting).
D. O. partial pressure	Pressure caused by the oxygen in a gas mixture or liquid.
D. O. saturation	Short name for the relative D. O. saturation. Note: The D. O. saturation value of air-saturated water and the D. O. saturation value of oxygen-saturated water are different.
Diaphragm	The junction is a porous body in the housing wall of reference electrodes or electrolyte bridges. It forms the electrical contact between two solutions and makes electrolyte exchange more difficult. The expression, junction, is also used for ground or junction-less transitions.
Electrode zero point	The zero point of a pH electrode is the pH value at which the electromotive force of the pH electrode at a specified temperature is zero. Normally, this is at 25 °C.
Electromotive force of an electrode	The electromotive force U of the electrode is the measurable electromotive force of an electrode in a solution. It equals the sum of all the galvanic voltages of the electrode. Its dependency on the pH results in the electrode function which is characterized by the parameters, slope and zero point.
Measured parameter	The measured parameter is the physical dimension determined by measuring, e. g. pH, conductivity or D. O. concentration.

Measured value	The measured value is the special value of a measured parameter to be determined. It is given as a combination of the numerical value and unit (e. g. 3 m; 0.5 s; 5.2 A; 373.15 K).
Measuring system	The measuring system comprises all the devices used for measuring, e. g. measuring instrument and probe. In addition, there is the cable and possibly an amplifier, terminal strip and armature.
Molality	Molality is the quantity (in Mol) of a dissolved substance in 1000 g solvent.
MultiCal®	WTW name stating that a measuring instrument provides several calibration procedures.
Offset potential	The measurable potential of a symmetrical electrode, the membrane of which is immersed in a solution with the pH of the nominal electrode zero point. The asymmetry is part of the offset potential.
ORP voltage	The ORP is caused by oxidizing or reducing substances dissolved in water if these substances become effective on an electrode surface (e. g. a gold or platinum surface).
OxiCal®	WTW name for a procedure to calibrate D. O. measuring systems in water vapor saturated air.
pH value	The pH is a measure of the acidic or basic effect of an aqueous solution. It corresponds to the negative decadic logarithm of the molal hydrogen ions activity divided by the unit of the molality. The practical pH value is the value of a pH measurement.
Potentiometry	Name of a measuring technique. The signal (depending on the measured parameter) of the electrode is the electrical potential. The electrical current remains constant.
Reference temperature	Fixed temperature value to compare temperature-dependent measured values. For conductivity measurements, the measured value is converted to a conductivity value at a reference temperature of 20 °C or 25 °C.
Reset	Restoring the original condition of all settings of a measuring system.
Resolution	Smallest difference between two measured values that can be displayed by a measuring instrument.
Salinity	The absolute salinity S_A of seawater corresponds to the relationship of the mass of dissolved salts to the mass of the solution (in g/Kg). In practice, this dimension cannot be measured directly. Therefore, the practical salinity is used for oceanographic monitoring. It is determined by measuring the electrical conductivity.
Salt content	General designation for the quantity of salt dissolved in water.
Slope	The slope of a linear calibration function.

Slope (relative)	Designation used by WTW in the D. O. measuring technique. It expresses the relationship of the slope value to the value of a theoretical reference probe of the same type of construction.
Standard solution	The standard solution is a solution where the measured value is known by definition. It is used to calibrate a measuring system.
Temperature coefficient	Value of the slope of a linear temperature function.
Temperature function	Name of a mathematical function expressing the temperature behavior of a test sample, a probe or part of a probe.
Test sample	Designation of the sample ready to be measured. Normally, a test sample is made by processing the original sample. The test sample and original sample are identical if the test sample was not processed.

A

Asymmetry 21
Auslieferungszustand 58
Authorized use 11
AutoCal DIN 21, 27
AutoCal TEC 21, 23
AutoRead 18
 D. O. 36

C

Calibration evaluation 22
 Dissolved oxygen 38
Calibration interval 41
Calibration procedures 21
Calibration protocol 21
 D. O. 37
ConCal 21, 30
Connecting a printer 53
Connecting the original plug-
in power supply unit 14

D

Data record 42
Data transmission (interval)
52
Display 9
Drift control 18
 D. O.
 36

E

Energy saving feature 15
Error messages 65

I

Initial commissioning 13
Initialization 58
Interval
 calibration 41
 data transmission 52
 save 44

K

Keys 8

M

Measuring mode when
switching on 15
Measuring the D. O. concen-
tration 34
Measuring the D. O. satura-
tion 35

O

Operating 10
Operating safety 12
ORP electrode 20
ORP voltage 20

P

Printing
 calibration protocol 21
Probe evaluation (D. O.) 38

R

Remote control 54
Replacing the batteries 61
Reset 58
Resetting the measuring pa-
rameters 59
RS232 interface 53

S

Safety 11
safety precautions 11
Salt content correction
 entering the salinity 40
 switching on 34
Save interval 44
Scope of delivery 13
Setting the baud rate 55
Setting the date 13, 57
Setting the time 13, 57
single-point calibration 21
Single-point calibration (pH)
21, 25, 28, 30
Slope 21, 37

T

Temperature sensor 34
 pH 17
Two-point calibration
 ConCal 30
two-point calibration 21
Two-point calibration (pH)
21, 25, 28, 30



Wissenschaftlich-Technische Werkstätten GmbH

Dr.-Karl-Slevogt-Straße 1
D-82362 Weilheim

Germany

Tel:	+49 (0) 881 183-0
	+49 (0) 881 183-100
Fax:	+49 (0) 881 183-420
E-Mail:	Info@WTW.com
Internet:	http://www.WTW.com