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**Mikroprozessor-
pH-Meter
CG 840
CG 840 B**

Operating Instructions

**Microprocessor
pH Meters
CG 840
CG 840 B**

Mode d'emploi

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Manual de instrucciones

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Operating Instructions Page 9 ... 16

Important notes: Before initial operation of the Microprocessor pH Meters CG 840 and CG 840 B please read and observe carefully the operating instructions. For safety reasons the Microprocessor pH Meters CG 840 and CG 840 B may only be used for the purposes described in these present operating instructions.

Please also observe the operating instructions for the units to be connected.

All specifications in this instruction manual are guidance values which are valid at the time of printing. However, for technical or commercial reasons or in the necessity to comply with the statutory stipulations of various countries, SCHOTT may perform additions to the Microprocessor pH Meters CG 840 and CG 840 B without changing the described properties.

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CONTENTS

PAGE

1	Front view of the Microprocessor pH Meters CG 840, CG 840 B	10
	Rear view of the Microprocessor pH Meters CG 840, CG 840 B	11
2	Functions of the Microprocessor pH Meters CG 840, CG 840 B	12
3	Commissioning	12
4	Calibrating the pH Meters	12
5	Single point calibration	13
6	pH measurement with combination pH electrode N 1042 A or N 2042 A (with integrated Pt 1000 temperature sensor)	13
7	pH measurement with combination pH electrode and separate Pt 1000 temperature sensor	13
8	pH measurement with manual temperature compensation	13
9	mV measurement	13
10	Temperature measurement	13
11	Chart strip recorder output	13
12	Reset function	14
13	Error messages	14
14	Battery change for the pH-Meter CG 840 B	14
15	Warning and safety instructions	14
16	Stored pH values of buffer solutions	15
17	Technical data	16

Important notes: Read these operating instructions carefully before using the unit. For safety reasons the unit may only be used for the purposes described in the operating instructions. This product is subject to continuous modification to adapt it to the state of the art. For this reason, despite great care, it may be possible that these operating instructions do not describe the properties of the unit in full. If you are in doubt please contact the Technical Application Department at our company.

1 Front view of the Microprocessor pH Meter CG 840


Fig. 1

1 Pilot light

2 Display

3 Keyboard

1 Rear view of the Microprocessor pH Meter CG 840


Fig. 2

- | | | |
|-----------------------------|--|--------------------------------------|
| 1 Connectors for electrodes | 2 Connector for temperature sensor | 3 Connector for chart strip recorder |
| 4 Mains connector | 5 Changeover switch for buffer solution sets | 6 Type specification plate |

Rear view of the Microprocessor pH Meter CG 840 B


Fig. 3

- | | | |
|---|--|--------------------------------------|
| 1 Connectors for electrodes | 2 Connector for temperature sensor | 3 Connector for chart strip recorder |
| 4 Battery compartment | 5 Changeover switch for buffer solution sets | 6 Type specification plate |
| 7 Low voltage power supply connection 7...15 V- | | |

2 Functions of the Microprocessor pH Meters CG 840 and CG 840 B

pH values, redox voltages and temperatures can be measured with the Microprocessor pH Meters CG 840 and CG 840 B. The temperature readings are shown in the display at the same time as the pH or mV readings.

The values of two sets of buffer solutions together with their temperature characteristics are stored in memory for matching to the electrode function ("calibration"). The calibration data for pH measurements are preserved in memory when the instruments are switched off. A chart recorder can be connected to the 4 mm sockets on the rear for the functions "pH" and "mV".

A clampstand Z 325 can be attached on the right or left underside for holding the electrode.

The pH Meter CG 840 is powered with AC mains voltage (90 to 250 V~). The pH Meter CG 840 B is powered with batteries or with a mains adapter constructed as mains plug (power supply voltage 7 to 15 V-).

3 Commissioning

The pH meter CG 840 is connected to the mains voltage by plugging the mains cable into the power connector (small equipment panel plug conforming to DIN 49 455) on the rear (Fig. 2, Pos. 4). The pH meter CG 840 B is delivered with batteries ready for operation. It also has a low voltage plug connector (2.1 mm diam.) for power supply with 7 to 15 V- (Fig. 3, Pos. 7) from a power unit contained in the mainsplug. When the instrument is powered with a plug power unit, switch-over from battery to mains is electronically automatic.

Connect the electrode to the sockets on the rear of the instrument, e.g. a Combination Electrode Type N 1042 A or N 2042 A with integrated temperature sensor (Pt 1000 resistance thermometer) (Fig. 2 + 3, Pos. 1 + 2).

The pH Meter CG 840 or CG 840 B is switched on by pressing one of the keys <ON/OFF>, <pH> or <mV>. A green light emitting diode (LED) on the CG 840 lights to indicate the operating state "On". The message <<WAIT>> appears for two seconds in the left top corner of the display until the instrument is ready for making measurements. The <<WAIT>> message also appears for two seconds when switching over between the measuring functions "pH" and "mV".

Before calibrating the instruments, select the desired set of buffer solutions by pressing the changeover switch for buffer sets (Fig. 2 + 3, Pos. 5). The selected values are shown in the digital display.

4 Calibrating the pH meter

The simplest procedure is to use a combination pH electrode with integrated temperature sensor (Pt 1000 resistance thermometer), e.g. type N 1042 A or N 2042 A. The buffer solution pH values stored in memory and valid at 25 °C according to DIN 19 266 and NBS are pH = 4.01, pH = 6.87 and pH = 9.18. The alternative integer buffer solution pH values are pH = 4.00, pH = 7.00 and pH = 10.00. The order in which the buffer solutions are called in the calibration procedure is arbitrary.

Dip the electrode into the first buffer solution which is to be used. Then press the corresponding buffer value key. The messages <<CAL>>, <<WAIT>>, <<P1>> and the corresponding pH value appear in the display. <<WAIT>> flashes when the buffer solution has been identified. <<WAIT>> disappears when the reading has stabilized. <<P2>> then flashes as prompt to continue with the next step of the calibration procedure.

Rinse the electrode with distilled water and repeat the procedure as described above with the second buffer solution. The message "P2" is displayed in this case. Temperature differences of the buffer solutions are standardized automatically by the pH Meters CG 840 and CG 840 B to 25 °C (take into account the settle time of the Pt 1000 resistance thermometer).

When the pH meter has taken over the calibration, the display shows <<CAL>> and <<OK>>, then the electrode zero point in <<pH>> and the slope in <<%>>. The pH meter is now ready for making measurements with this electrode. To display the slope and zero point again, press one of the buffer value keys and then the <pH> key.

Switching over from one set of buffer solutions to the other the calibration parameters are still stored.

5 Single point calibration

A simplified calibration procedure with only one buffer solution is possible too. For this purpose, proceed as described in 4. When the flashing indication $\langle \langle P2 \rangle \rangle$ appears in the display, terminate the calibration procedure by pressing the $\langle \text{pH} \rangle$ key. Only the zero point correction is made in this single point calibration procedure. The slope value is retained from the most recent previous calibration. If a reset (see 12) has been made, the slope setting is 98%.

6 pH measurement with combination pH electrode N 1042 A or N 2042 A (with integrated Pt 1000 temperature sensor)

Press the $\langle \text{pH} \rangle$ key to switch to the "pH" function. The function indications $\langle \langle \text{pH} \rangle \rangle$ and $\langle \langle ^\circ\text{C auto} \rangle \rangle$ are shown in the display (see Fig. 1, Pos. 2).

Dip the electrode into the sample solution. After waiting for a suitable time to establish temperature equilibrium, take the temperature compensated pH reading.

7 pH measurement with combination pH electrode and separate Pt 1000 temperature sensor

pH measurement with automatic temperature compensation is also possible with a combination pH electrode and a separate Pt 1000 temperature sensor. In this configuration, connect one plug of the temperature sensor to the Pt 1000 socket and the other plug to the reference electrode socket (Fig. 2 + 3, Pos. 2).

8 pH measurement with manual temperature compensation

When operating with an electrode without integrated temperature sensor (e.g. Combination Electrode N 62), the manual temperature setting is displayed together with the message $\langle \langle ^\circ\text{C man} \rangle \rangle$.

This temperature can be set to the sample temperature with the keys $\langle ^\circ\text{C} \wedge \rangle$ and $\langle ^\circ\text{C} \vee \rangle$. Manual temperature setting is not possible in the "mV" function and no temperature value is displayed.

9 mV measurement

Press the $\langle \text{mV} \rangle$ key to switch the instrument to the "mV" measuring function. The function display $\langle \langle \text{mV} \rangle \rangle$ lights.

Dip the electrode combination (e.g. Combination Pt Electrode Type Pt 62 for redox voltage measurements) into the sample solution. The mV readings appear directly in the display.

The temperature in $^\circ\text{C}$ is displayed only when a Pt 1000 resistance thermometer is connected (Fig. 2 + 3, Pos. 2).

10 Temperature measurement

Dip the Pt 1000 temperature sensor or the Combination Electrode Type N 1042 A or N 2042 A (with integrated temperature sensor) into the sample solution.

The temperature reading is shown in the display. The pH or mV reading shown simultaneously is without significance if only the temperature sensor is connected.

11 Chart strip recorder output

The instrument is equipped with a chart strip recorder output (Fig. 2 + 3, Pos. 3) with $R_i = 2 \text{ k}\Omega$ source resistance, for connecting a servo chart strip recorder with floating input:

pH: 118 mV/ Δ pH at 25 $^\circ\text{C}$, 100 % slope
U: 1 mV/mV

12 Reset function

If the pH meter gives incorrect readings for no apparent reason, a data reset must be made. For this purpose, switch off the pH meter but do not disconnect it from the mains voltage (CG 840 B: Leave the batteries in the instrument). Then press the <CAL 4>-key and the <ON/OFF>-key at the same time. Only the message <<RESET>> appears in the display. The instrument is ready for making measurements after pressing the <ON/OFF>-key again. When the pH meter has been switched on in this manner, the slope is set to 98% and the zero point to pH = 7.00. Recalibration with the electrode is then necessary if the electrode function differs from these values (see "calibration").

13 Error messages

- Err1: Measuring range overshoot. This error display disappears automatically when the reading returns to a value lying within the measuring range.
- Err2: The settle time of the electrode is too long (> 1 minute). This error display can appear during calibration. It usually means that the electrode has not been watered or that it is too old. "Err2" is also displayed when no electrode is connected.
- Err3: A buffer solution is being used which does not correspond to the <CAL> key which has been pressed. To cure this error status, press the correct key.
- Err4: The temperature during calibration lies below 0 °C or above +100 °C. In this case, correct the temperature of the buffer solutions or the manual temperature setting.
- Err%: The actual slope of the electrode is less than 85% or more than 105% of the theoretical value.
- ErrΔpH: The zero point of the electrode lies below pH = 5.5 or above pH = 8.0 .

Attention: Error messages Err1 and ERR3 may also be caused by:

- 1) drifting of the electrode tension too strong (ageing of electrode)
- 2) exchange of an electrode for a new one with different slope and different zeropoint
- 3) wrong sequence of single operating steps

In all these cases the RESET FUNCTION (see point 12) helps to make the pH meters CG 840 and CG 840 B function again.

14 Battery change for the pH Meter CG 840 B

A red light emitting diode (LED) in the left part of the display area indicates when the battery voltage is low.

To replace the batteries, take off the cover of the battery compartment at the rear of the instrument. Battery replacement will become necessary only after approx. 2 years, because the operating life of the batteries is very long (approx. 3500 hours). The calibration data are preserved for approx. 2 minutes without batteries. If the battery replacement takes longer than this or if the batteries were completely discharged, a reset must be carried out (see Paragraph 12) after battery replacement. A battery set for this instrument consists of four IEC-LR6 Mignon cells of **equal** type (alkaline manganese type). The used-up batteries have to be disposed off in accordance with environmental protection regulations.

15 Warning and safety instructions

The Microprocessor pH Meters CG 840 and CG 840 B may be used only for the purposes described in this operating instructions manual.

Safety considerations require that the Microprocessor pH Meters CG 840 and CG 840 B shall be opened only by authorized persons, e. g. work on the electrical circuitry may be carried out only by properly trained persons. Warranty becomes void in the case of unauthorized tampering inside the instrument or in the case of negligent or deliberate damage. For safety reasons, the instruments may be utilized only for the purposes described in this operating instructions manual.

These instruments conform to the protection class II. They have been manufactured and tested in accordance with DIN 57 411 Part 1/VDE 0411, Part 1, protection measures for electronic measuring equipment, and have left the factory in a safe and properly functioning condition. In order to preserve this condition and ensure safe operation, the user must observe the instructions and warnings contained in this operating instructions manual.

Before switching on, make sure that the mains voltage setting on the instruments agrees with the nominal value of the available mains voltage. The voltage setting is specified on the type specification plate.

The instruments must not be operated or stored in rooms in which water condensation on the equipment units takes place frequently.

Only equal leakproof batteries (4 x Mignon IEC-LR 6) may be used in the pH Meter CG 840 B. Always replace all four battery cells at the same time with four new cells of the same type. Do not use accumulator cells, because the instrument could be destroyed in the case of deep discharge of accumulator cells.

Repairs may be carried out only by properly trained persons. Warranty becomes void in the case of improper tampering inside or damage to the unit.

The buffer solutions used for calibration must correspond exactly to the pH values specified on the keys.

Please also observe the operating instructions for the connected devices, e. g. for the electrodes.

Declaration by the manufacturer

We herewith certify that the Microprocessor pH Meters CG 840 and CG 840 B are equipped with radio frequency interference suppression in accordance with the stipulations of the Official Bulletin Amtsbl Vfg 1046/1984. Commercial distribution of these units have been reported to the German Federal Post Office which has been given the right to check the serial production for compliance with the regulations.

16 Stored pH values of buffer solutions depending on temperature

Value pH	4,00	4,01	6,87	7,00	9,18	10,00
Buffer solution	Kalium-hydrogen-phthalat	Kalium-hydrogen-phthalat	Phosphat	Phosphat-citrat	Borax	Glykokoll-NaOH
Temperature °C						
15	3,99	4,00	6,90	7,03	9,28	10,29
25	4,00	4,01	6,87	7,00	9,18	10,00
35	4,01	4,02	6,84	6,98	9,10	9,76
45	4,03	4,04	6,83	6,97	9,04	9,53
55	4,06	4,07	6,83	6,97	8,99	9,32
65	4,10	4,11	6,84	6,97	8,94	9,11
75	4,14	4,15	6,85	6,98	8,90	8,92

17 Technical Data

Numerical display:	3½digit LCD display for pH, mV; H = 18 mm 2½digit LCD display for temperature; H = 7 mm
Segment displays:	pH, mV, WAIT, Δ, RESET, °C auto, %, °C man, CAL, buffer sets
Pilot light CG 840:	green LED
Battery low indicator CG 840 B:	red LED
Measuring ranges (1 digit =) :	pH: 0.00... 14.00 (pH = 0,01) U: - 1000 ... + 1700 mV (1 mV) T: - 199 ... + 199 °C (0.1 K)
Tolerance:	$\Delta \text{pH} \leq 0.01 \pm 1 \text{ digit}$ $\Delta U \leq 1 \text{ mV} \pm 1 \text{ digit}$ $\Delta T \leq 0.2 \text{ K} (-150...+150 \text{ °C}) \pm 1 \text{ digit}$ (plus the tolerance of the Pt 1000)
Input resistance (20 °C):	$\geq 10^{12} \Omega$
Offset current (20 °C):	$\leq 10^{-12} \text{ A}$
Temperature compensation:	manual: - 1 ... + 100 °C automatic: - 199 ... + 199 °C (with Pt 1000)
Slope matching:	85 ... 105 %
Zero point matching:	pH = 5.5 ... 8
Calibration, automatic:	with 2 buffer sets (switch selected) 1: pH = 4.00/7.00/10.00 2: pH = 4.01/6.87/ 9.18 (DIN 19 266 + NBS)
Electrode connector:	Socket according to DIN 19 262, 2 and 4 mm sockets for reference electrode
Connector for resistance thermometer (Pt 1000):	2 and 4 mm sockets
Chart recorder output:	4 mm sockets $R_i = 2 \text{ k}\Omega$ pH: 118 mV/Δ pH at 25 °C and 100 % slope U: 1 mV/mV for servo chart recorders with floating input
Power supply CG 840:	Mains connection: 90...250 V~, 50...60 Hz, 2,5 VA, protection class II Connection: Small instruments plug, according to DIN 49 455 (built-in small instruments plug, pin distance 8.6 mm)
Power supply CG 840 B:	Batteries, 4 Mignon cells IEC-LR6, external access for replacement, or mains power supply via transformer in mains plug (7 to 15 V-) Connection: Low voltage built-in plug, pluspol on the pin contact (internal contact = ø 2,1 mm) USA/Japan
Ambient temperature:	0 ... 50 °C
Case:	ABS table cabinet, shielded, sprayproof
Dimensions:	240 x 80 x 180 mm (W x H x D)
Weight:	approx. 0.8 kg

All specifications in this instructions manual are guidance values which are valid at the time of printing. However, SCHOTT-GERÄTE may make changes at any time for technical or commercial reasons or in the case of necessity to comply with the statutory stipulations of various countries.

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SCHOTT-GERÄTE GmbH